CALCIUM ELECTRODE

ION-SELECTIVE ELECTRODE FOR MEASUREMENT OF CALCIUM IONS

CALCIUM ELECTRODE

Cat. No E-712

PACKAGE

1 CALCIUM ELECTRODE

INTENDED USE

The CALCIUM ELECTRODE is a membrane electrode designed for in vitro calcium ions measurements in blood, serum, plasma samples. Electrode is intended for use in electrolyte analysers based on ISE system.

STABILITY AND STORAGE

- Store at room temperature (18-25°C).
- The electrode must be installed before the expiry date indicated on the product label, but can remain in the instrument after the date.

INSTALLATION

1. Check that the electrode is $\frac{3}{4}$ filled with fill solution.
2. Check for air bubbles in the fill chamber near sample path. If there are any, tap the electrode body to dislodge air bubbles.
3. Verify that the O-rings are in place on the electrode body.
4. Carefully install electrode in your system.
5. Allow the electrode to equilibrate for at least 15 minutes.
6. Perform a two-point calibration.
7. Analyse a control material sample.

NOTES

- Product for in vitro diagnostic use only.
- Electrodes should be used by suitably qualified laboratory personnel only in accordance with intended purpose.

WASTE MANAGEMENT

Please refer to local legal requirements.

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PZ CORMAY S.A.