ISE LINEARITY CONTROL

CONTROL KIT FOR CALIBRATION AND LINEARITY OF ELECTROLYTES DETERMINATIONS

ISE LINEARITY CONTROL

LEVEL 1	4 x 1.8 ml
LEVEL 2	4 x 1.8 ml
LEVEL 3	4 x 1.8 ml
LEVEL 4	4 x 1.8 ml
LEVEL 5	4 x 1.8 ml

Cat. No E-720

INTENDED USE

ISE LINEARITY CONTROL kit is designed to control the calibration and linearity of electrolytes determinations in ISE electrolyte analysers. The controls are buffered solution of electrolytes: sodium, potassium, chloride, lithium and calcium provided in five levels. The controls contain no human based materials.

STABILITY AND STORAGE

- Unopened controls are stable at 2-8°C until expiry date printed on the product label.
- Avoid freezing and exposure to temperatures greater than 30°C.
- Controls may also be stored at room temperature (up to 25°C) for 6 months, provided the labeled expiration date is not exceeded.
- Do not store once opened ampoules.

PREPARATION PROCEDURE

- Ampoules stored at 2-8°C should be left at least 4 hours to equilibrate to 20-25°C prior to testing.
- It is recommended to proceed with control materials in the same manner as patient samples following the instrument instructions for analyser regarding the use of control materials.
- The control procedure should be carried out from the lowest to the highest level of control for three ampoules of each level. The fourth ampoule of each level is included in the event of accidental damage or sampling error.
- Perform a two-point calibration on your instrument before beginning the calibration verification procedure.

PERFORMANCE

1. Before use mix the controls: held the ampoule at the top and bottom (with forefinger and thumb) and shake 15-20 times (about 10 seconds) to mix the solution. Tap the ampoule to restore the liquid to the bottom on the ampoule.
2. Carefully open the ampoule by snapping off the tip at the score.
3. Immediately introduce the liquid from the ampoule to the analyser, following the instrument manufacturer's instructions for sampling a control material.
4. Record all values for each analyte.
5. Repeat steps 1 through 4 for the remaining 2 ampoules of level 1 (ISE Linearity Control Level 1).
6. Test levels 2, 3, 4, and 5 the same way. Record all values.
7. Calculate the mean value for each test analyte and compare your mean to the range on the Assigned Values on the Value Sheet. Graph the linearity of tests results.

NOTES

- Product for in vitro diagnostic use only.
- Controls should be used by suitably qualified laboratory personnel only in accordance with intended purpose.
- This control are not intended for use as a calibrator.
- Do not use control after expiry date!
- Controls are not a blood-based material, it may not detect certain malfunctions, which would affect the testing of blood. The ISE Linearity Control is sensitive to many instrument related factors that could affect analytical results.

ASSIGNED VALUES

Assigned Values are calculated from results of multiple determinations performed on randomly selected samples from each lot. It is recommended that each laboratory establish its own means and acceptable ranges. The mean value established should fall within the Assigned Values shown on the Value Sheet.

WASTE MANAGEMENT

Please refer to local legal requirements.

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ISE LINEARITY CONTROL

Nr kat. E-720

CONTROL 1	4 x 1.8 ml
CONTROL 2	4 x 1.8 ml
CONTROL 3	4 x 1.8 ml
CONTROL 4	4 x 1.8 ml
CONTROL 5	4 x 1.8 ml

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PZ CORMAY S.A.