## GGT

| Nazwa zestawu       | (PL)  | Nr kat. |
|---------------------|-------|---------|
| Liquick Cor-GGT 30  | 1-226 |         |
| Liquick Cor-GGT 60  | 1-224 |         |
| Liquick Cor-GGT 120 | 3-335 |         |
| HC-GGT              | 4-524 |         |
| OS-GGT              | 9-423 |         |
| B50-GGT             | 5-530 |         |

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### LITERATURA

- Szasz G., Weimann G., Suhler F., Wahlefeld A.W., Persijn J.P.: Z. Klin. Chem. Klin. Biochem., 12 (1974): 228.
- Persijn J.P., van der Silk W.: J. Clin. Chem. Clin. Biochem. 14 (1976): 421-427
- Burtis C.A., Ashwood E.R., ed. Tietz Textbook of Clinical Chemistry, 2nd ed. Philadelphia, PA: WB Saunders, (1994), 850-1.
- Tietz N.W., ed. Clinical Guide to Laboratory Tests, 3rd ed. Philadelphia, PA (1995), WB Saunders: 286
- Kaplan L.A., Pesce A.J.: Clinical Chemistry. Theory, analysis and correlation 3<sup>rd</sup> Ed. (1996), the C.V. Mosby Company, St. Louis: 1072.
- Shaw L.M., Strömme J.H., London J.L., Theodorsen L. International Federation Clinical Chemistry, (IFCC), Scientific Committee, Analytical Section: IFCC method for gamma-glutamyltransferase J Clin Chem Clin Biochem. Vol 21(1983): 633-46
- Schumann G.I et al. International Federation of Clinical Chemistry and Laboratory Medicine. IFCC primary reference procedures for the measurement of catalytic activity concentrations of enzymes at 37 degrees C. International Federation of Clinical Chemistry and Laboratory Medicine. Part 6. Reference procedure for the measurement of gamma-glutamyltransferase. Clin Chem Lab Med. (2002), 40(7): 734-8.
- Schumann G.I., Klauke R. New IFCC reference procedures for the determination of catalytic activity concentrations of five enzymes in serum: preliminary upper reference limits obtained in hospitalized subjects. Clin Chim Acta 327 (2003): 69-79

Data wydania: 10. 2023.

## GGT

| Kit name            | (EN)<br>Cat. No |
|---------------------|-----------------|
| Liquick Cor-GGT 30  | 1-226           |
| Liquick Cor-GGT 60  | 1-224           |
| Liquick Cor-GGT 120 | 3-335           |
| HC-GGT              | 4-524           |
| OS-GGT              | 9-423           |
| B50-GGT             | 5-530           |

### INTENDED USE

Diagnostic kit for determination of  $\gamma$ -glutamyltransferase activity used both for manual assay and in several automatic analysers.


The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

### INTRODUCTION

$\gamma$ -Glutamyltransferase (GGT, GGTP) is a membrane localized enzyme that catalyzes the transfer of glutamyl groups from glutathione to amino acids or peptides. Large GGT amounts are present in secretory organs: kidney, liver, bile duct, pancreas. Although the GGT activity is highest in renal tissue, serum GGT is generally elevated as a result of liver disease. Since alcohol induces GGT production, measurement of GGT activity is used for monitoring of abstinence in withdrawal treatment.

### METHOD PRINCIPLE

IFCC method. Kinetic method with L- $\gamma$ -glutamyl-3-carboxy-4-nitroanilide.

The rate of 5-amino-2-nitrobenzoate creation measured colorimetrically is directly proportional to  $\gamma$ -glutamyltransferase activity.

### REAGENTS

| Package   | Liquick<br>Cor-GGT<br>30 | Liquick<br>Cor-GGT<br>60 | Liquick<br>Cor-GGT<br>120 |
|-----------|--------------------------|--------------------------|---------------------------|
| 1-REAGENT | 5 x 24 ml                | 5 x 48 ml                | 5 x 96 ml                 |
| 2-REAGENT | 1 x 30 ml                | 1 x 60 ml                | 1 x 120 ml                |
|           | HC-GGT                   | OS-GGT                   | B50-GGT                   |
| 1-REAGENT | 6 x 76 ml                | 4 x 53.5 ml              | 3 x 58.5 ml               |
| 2-REAGENT | 6 x 19.5 ml              | 4 x 15.5 ml              | 3 x 17.5 ml               |

The reagents when stored at 2-8°C are stable up to expiry date printed on the package. The reagents are stable for 12 weeks on board the analyser at 2-10°C.

### Working reagent preparation and stability

Assay can be performed with use of separate 1-REAGENT and 2-REAGENT reagents or with use of working reagent. For working reagent preparation mix gently 4 parts of 1-REAGENT with 1 part of 2-REAGENT. Avoid foaming.

Stability of working reagent in darkness: 4 weeks at 2-8°C  
5 days at 15-25°C

### Concentrations in the test

|                                                |            |
|------------------------------------------------|------------|
| Tris (pH 8.25)                                 | 100 mmol/l |
| glycylglycine                                  | 100 mmol/l |
| L- $\gamma$ -glutamyl-3-carboxy-4-nitroanilide | 4 mmol/l   |

### Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.
- 1-REAGENT and 2-REAGENT meeting the criteria for classification in accordance with Regulation (EC) No 1272/2008.

### Ingredients:

1-REAGENT and 2-REAGENT contain reaction mass of 5-chloro-2-methyl-2H-isothiazol-3-one and 2-methyl-2H-isothiazol-3-one (3:1).

### Warning

H317 May cause an allergic skin reaction.  
H412 Harmful to aquatic life with long lasting effects.

P280 Wear protective gloves, protective clothing, eye protection or face protection.

P302 + P352 IF ON SKIN: Wash with plenty of soap and water.

P273 Avoid release to the environment.

### ADDITIONAL EQUIPMENT

- automatic analyzer or photometer able to read at 405 nm;
- thermostat at 25°C, 30°C or 37°C;
- general laboratory equipment;

### SPECIMEN

Serum, EDTA plasma free from haemolysis.  
Do not use citrate, oxalate and fluoride as anticoagulants because of GGT activity inhibition!

Heparin causes turbidity in the reaction mixture!  
GGT activity remains stable in specimen up to 2 days at 15-25°C or 1 week at 2-8°C or 1 month at -25°C. Freezing of sample causes a loss of enzyme activity. Frozen specimens should be thawed and kept at room temperature for 18 to 24 hours before measurement to achieve full enzyme reactivation.  
Nevertheless it is recommended to perform the assay with freshly collected samples!

### PROCEDURE

Applications for analyzers are available on request.

### Manual procedure

|             |                    |
|-------------|--------------------|
| wavelength  | 405 nm             |
| temperature | 25°C / 30°C / 37°C |
| cuvette     | 1 cm               |

### Sample Start method

Pipette into the cuvette:

|                 |              |
|-----------------|--------------|
| working reagent | 1000 $\mu$ l |
|-----------------|--------------|

Bring up to the temperature of determination. Then add:

|        |             |
|--------|-------------|
| sample | 100 $\mu$ l |
|--------|-------------|

Mix and incubate at adequate temperature. After about 1 min. read the absorbance against air or water. Repeat the reading after exactly 1, 2 and 3 minutes. Calculate the mean absorbance change per minute ( $\Delta A/\text{min.}$ ).

### Calculation

GGT activity [U/l] =  $\Delta A/\text{min.} \times 1511$

1 U/l = 0.0167  $\mu$ kat/l

### Reagent Start method

The determination can be also performed with use of separate 1-REAGENT and 2-REAGENT reagents.

Pipette into the cuvette:

|           |              |
|-----------|--------------|
| 1-REAGENT | 1000 $\mu$ l |
|-----------|--------------|

Bring up to the temperature of determination. Then add:

|        |             |
|--------|-------------|
| sample | 100 $\mu$ l |
|--------|-------------|

Mix well, incubate for 1 min. Then add:

|           |             |
|-----------|-------------|
| 2-REAGENT | 250 $\mu$ l |
|-----------|-------------|

Mix well; perform measurement as described for Sample Start method.

### Calculation

GGT activity [U/l] =  $\Delta A/\text{min.} \times 1810$

1 U/l = 0.0167  $\mu$ kat/l

### REFERENCE VALUES <sup>7</sup>

| serum / plasma | 37°C                         |
|----------------|------------------------------|
| female         | < 38 U/l (0.633 $\mu$ kat/l) |
| male           | < 55 U/l (0.917 $\mu$ kat/l) |

It is recommended for each laboratory to establish its own reference ranges for local population.

### QUALITY CONTROL

For internal quality control it is recommended to use the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173) with each batch of samples.

For the calibration of automatic analysers systems the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) or LEVEL 2 (Cat. No 5-175; 5-177) is recommended.

The calibration curve should be prepared every 12 weeks, with change of reagent lot number or as required e.g. quality control findings outside the specified range.

### PERFORMANCE CHARACTERISTICS

These metrological characteristics have been obtained using the automatic analyser Biolis 24i Premium. Results may vary if a different instrument or a manual procedure is used.

- Sensitivity:** 11.2 U/l (0.187  $\mu$ kat/l).
- Linearity:** up to 580 U/l (9.667  $\mu$ kat/l).
- Specificity / Interferences**  
Haemoglobin up to 0.16 g/dl, ascorbate up to 62 mg/l, bilirubin up to 20 mg/dl and triglycerides up to 1000 mg/dl do not interfere with the test.

### Precision

| Repeatability (run to run)<br>n = 20 | Mean<br>[U/l] | SD<br>[U/l] | CV<br>[%] |
|--------------------------------------|---------------|-------------|-----------|
| level 1                              | 29.06         | 1.44        | 4.97      |
| level 2                              | 134.05        | 1.86        | 1.39      |

| Reproducibility (day to day)<br>n = 80 | Mean<br>[U/l] | SD<br>[U/l] | CV<br>[%] |
|----------------------------------------|---------------|-------------|-----------|
| level 1                                | 29.05         | 2.14        | 7.36      |
| level 2                                | 132.33        | 3.50        | 2.64      |

### Method comparison

A comparison between GGT values determined at Biolis 24i Premium (y) and at ADVIA 1650 (x) using 30 samples gave following results:

$y = 1.0181x - 2.5022$  U/l;

$R = 0.9995$  (R – correlation coefficient)

### WASTE MANAGEMENT

Please refer to local legal requirements.

### LITERATURE

- Szasz G., Weimann G. Suhler F., Wahlefrld A.W., Persijn J.P.: Z. Klin. Chem. Biochem., 12 (1974): 228.
- Persijn J.P., van der Silk W.: J. Clin. Chem. Clin. Biochem. 14 (1976): 421-427
- Burtis C.A., Ashwood E.R., ed. Tietz Textbook of Clinical Chemistry, 2nd ed. Philadelphia, PA: WB Saunders, (1994), 850-1.
- Tietz N.W., ed. Clinical Guide to Laboratory Tests, 3rd ed. Philadelphia, PA (1995), WB Saunders: 286
- Kaplan L.A., Pesce A.J.: Clinical Chemistry. Theory, analysis and correlation 3<sup>rd</sup> Ed. (1996), the C.V. Mosby Company, St. Louis: 1072.
- Shaw LM, Strømme JH, London JL, Theodorsen L. International Federation Clinical Chemistry, (IFCC), Scientific Committee, Analytical Section: IFCC method for gamma-glutamyltransferase J Clin ChemClin Biochem. Vol 21(1983): 633-46
- Schumann Gl et al. International Federation of Clinical Chemistry and Laboratory Medicine. IFCC primary reference procedures for the measurement of catalytic activity concentrations of enzymes at 37 degrees C. International Federation of Clinical Chemistry and Laboratory Medicine. Part 6. Reference procedure for the measurement of gamma-glutamyltransferase. Clin Chem Lab Med. (2002), 40(7): 734-8.
- Schumann Gl, Klauke R. New IFCC reference procedures for the determination of catalytic activity concentrations of five enzymes in serum: preliminary upper reference limits obtained in hospitalized subjects. Clin Chim Acta 327 (2003): 69-79

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- Szasz G., Weimann G. Suhler F., Wahlefrld A.W., Persijn J.P.: Z. Klin. Chem. Klin. Biochem., 12 (1974): 228.
- Persijn J.P., van der Silk W.: J. Clin. Chem. Clin. Biochem. 14 (1976): 421-427
- Burtis C.A., Ashwood E.R., ed. Tietz Textbook of Clinical Chemistry, 2nd ed. Philadelphia, PA: WB Saunders, (1994), 850-1.
- Tietz N.W., ed. Clinical Guide to Laboratory Tests, 3rd ed., Philadelphia, PA (1995), WB Saunders: 286
- Kaplan L.A., Pesce A.J.: Clinical Chemistry. Theory, analysis and correlation 3<sup>rd</sup> Ed. (1996), the C.V. Mosby Company, St. Louis: 1072.
- Shaw LM, Strømme JH, London JL, Theodorsen L. International Federation Clinical Chemistry, (IFCC), Scientific Committee, Analytical Section: IFCC method for gamma-glutamyltransferase J Clin ChemClin Biochem. Vol 21(1983): 633-46
- Schumann GI et al. International Federation of Clinical Chemistry and Laboratory Medicine. IFCC primary reference procedures for the measurement of catalytic activity concentrations of enzymes at 37 degrees C. International Federation of Clinical Chemistry and Laboratory Medicine. Part 6. Reference procedure for the measurement of gamma-glutamyltransferase. Clin Chem Lab Med. (2002), 40(7): 734-8.
- Schumann GI, Klauke R. New IFCC reference procedures for the determination of catalytic activity concentrations of five enzymes in serum: preliminary upper reference limits obtained in hospitalized subjects. Clin Chim Acta 372 (2003): 69-79

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