FOR INFORMATION ONLY.
WHEN PERFORMING
THE ASSAY ALWAYS REFER
TO PACKAGE ALWAYS INSERT
SUPPLIED WITH THE KIT

CanChek Tumor Marker Control Serum Levels 1 and 2

REF

107-20

<div style="display: flex; align-items: center;"> <div style="border: 1px solid black; padding: 2px; margin-right: 5px;">LOT</div> <div style="margin-right: 20px;">Kit lotnr</div> <div style="margin-right: 20px;"> </div> <div>Exp.date</div> </div>		<div style="display: flex; align-items: center;"> <div style="border: 1px solid black; padding: 2px; margin-right: 20px;">IVD</div> <div style="font-size: 2em;">CE</div> </div>		
		Range/Bereich/Plages de valeurs/ Intervallo/Rango/Limites/Område:		
		SCC (µg/L)	NSE (µg/L)	S100B (ng/L)
<div style="display: flex; align-items: center;"> <div style="border: 1px solid black; padding: 2px; margin-right: 5px;">CONT</div> </div>				
<div style="display: flex; align-items: center;"> <div style="width: 150px;">3 x 0.75 mL Level 1</div> <div style="margin-left: 20px;"> <div style="border: 1px solid black; padding: 2px; margin-right: 5px;">LOT</div> </div> <div>Lot 1</div> </div>		Gräns SCC,1	Gräns NSE, 1	Gräns S100B, 1
<div style="display: flex; align-items: center;"> <div style="width: 150px;">3 x 0.75 mL Level 2</div> <div style="margin-left: 20px;"> <div style="border: 1px solid black; padding: 2px; margin-right: 5px;">LOT</div> </div> <div>Lot 2</div> </div>		Gräns SCC, 2	Gräns NSE, 2	Gräns S100B, 2

GB EXPLANATION OF SYMBOLS
DE BEDEUTUNG DER SYMBOLE
ES EXPLICACIÓN DE SÍMBOLOS
IT SIGNIFICATO DEI SIMBOLI
FR EXPLICATION DES SYMBOLES
DK SYMBOLFORKLARING
GR ΕΠΕΞΗΓΗΣΗ ΤΩΝ ΣΥΜΒΟΛΩΝ
SE SYMBOLFÖRKLARING

Use By/Verwendbar bis/
Fecha de caducidad/Utilizzare entro/
Utiliser jusque/Holdbar til/
Ημερομηνία λήξης/Bäst före datum

Batch code/Chargenbezeichnung/
Codigo de lote/Codice del lotto/Code du lot/
Lotnummer/Αριθμός Παρτίδας/Lotnummer

Date of manufacture/Herstellungsdatum/
Fecha de fabricación/Data di fabbricazione/
Date de fabrication/Produktionsdato/
Ημερομηνία Παραγωγής/Tillverkningsdatum

REF

Catalogue number/Bestellnummer/
Número de catálogo/Numero di catalogo/
Référence du catalogue/Katalognummer/
Αριθμός καταλόγου/Produktnummer

Manufacturer/Hersteller/Fabricante/Fabbricante/
Fabricant/Producent/τασκευαστής/Tillverkare

In Vitro Diagnostic Medical Device/
In Vitro Diagnostikum/Producto
sanitario para diagnóstico in vitro/
Dispositivo medico diagnostico in vitro/
Dispositif médical de diagnostic in
vitro/Medicinsk udstyr til in vitro-
diagnostik/In Vitro Διαγνωστικό
Ιατροτεχνολογικό προϊόν/
Endast för in vitro-diagnostik

Temperature limitation/
Zulässiger Temperaturbereich/
Limite de temperatura/
Limiti di temperatura/
Limites de température/
Temperaturbegrænsning/
Περιορισμοί θερμοκρασίας/
Temperaturgräns

Consult Instructions for Use/
Gebrauchsanweisung beachten/
Consulte las instrucciones de uso/
Consultare le istruzioni per l'uso/
Consulter les instructions d'utilisation/
Se brugsanvisning/
Συμβουλευτείτε τις οδηγίες χρήσης/
Se bruksanvisning

Biological risks/
Biogefährdung/Riesgo biológico/
Rischio biologico/Risques
biologiques/Biologisk fare/
Βιολογικοί κίνδυνοι/Biologisk risk

Contents of kit/Inhalt/Contenido/
Contenido/Contenu/Inndhold/
ανιδραστήρια/Kit innehåll

INTENDED USE

CanChek Tumor Marker Control Serum is intended for use as an assayed quality control serum to monitor the accuracy of laboratory testing for Neuron Specific Enolase (NSE), S100B protein and Squamous Cell Carcinoma Antigen (SCCA).

SUMMARY AND PRINCIPLE

The use of a quality control serum is indicated to verify the performance characteristics of an in-vitro diagnostic laboratory test and is an integral part of good laboratory practice. The CanChek controls are available in two levels in order to provide performance monitoring over the clinical range.

REAGENTS

CanChek Tumor Marker Control Serum is prepared from human serum with added human proteins, excipient and preservative. To be reconstituted with distilled water before use.

STORAGE AND STABILITY

CanChek Tumor Marker Control Serum is stable until the expiry date stated on the vials when stored unopened at 2–8°C. Once reconstituted, all analytes are stable for 2 weeks when stored tightly capped at 2–8°C or 3 months when stored at –20°C or below.

PREPARATION OF REAGENT

Add exactly 0.75 mL of distilled water to each vial. Replace the stopper and allow standing for 15 minutes, swirling occasionally. Gently invert the vial several times to ensure homogeneity before use.

WARNINGS AND PRECAUTIONS**For in vitro diagnostic use.**

- For professional use only.
- Please refer to the US Department of Health and Human Services (Bethesda, Md., US) publication No. (CDC) 88-8395 on laboratory safety or any other local or national regulation.
- Follow local guidelines for disposal of all waste material.

Caution

Each donor unit used in the preparation of human source reagent has been tested and found to be Non-Reactive for HIV-1/2 Antibody, HCV Antibody and Hepatitis B Surface Antigen (HBsAg). Since no method can completely rule out the presence of blood borne diseases, the handling and disposal of human source reagents from this product should be made as if they were potentially infectious.

VALUE ASSIGNMENT

The assigned value ranges printed in this insert were derived from replicate analyses using CanAg EIA kits and are valid only for this lot of CanChek Tumor Marker Control. Since no common reference material is available for NSE, S100B or SCC antigen, the assigned value ranges are traceable to a set of in-house reference standards.

LIMITATIONS

CanChek Tumor Marker Control Serum should not be used as a standard.

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CanAg® is a registered trademark of Fujirebio Diagnostics AB
Fujirebio Diagnostics AB