

# IMMUNOLOGY CONTROL 1

REF 1578-1195-04      Packaging: 4 x 3 mL



## INTENDED USE

Material for internal quality control of the listed special protein assays using the related Medicon reagents with Diatron Pictus® P700 and P500 analyzers. For in vitro diagnostic use only by trained laboratory professionals.

## COMPOSITION

Liquid, human serum-based control material with additives and preservative. The concentrations/activities of the proteins are LOT depended.

### ⚠️ WARNINGS – PRECAUTIONS

- For in vitro use only.
- Every donor used for the preparation of the material was tested and found negative for HbsAg, anti-HCV και anti-HIV 1 and 2 with FDA approved methods. As there is no known test method that can offer complete assurance that products derived from human blood will not transmit infectious agents this product should be handled as a potentially infectious material.
- Contains sodium azide (NaN<sub>3</sub>) <0.1% as preservative. Avoid contact with skin, eyes and mucous membranes.
- Any serious incident that may occur in relation to this device must be reported by the user to the manufacturer and the competent authority of the country in which the user and/or the patient is established!

### ⚠️ PREPARATION

Material ready to use. Swirl the vials gently to acquire homogenous mix. Avoid foaming. Replace cap immediately after use and store at 2-8°C. Unsuitable storage, handling, or errors during the analytical procedure may give erroneous results.

## MATERIALS REQUIRED BUT NOT PROVIDED WITH THE KIT

- Special proteins reagents
- Automated biochemical analyzer
- Common laboratory equipment.

### ⚠️ STORAGE - STABILITY

The material is stable, unopened up to the expiry date stated on the label when stored at 2-8°C. Once opened the material is stable for 1 month when stored tightly capped at 2-8°C.

### ⚠️ DETERIORATION

**The material should not be used:**

- After prolonged exposure to high temperature
- After the expiry date or when microbial growth is evident.

## TEST PROCEDURE

Refer to the user's manual of the analyzer for calibration and quality control process.

Each laboratory should establish its own control frequency, however good laboratory practice suggests that controls be tested each date patient samples are tested and each time calibration is performed.

The results obtained by any individual laboratory may vary from the given mean value but should fall within the corresponding acceptable ranges given in the enclosed table.

If any trends or sudden shifts in values are detected, review all operating parameters. Each laboratory should establish guidelines for corrective actions to be taken if controls do not recover within the specified limits. Make sure the LOT on the vial is the same as on the value sheet accompanying the material.

### ⚠️ WASTE DISPOSAL

Dispose of all material according to local guidelines for potentially infectious materials. Sodium azide forms explosive compounds with lead or copper. Flush waste pipes with plenty of water after disposal of undiluted material to avoid sodium azide build up in the drain. Material safety data sheet is available by MEDICON on request.

## SPECIAL PROTEINS ASSAYS ASSIGNED VALUES are lot number specific

Please refer to the value sheet for the specific lot available at <https://medicondoc.com>.

Contents:

|                       |                |                     |
|-----------------------|----------------|---------------------|
| Albumin               | Ceruloplasmin  | Immunoglobulin E    |
| A1-acid glycoprotein  | Complement 3   | Kappa light chains  |
| A1-antitrypsin        | Complement 4   | Lambda light chains |
| A2-macroglobulin      | Haptoglobin    | Prealbumin          |
| Antistreptolysin-O    | Immunoglobulin | Rheumatoid factor   |
| B2-microglobulin      | Immunoglobulin | Transferrin         |
| C-reactive protein    | Immunoglobulin | Ferritin            |
| C-reactive protein La |                |                     |

## SYMBOLS



Manufacturer



In vitro diagnostic medical device



Temperature Limit



Catalogue Number



Caution



Biological Risks