

**AccuDiag™
HCV Ab ELISA**

Cat# 1707-12



IgG ANTIBODIES TO HEPATITIS C VIRUS ELISA KIT

Two-Step Incubation, Indirect Principle

Test	HCV Ab ELISA
Method	Enzyme Linked Immunosorbent Assay
Principle	Solid Phase Indirect ELISA
Detection Range	Qualitative: Positive & Negative Control
Sample	10ul
Specificity	99.55%
Sensitivity	100%
Total Time	~ 75 min
Shelf Life	12-14 Months from the manufacturing date

** Laboratory results can never be the only base of a medical report. The patient history and further tests have to be taken into account.*

INTENDED USE

The purpose of the HCV Ab ELISA Test is for clinical lab diagnosis of patients who are suspected of having a hepatitis C virus infection, and for blood donor screening. This HCV Ab ELISA test is an enzyme-linked immunosorbent assay for in vitro qualitative identification of IgG antibodies to hepatitis C virus in human serum/plasma.

SUMMARY AND EXPLANATION

For more than 50% of patients, HCV develops into chronic hepatitis and has become the principal cause of cirrhosis of the liver and hepatocellular carcinomas. HCV is now accepted as the major agent for non-A and non-B hepatitis transfusion infection. HCV is defined as an envelope, single stranded positive sense RNA (9.5kb) affiliated with the Flaviviridae family. Classification of hepatitis C virus is established through six major genotypes and series of subtypes of HCV.

Blood donor transmissions of hepatitis C virus have notably decreased since screening started in 1990. Three generations of HCV ELISA tests have been established and each generation has resulted in improvement in sensitivity of detecting anti-HCV. 1st generation - limited sensitivity and specificity, using recombinant proteins complementary to NS4 region (c100-3). 2nd generation - noticeable improvement in sensitivity and specificity employing recombinant/synthetic antigens from Core (c22) and nonstructural regions NS3 (c33c, c100-3) and NS4 (c100-3, c200). A third generation was established because antibodies to NS5 non-structural protein of HCV were being developed by infected individuals. Thus, in addition to NS3 (c200), NS4 (c200), and Core (c22), the 3rd generation tests included antigens from NS5 region. This 3rd generation of tests have improved sensitivity and have lessened the time period


(60 days) between infection and the subsequent appearance of antibodies to HCV.

TEST PRINCIPLE

The principle of the HCV ELISA test involves a two-step incubation procedure in which solid phase, indirect ELISA technique for HCV antibodies is established. In this third generation HCV ELISA test, recombinant, highly immunoreactive antigens corresponding to the core and non-structural regions of HCV, are pre-coated on the polystyrene microwell strips.

The test begins in the first incubation stage where anti-HCV specific antibodies (if present) are bound to the solid phase pre-coated HCV antigens. Next, it is important to wash the wells to remove unbound serum proteins. Added after this are rabbit anti-human IgG antibodies (anti-IgG) which are conjugated to the enzyme horseradish peroxidase (HRP-Conjugate). The second incubation stage involves these HRP-conjugated antibodies, which will be bound to any antigen-antibody (IgG) complexes previously formed. Next, wash the wells to remove unbound HRP-conjugate. Chromogen solutions containing Tetramethylbenzidine (TMB) and urea peroxide are then added to the wells. A blue-colored product develops in the presence of the antigen-antibody-anti-IgG (HRP) immunocomplex and when the colorless chromogens are hydrolyzed by the bound HRP conjugate. After stopping the reaction with sulfuric acid, the blue color turns yellow. The color intensity can be gauged proportionally to the amount of antibody captured in the wells, and to the amount in the sample, respectively. The wells remain colorless if the anti-HCV result is negative.

Assay principle scheme: Indirect ELISA



Incubation1	Incubation 2	Immobilized Complex	Coloring	results
30min.	30min.		15min	

Ag(p)–pre-coated HCV antigens(core, NS3/4,NS5);
Ab(s)–HCV antibodies in sample (IgG);
ENZ–HRP conjugated rabbit anti-human IgG;

SPECIMEN COLLECTION AND PREPARATION

- Sample Collection:** Either fresh serum or plasma samples can be used for this assay. Blood collected by venipuncture should be allowed to clot naturally and completely – the serum/plasma must be separated from the clot as early as possible as to avoid hemolysis of the RBC. Care should be taken to ensure that the serum samples are clear and not contaminated by microorganisms. Any visible particulate matters in the sample should be removed by centrifugation at 3000 RPM for at least 20 minutes at room temperature, or by filtration on 0.22u filters. Plasma samples collected into EDTA, sodium citrate or heparin may be tested, but highly lipemic, icteric, or hemolized samples should not be used as they could give erroneous results in the assay. Do not heat inactivate samples. This can cause deterioration of the target proteins in the sample.
- Transportation and Storage:** Store samples at 2-8 °C. Samples not required for assaying within 3 days should be stored frozen (-20 °C or lower). Multiple freeze-thaw cycles should be avoided. For shipment, samples should be packaged and labeled in accordance with the existing local and international regulations for transport of clinical samples and ethological agents.

MATERIALS AND COMPONENTS

Materials provided with the test kits

- Microwell Plate**

Blank microwell strips fixed on white strip holder.
 The plate is sealed in aluminum pouch with desiccant.
8x12/12x8-well strips per plate.

Each well contains recombinant HCV antigens.
 The microwell strips can be broken to be used separately.
 Place unused wells or strips in the plastic sealable storage bag together with the desiccant and return to 2~8 °C

2. Negative Control

Blue-colored liquid filled in a vial with green screw cap.
 0.2ml per vial.
 Protein-stabilized buffer tested non-reactive for HCV antibodies.
 Preservatives: 0.1% ProClin 300.
 Ready to use as supplied.
 Once open, stable for one month at 2-8°C.

3. Positive Control

Red-colored liquid filled in a vial with red screw cap.
 0.2ml per vial.
 Anti-HCV antibodies diluted in protein-stabilized buffer.
 Preservatives: 0.1% ProClin 300.
 Ready to use as supplied. Once open, stable for one month at 2-8°C.

4. Specimen Diluent

Blue liquid filled in a white vial with blue screw cap.
 13ml per vial.
 Protein-stabilized buffer, casein, and sucrose solution.
 Ready to use as supplied. Once open, stable for one month at 2-8°C.

5. HRP-Conjugated Reagent

Red-colored liquid filled in a white vial with red / orange screw cap.
 13ml per vial.
 Horseradish peroxidase-conjugated rabbit anti-human IgG antibodies.
 Ready to use as supplied. Once open, stable for one month at 2-8°C

6. Wash Buffer

Colorless liquid filled in a clear bottle with white screw cap.
 50ml per bottle.
 PH 7.4, 20 × PBS (Contains Tween-20 as a detergent).
 The concentration must be diluted **1 to 20** with distilled/deionized water before use. Once diluted, stable for one week at room temperature or for two weeks at 2-8°C.

7. Chromogen Solution A

Colorless liquid filled in a white vial with green screw cap.
 8ml per vial.
 Urea peroxide solution.
 Ready to use as supplied.
 Once open, stable for one month at 2-8°C

8. Chromogen Solution B

Colorless liquid filled in a Brown vial with brown/ black screw cap.
 8 ml per vial.
 TMB solution (Tetramethyl-benzidine dissolved in citric acid).
 Ready to use as supplied.
 Once open, stable for one month at 2-8°C.

9. Stop Solution

Colorless liquid filled in a white vial with yellow/white screw cap.
 8 ml per vial.
 Diluted sulfuric acid solution (2.0M H₂SO₄)

10. Plastic Sealable Bag

For enclosing the strips not in use.

11. Cardboard Plate Cover

To cover the plates during incubation and prevent evaporation or contamination of the wells.
 2 Sheets

12. Package Insert

4. Disposable V-shaped troughs.
5. Dispensing system and/or pipette (single or multichannel), disposable pipette tips
6. Absorbent tissue or clean towel.
7. Dry incubator or water bath, 37±0.5°C.
8. Microshaker for dissolving and mixing conjugate with samples.
9. Microwell plate reader, single wavelength 450nm or dual wavelength 450nm and 630nm.
10. Microwell aspiration/wash system.

SPECIAL INSTRUCTIONS FOR WASHING

1. A good washing procedure is essential to obtain correct and precise analytical data.
2. It is therefore recommended to use a good quality ELISA microplate washer, maintained at the best level of washing performances. In general, no less than 5 automatic washing cycles with dispensing of 350-400µl/well, are sufficient to avoid false positive reactions and high background (all wells turn yellow).
3. To avoid cross-contaminations of the plate with sample or HRP-conjugate, after incubation do not discard the content of the wells, but allow the plate washer to aspirate it automatically.
4. It is recommended that the washing system should be calibrated on the kit itself in order to match the declared analytical performances. Assure that the microplate washer's liquid dispensing channels are not blocked or contaminated, and sufficient volume of Wash buffer is dispensed each time into the wells.
5. In case of manual washing, we suggest to perform at least 5cycles, dispensing 350-400µl/well and aspirating the liquid for 5times. If poor results (high background) are observed, increase the washing cycles or soaking time per well.
6. In any case, the liquid aspirated out the strips should be treated with a sodium hypochlorite solution (final concentration of 2.5%) for 24 hours, before liquids are disposed in an appropriate way.
7. The concentrated Washing solution should be diluted 1 to 20 before use. For one plate, mix 50 ml of the concentrate with 950ml of water for a final volume of 1000ml diluted Wash Buffer. If less than a whole plate is used, prepare the proportional volume of solution.

ASSAY PROCEDURE

1. **Reagents preparation:** Allow the reagents and samples to reach room temperature (**18-30°C**) for at least 15-30minutes. Check the Wash buffer concentrate for the presence of salt crystals. If crystals have formed in the solution, resolubilize by warming at 37° until crystals dissolve. Dilute the stock Wash Buffer **1 to 20** with distilled or deionized water. Use only clean vessels to dilute the buffer.
2. **Numbering Wells:** Set the strips needed in strip-holder. and number sufficient number of wells including three for the Negative control (e.g. **B1, C1, D1**), two for the Positive control (e.g. **E1, F1**) and one Blank (e.g. **A1**, neither samples nor HRP-Conjugate should be added into the Blank well). If the results will be determined by using dual wavelength plate reader, the requirement for use of Blank well could be omitted. Use only number of strips required for the test.
3. **Adding Diluent:** Add **100 µl** Specimen Diluent into each well except the blank.
4. **Adding Sample:** Add **10µl** of Positive control, Negative control, and Specimen into their respective wells. Note: Use a separate disposal pipette tip for each specimen, Negative and Positive Control as to avoid cross-contamination. Mix by tapping the plate gently.
5. **Incubating (1):** Cover the plate with the plate cover and incubate for **30minutes at 37°C**. It is recommended to use thermostat-controlled water tank to assure the temperature stability and humidity during the incubation. If dry incubator is used, do not open the door frequently.
6. **Washing 1:** After the end of the incubation, remove and discard the plate cover. Wash each well **5times** with diluted Wash buffer. Each time, allow

Materials required but not provided

1. Freshly distilled or deionized water.
2. Disposable gloves and timer.
3. Appropriate waste containers for potentially contaminated materials.

the microwells to soak for 30- 60seconds. After the final washing cycle, turn the strips plate onto blotting paper or clean towel, and tap it to remove any remainders

7. **Adding HRP Conjugate:** Add 100µl HRP-Conjugate to each well except the Blank.
8. **HRP-Conjugate Incubation (2):** Cover the plate with the plate cover and incubate for 30 minutes at 37°C.
9. **Washing (2):** At the end of the incubation, remove and discard the plate cover. Wash each well 5times with diluted Wash buffer as in Step 6.
10. **Coloring:** Dispense 50µl of Chromogen A and 50µl Chromogen B solution into each well including the **Blank**, and mix by tapping the plate gently. Incubate the plate at 37° for 15minutes avoiding light. The enzymatic reaction between the Chromogen A/B solutions produces blue color in Positive control and anti-HCV positive sample wells.
11. **Stopping Reaction:** Using a multichannel pipette or manually, add 50µl Stop Solution into each well and mix gently. Intensive yellow color develops in Positive control and anti-HCV positive sample wells.
12. **Measuring the Absorbance:** Calibrate the plate reader with the Blank well and read the absorbance at 450nm.If a dual filter instrument is used, set the reference wavelength at 630nm. Calculate the Cut-off value and evaluate the results (**Note:** read the absorbance within 5 minutes after stopping the reaction).

RESULTS

Each microplate should be considered separately when calculating and interpreting results of the assay, regardless of the number of plates concurrently processed. The results are calculated by relating each sample's optical density (OD) value to the Cut-off value (C.O.) of the plate. If the Cut-off reading is based on single filter plate reader, the results should be calculated by subtracting the Blank well OD value from the print report values of samples and controls. In case the reading is based on dual filter plate reader, do not subtract the Blank well OD from the print report values of samples and controls.

1. **Calculation of Cut-off value:** $\text{Cut-off value (C.O.)} = *Nc + 0.12$

*Nc = the mean absorbance value for three negative controls.

Important: If the mean OD value of the negative control is lower than 0.02, take it as 0.02 .If higher than 0.02 see the Quality control range.

Example:

1. Calculation of Nc:

Well No	B1	C1	D1
Negative controls OD value	0.02	0.012	0.016

Nc= 0.016 (the Nc value is lower than 0.02 so take it as 0.02)

2. Calculation of Cut-off value: $\text{Cut-off (C.O.)} = 0.02 + 0.12 = 0.140$

If one of the Negative control values does not meet the Quality control range specifications, it should be discarded and the mean value is calculated again using the remaining two values. If more than one negative control OD value does not meet the Quality control range specifications, the test is invalid and must be repeated.

2. Quality control range:

The test results are valid if the Quality Control criteria are verified. It is recommended that each laboratory must establish appropriate quality control system with quality control material similar to or identical with the patient sample being analyzed

1. The OD value of the Blank well, which contains only Chromogens and Stop solution, is less than 0.080 at 450 nm.
2. The OD value of the Positive control must be equal to or greater than 0.800 at 450/630nm, or at 450nm after blanking.
3. The OD value of the Negative control must be less than 0.100 at 450/630nm or at 450nm after blanking.

3. Interpretations of the results:

(S = the individual optical density (OD) of each specimen)

Negative Results (S/C.O. <1): samples giving an optical density greater than, or equal to the Cut-off value are considered initially reactive, which indicates that antibodies to hepatitis C virus have probably been detected using this anti-HCV ELISA kit. Retesting in duplicates of any initially reactive sample is recommended. Repeatedly reactive samples could be considered positive for antibodies to HCV and therefore the patient is probably infected with hepatitis C virus. Blood unit positive for HCV antibodies should be immediately discarded.

Positive Results (S/C.O. ≥1): samples giving an optical density greater than, or equal to the Cut-off value are considered initially reactive, which indicates that antibodies to hepatitis C virus have probably been detected using this anti-HCV ELISA kit. Retesting in duplicates of any initially reactive sample is recommended. Repeatedly reactive samples could be considered positive for antibodies to HCV and therefore the patient is probably infected with hepatitis C virus. Blood unit positive for HCV antibodies should be immediately discarded.

Borderline: Samples with optical density O.D. ≤ Cut-off ×2 are considered borderline and retesting of those samples in duplicates is recommended. Repeatedly positive samples could be considered positive for hepatitis C virus infection.

Follow-up and supplementary testing of any anti-HCV positive samples with other analytical system (e.g. RIBA, WB) is required to confirm the diagnosis

PERFORMANCE CHARACTERISTICS

Clinical Specificity: A blood donor population of 2948 individuals was tested with 3 different kits from different manufacturers. The specificity of this anti-HCV ELISA kit was 99.55%.

Manufactures	-	+	False pos.	Specificity
Manufacturer 1*	2896	38	14	99.52
Manufacturer 2*	2895	38	15	99.48
anti-HCV ELISA	2897	38	13	99.55

Clinical Sensitivity: Among 480 clinical hepatitis C patients confirmed positive by RIBA 3.0, 480 were positive when tested with this anti-HCV EISA kit. The sensitivity was 100%.

Analytical Specificity:

1. No cross reactivity observed with samples from patients infected with HAV, HIV, HBV, CMV, and TP.
No interference from rheumatoid factors up to 2000U/ml
2. This assay performance characteristics are unaffected from elevated concentrations of bilirubin, hemoglobin, and triolein.
3. Frozen specimens have been tested to check for interferences due to collection and storage.

Performance on BBI low and mixed titers anti-HCV panels:
**BBI PHV105 ANTI-HCV LOW TITER
 PERFORMANCE PANEL**

No.	BBI Ref. Data			Anti-HCV ELISA		
	Western Blot	EIA Result, S/CO			Lot 1	Lot 2
	RIBA 3.0	EIA 2.0	EIA 3.0			
1	POS	2.2	>4.7	6.6	4.5	
2	POS	>4.9	4.3	6.3	14.9	
3	POS	1.3	4.5	6.7	5.4	
4	POS	2.7	>4.7	4.2	4.7	
5	POS	2.1	4.4	6.6	5.5	
6	POS	2.1	3.3	4.8	6.0	
7	IND	2.7	1.7	4.0	2.8	
8	IND	3.8	2.9	3.6	3.3	
9	IND	1.1	1.7	2.5	2.5	
10	NEG	0.1	0.0	0.0	0.0	
11	IND	2.3	>4.7	3.6	3.0	
12	POS	2.5	4.4	2.4	2.0	
13	IND	2.2	3.5	1.4	1.6	
14	IND	2.7	3.7	5.6	4.2	
15	POS	2.8	1.8	2.3	2.8	

**BBI PHV205 ANTI-HCV MIXED TITER PERFORMANCE
 PANEL**

No.	BBI Ref. Data			Anti-HCV ELISA		
	Western Blot	EIA Result, S/CO			Lot 1	Lot 2
	RIBA 3.0	EIA 3.0	EIA 3.0			
1	POS	6.8	>4.9	7.2	6.4	
2	NEG	0.1	0.1	0.0	0.0	
3	POS	6.8	>4.9	7.8	11.3	
4	POS	6.6	>4.9	8.0	10.9	
5	IND	3.3	2.1	7.2	3.6	
6	IND	>4.9	1.5	5.5	4.2	
7	POS	1.5	2.4	5.4	4.2	
8	IND	4.0	2.3	5.5	4.1	
9	POS	>4.9	2.2	3.3	2.1	
10	POS	3.9	4.5	7.6	6.8	
11	POS	>4.9	>4.9	7.0	10.4	
12	POS	>4.9	2.9	4.8	3.7	
13	POS	>4.9	>4.9	9.4	9.8	
14	POS	>4.9	>4.9	8.0	7.7	
15	POS	>4.9	4.4	8.7	7.3	
16	POS	>4.9	>4.9	7.7	8.6	
17	POS	>4.9	>4.9	8.8	12.0	
18	POS	>4.9	>4.9	6.3	12.7	
19	POS	>4.9	>4.9	6.1	11.3	
20	POS	>4.9	>4.9	6.0	21.4	
21	POS	>4.9	>4.9	5.9	21.4	
22	POS	>4.9	>4.9	6.1	21.4	
23	POS	>4.9	>4.9	6.4	16.9	
24	POS	> 4.9	> 4.9	6.5	21.4	
25	NEG	0.0	NEG	0.0	0.0	

HCV genotype antibody testing (reference standards)

GENOTYPE	SAMPLES	POSITIVE
1a-b	15	15
2a-b	13	13
3a-b	10	10
4h	6	6
5	12	12
6	18	18
TOTAL	74	74

Reproducibility		Within run		Between run	
Sample	No	Mean OD	CV%	Mean OD	CV%
Weak positive	10	0.436	9.1%	0.401	9.5%
Moderate positive	10	0.946	7.0%	0.916	7.5%
Strong positive 1	10	1.917	4.4%	1.895	4.2%
Strong positive 2	10	2.372	3.8%	2.309	4.0%

LIMITATIONS OF PROCEDURE

1. Non-repeatable positive result may occur due to the general biological characteristics of ELISA assays. The assay is design to achieve very high performance characteristics of sensitivity and specificity and the “indirect model” minimizes the unspecific reactions, which can occur due to interference between unknown meters in sample and the rabbit anti-human IgG used as a conjugate. Antibodies may be undetectable during the early stages of the disease and in some immunosuppressed individuals.
2. Positive results must be confirmed with another available method. Any positive result must be interpreted together with the patient clinical information and other laboratory results.
3. Common sources for mistakes: kits beyond the expiry date, bad washing procedures, contaminated reagents, incorrect assay procedure steps, insufficient aspiration during washing, failure to add samples or reagents, equipment, timing, volumes, sample nature and quality.
4. If, after retesting of the initially reactive samples, the assay results are negative, these samples should be considered as non-repeatable (false positive) and interpreted as negative. As with many very sensitive ELISA assays, false positive results can occur due to the several reasons, most of which are related but not limited to inadequate washing step.
5. The prevalence of the marker will affect the assay’s predictive values.
6. This kit is intended ONLY for testing of individual serum or plasma samples. Do not use it for testing of cadaver samples, saliva, urine or other body fluids, or pooled (mixed) blood.

INDICATIONS OF INSTABILITY OR DETERIORATION OF THE REAGENTS

1. Values of the Positive or Negative controls, which are out of the indicated Quality control range, are indicator of possible deterioration of the reagents and/or operator or equipment errors. In such case, the results should be considered as invalid and the samples must be retested. In case of constant erroneous results classified as due to deterioration or instability of the reagents, immediately substitute the reagents with new ones.
2. If after mixing of the Chromogen A and B solutions into the wells, the, the color of the mixture turns blue within few minutes, do not continue carrying out the testing and replace the reagents with fresh ones.

VALIDITY

Please do not use this kit beyond the expiry date indicated on the kit box and reagent labels.

PRECAUTIONS

This kit is intended **FOR IN VITRO USE ONLY**



FOR PROFESSIONAL USE ONLY

The ELISA assay is a time and temperature sensitive method. To avoid incorrect result, strictly follow the test procedure steps and do not modify them.

1. Do not exchange reagents from different lots, or use reagents from other commercially available kits. The components of the kit are precisely matched as to achieve optimal performance during testing.
2. Make sure that all reagents are within the validity indicated on the kit box and are of the same lot. Never use reagents beyond the expiry date stated on reagents labels or on the kit box.
3. **CAUTION - CRITICAL STEP:** Allow the reagents and samples to stabilize at room temperature (18-30°) before use. Shake reagent gently before, and return to 2-8° immediately after use.
4. Use only sufficient volume of sample as indicated in the procedure steps. Failure to do so, may cause in low sensitivity of the assay.
5. Do not touch the bottom exterior of the wells; fingerprints or scratches may interfere with microwell reading.
6. When reading the results, ensure that the plate bottom is dry and there are no air-bubbles inside the wells.
7. Never allow the microplate wells to dry after the washing step. Immediately proceed to the next step. Avoid the formation of air-bubbles when adding the reagents.
8. Avoid assay steps long time interruptions. Assure same working conditions for all wells.
9. Calibrate the pipette frequently to assure the accuracy of samples/reagents dispensing. Always use different disposal pipette tips for each specimen and reagents as to avoid cross-contaminations. Never pipette solutions by mouth.
10. The use of automatic pipettes is recommended.
11. Assure that the incubation temperature is 37° inside the incubator.
12. When adding samples, avoid touching the well's bottom with the pipette tip.
13. When reading the results with a plate reader, it is recommended to determine the absorbance at 450nm or at 450nm with reference at 630nm.
14. All specimens from human origin should be considered as potentially infectious.
15. Materials from human origin may have been used in the kit. These materials have been tested with tests kits with accepted performance and found negative for antibodies to HIV 1/2, HCV, TP and HBsAg. However, there is no analytical method that can assure that infectious agents in the specimens or reagents are completely absent. Therefore, handle reagents and specimens with extreme caution as if capable of transmitting infectious diseases. Strict adherence to GLP (Good Laboratory Practice) regulations can ensure the personal safety. Never eat, drink, smoke, or apply cosmetics in the assay laboratory.
16. Bovine derived sera may have been used in this kit. Bovine serum albumin (BSA) and fetal calf sera (FCS) are derived from animals from BSE/TSE free-geographical areas.
17. The pipette tips, vials, strips and sample containers should be collected and autoclaved for 1hour at 121°C or treated with 10% sodium hypochlorite for 30minutes to decontaminate before any further steps for disposal.
18. The Stop solution (2M H₂SO₄) is a strong acid. Corrosive. Use it with appropriate care. Wipe up spills immediately or wash with water if come into contact with the skin or eyes. ProClin 300 used as a preservative can cause sensation of the skin.
19. The enzymatic activity of the HRP-conjugate might be affected from dust, reactive chemical, and substances like sodium hypochlorite, acids, alkalins etc. Do not perform the assay in the presence of such substances.
20. Materials Safety Data Sheet (MSDS) available upon request.
21. If using fully automated microplate processing system, during incubation, do not cover the plates with the plate cover. The tapping out of the remainders inside the plate after washing, can also be omitted.

STORAGE

The components of the kit will remain stable through the expiration date indicated on the label and package when stored between 2-8 °C, **do not freeze**. To assure maximum performance of this anti-HCV ELISA kit, during storage protect the reagents from contamination with microorganism or chemicals.

<p>ISO 13485 ISO 9001</p> 	
 <p>Diagnostic Automation/ Cortez Diagnostics, Inc. 23961 Craftsman Road, Suite E/F, Calabasas, California 91302 USA</p>	
<p>Date Adopted 2008-01-05</p>	<p>Reference No. 1707-12 AccuDiag™-HCV Ab -2012</p>
<p>Revision B Date: 07-25-2013</p>	