

8

9

# **AccuDiag**<sup>TM</sup> Human Allergen Specific IgG **ELISA Kit**

# Cat# 1809-9

See external Label	2°C 8°C	$\sum$	192 Tests
			1/2 1050

Test	Human Allergen Specific IgG	
Method	Enzyme Linked Immunosorbent Assay	
Principle	Sandwich Complex	
Detection Range	Class I - Class IV	
Sample	10 µL serum/ plasma	
Specificity	N/A	
Sensitivity	N/A	
Total Time	~105 min.	
Shelf Life	12-14 Months from the manufacturing date	

## INTENDED USE

To quantitate total human Immunoglobulin G (IgG) Specific for Common Allergens

# TEST PRINCIPLE

Solid phase capture sandwich ELISA assay using a microwell format.

# Patient and Standard Dilutions:

Dilute each serum or plasma specimen to be tested 1:100 in the Specimen Diluent provided e.g. 10ul of specimen into 990ul of Specimen Diluent.

# MATERIALS AND COMPONENTS

#### Materials provided with the test kits

- 1. 2-Allergen coated microwell strips 12x8 with plastic frame
- 2. 2-HRP conjugated goat anti-human IgG-12mL
- 3. 2-TMB/peroxide substrate color developer -12mL
- 4. IgG Specimen Diluent 1x 60 mL
- 5. 2-Sulfuric acid termination reagent (0.5N) -12mL
- 6. 15 X Wash buffer concentrate 60mL

# ASSAY PROCEDURE

\*Allow each reagent to reach room temperature before use

- 1. Add 100uL of *diluted* specimen to each microwell
- 2. Incubate at room temperature for 45 minutes
- 3. Decant and wash each microwell four times with wash buffer (dilute buffer 1:15 with reagent grade water)
- 4. Add 100uL of HRP conjugated goat anti-human IgG to each well
- 5. Incubate at room temperature for 45 minutes
- 6. Decant and wash as in step 3

Diagnostic Automation/Cortez Diagnostics, Inc.

23961 Craftsman Road, Suite D/E/F, Calabasas, CA 91302 USA Phone: 818.591.3030 Fax 818.591.8383 Email:onestep@rapidtest.com Website: www.rapidtest.com

Zero the microwell reader at 450nm using the specimen diluent zero control well 10. Determine the optical density (O.D.) of the remaining wells

Terminate the reaction with 100uL of 0.5N sulfuric acid

11. Optical Densities are interpolated into Classes:

## RESULTS

minutes

O.D.@450nm	Class
0.000-0.250	Class O (Negative)
0.250-0.350	Class I (Equivocal)
0.350-0.450	Class II
0.450-0.550	Class III
>0.550	Class IV

## LIMITATIONS OF PROCEDURE

No single assay should be used as the only basis for arriving at a diagnostic conclusion. For research use only.

#### Shelf Life

The expiration date for the package and each component is stated on the label(s). Store components at 2-8°C.

