LDL DIRECT

| | (PL) | CORMAY | HC- | OS- | B50- |
|-----------------------|---------|------------|------------|-----------|-----------|
| Nazwa zestawu | Nr kat. | LDL DIRECT | LDL | LDL | LDL |
| CORMAY LDL DIRECT | 2-180 | | | | |
| CORMAY LDL DIRECT 30 | 2-191 | 1-REAGENT | 2 x 37 ml | 2 x 34 ml | 1 x 48 ml |
| CORMAY LDL DIRECT 60 | 2-192 | | | | |
| CORMAY LDL DIRECT 120 | 2-193 | 2-REAGENT | 1 x 100 ml | 2 x 12 ml | 1 x 18 ml |
| HC- LDL DIRECT | 4-569 | | | | |
| OS- LDL DIRECT | 9-412 | | | | |
| B50- LDL DIRECT | 5-509 | | | | |

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The reagent is stable up to the kit expiry date printed on the package when stored at 2-8°C. The reagents are stable for 12 weeks on board the analyser at 2-10°C.

Concentrations in the reagent

| 1-REAGENT | 2-REAGENT |
|---|----------------|
| Buffer | |
| Detergent 1 | < 1.0 % |
| Cholesterol esterase (<i>Pseudomonas</i> spp.) | < 1500 U/l |
| Cholesterol oxidase (<i>Cellulomonas</i> spp.) | < 1500 U/l |
| Peroxidase (horseradish) | < 1300 ppg U/l |
| 4-aminoantipyrine | < 0.1 % |
| Ascorbic Acid Oxidase (<i>Curcubita</i> spp.) | < 3000 U/l |
| Preservative | |
| Detergent 2 | < 1.0 % |
| N,N-bis(sulfobutyl)-toluidine, disodium (DSBmT) | < 1.0 mM |
| Preservative | |

Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Do not freeze reagents.
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.
- 1-REAGENT and 2-REAGENT meeting the criteria for classification in accordance with Regulation (EC) No 1272/2008.

Ingredients:

1-REAGENT and 2-REAGENT contain reaction mass of 5-chloro-2-methyl-2H-isothiazol-3-one and 2-methyl-2H-isothiazol-3-one (3:1).

Warning

H317 May cause an allergic skin reaction.
P280 Wear protective gloves, protective clothing, eye protection or face protection.
P302 + P352 IF ON SKIN: Wash with plenty of soap and water.

SPECIMEN

Serum, sodium heparinized or EDTA plasma. Anticoagulants containing citrate should not be used. **Serum:** Collect whole blood by venipuncture and allow to clot. Centrifuge and remove the serum as soon as possible after collection (within 3 hours).

Plasma: Specimens may be collected in EDTA or sodium heparin. Centrifuge and remove the plasma as soon as possible after collection (within 3 hours).

If not analysed promptly, specimens may be stored at 2-8°C for up to 5 days. If specimens need to be stored for more than 5 days, they may be preserved at -80°C. Samples may be frozen once.

Nevertheless it is recommended to perform the assay with freshly collected samples.

ADDITIONAL EQUIPMENT

- automated analyser or photometer able to read at 630 nm;
- thermostat at 37°C;
- general laboratory equipment;

PROCEDURE

The reagents are ready to use.

Applications for automatic analysers are available on request.

Manual procedure

| | |
|-------------|--------|
| wavelength | 630 nm |
| temperature | 37°C |
| cuvette | 1 cm |

Pipette into the cuvettes:

| | test (T) | standard (S) |
|-----------|-------------|-----------------|
| 1-REAGENT | 900 µl | 900 µl |

Bring up to the temperature of determination (37 °C). Then add

| | | |
|------------|------|------|
| calibrator | - | 9 µl |
| sample | 9 µl | - |

Mix well and after 5 min. of incubation at 37°C add:

| | | |
|-----------|--------|--------|
| 2-REAGENT | 300 µl | 300 µl |
|-----------|--------|--------|

Mix well and after exactly 5 min. of incubation read the absorbance of standard A(S) and test A(T) against air or water.

Calculation

$$\text{LDL cholesterol concentration} = \frac{A(T)}{A(S)} \times \text{calibrator concentration}$$

REFERENCE VALUES ⁷

| Adults | | |
|-----------------------|---------|-----------|
| NCEP* Classification: | mg/dl | mmol/l |
| Optional | < 100 | < 2.59 |
| Near optional | < 130 | < 3.37 |
| Bordeline high | 130-159 | 3.37-4.12 |
| High | 160-189 | 4.14-4.90 |
| Very high | ≥ 190 | ≥ 4.92 |

* National Cholesterol Education Program

As LDL cholesterol is affected by a number of factors such as smoking, exercise, hormones, age and sex, each laboratory should establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use CORMAY LIPID CONTROL 1 (Cat. No 5-179) and CORMAY LIPID CONTROL 2 (Cat. No 5-180) or CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173) with each batch of samples.

For the calibration the CORMAY HDL/LDL CALIBRATOR (Cat. No 5-178) is recommended.

The calibration curve should be prepared every 12, with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using Multi+ for manual assay and automatic analyser Biolis 24i Premium. Results may vary if a different instrument is used.

Sensitivity:

5.5 mg/dl (0.14 mmol/l) – Multi+
4.2 mg/dl (0.11 mmol/l) – Biolis 24i Premium

Linearity:

up to 316 mg/dl (8,18 mmol/l) – Multi+
up to 700 mg/dl (18,13 mmol/l) - Biolis 24i Premium

For higher concentrations dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

Specificity / Interferences

Triglycerides up to 1293 mg/dl, bilirubin conjugated up to 20 mg/dl, bilirubin total up to 20 mg/dl, haemoglobin up to 0.5 g/dl, ascorbic acid up to 500 mg/l and gamma-globulins up to 5000 mg/dl do not interfere with the test.

N-acetyl-p-benzoquinone imine (NAPQI), a metabolite of paracetamol (acetaminophen), may cause falsely low results for patients with a toxic level of paracetamol.

Precision (Multi+)

| Repeatability (run to run) n = 10 | Mean [mg/dl] | SD [mg/dl] | CV [%] |
|--------------------------------------|-----------------|---------------|-----------|
| level 1 | 62.72 | 1.22 | 1.94 |
| level 2 | 158.77 | 0.69 | 0.44 |

Precision (Biolis 24i Premium)

| Repeatability (run to run) n = 10 | Mean [mg/dl] | SD [mg/dl] | CV [%] |
|--|-----------------|---------------|-----------|
| level 1 | 135.60 | 4.35 | 3.21 |
| level 2 | 187.07 | 5.19 | 2.77 |
| Reproducibility (day to day) n = 20 | Mean [mg/dl] | SD [mg/dl] | CV [%] |
| level 1 | 135.06 | 3.67 | 2.71 |
| level 2 | 187.24 | 5.62 | 3.00 |

Method comparison

A comparison between LDL cholesterol values determined at Multi+ (y) and at Hitachi 912 (x) using 48 samples gave following results:

$$y = 0.8968x + 11.543 \text{ mg/dl};$$

$$R = 0.982 \quad (R - \text{correlation coefficient})$$

A comparison between LDL cholesterol values determined at Biolis 24i Premium (y) and at Hitachi 912 (x) using 29 samples gave following results:

$$y = 1.0365x + 3.7224 \text{ mg/dl};$$

$$R = 0.996 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT

Please refer to local legal requirements.

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