## CHOL

Nazwa zestawu	(PL) Nr kat.
Liquick Cor-CHOL 30	2-211
Liquick Cor-CHOL 60	2-204
Liquick Cor-CHOL 120	2-205
HC-CHOL	4-504
OS-CHOL	9-402
B50-CHOL	5-507

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### LITERATURA

- Siedel J., Hägele E.O., Zigenhorn J., Wahlefeld A.W.: Clin. Chem. 29, 6 (1983).
- Tel R.M., Berends G.T.: J Clin.Chem. Clin. Biochem. 18, 10 (1980).
- Rautela G.S., Liedtke R.J.: Clin. Chem. 24, 1 (1978)
- Schettler G., Nussel E.: Arbeitsmed. Sozialmed. Praventivmed. 10, 25 (1975).
- Richmond W.: Clin. Chem. 19, 1350 (1973).
- Roeschlau P., Bernt E., Gruber W.: J. Clin. Chem. Biochem. 12, 403 (1974).
- Tietz N.W., ed. Clinical Guide to Laboratory Tests, 3rd ed. Philadelphia, PA: WB Saunders, 130 (1995).
- Trinder P.: Ann. Clin. Biochem. 6, 24 (1969).
- Dembińska-Kieć A., Naskalski J.W., Solnica B.: Diagnostyka laboratoryjna z elementami biochemii klinicznej, Wyd. IV, str. 670, 875 (2017).
- Cybulska B, Szostak WB, Filipiak KJ, et al. Polish Forum for Prevention Guidelines on Dyslipidaemia: update 2016. Kardiol Pol. 2017; 75(2): 187–190, doi: 10.5603/KP.2017.0031.
- Piepoli MF, Hoes AW, Agewall S, et al. 2016 European Guidelines on cardiovascular disease prevention in clinical practice: The Sixth Joint Task Force of the European Society of Cardiology and Other Societies on Cardiovascular Disease Prevention in Clinical Practice (constituted by representatives of 10 societies and by invited experts) Developed with the special contribution of the European Association for Cardiovascular Prevention & Rehabilitation (EACPR). Eur Heart J. 2016; 37(29): 2315–2381, doi: 10.1093/eurheartj/ehw106, indexed in Pubmed: 27222591.
- NCEP Expert Panel. Arch Inter Med (148), 36-69, (1988).
- Jacobso D.S., DeMott W.R., Grady H.J., et. al., ed., Laboratory Tests Handbook, 4th ed., Hudson, Lexi-Comp, 143, (1996).

Data wydania: 06. 2023.

## CHOL

<b>Kit name</b>	<b>(EN)</b>	<b>Cat. No</b>
Liquick Cor-CHOL 30		2-211
Liquick Cor-CHOL 60		2-204
Liquick Cor-CHOL 120		2-205
HC-CHOL		4-504
OS-CHOL		9-402
B50-CHOL		5-507

### INTENDED USE

Diagnostic kit for determination of total cholesterol concentration intended to use for manual assay and in several automatic analyzers.



The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

### INTRODUCTION

Cholesterol is essential structural component of cell membranes and precursor of bile acids and all steroids hormones. This is why cholesterol has enormous significance for organism normal functioning. But there is also well established association between blood cholesterol concentration and coronary heart disease. Measurement of cholesterol serum level is valuable in prevention and monitoring cardiovascular disease. This determination is useful also for evaluation of intestine absorption, liver and gallbladder function.

### METHOD PRINCIPLE

Enzymatic, colorimetric method with cholesterol esterase and cholesterol oxidase.

The colour intensity is proportional to the cholesterol concentration.

### REAGENTS

<b>Package</b>	<b>Liquick Cor-CHOL 30</b>	<b>Liquick Cor-CHOL 60</b>	<b>Liquick Cor-CHOL 120</b>
1-REAGENT	6 x 30 ml	6 x 60 ml	6 x 120 ml
2-STANDARD	1 x 2 ml	1 x 2 ml	-
	<b>HC-CHOL</b>	<b>OS-CHOL</b>	<b>B50-CHOL</b>
1-REAGENT	6 x 96.5 ml	2 x 56 ml	4 x 58.5 ml

2-STANDARD is a standard solution of cholesterol with a concentration within the range 180 - 220 mg/dl (4.66 – 5.70 mmol/l). Exact concentration is printed on the label of the each vial.

### Working reagent preparation and stability

CHOL

The reagent is ready to use.  
The reagent when stored at 2-8°C is stable up to expiry date printed on the package. The reagents are stable for 12 weeks on board the analyzer at 2-10°C.

### Concentrations in the test

Good's buffer (pH 6.4)	< 120 mmol/l
phenol	< 6 mmol/l
4-aminoantipyrine	< 0.4 mmol/l
cholesterol esterase (CHE)	< 4 µkat/l
cholesterol oxidase (CHO)	< 5 µkat/l
peroxidase (POD)	< 24 µkat/l
preservatives, detergents, stabilizers	

### Warnings and notes

- Protect from direct sunlight and avoid contamination!
- The reagents are usable when the absorbance of the working reagent is less than 0.100 (read against distilled water, wavelength λ=500 nm).
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.
- EUH210 Safety data sheet available on request.
- 1-REAGENT meeting the criteria for classification in accordance with Regulation (EC) No 1272/2008.

### Warning

H319 Causes serious eye irritation.  
P280 Wear protective gloves, protective clothing, eye protection or face protection.  
P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

### ADDITIONAL EQUIPMENT

- automatic analyzer or photometer able to read at 500 nm (Hg 546 nm);
- thermostat at 37°C;
- general laboratory equipment;

### SPECIMEN

Serum, EDTA or heparinized plasma (recommended: heparine lithium, sodium or ammonium salt) free from hemolysis. Blood should be collected only if the patient has been fasting for minimum of 12 hours. Before blood collection patient should stay in rest position for about 30 minutes. Venous blood is recommended for cholesterol measurement.

Plasma cholesterol values have been reported to be 3% to 5% lower than serum cholesterol values. Serum should be separated from red blood cells as soon as possible after blood collection.

Serum and plasma can be stored up to 3 days at 2-8°C or 6 months at -20°C. Nevertheless it is recommended to perform the assay with freshly collected samples!

### PROCEDURE

Applications for analysers are available on request.

### Manual procedure

wavelength	500 nm (Hg 546 nm)
temperature	20-25°C / 37°C
cuvette	1 cm

Pipette into the cuvettes:

	reagent blank (RB)	test (T)	standard (S)
1-REAGENT	1000 µl	1000 µl	1000 µl
Bring up to the temperature of determination. Then add:			
standard / calibrator	-	-	10 µl
sample	-	10 µl	-

Mix well, incubate for 5 min. at 37°C or 10 min. at 20-25°C. Read the absorbance of the test A(T) and standard A(S) against reagent blank (RB).

### Calculation

$$\text{cholesterol concentration} = \frac{A(T)}{A(S)} \times \text{standard / calibrator concentration}$$

### REFERENCE VALUES

serum / plasma	mg/dl	mmol/l
children <sup>9</sup>	< 1 year	66.1 – 228.5
	1 – 19 years	111.4 – 202.2
adults <sup>9,10,11</sup>		< 190
		< 5.00

It is recommended for each laboratory to establish its own reference ranges for local population.

### QUALITY CONTROL

For internal quality control it is recommended to use the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173) with each batch of samples. For the calibration is also recommended the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176), LEVEL 2 (Cat. No 5-175; 5-177) or CHOLESTEROL STANDARD 200 (Cat. No 5-118).

The calibration curve should be prepared every 12 weeks, with change of reagent lot number or as required e.g. quality control findings outside the specified range.

### PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analyzer Biolis 24i Premium. Results may vary if a different instrument or a manual procedure is used.

#### Sensitivity

13 mg/dl (0.337 mmol/l)

#### Linearity

up to 730 mg/dl (18.9 mmol/l)

#### Specificity / Interferences

Haemoglobin up to 0.31 g/dl, bilirubin up to 20 mg/dl, ascorbate up to 62 mg/l and triglycerides up to 1000 mg/dl do not interfere with the test.

N-acetyl-p-benzoquinone imine (NAPQI), a metabolite of paracetamol (acetaminophen), may cause falsely low results for patients with a toxic level of paracetamol.

#### Precision

Repeatability (run to run) n = 20	Mean [mg/dl]	SD [mg/dl]	CV [%]
level 1	254.84	2.62	1.03
level 2	106.28	1.11	1.05
Reproducibility (day to day) n = 20	Mean [mg/dl]	SD [mg/dl]	CV [%]
level 1	256.65	3.41	1.33
level 2	110.86	2.61	2.35

#### Method comparison

CHOL

A comparison between cholesterol values determined at **Biolis 24i Premium** (y) and at **ADVIA 1650** (x) using 100 samples gave following results:

$$y = 1.0132x + 2.995 \text{ mg/dl};$$

$$R = 0.9964 \quad (R - \text{correlation coefficient})$$

### TRACEABILITY

CHOLESTEROL STANDARD 200 are traceable to the SRM 1951C reference material.

### WASTE MANAGEMENT

Please refer to local legal requirements.

### LITERATURE

- Siedel J., Hägele E.O., Zigenhorn J., Wahlefeld A.W.: Clin. Chem. 29, 6 (1983).
- Tel R.M., Berends G.T.: J Clin.Chem. Clin. Biochem. 18, 10 (1980).
- Rautela G.S., Liedtke R.J.: Clin. Chem. 24, 1 (1978)
- Schettler G., Nussel E.: Arbeitsmed. Sozialmed. Präventivmed. 10, 25 (1975).
- Richmond W.: Clin. Chem. 19, 1350 (1973).
- Roeschlau P., Bernt E., Gruber W.: J. Clin. Chem. Biochem. 12, 403 (1974).
- Tietz N.W., ed. Clinical Guide to Laboratory Tests, 3rd ed. Philadelphia, PA: WB Saunders, 130 (1995).
- Trinder P.: Ann. Clin. Biochem. 6, 24 (1969).
- Dembińska-Kieć A., Naskalski J.W., Solnica B.: Diagnostyka laboratoryjna z elementami biochemii klinicznej, Wyd. IV, str. 670, 875 (2017).
- Cybulska B, Szostak WB, Filipiak KJ, et al. Polish Forum for Prevention Guidelines on Dyslipidaemia: update 2016. Kardiol Pol. 2017; 75(2): 187–190, doi: 10.5603/KP.2017.0031.
- Piepoli MF, Hoes AW, Agewall S, et al. 2016 European Guidelines on cardiovascular disease prevention in clinical practice: The Sixth Joint Task Force of the European Society of Cardiology and Other Societies on Cardiovascular Disease Prevention in Clinical Practice (constituted by representatives of 10 societies and by invited experts) Developed with the special contribution of the European Association for Cardiovascular Prevention & Rehabilitation (EACPR). Eur Heart J. 2016; 37(29): 2315–2381, doi: 10.1093/eurheartj/ehw106, indexed in Pubmed: 27222591.
- NCEP Expert Panel. Arch Inter Med (148), 36-69, (1988).
- Jacobso D.S., DeMott W.R., Grady H.J., et. al., ed., Laboratory Tests Handbook, 4th ed., Hudson, Lexi-Comp, 143, (1996).

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1. Siedel J., Hägele E.O., Zigenhorn J., Wahlefeld A.W.: Clin. Chem. 29, 6 (1983).
2. Tel R.M., Berends G.T.: J Clin.Chem. Clin. Biochem. 18, 10 (1980).
3. Rautela G.S., Liedtke R.J.: Clin. Chem. 24, 1 (1978)
4. Schettler G., Nussel E.: Arbeitsmed. Sozialmed. Praventivmed. 10, 25 (1975).
5. Richmond W.: Clin. Chem. 19, 1350 (1973).
6. Roeschlau P., Bernt E., Gruber W.: J. Clin. Chem. Biochem. 12, 403 (1974).
7. Tietz N.W., ed. Clinical Guide to Laboratory Tests, 3rd ed. Philadelphia, PA: WB Saunders, 130 (1995).
8. Trinder P.: Ann. Clin. Biochem. 6, 24 (1969).
9. Dembińska-Kieć A., Naskalski J.W., Solnica B.: Diagnostyka laboratoryjna z elementami biochemii klinicznej, Wyd. IV, str. 670, 875 (2017).
10. Cybulska B, Szostak WB, Filipiak KJ, et al. Polish Forum for Prevention Guidelines on Dyslipidaemia: update 2016. Kardiol Pol. 2017; 75(2): 187–190, doi: 10.5603/KP.2017.0031.
11. Piepoli MF, Hoes AW, Agewall S, et al. 2016 European Guidelines on cardiovascular disease prevention in clinical practice: The Sixth Joint Task Force of the European Society of Cardiology and Other Societies on Cardiovascular Disease Prevention in Clinical Practice (constituted by representatives of 10 societies and by invited experts) Developed with the special contribution of the European Association for Cardiovascular Prevention & Rehabilitation (EACPR). Eur Heart J. 2016; 37(29): 2315–2381, doi: 10.1093/eurheartj/ehw106, indexed in Pubmed: 27222591.
12. NCEP Expert Panel. Arch Inter Med (148), 36-69, (1988).
13. Jacobso D.S., DeMott W.R., Grady H.J., et. al., ed., Laboratory Tests Handbook, 4th ed., Hudson, Lexi-Comp, 143, (1996).

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