GLUCOSE

Nazwa zestawu	(PL)	Nr kat.
Liquick Cor-GLUCOSE 30	2-219	2-219
Liquick Cor-GLUCOSE 60	2-201	2-201
Liquick Cor-GLUCOSE 120	2-202	2-202
HC- GLUCOSE	4-501	4-501
OS- GLUCOSE	9-401	9-401
B50- GLUCOSE	5-510	5-510

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Data wydania: 06. 2023.

GLUCOSE

Kit name	(EN) Cat. No
Liquick Cor-GLUCOSE 30	2-219
Liquick Cor-GLUCOSE 60	2-201
Liquick Cor-GLUCOSE 120	2-202
HC- GLUCOSE	4-501
OS- GLUCOSE	9-401
B50- GLUCOSE	5-510

INTENDED USE



GLUCOSE reagent is intended to determine quantitatively glucose level in plasma, serum and cerebrospinal fluid. It is intended to diagnosis, monitoring and as an aid in the diagnosis of clinical conditions connected with abnormal glucose level. GLUCOSE reagent is intended to use to use on automatic analysers and manually. 2-STANDARD is intended to calibrate manual determinations of glucose, performed with Liquick Cor-GLUCOSE 30 and Liquick Cor-GLUCOSE 60 kit. It is only for *in vitro* diagnostics, for healthcare professional users.

SUMMARY ^{1,2}

Determinations of glucose level in serum, plasma and cerebrospinal fluid are mainly used in diagnosis and treatment monitoring of diabetes mellitus. Measurement of glucose concentration is also used as an aid in diagnosis of clinical conditions associated with hypoglycemia and hyperglycemia. Increased glucose levels can be observed in transient hyperglycemia, acute stress reaction, Cushing's syndrome, glucocorticoid therapy, hyperthyroidism. Decreased concentrations of glucose can be associated with neonatal hypoglycemia or hypothyroidism.

METHOD PRINCIPLE ³

Colorimetric, enzymatic method with glucose oxidase.

(red colour)

The colour intensity is proportional to the glucose concentration.

REAGENTS

Package	Liquick Cor- GLUCOSE 30	Liquick Cor- GLUCOSE 60	Liquick Cor- GLUCOSE 120
1-REAGENT	6 x 30 ml	6 x 60 ml	6 x 120 ml
2-STANDARD	1 x 2 ml	1 x 2 ml	-
	HC- GLUCOSE	OS- GLUCOSE	B50- GLUCOSE
1-REAGENT	6 x 96.5 ml	4 x 53.5 ml	4 x 58 ml

2-STANDARD is a standard solution of glucose with a concentration within the range 90 - 110 mg/dl (5.0 - 6.1 mmol/l). The exact concentration is printed on the label of each vial.

CONCENTRATION OF THE ACTIVE INGREDIENTS IN THE REAGENT

1-REAGENT	
phenol	5 mmol/l
glucose oxidase (GOD)	333.3 μ kat/l
peroxidase (POD)	38.33 μ kat/l
4- aminoantipyrine (4-AA)	0.75 mmol/l
phosphate buffer	
metal ion chelator	
stabilizer	
preservatives	
2-STANDARD	
glucose	100 mg/dl
acetate buffer	
preservative	

REAGENT STABILITY

The reagent when stored at 2-8°C is stable up to expiry date printed on the package. The reagents are stable for 12 weeks on board the analyser at 2-10°C (Biolis 30i).

WARNINGS AND NOTES

- Protect from direct sunlight and avoid contamination!
- The reagents are usable when the absorbance is less than 0.300 (read against distilled water, wavelength $\lambda=500$ nm, cuvette l=1 cm, at temp. 25°C).
- Do not use reagent beyond expiry date printed on the package.
- Do not mix reagents from different kits or lots.
- Use personal protective equipment to prevent contact with samples, reagents and controls.
- EUH208 Contains reaction mass of 5-chloro-2-methyl-2H-isothiazol-3-one and 2-methyl-2H-isothiazol-3-one (3:1). May produce an allergic reaction.
- EUH210 Safety data sheet available on request.

EQUIPMENT REQUIRED BUT NOT PROVIDED

- automatic analyzer or photometer able to read at 500 nm (Hg 546 nm);
- thermostat at 37°C;
- general laboratory equipment;

SPECIMEN ^{4,5,6}

EDTA or heparinized plasma in tubes containing sodium fluoride or sodium iodoacetate additive/ serum, free from hemolysis, cerebrospinal fluid.

Plasma / Serum. Serum and plasma specimens should be separated from cells within 30 minutes after collection. Plasma specimen which is not assayed immediately after collection should be kept in tubes containing sodium fluoride or sodium iodoacetate. These compounds adding prevent glycolysis and stabilize glucose level. Serum and plasma can be stored up to 8 hours at 25°C or 72 hours at 4 °C.

Plasma is the specimen recommended for the glucose determination in the blood.

Cerebrospinal fluid. Glucose concentration in cerebrospinal fluid should be measured directly after specimen collection. Cerebrospinal fluid must be analysed simultaneously with a blood sample.

After centrifuge CSF sample can be stored up to 24 hours at 4°C.

Nevertheless it is recommended to perform the assay with freshly collected samples!

Follow tube manufacturers' instructions carefully when using collection tubes.


Human-origin material should be handled as potentially infectious. Standard precautions in normal laboratory work are required.

ASSAY PROCEDURE

Reagents are ready to use.

Carry-over effect may occur when the reagent is used on automated chemistry analyzers. It is based on interference between specific test combinations and results in either underestimating or inflating patient sample results. To minimize this effect, the preventive actions should be taken, such as ordering an extra washing program or performing tests in a separate order.

General procedure for automatic determinations:

Main Wavelength: 505 nm

2nd wavelength: 700 nm

Method: End Point

Calibration Type: Linear

Direction: Increasing

R1 Volume: 250 μ l

Sample Volume: 3 μ l

Temperature: 37°C

Application guides for analysers are available on request.

The application guide should be verified and validated by the user prior to testing patient samples.

Manual procedure

wavelength	500 nm (Hg 546 nm)
temperature	37°C
cuvette	1 cm

Pipette into the cuvettes:

	reagent blank (RB)	test (T)	standard (S)
1-REAGENT	1000 μ l	1000 μ l	1000 μ l

Bring up to the temperature of determination. Then add:

standard / calibrator	-	-	10 μ l
sample	-	10 μ l	-

Mix well, incubate for 5 min. at 37°C. Read the absorbance of the test A(T) and standard A(S) against reagent blank A(RB).

Calculation

$$\text{glucose concentration} = \frac{A(T)}{A(S)} \times \text{standard / calibrator concentration}$$

QUALITY CONTROL

For internal quality control it is recommended to use the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173) with each batch of samples.

For calibration when using the manual method the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176), LEVEL 2 (Cat. No 5-175; 5-177), GLUCOSE STANDARD 100 (Cat. No 5-121) or 2-STANDARD are recommended.

For the calibration of automatic analysers systems the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) is recommended. Deionised water should be used as a calibrator 0.

Calibration stability depends on the type of analyser used for analysis. The calibration curve should be prepared every 10 weeks (Biolis 30i).

Calibration is recommended in the following cases:

- after each change of lot,
- after instrument service,
- if controls lie outside the expected range,
- each time a new reagent kit is used.

If the quality control results do not fall within the expected values or within the range determined in the laboratory, despite a successful calibration procedure, do not report results. In this case, please take the following actions:

- verify reagents are not out of expiration date.
- verify that the required maintenance has been carried out.
- verify that the procedure has been performed in accordance with the instructions for use.
- contact the Service Department or distributor for assistance.

REFERENCE VALUES ^{4,7,8,9}

	mg/dl	mmol/l
plasma, serum ^{7,8,9}	70 - 99	3.9 - 5.5
cerebrospinal fluid ⁴	40 - 70	2.2 - 3.9

It is recommended for each laboratory to establish its own reference ranges for local population.

It is recommended for each laboratory to establish its own reference ranges for local population and method principle.

The assay results should be used in conjunction with other data, such as symptoms, results of other tests and clinical history to make clinical decisions. It is not recommended to make clinical diagnosis based on a single result.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using a Multi+ analyser for manual assay and automatic analysers Biolis 30i and/or BS-400. Results may vary if a different instrument or a manual procedure is used.

LoB (Limit of Blank):

0.0 mg/dl (0.0 mmol/l) - BS-400

LoD (Limit of Detection):

0.3 mg/dl (0.017 mmol/l) - BS-400

LoQ (Limit of Quantitation):

12.0 mg/dl (0.67 mmol/l) - Multi+
6.0 mg/dl (0.33 mmol/l) - Biolis 30i

- **Linearity:**
up to 350 mg/dl (19.4 mmol/l) – Multi+
up to 600 mg/dl (33.3 mmol/l) – Biolis 30i

For higher concentration dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

- **Measurement range:**

Multi+
12.0 mg/dl (0.67 mmol/l) - 350 mg/dl (19.4 mmol/l)

Biolis 30i
6.0 mg/dl (0.33 mmol/l) - 600 mg/dl (33.3 mmol/l)

- **Specificity / Interferences**

a) samples containing low glucose concentration

Haemoglobin up to 0.63 g/dl, ascorbate up to 62 mg/l, bilirubin up to 20 mg/dl and triglycerides up to 500 mg/dl do not interfere with the test.

b) samples containing high glucose concentration

Haemoglobin up to 2.50 g/dl, ascorbate up to 62 mg/l, bilirubin up to 20 mg/dl and triglycerides up to 1000 mg/dl do not interfere with the test.

- **Precision (Multi+)**

Repeatability (run to run) n = 20	Mean [mg/dl]	SD [mg/dl]	CV [%]
level 1	94.4	1.70	1.8
level 2	301.6	6.10	2.0

- **Precision (Biolis 30i)**

Repeatability (run to run) n = 20	Mean [mg/dl]	SD [mg/dl]	CV [%]
level 1	86.6	0.45	0.5
level 2	293.6	2.38	0.8
Reproducibility (day to day) n = 80	Mean [mg/dl]	SD [mg/dl]	CV [%]
level 1	87.8	2.23	2.5
level 2	291.6	5.32	1.8

- **Method comparison**

A comparison between glucose values determined at **Multi+** (y) and at **BECKMAN COULTER AU680** (x) using 21 plasma samples gave following results:

$$y = 0.9547x + 4.2809 \text{ mg/dl};$$

$$R = 0.999 \quad (R - \text{correlation coefficient})$$

A comparison between glucose values determined at **Biolis 30i** (y) and at **SIEMENS ADVIA 1800** (x) using 59 plasma samples gave following results:

$$y = 1.0133x + 0.1086 \text{ mg/dl};$$

$$R = 0.999 \quad (R - \text{correlation coefficient})$$

A comparison between glucose values determined at **Biolis 30i** (y) and at **SIEMENS ADVIA 1800** (x) using 56 serum samples gave following results:

$$y = 1.0559x - 1.8604 \text{ mg/dl};$$

$$R = 1.000 \quad (R - \text{correlation coefficient})$$

A comparison between glucose values determined at **Biolis 30i** (y) and at **BECKMAN COULTER AU680** (x) using 30 PMR samples gave following results:

$$y = 1.02x + 0.4699 \text{ mg/dl};$$

$$R = 0.992 \quad (R - \text{correlation coefficient})$$

TRACEABILITY

2-STANDARD and GLUCOSE STANDARD 100 are traceable to the SRM 965B reference material.

WASTE MANAGEMENT¹⁰

After use, the reagents should be handled as potentially infectious and disposed of in accordance with local legal requirements.

- Reagents residues:
1-REAGENT - 18 01 07
2-STANDARD - 18 01 07
- Empty packages:
1-REAGENT - 15 01 02
2-STANDARD – 15 01 07
- Wastewater from the analyzer:
1-REAGENT - 18 01 03*
2-STANDARD/ 1-REAGENT - 18 01 03*

INCIDENTS¹¹

Any serious incident that has occurred in relations to the device shall be reported to the manufacturer (website address: incidents@cornmay.pl) and the competent authority of the Member State in which the user and/or the patient is established.

Serious incident means any incident that directly or indirectly led, might have led or might lead to any of the following:

- the death of a patient, user or other person,
- the temporary or permanent serious deterioration of a patient's, user's or other person's state of health,
- a serious public health threat.

LITERATURE

1. Padhi, S., Nayak, A. K., Behera, A. Type II diabetes mellitus: a review on recent drug based therapeutics. Biomedicine & Pharmacotherapy, 131, 110708 (2020).
2. Pisarczyk-Wiza D., Zozulińska-Ziólkiewicz D. Glikokortykosteroidy a zaburzenia metabolizmu glukozy. Diabetologia Kliniczna 4,3, 110-116 (2015).
3. Barham P., Trinder P. An improved colour reagent for the determination of blood glucose by the oxidase system: Analyst 97, 142-145 (1972).
4. Alan H.B. Wu: Tietz Clinical Guide to Laboratory Tests, 4th ed. WB Saunders, 444-450 (2006).
5. Rifai N., Horvath A.R., Wittwer C., ed. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics 6th ed., Elsevier, St. Louis, USA, 528 (2018).
6. Dujmovic and F. Deisenhammer, Stability of cerebrospinal fluid/serum glucose ratio and cerebrospinal fluid lactate concentrations over 24 h: analysis of repeated measurements, Clinical Chemistry and Laboratory Medicine 48, 2, 209-212 (2010).
7. 2021 Guidelines on the management of patients with diabetes. A position of Diabetes Poland
8. Sacks, David B., et al. Guidelines and recommendations for laboratory analysis in the diagnosis and management of diabetes mellitus. Clinical Chemistry 48, 3, 436-472 (2002).
9. Miles RR, Roberts RF, Putnam AR, Roberts WL. Comparison of serum and heparinized plasma samples for measurement of chemistry analytes. Clinical Chemistry 50, 1704-1706 (2004).
10. European Commission notice on technical guidance on the classification of waste (2018/C 124/01) of 9 April 2018.

11. Regulation (EU) 2017/746 of the European Parliament and of the Council of 5 April 2017 on in vitro diagnostic medical devices.

LIST OF CHANGES

Previous version: 04	Current version: 05
Sections updated: <i>REAGENTS; WARNING AND NOTES.</i>	

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- Padhi, S., Nayak, A. K., Behera, A. Type II diabetes mellitus: a review on recent drug based therapeutics. Biomedicine & Pharmacotherapy, 131, 110708 (2020).
- Pisarczyk-Wiza D., Zozulińska-Ziółkiewicz D. Glikokortykosteroidy a zaburzenia metabolizmu glukozy. Diabetologia Kliniczna 4,3, 110-116 (2015).
- Barham P., Trinder P. An improved colour reagent for the determination of blood glucose by the oxidase system: Analyst 97, 142-145 (1972).
- Alan H.B. Wu: Tietz Clinical Guide to Laboratory Tests, 4th ed. WB Saunders, 444-450 (2006).
- Rifai N., Horvath A.R., Wittwer C., ed. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics 6th ed., Elsevier, St. Louis, USA, 528 (2018).
- Dujmovic and F. Deisenhammer, Stability of cerebrospinal fluid/serum glucose ratio and cerebrospinal fluid lactate concentrations over 24 h: analysis of repeated measurements, Clinical Chemistry and Laboratory Medicine 48, 2, 209-212 (2010).
- 2021 Guidelines on the management of patients with diabetes. A position of Diabetes Poland
- Sacks, David B., et al. Guidelines and recommendations for laboratory analysis in the diagnosis and management of diabetes mellitus. Clinical Chemistry 48, 3, 436-472 (2002).
- Miles RR, Roberts RF, Putnam AR, Roberts WL. Comparison of serum and heparinized plasma samples for measurement of chemistry analytes. Clinical Chemistry 50, 1704-1706 (2004).
- European Commission notice on technical guidance on the classification of waste (2018/C 124/01) of 9 April 2018.
- Regulation (EU) 2017/746 of the European Parliament and of the Council of 5 April 2017 on in vitro diagnostic medical devices.

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