## CREATININE

Nazwa zestawu	(PL)
Liquick Cor-CREATININE 30	2-232
Liquick Cor-CREATININE 60	2-233
Liquick Cor-CREATININE 120	3-332
HC- CREATININE	4-533
OS- CREATININE	9-408
B50- CREATININE	5-513

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### INTENDED USE

Diagnostic kit for determination of creatinine concentration used both for manual assay and in several automatic analysers. The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

### INTRODUCTION

Creatinine is a product of creatine nonenzymatic dehydration in skeletal muscle. The amount of creatinine generated and excreted by kidney is proportional to muscle mass and usually is higher in men than women. Daily creatinine generation remains fairly constant, with the exception of crushing injury or degenerative diseases that cause massive damage to muscle. Creatinine blood and urine level depends on glomerular filtration so creatinine clearance is excellent index of renal function.

### METHOD PRINCIPLE

Modified Jaffe's method, without deproteinization. In alkaline solution picrate reacts with creatinine to form a yellow-red 2,4,6-trinitrocyclohexadienate. The colour intensity is proportional to the creatinine concentration.

### REAGENTS

Package	Liquick Cor - CREATININE 30	Liquick Cor - CREATININE 60	Liquick Cor - CREATININE 120
1-REAGENT	5 x 24 ml	5 x 48 ml	5 x 96 ml
2-REAGENT	1 x 30 ml	1 x 60 ml	1 x 120 ml

	HC - CREATININE	OS - CREATININE	B50 - CREATININE
1-REAGENT	6 x 76 ml	2 x 56 ml	4 x 58.5 ml
2-REAGENT	6 x 19.5 ml	2 x 18.5 ml	4 x 17 ml

The reagents, stored at 15-25°C are stable up to expiry date printed on the package. The reagents stored on board of the analyser at 2-10°C are stable for 2 weeks.

### Working reagent preparation and stability

Assay can be performed with use of separate 1-REAGENT and 2-REAGENT reagents or with use of working reagent. For working reagent preparation mix gently 4 parts of 1-REAGENT with 1 part of 2-REAGENT. Avoid foaming.

Stability of working reagent: 4 weeks at 2-8°C  
7 days at 15-25°C

Working reagent should be stored in tightly closed vials. If stored in open vial, it is stable for 1 day at 15-25°C!

### Concentrations in the reagent

1-REAGENT sodium hydroxide  $\leq 450$  mmol/l

carbonate buffer  $\leq 150$  mmol/l  
2-REAGENT  
picric acid  $\leq 38.8$  mmol/l

### Warnings and notes

- Protect from direct sunlight and avoid contamination!
- The reagents are usable when the absorbance of the working reagent is less than 0.225 (read against distilled water, wavelength  $\lambda = 500$  nm, cuvette l = 1 cm, at temp. 25°C).
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.
- 1-REAGENT and 2-REAGENT meeting the criteria for classification in accordance with Regulation (EC) No 1272/2008.

### Ingredients:

1-REAGENT contains sodium hydroxide.

2-REAGENT contains picric acid.

### Danger

H314 Causes severe skin burns and eye damage.  
P280 Wear protective gloves, protective clothing, eye protection or face protection.  
P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes.

Remove contact lenses, if present and easy to do. Continue rinsing.

P310 Immediately call a POISON CENTER or doctor.

### ADDITIONAL EQUIPMENT

- automatic analyzer or photometer able to read at 500 nm (492 nm);
- thermostat at 37°C;
- general laboratory equipment;

### SPECIMEN

Serum, EDTA or heparinized plasma free from hemolysis, 24-hours urine, without preservatives<sup>9</sup>.

**Urine preparation:** before analysis urine sample should be diluted 100-fold with 0.9% NaCl, and the results multiplied by 100. Mix well probes before measurement. Specimen can be stored up to 7 days at 2-8°C. For longer storage samples should be frozen at -20°C.

Nevertheless it is recommended to perform the assay with freshly collected samples!

### PROCEDURE

Applications for analyzers are available on request.

### Manual procedure

wavelength 500 nm (492 nm)  
temperature 37°C  
cuvette 1 cm  
reaction type Fixed time  
blank deionized water

### Sample Start method

Pipette into the cuvettes:

	test (T)	standard (S)
working reagent	1000 $\mu$ l	1000 $\mu$ l

Bring up to the temperature of determination. Then add:

calibrator	-	120 $\mu$ l
sample	120 $\mu$ l	-

Mix well and after exactly 70 sec read absorbance A1 of the test (T) and standard (S) against air. After next 40 sec repeat absorbance reading (A2) and calculate  $\Delta A$  (A2 - A1) for the test and standard.

### Calculation

$$\text{creatinine concentration} = \frac{\Delta A(T)}{\Delta A(S)} \times \text{calibrator concentration}$$

### REFERENCE VALUES<sup>7</sup>

serum / plasma	mg/dl	$\mu$ mol/l
female	0.6 - 1.1	53 - 97
male	0.7 - 1.3	62 - 115
24-hours urine	mg/kg/24h	$\mu$ mol/kg/24h
female	11 - 20	97 - 177
male	14 - 26	124 - 230

It is recommended for each laboratory to establish its own reference ranges for local population.

### QUALITY CONTROL

For internal quality control it is recommended to use, with each batch of samples, the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173) for determination in serum or CORMAY URINE CONTROL LEVEL 1 (Cat. No 5-161) or LEVEL 2 (Cat. No 5-162) for determination in urine.

For calibration of the automatic analysers systems or when using the manual method CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174 and 5-176) and LEVEL 2 (Cat. No 5-175 and 5-177) are recommended.

The calibration curve should be prepared every week, with change of reagent lot number or as required e.g. quality control findings outside the specified range.

### PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analyser Biolis 30i. Results may vary if a different instrument or a manual procedure is used.

- LoB (Limit of Blank):** 0.01 mg/dl (0.884  $\mu$ mol/l)
- LoD (Limit of Detection):** 0.04 mg/dl (3.54  $\mu$ mol/l)
- LoQ (Limit of Quantitation):** 0.6 mg/dl (53.04  $\mu$ mol/l)
- Linearity:** up to 17.5 mg/dl (1547  $\mu$ mol/l)

For higher concentration dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

### Specificity / Interferences

Haemoglobin up to 2.5 g/dl, triglycerides up to 500 mg/dl, ascorbate up to 62 mg/l and bilirubin up to 20 mg/dl do not interfere with the test.

### Precision

Repeatability (run to run) n = 20	Mean	SD	CV
	[mg/dl]	[mg/dl]	[%]
level 1	2.03	0.03	1.51
level 2	4.89	0.04	0.89
Reproducibility (day to day) n = 80	Mean	SD	CV
	[mg/dl]	[mg/dl]	[%]
level 1	2.1	0.07	3.6
level 2	5.0	0.14	2.7

### Method comparison

A comparison between creatinine values determined at **Biolis 30i** (y) and at **BECKMAN COULTER AU680** (x) using 62 serum samples gave following results:

$$y = 0.948x + 0.0411 \text{ mg/dl;}$$

$$R = 0.998 \quad (R - \text{correlation coefficient})$$

A comparison between creatinine values determined at **Biolis 30i** (y) and at **BECKMAN COULTER AU680** (x) using 61 urine samples gave following results:

$$y = 0.9441x - 0.0069 \text{ mg/dl;}$$

$$R = 0.999 \quad (R - \text{correlation coefficient})$$

### WASTE MANAGEMENT

Please refer to local legal requirements.

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