## BIL DIRECT

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- Dhungana N., Morris C., Krasowski M.D. Operational impact of using a vanadate oxidase method for direct bilirubin measurements at an academic medical center clinical laboratory. Practical laboratory medicine 8: 77-85 (2017).
- Tietz N.W., ed. Clinical Guide to Laboratory Test. 3rd ed. Philadelphia: WB Saunders 518-525 (1995).
- Kaplan L.A., Pesce A.J., Clinical Chemistry: Theory, analysis, and correlation, 2nd ed., The C.V. Mosby Company, St. Louis, USA, 1013 (1989).
- World Health Organization (2002). Use of Anticoagulants in Diagnostic Laboratory Investigations (Stability of Blood, Plasma and serum samples); WHO/DIL/LAB/99.1 Rev.2: 1-64.
- Zawiadomienie Komisji Europejskiej dotyczące wytycznych technicznych w sprawie klasyfikacji odpadów 2018/C 124/01 z dnia 9 kwietnia 2018 r.
- Rozporządzenie Parlamentu Europejskiego i Rady (UE) 2017/746 z dnia 5 kwietnia 2017 r. w sprawie wyrobów medycznych do diagnostyki in vitro.

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## BIL DIRECT

Kit name	(EN)	Cat. No	stabilizers
Liquick Cor-BIL DIRECT 30		2-247	detergent
Liquick Cor-BIL DIRECT 60		2-248	pH adjuster
Liquick Cor-BIL DIRECT 120		3-329	preservative
HC-BIL DIRECT		4-548	<b>2-BIL DIRECT/ 2-REAGENT</b>
OS-BIL DIRECT		9-406	sodium metavanadate 4.1 mmol/l
B50-BIL DIRECT		5-506	phosphate buffer
			pH adjuster
			stabilizer

## INTENDED USE

BIL DIRECT reagent is intended to determine quantitatively direct bilirubin level in serum. It is intended to diagnosis, monitoring and as an aid to diagnosis of clinical conditions connected with abnormal bilirubin level. BIL DIRECT reagent is intended to use on automatic analysers. It is only for *in vitro* diagnostics, for healthcare professional users.

## SUMMARY<sup>1,2,3</sup>

Determinations of direct bilirubin level in serum are mainly used in diagnosis of liver fuctions. Measurement of direct bilirubin concentration is used as an aid in diagnosis of clinical conditions associated with hyperbilirubinemia. Increased direct bilirubin levels can be observed in mechanical jaundice, Dubin-Jonson syndrome, bile ducts diseases, gallbladder diseases. Direct bilirubin is also the marker of perforation in acute sigmoid diverticulitis.

## METHOD PRINCIPLE<sup>4,5</sup>

Method is based on chemical oxidation, utilizing vanadate as an oxidizing agent.

In the presence of detergent and and vanadate in an acidic solution, conjugated (direct) bilirubin is oxidized to produce biliverdin.

This oxidation reaction causes change of the yellow colour, which is specific to bilirubin to the green colour typical for biliverdin. Therefore, the direct bilirubin concentration in the sample can be obtained by measuring the absorbance before and after the vanadate oxidation.

## REAGENTS

	Liquick Cor-BIL DIRECT 30	Liquick Cor-BIL DIRECT 60	Liquick Cor-BIL DIRECT 120
1-BIL DIRECT	4 x 27 ml	4 x 54 ml	5 x 96 ml
2-BIL DIRECT	1 x 27 ml	1 x 54 ml	1 x 120 ml

	HC-BIL DIRECT	OS-BIL DIRECT	B50-BIL DIRECT
1-REAGENT	6 x 70 ml	6 x 49 ml	4 x 58.5 ml
2-REAGENT	6 x 18 ml	6 x 14 ml	4 x 17.5 ml

## CONCENTRATION OF THE ACTIVE INGREDIENTS

### 1-BIL DIRECT/ 1-REAGENT

citrate buffer  
hydroxylammonium chloride

## REAGENT STABILITY

Reagents stored at 2-8°C are stable until expiry date printed on the package. The reagents are stable for 12 weeks on board the analyser at 2-10°C (Biolis 30i).

## WARNINGS AND NOTES

- Protect from direct sunlight and avoid contamination!
- Do not freeze reagents.
- Do not interchange caps.
- Reagent bottles should be shaken before use by gently inverting several times.
- The appearance of turbidity or control sera values outside the manufacturer's acceptable range may indicate of the reagents instability.
- Lack of significant changes in the color of the reaction mixture at the samples with low bilirubin concentration does not indicate the assay malfunction.
- Do not use the reagent beyond expiry date printed on the package.
- Do not mix reagents from different kits or lots.
- Use personal protective equipment to prevent contact with samples, reagents and controls.
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.
- 1-BIL DIRECT/ 1-REAGENT contains hydroxylammonium chloride. May produce an allergic reaction (EUH208).

## SPECIMEN<sup>6,7,8</sup>

Serum free from hemolysis. Serum should be separated from red blood cells as soon as possible after blood collection.

Lipemic specimens may show falsely decreased bilirubin concentration thus fasting specimen is recommended.

Because bilirubin is photooxidized when exposed to light, specimen should be protected from direct exposure to either artificial light or sunlight. Therefore it is essential to store specimens in the dark! Serum can be stored at 4-8°C or at room temperature (20-25°C) up to 2 days.

Nevertheless it is recommended to perform the assay with freshly collected samples!

Follow tube manufacturers' instructions carefully when using collection tubes .

Human-origin material should be handled as potentially infectious. Standard precautions in normal laboratory work are required.


## ASSAY PROCEDURE

Reagents are ready to use.

Carry-over effect may occur when the reagent is used on automated chemistry analyzers. It is based on interference between specific test combinations and results in either underestimating or inflating patient sample results. To minimize this effect, the preventive actions should be taken, such as ordering an extra washing program or performing tests in a separate order.

### General procedure for automatic determinations:

R1: 240 µl

**Main Wavelength:** 450 nm

**2nd wavelength:** 546 nm

**Method:** End Point

**Calibration Type:** Linear

**Direction:** Decreasing

**R1 Volume:** 240 µl

**Sample Volume:** 10 µl

**R2 Volume:** 60 µl

**Temperature:** 37°C

Application guides for analysers are available on request. The application guide should be verified and validated by the user prior to testing patient samples.

## QUALITY CONTROL

For internal quality control it is recommended to use with each batch of samples, the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173). For the calibration of automatic analysers systems the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) are recommended.

Calibration stability depends on type of analyser used for analysis. The calibration curve should be prepared every 10 weeks (Biolis 30i).

Calibration is recommended in the following cases:

- after each change of lot,
- after instrument service,
- if controls lie outside the expected range,
- each time a new reagent kit is used.

If the quality control results do not fall within the expected values or within the range determined in the laboratory, despite a successful calibration procedure, do not report results. In this case, please take the following actions:

- verify reagents are not out of expiration date.
- verify that the required maintenance has been carried out.
- verify that the procedure has been performed in accordance with the instructions for use.
- contact the Service Department or distributor for assistance.

## REFERENCE VALUES<sup>5</sup>

serum (adults)	0.1 – 0.3 mg/dl 1.7 – 5.1 µmol/l
----------------	-------------------------------------

It is recommended for each laboratory to establish its own reference ranges for local population and method principle. The assay results should be used in conjunction with other data, such as symptoms, results of other tests and clinical history to make clinical decisions. It is not recommended to make clinical diagnosis based on a single result.

## PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analyser Biolis 30i or BS-400. Results may vary if a different instrument is used.

- LoB (Limit of Blank):**  
0.00 mg/dl (0.00 µmol/l) – BS-400
- LoD (Limit of Detection):**  
0.01 mg/dl (0.171 µmol/l) – BS-400
- LoQ (Limit of Quantitation):**  
0.04 mg/dl (0.684 µmol/l) – Biolis 30i
- Linearity:**  
up to 40 mg/dl (684 µmol/l) – Biolis 30i
- Measurement range:**  
0.04 mg/dl - 40 mg/dl (0.684 µmol/l - 684 µmol/l) – Biolis 30i
- Specificity / Interferences**  
Ascorbic acid up to 62 mg/l and triglycerides up to 550 mg/dl (in samples containing low concentration of direct bilirubin) and up to 900 mg/dl (in samples containing high concentration of direct bilirubin) do not interfere with the test. Haemoglobin interferes even in small amount with the determination.

### Precision (Biolis 30i)

Repeatability (run to run) n = 20	Mean [mg/dl]	SD [mg/dl]	CV [%]
level 1	0.35	0.01	2.9
level 2	1.00	0.01	1.0
Reproducibility (day to day) n = 80	Mean [mg/dl]	SD [mg/dl]	CV [%]
level 1	0.36	0.01	3.2
level 2	1.02	0.07	6.7

### Method comparison

A comparison between direct bilirubin values determined at **BIOLIS 30i** (y) and at **BECKMAN COULTER AU680** (x) using 120 serum samples gave following results:

$$y = 1.0564 x - 0.0196 \text{ mg/dl};$$

$$R = 0.999 \quad (R - \text{correlation coefficient})$$

## WASTE MANAGEMENT<sup>9</sup>

After use, the reagents should be handled as potentially infectious and disposed of in accordance with local legal requirements.

- Residual reagents:  
1-BIL DIRECT/ 1-REAGENT - 18 01 07  
2- BIL DIRECT/ 2-REAGENT - 18 01 07

- Empty packages:  
1-BIL DIRECT/ 1-REAGENT - 15 01 02  
2- BIL DIRECT/ 2-REAGENT – 15 01 02
- Wastewater from the analyzer:  
1-BIL DIRECT/ 1-REAGENT - 18 01 03\*  
2- BIL DIRECT/ 2-REAGENT - 18 01 03\*

## INCIDENTS<sup>10</sup>

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer (website address: incidents@cormay.pl) and the competent authority of the Member State in which the user and/or the patient is established.

**Serious incident** means any incident that directly or indirectly led, might have led or might lead to any of the following:

- the death of a patient, user or other person,
- the temporary or permanent serious deterioration of a patient's, user's or other person's state of health,
- a serious public health threat.

## LITERATURE

- Rifai N., Horvath A.R., Wittwer C., ed. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics 6th ed., St. Louis, Missouri, Elsevier, 767-771 (2018).
- Käser S.A., Fankhauser G., Glauser P.M. et al. Diagnostic Value of Inflammation Markers in Predicting Perforation in Acute Sigmoid Diverticulitis. World J Surg 34, 2717–2722 (2010).
- Pagana K, Pagana T, Pagana T. Mosby's Diagnostic & Laboratory Test Reference 14th ed., St. Louis, Missouri, Elsevier, 137-139 (2019).
- Tokuda K., Tanimoto K. A new method of measuring bilirubin in serum by vanadic acid. Jpn J Clin Chem (in Japanese with English abstract) 22: 116–22 (1993).
- Dhungana N., Morris C., Krasowski M.D. Operational impact of using a vanadate oxidase method for direct bilirubin measurements at an academic medical center clinical laboratory. Practical laboratory medicine 8: 77-85 (2017).
- Tietz N.W., ed. Clinical Guide to Laboratory Test. 3rd ed. Philadelphia: WB Saunders 518-525 (1995).
- Kaplan L.A., Pesce A.J., Clinical Chemistry: Theory, analysis, and correlation, 2nd ed., The C.V. Mosby Company, St. Louis, USA, 1013 (1989).
- World Health Organization (2002). Use of Anticoagulants in Diagnostic Laboratory Investigations (Stability of Blood, Plasma and serum samples); WHO/DIL/LAB/99.1 Rev.2: 1-64.
- European Commission notice on technical guidance on the classification of waste (2018/C 124/01) of 9 April 2018.
- Regulation (EU) 2017/746 of the European Parliament and of the Council of 5 April 2017 on in vitro diagnostic medical devices.

## LIST OF CHANGES

Previous version: 03	Current version: 04
Sections updated: <i>INTENDED USE; SUMMARY; CONCENTRATION OF THE ACTIVE INGREDIENTS; WARNING AND NOTES; SPECIMEN; ASSAY PROCEDURE; QUALITY CONTROL; REFERENCE VALUES; PERFORMANCE CHARACTERISTICS; WASTE MANAGEMENT; LITERATURE</i>	
Sections removed: <i>ADDITIONAL EQUIPMENT</i>	
Sections added: <i>REAGENT STABILITY, INCIDENTS</i>	

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4. Tokuda K., Tanimoto K. A new method of measuring bilirubin in serum by vanadic acid. Jpn J Clin Chem (in Japanese with English abstract) 22: 116–22 (1993).
5. Dhungana N., Morris C., Krasowski M.D. Operational impact of using a vanadate oxidase method for direct bilirubin measurements at an academic medical center clinical laboratory. Practical laboratory medicine 8: 77-85 (2017).
6. Tietz N.W., ed. Clinical Guide to Laboratory Test. 3rd ed. Philadelphia: WB Saunders 518-525 (1995).
7. Kaplan L.A., Pesce A.J., Clinical Chemistry: Theory, analysis, and correlation, 2nd ed., The C.V. Mosby Company, St. Louis, USA, 1013 (1989).
8. World Health Organization (2002). Use of Anticoagulants in Diagnostic Laboratory Investigations (Stability of Blood, Plasma and serum samples); WHO/DIL/LAB/99.1 Rev.2: 1-64.
9. European Commission notice on technical guidance on the classification of waste (2018/C 124/01) of 9 April 2018.
10. Regulation (EU) 2017/746 of the European Parliament and of the Council of 5 April 2017 on in vitro diagnostic medical devices.

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