

# Liquick Cor-TG mono

## DIAGNOSTIC KIT FOR DETERMINATION OF TRIGLYCERIDES CONCENTRATION





Kit name	Cat. No
Liquick Cor-TG mono mini	2-281
Liquick Cor-TG mono 30	2-282
Liquick Cor-TG mono 60	2-283
Liquick Cor-TG mono 120	2-284

### INTRODUCTION

Triglycerides are built of glycerol molecule esterified with three fatty acids molecules. Triglycerides are delivered with food or are synthesized endogenously in liver. Triglycerides stored in adipose tissue constitute a reserve of energy. Elevated triglycerides serum level is a risk factor of atherosclerosis. Triglycerides measurement is useful for hyperlipidemia diagnosis and treatment or for estimation of atherosclerosis progression.

### METHOD PRINCIPLE

Colorimetric, enzymatic method with glycerophosphate oxidase.



The colour intensity is proportional to the triglycerides concentration.

### REAGENTS

#### Package

	Liquick Cor-TG mono mini	Liquick Cor-TG mono 30	Liquick Cor-TG mono 60	Liquick Cor-TG mono 120
1-TG mono	2 x 30 ml	6 x 30 ml	6 x 60 ml	6 x 120 ml
2-STANDARD	1 x 1 ml	1 x 2 ml	-	-

2-STANDARD is triglycerides standard solution: 2.5 mmol/l (220 mg/dl).

The reagent when stored at 2-8°C are stable up to expiry date printed on the package. The reagents are stable for 8 weeks on board the analyser at 2-10°C. Protect from light and avoid contamination!

#### Concentrations in the test

buffer TRIS (PH 8.0)	200 mmol/l
4-aminoantipyrine (4-AA)	< 0.4 mmol/l
ATP	< 1.5 mmol/l
Mg <sup>2+</sup>	< 1.6 mmol/l
4-chlorophenol	< 2.5 mmol/l
chlorophenicol	1.6 mmol/l
potassium hexacyanoferrate (II)	< 1 mmol/l
FAD-2Na	< 1 mmol/l
glycerol kinase (GK)	~2500 U/l
glycerol phosphate oxidase (GPO)	~2500 U/l
peroxidase (POD)	~1900 U/l
lipoprotein lipase (LPL)	~2000 U/l
detergents, preservatives	

#### Warnings and notes

- Product for in vitro diagnostic use only.
- The reagents should be used by suitably qualified laboratory personnel only in accordance with intended purpose.
- The reagents contain <0.1% sodium azide as a preservative. Avoid contact with skin and mucous membranes.

### ADDITIONAL EQUIPMENT

- automatic analyser or photometer able to read at 550 nm;
- thermostat at 37°C;
- general laboratory equipment;

### SPECIMEN

Serum, EDTA or heparinized plasma (recommended: heparine lithium, sodium or ammonium salt) free from hemolysis.

Blood should be collected only if the patient has been fasting for minimum of 12 hours. Before blood collection patient should stay in rest position for about 30 minutes. Venous blood is recommended for triglycerides measurement.

Plasma triglycerides values have been reported to be 2% to 4% lower than serum triglycerides values.

Serum and plasma can be stored up to 3 days at 2-8°C or up to 3 months at -20°C. Nevertheless it is recommended to perform the assay with freshly collected samples!

### PROCEDURE

These reagents may be used both for manual assay and in several automatic analysers. Applications for them are available on request.

#### Manual procedure

wavelength	550 nm (546 nm)
temperature	37°C
cuvette	1 cm

Pipette into the cuvettes:

	reagent blank (RB)	test (T)	standard (S)
1-TG mono	1000 µl	1000 µl	1000 µl

Bring up to the temperature of determination. Then add:

standard / calibrator	-	-	10 µl
sample	-	10 µl	-

Mix well, incubate for 5 min. at the temperature of determination. Read the absorbance of standard (S) and test (T) against reagent blank (RB).

#### Calculation

$$\text{triglycerides concentration} = \frac{A(T)}{A(S)} \times \text{standard / calibrator concentration}$$

### REFERENCE VALUES <sup>7</sup>

serum, plasma	< 150 mg/dl < 1.7 mmol/l
---------------	-----------------------------

It is recommended for each laboratory to establish its own reference ranges for local population.

### QUALITY CONTROL

For internal quality control it is recommended to use the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173) with each batch of samples.

For the calibration is also recommended the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176), LEVEL 2 (Cat. No 5-175; 5-177) or TRIGLYCERIDES STANDARD 220 (Cat. No 5-130).

The calibration curve should be prepared every 8 weeks, with change of reagent lot number or as required e.g. quality control findings outside the specified range.

### PERFORMANCE CHARACTERISTICS

These metrological characteristics have been obtained using automatic analyser Prestige 24i. Results may vary if a different instrument or a manual procedure is used.

- **Sensitivity:** 5.3 mg/dl (0.06 mmol/l).
- **Linearity:** up to 1100 mg/dl (12.43 mmol/l).  
For higher triglycerides concentrations dilute the sample with 0.9% NaCl in the ratio of 1 to 4 and repeat the assay. Multiply the result by 5.
- **Specificity / Interferences**  
Haemoglobin up to 0.31 g/dl, bilirubin up to 8.6 mg/dl and ascorbate up to 31 mg/l do not interfere with the test.

▪ **Precision**

Repeatability (run to run) n = 20	Mean [mg/dl]	SD [mg/dl]	CV [%]
level 1	88.86	0.55	0.61
level 2	170.10	1.49	0.87

Reproducibility (day to day) n = 80	Mean [mg/dl]	SD [mg/dl]	CV [%]
level 1	83.51	2.43	2.91
level 2	166.08	4.19	2.52

▪ **Method comparison**

A comparison between triglycerides values determined on Prestige 24i (y) and on COBAS INTEGRA 400 (x) using 96 samples gave following results:

$$y = 0.9785 x + 4.1353 \text{ mg/dl};$$

$$R = 0.9946 \quad (R - \text{correlation coefficient})$$

**TRACEABILITY**

TRIGLYCERIDES STANDARD 220 is traceable to the SRM 1951B reference material.

**WASTE MANAGEMENT**

Please refer to local legal requirements.

**LITERATURE**

1. Jacobs N.J., Van Denmark P.: J. Arch. Biochem. Biophys. 88, 250-255 (1960).
2. Kodischek L.K., Umbreit W.W.: J. Bacteriol. 98, 1063-1068 (1969).
3. Trinder P.: Ann. Clin. Biochem. 6, 24-27 (1969).
4. Schettler G., Nussel E.: Arb. Med. Soz. Med. Prav. Med. 10, 25 (1975).
5. Tietz N.W., ed. Clinical Guide to Laboratory Tests, 3rd ed. Philadelphia, PA: WB Saunders, 610, (1995).
6. Burtis C.A., Ashwood E.R., ed. Tietz Textbook of Clinical Chemistry, 2nd ed. Philadelphia, PA: WB Saunders, 2209, (1994).
7. Alan H.B. Wu. editor. Tietz Clinical Guide to Laboratory Tests, 4<sup>th</sup> ed. St. Louis: W.B Saunders Company; 2006, p.1074.

**Date of issue:** 01. 2018.

**MANUFACTURER**

**PZ CORMAY S.A.**  
ul. Wiosenna 22,  
05-092 Łomianki, POLAND  
tel.: +48 (0) 22 751 79 10  
fax: +48 (0) 22 751 79 14  
<http://www.cormay.pl>

01/18/01/18