Liquick Cor - dTIBC

Nazwa zestawu	(PL) Nr kat.
Liquick Cor-dTIBC 500	2-306
Liquick Cor-dTIBC "bulk"	2-307

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INTENDED USE

Diagnostic kit for determination of total iron binding capacity intended to use in several automatic analyzers.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Iron is transported as Fe (III) bound to the apotransferrin. The apotransferrin-Fe (III) complex is called transferrin. Normally only about one third of the iron binding sites of transferrin are occupied by Fe (III). TIBC (Total iron-binding capacity) is a measurement for the maximum iron concentration that transferrin can bind.


Serum TIBC levels vary in disorders of iron metabolism, e.g. elevated TIBC level when iron levels are low is useful in the diagnosis and monitoring of iron deficiency anaemia and also anaemia in late pregnancy. Decreased level of TIBC may indicate disorders such as anaemia associated with chronic inflammatory disorders, haemochromatosis, other anaemias, thalassemia, cirrhosis and malignancies.

METHOD PRINCIPLE


Colorimetric method, direct, with chromazurol B.

The assay consists of 2 reaction steps:

1. The sample is added to 1-Reagent acidic buffer which contains excess iron and iron-binding dye (Chromazurol B). The low pH causes transferrin to release all iron atoms which then are bound by the iron-binding dye along with the excess iron, producing a coloured dye-iron complex.

2. Neutral buffer 2-Reagent addition results in pH change and causes a large increase in affinity of transferrin for iron. Then serum transferrin rapidly binding iron ions separating it from the dye-iron complex. The observed decrease in the absorbance of a coloured dye-iron complex is directly proportional to the total iron-binding capacity TIBC in the serum sample.

REAGENTS

Package	Liquick Cor-dTIBC 500	Liquick Cor-dTIBC "bulk"
1-dTIBC	3 x 400 ml	--*
2-dTIBC	1 x 300 ml	--*

*reagent volume is printed on the label.

The reagents when stored at 2-8°C are stable up to expiry date printed on the package. The reagents are stable for 11 weeks on board the analyser.

Concentrations in the test

1-Reagent		
sodium acetate		200 mmol/l
iron chloride		0.02 mmol/l
chromazurol B		< 0.25 mmol/l
hexadecyltrimethylammonium bromide (CTAB)		< 3.0 mmol/l
preservatives, conservants		
2-Reagent		
sodium bicarbonate		325 mmol/l
buffer, preservatives, conservants		

Warnings and notes

- Do not freeze the reagents.
- Protect from light and air, avoid contamination!
- Contaminated glassware is the greatest source of error. The use of disposable plastic ware is recommended.
- Because of the great quantities of iron in the environment, scrupulous care is necessary to ensure that glassware, water and reagents do not become contaminated with iron.
- For internal quality control, non-frozen control serum is recommended.
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.

SPECIMEN

Serum free from haemolysis. Presence of fibrin in the serum sample caused by incomplete clot formation causes non-specific interference with TIBC determination.

Do not use plasma samples as anticoagulants affect TIBC test results.

Serum can be stored up to 3 days at 4°C or up to 6 months at -20°C or also at -70°C indefinitely.

Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

Applications for them are available on request.

1-Reagent and 2-Reagent are ready to use.

For reagent blank deionised water is recommended.

REFERENCE VALUES 2

serum	µg/dl	µmol/l
adults	250 – 425	44,8 – 71,6

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173) with each batch of samples.

For the calibration of automatic analysers systems the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) or CORMAY MULTICALIBRATOR LEVEL 2 (Cat. No 5-175; 5-177) is recommended depending on the calibrator lot number.

Calibrator and deionised water should be used for calibration. The calibration curve should be prepared every 11 weeks, with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analyzer BS-400. Results may vary if a different instrument is used.

Sensitivity

100 µg/dl (17.90 µmol/l).

Linearity

up to 600 µg/dl (107.4 µmol/l).

For higher concentration dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

Specificity / Interferences

Haemoglobin up to 1.25 g/dl, ascorbate up to 31 mg/l, bilirubin up to 20 mg/dl, triglycerides up to 1000 mg/dl, copper up to 3 mg/dl, zinc up to 250 µg/dl, iron up to 357,5 µg/dl, Desferal up to 11.5 µg/ml, Cuprimine up to 250 µg/dl and Iron Dextran (Imferon) up to 1430 µg/ml do not interfere with the test.

Precision

Repeatability (run to run) n = 10	Mean [µg/dl]	SD [µg/dl]	CV [%]
level 1	322	1.2	0.4
level 2	223	0.8	0.4
Reproducibility (day to day) n = 10	Mean [µg/dl]	SD [µg/dl]	CV [%]
level 1	357	5.0	1.4
level 2	245	2.4	1.00

Method comparison

A comparison between TIBC values determined at **BS-400** (y) and at **ADVIA 1650** (x) using 149 samples gave following results:

$$y = 0.9628x + 6.373 \mu\text{g/dl};$$

$$R = 0.960 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT

Please refer to local legal requirements.

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