# Vitamin D CONTROL

REF

2130152 - Tri-level Kit: 6 x 2 mL  
2130159 - Level 1 Kit: 6 x 2 mL  
2130160 - Level 2 Kit: 6 x 2 mL  
2130161 - Level 3 Kit: 6 x 2 mL

IVD

EN EXPLANATION OF SYMBOLS  
CS VÝZNAM SYMBOLŮ  
DA SYMBOLFORKLARING  
DE ERKLÄRUNG DER SYMBOLE  
EL ΕΠΕΞΗΓΗΣΗ ΤΩΝ ΣΥΜΒΟΛΩΝ  
ES SIGNIFICADO DE LOS SÍMBOLOS  
ET SÜMBOLITE SELGITUS  
FR EXPLICATION DES SYMBOLES  
HU JELMAGYARÁZAT  
IT SPIEGAZIONE DEI SIMBOLI  
JA 記号の説明  
LV SIMBOLU SKAIDROJUMS  
LT ŽYMĖJIMŲ PAAIŠKINIMAS  
NL VERKLARING DER SYMBOLEN  
PL OBJAŚNIENIE SYMBOLI  
PT EXPLICAÇÃO DOS SÍMBOLOS  
RO EXPLICAREA SIMBOLURILOR  
SK VÝZNAM SYMBOLOV  
SL RAZLAGA SIMBOLOV  
SV SYMBOLFÖRKLARING  
ZH 符号注解

LOT

Batch code/Číslo šarže/  
Lotnummer/  
Chargenbezeichnung/Αριθμός  
Παρτίδας/Código de lote/Partii  
kood/Code du lot/Sarzzszám/  
Codice del lotto/バッチ コード/  
Partijas kods/Partijos kods/Lot  
number/Kod partii/Código do  
lote/Lot nr./Číslo šarže/Številka  
serije/Satsnummer/批号

Date of manufacture/Datum  
výroby/Produktionsdato/  
Herstellungsdatum/  
Ημερομηνία παραγωγής/Fecha  
de fabricación/Tootmiskuupäev/  
Date de fabrication/Gyártás  
időpontja/Data di fabbricazione/  
製造日/Izgatavošanas datums/  
Pagaminimo data/Productiedatum/  
Data produkcyj/Data de fabrico/  
Data de fabricație/Dátum  
výroby/Datum izdelave/  
Tillverkningsdatum/生产日期

Use By/Použitelné do/  
Holdbar til/Verwendbar bis/  
Ημερομηνία λήξης/Fecha de  
caducidad/Kasutada enne/Utiliser  
jusque/Felhasználható/Utilizaze  
entro/使用期限/Izlietot līdz/  
Sunaudoti iki/Houdbaar tot/Użytyć  
przed/Prazo de validade/Data  
expirării/Použiteneľ do/Uporabno  
do/Använd före/此日期前使用

**REF**

Catalogue number/Katalogové číslo/  
Katalognummer/Bestellnummer/  
Αριθμός καταλόγου/Número de  
catálogo/Katalogi number/Référence  
du catalogue/Katalógusszám/Numero  
di catalogo/カタログ番号/Kataloga  
numurs/Katalogo numeris/Catalogus  
nummer/Numer katalogowy/Referência  
de catálogo/Număr de catalog/  
Katalógové číslo/Kataloška številka/  
Katalognummer/目録番号

Manufacturer/Výrobce/Producent/  
Hersteller/Κατασκευαστής/  
Fabricante/Tootja/Fabricant/Gyártó/  
Fabbicante/製造業者/Izgatavotājs/  
Gamintojas/Fabrikant/Producent/  
Fabricante/Producător/Výrobca/  
Izdelovalec/Tilverkare/生产商

**IVD**

In Vitro Diagnostic Medical Device/  
In Vitro diagnostický zdravotnícký  
prostředek/Medicinsk udstyr  
til in vitro-diagnostik/In-vitro-  
Diagnostikum/In Vitro Διαγνωστικό  
Ιατροτεχνολογικό προϊόν/Producto  
sanitario para diagnóstico in vitro/  
In vitro diagnostiline meditsiiniseade/  
Dispositif médical de diagnostic in  
vitro/In vitro diagnosztikum/Dispositivo  
medico-diagnostico in vitro/体外診断  
医療機器/Medicinas ierice, ko lieto in  
vitro diagnostikā/In vitro diagnostikos  
medicinos prietaisais/Medisch  
hulpmiddel voor in-vitro diagnostiek/  
Wyrób do diagnostyki In Vitro/  
Dispositivo médico para diagnóstico in  
vitro/Dispozitiv medical de diagnostic  
in-vitro/Zdravotnícka pomôcka in  
vitro/In vitro diagnostični pripomoček/  
Medicintekniska produktier för in vitro  
diagnostik/体外诊断医疗器械

Temperature limitation/  
Templotni rozmezi od do/  
Temperaturbegrænsning/  
Temperaturbegrenzung/  
Περιορισμοί θερμοκρασίας/  
Limite de temperatura/  
Temperatuuriirangud/  
Limites de température/  
Hömerséklettartomány/  
Limiti di temperatura/温度制限/  
Temperatūras ierobežojumi/  
Temperatūras apribojimas/  
Temperatuurlimiet/  
Przestrzegać zakresu temperatury/  
Limites de temperatura/  
Limite de temperatură/  
Teplotné rozmedzie od do/  
Omejitve temperature/  
Temperaturbegrænsning/  
温度限制

**EC REP**

Authorized representative/Pověřený  
zástupce/Autoriseret representant/  
Bevollmächtigter Vertreter/  
Εξουσιοδοτημένος αντιπρόσωπος/  
Representante autorizado/  
Volitatud esindaja/Márkakereskedő/  
Rappresentante autorizzato/  
委任された代表/Pilnvarots pārstāvis/  
Igalilotsis atstovas/Geautoriseerd  
vertegenwoordiger/Autoryzowany  
przedstawiciel/Representante autorizado/  
Representant autorizat/Autorizovaný  
zástupca/Пооблашчени представник/  
Auktoriserad representant/授权代表

Consult Instructions for Use/Viz  
návod k použití/Se brugsanvisning/  
Gebrauchsanweisung beachten/  
Συμβουλευτείτε τις οδηγίες χρήσης/  
Consulte las instrucciones de uso/  
Lugege kasutusjuhendit/Consulter les  
instructions d'utilisation/ Nézze meg  
a Használati utasítást/Consultare  
le istruzioni per l'uso/使用説明書  
を参照/lepazīstieties ar lietošanas  
instrukciju/Žr. naudojimo instrukcija/  
Raadpleeg de gebruiksaanwijzing/  
Sprawdź w instrukcji obsługi/Consulte  
as instruções de utilização/Consultați  
instrucțiunile de utilizare/Vid' návod  
na použitie/Oglejte si navodila za  
uporabo/Se brugsanvisningen/  
请参考“使用说明”

Biological risks/Biologicky nebezpečné/  
Biologisk fare/Biogefährdung/Biολογικοί  
κίνδυνοι/Riesgo biológico/Bioloogiline  
oht/Risques biologiques/Biológiai  
kockázat/Rischio biologico/  
生物学的 リスク/Biologiskie riski/  
Biologinis pavojus/Biologiskie risko/  
Ryzyko biologiczne/Risco biológico/  
Risc biologic/Biološka tveganja/Biološka  
tveganja/Biologisk risk/生物风险

**CONT**

Contents of kit/Obsah sady/Kittets  
indhold/Inhalt des Kits/Περιεχόμενα  
του κιτ/Contenido del kit/Komplekti  
sisu/Contenu du kit/A készlet  
tartalma/Contenido del kit/キットの  
内容/Komplekta saturs/Komplekto  
turinys/Inhoud van de set/Zawartość  
zestawu/ Conteúdo do kit/Conținutul  
kitului/Obsah súpravy/Vsebina  
kompleta/Satsinnehäll/试剂盒内容

Reconstitute with/Rozedzte pomocí/  
Rekonstituertes med/Rekonstituieren mit/  
Αναστάση με/Reconstituir con/Lahusti  
Reconstituer avec/Feloldáshoz/Ricostituire  
con/次により調製する/Atškaidīt ar/Ištirpinti  
naudojant/Reconstituie met/Odworzyć za  
pomocą/Reconstituire com/Reconstituire cu/  
Rozriedzte pomocou/Rekonstituiraite z/s/  
Rekonstituera med/复溶剂

EN

CS

DA

DE

EL

ES

ET

FR

HU

IT

JA

LT

LV

NL

PL

PT

RO

SK

SL

SV

ZH

## INTENDED USE

For In Vitro Diagnostic Use Only.

Fujirebio Diagnostics Vitamin D Control is intended for use as an assayed control serum to monitor the precision of laboratory testing procedures for the analysis of Vitamin D.

## SUMMARY AND PRINCIPLE

This quality control product can be used as an objective judgment of the laboratory's procedures and personnel techniques. It is a valuable tool to assess good laboratory practices. Three levels of control are available to compare observations with expected ranges therefore assuring consistent performance of the testing system within the clinical range.

## REAGENTS

This product is prepared from human serum, protein (bovine), purified biochemical materials, and chemicals. Contains Proclin 300® / Gentamicin as preservatives. Level 1 of this control contains 25(OH) Vitamin D<sub>3</sub> only. Level 2 and Level 3 contain both 25(OH) Vitamin D<sub>2</sub> and 25(OH) Vitamin D<sub>3</sub>. Each level of control contains clinically relevant concentrations of Vitamin D (1). The control is provided in lyophilized form.

## STORAGE AND STABILITY

**Unopened:** This product is stable until the expiration date when stored un-opened at 2-8°C.

**Reconstituted:** Once the control is reconstituted, the product is stable for 14 days when stored tightly capped at 2-8°C. The product is stable for 60 days when stored at ≤ -10°C. Controls may be frozen and thawed repeatedly for up to 9 cycles. Do not store at ambient temperature.

## RECONSTITUTION

Using a calibrated pipette, reconstitute each vial with 2.0 mL of distilled or deionized water. Replace the stopper and lightly tap the vial 5-6 times to loosen the pellet and allow the controls to stand at room temperature, 15-30°C, for at least 30 minutes. Invert three times and swirl before use to ensure homogeneity.

## PROCEDURE

Follow all instructions provided with the assay or instrument of interest while running the control. The control is to be run as you would a patient sample. Dispose of any discarded materials in accordance with the requirements of your local waste management authorities. Quality control materials should be used in accordance with local, state, and/or federal regulations or accreditation requirements. In the event of damage to packaging, contact your local distributor or Fujirebio Diagnostics at the numbers listed on the last page of this package insert.

## WARNING

Human source material. Treat as potentially infectious. Each human donor unit used to manufacture this control was tested by FDA accepted methods and found nonreactive for Hepatitis B Surface Antigen (HBsAg), antibody to Hepatitis C (HCV), antibody to HIV-1/HIV-2 and antibody to Treponema Pallidum (Syphilis). This product may also contain other human pathogens for which there are no approved tests. All human source material should be considered potentially infectious and handled with the same precautions used with patient specimens.

Never mouth pipette. To avoid cross-contaminations, ensure that the stoppers and caps are not switched during use.

## LIMITATIONS

1. Do not use this product beyond the expiration date.
2. Discard the vial if there are indications of microbial contamination or cloudiness in the reconstituted control.
3. This product is not to be used as a standard.

## ASSIGNMENT OF VALUES

The assigned ranges shown have been established for the assays listed in the Assigned Values Sheet only. The ranges listed resulted from replicate analyses and are specific for each lot of product. The tests listed were performed by the reagent manufacturer and/or independent laboratories using manufacturer supported reagents (available at the time of assay) and a sampling of this lot of control. Your laboratory means should fall within the printed acceptable range; however, laboratory means may vary from the listed values during the life of this control. Variations between laboratories and variations over time may be caused by differences in laboratory technique, instrumentation and reagents, or by modifications made to manufacturer's test methods. It is recommended that each laboratory establish its own means and acceptable ranges and use the printed ranges as guidance only.

To convert ng/mL to nmol/L use a conversion factor of 2.496 (1).

## SPECIFIC PERFORMANCE CHARACTERISTICS

This is a lyophilized product manufactured under strict quality standards. This control must be stored and handled as directed to obtain expected values.

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**Fujirebio Diagnostics, Inc.**
**Fujirebio Diagnostics, AB**