ACCENT-200 ACID WASHING SOLUTION

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Data wydania: 03. 2022.

ACCENT-200 ACID WASHING SOLUTION

Reagent name	Variant	Cat.No
ACCENT-200 ACID WASHING SOLUTION	3 x 35 ml	3-109

INTENDED USE

ACCENT-200 ACID WASHING SOLUTION is intended to use for cleaning the reagent and sample probes of the biochemistry analyzers: ACCENT-200 (II GEN) / BS-200, ACCENT-220S / BS-180 and BS-120 during routine analyzer operation and to eliminate carryover effect. It is only for *in vitro* diagnostics, for healthcare professional users.

OPERATION PRINCIPLE

ACCENT-200 ACID WASHING SOLUTION is a washing solution placed on board the analyzer. Removal of biological material and reagents residues, especially protein-containing materials and microbial deposits from the reagent and sample probes takes place due to the cleaning properties of the reagent. ACCENT-200 ACID WASHING SOLUTION contains hydrochloric acid that dissolves and removes contaminations.

CONCENTRATION OF THE ACTIVE INGREDIENTS IN THE REAGENT

hydrochloric acid 278 mmol/l
preservatives

REAGENT STABILITY

The reagent must be used before the expiry date indicated on the label.

The reagent stored on board of the analyzer at 2-10°C is stable for 30 days.

WARNING AND NOTES

- Protect from direct sunlight!
- Do not use the reagent beyond expiry date printed on the label.
- Do not combine the remaining reagent with the reagent from a new package.
- Use personal protective equipment to prevent contact with samples, reagents and controls.
- Avoid discharge to sewage system or to environment.
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.
- The product meets the criteria for classification in accordance with Regulation (EC) No 1272/2008.

Ingredients:

ACCENT-200 ACID WASHING SOLUTION contains hydrochloric acid.

H314 Causes severe skin burns and eye damage.
H318 Causes serious eye damage.
P280 Wear protective gloves, protective clothing, eye protection or face protection.

P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P310 Immediately call a POISON CENTER or doctor.

ACCENT-200 ACID WASHING SOLUTION

ASSAY PROCEDURE

ACCENT-200 ACID WASHING SOLUTION is ready to use. Use according to the analyzer's operation manual and the guidelines contained in the instruction: 51_03_24_001_ACCENT-200_CARRYOVER.

WASTE MANAGEMENT¹

After use, the reagents should be handled as potentially infectious and disposed of in accordance with local legal requirements

- Reagents residues: 18 01 06*
- Empty packages: 15 01 10*
- Wastewater from the analyzer: 18 01 03*

INCIDENTS²

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer (website address: incidents@cormay.pl) and the competent authority of the Member State in which the user and/or the patient is established.

Serious incident means any incident that directly or indirectly led, might have led or might lead to any of the following:

- the death of a patient, user or other person,
- the temporary or permanent serious deterioration of a patient's, user's or other person's state of health,
- a serious public health threat.

LITERATURE

- European Commission notice on technical guidance on the classification of waste (2018/C 124/01) of 9 April 2018.
- Regulation (EU) 2017/746 of the European Parliament and of the Council of 5 April 2017 on *in vitro* diagnostic medical devices.

LIST OF CHANGES

Previous version: 02	Current version: 03
<i>Sections updated: INTENDED USE, OPERATION PRINCIPLE, CONCENTRATION OF THE ACTIVE INGREDIENTS IN THE REAGENT, WARNING AND NOTES, WASTE MANAGEMENT</i>	
<i>Sections added: REAGENT STABILITY, INCIDENTS, LITERATURE</i>	

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