A-400 ACID WASHING SOLUTION

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A-400 ACID WASHING SOLUTION

Reagent name	Variant	Cat.No
A-400 ACID WASHING SOLUTION	7 x 40 ml	3-111

INTENDED USE

A-400 ACID WASHING SOLUTION is intended to use for cleaning the reagent and sample probes of the biochemistry analyzers BS-400 and BS-480 during routine analyzer operation and to eliminate carryover effect. It is only for *in vitro* diagnostics, for healthcare professional users.

OPERATION PRINCIPLE

A-400 ACID WASHING SOLUTION is a washing solution placed on board the analyzer. Removal of biological material and reagents residues, especially protein-containing materials and microbial deposits from the reagent and sample probes takes place due to the cleaning properties of the reagent. A-400 ACID WASHING SOLUTION contains hydrochloric acid that dissolves and removes contaminations.

CONCENTRATION OF THE ACTIVE INGREDIENTS IN THE REAGENT

hydrochloric acid 278 mmol/l
preservatives

REAGENT STABILITY

The reagent must be used before the expiry date indicated on the label.

The reagent stored on board of the analyzer at 2-10°C is stable for 30 days.

WARNING AND NOTES

- Protect from direct sunlight!
- Do not use the reagent beyond expiry date printed on the label.
- Do not combine the remaining reagent with the reagent from a new package.
- Use personal protective equipment to prevent contact with samples, reagents and controls.
- Avoid discharge to sewage system or to environment.
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.
- The product meets the criteria for classification in accordance with Regulation (EC) No 1272/2008.

Ingredients:

A-400 ACID WASHING SOLUTION contains hydrochloric acid.

Danger

- H314 Causes severe skin burns and eye damage.
- H318 Causes serious eye damage.
- P280 Wear protective gloves, protective clothing, eye protection or face protection.

P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P310 Immediately call a POISON CENTER or doctor.

ASSAY PROCEDURE

A-400 ACID WASHING SOLUTION is ready to use.

Use according to the analyzer's operation manual and the guidelines contained in the instructions:

- 51_03_24_002_BS-400_CARRYOVER
- 51_03_24_003_BS-480_CARRYOVER

WASTE MANAGEMENT¹

After use, the reagents should be handled as potentially infectious and disposed of in accordance with local legal requirements

- Reagents residues: 18 01 06*
- Empty packages: 15 01 10*
- Wastewater from the analyzer: 18 01 03*

INCIDENTS²

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer (website address: incidents@cormay.pl) and the competent authority of the Member State in which the user and/or the patient is established.

Serious incident means any incident that directly or indirectly led, might have led or might lead to any of the following:

- the death of a patient, user or other person,
- the temporary or permanent serious deterioration of a patient's, user's or other person's state of health,
- a serious public health threat.

LITERATURE

- European Commission notice on technical guidance on the classification of waste (2018/C 124/01) of 9 April 2018.
- Regulation (EU) 2017/746 of the European Parliament and of the Council of 5 April 2017 on *in vitro* diagnostic medical devices.

LIST OF CHANGES

Previous version: 03	Current version: 04
<i>Sections updated: INTENDED USE, OPERATION PRINCIPLE, CONCENTRATION OF THE ACTIVE INGREDIENTS IN THE REAGENT, WARNING AND NOTES, WASTE MANAGEMENT</i>	
<i>Sections added: REAGENT STABILITY, INCIDENTS, LITERATURE</i>	

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