## FERRUM

Nazwa zestawu	(PL)	Nr kat.
Liquick Cor-FERRUM 30	3-257	
Liquick Cor-FERRUM 60	3-258	
Liquick Cor-FERRUM 120	3-334	
HC-FERRUM	4-558	
OS-FERRUM	9-416	
B50-FERRUM	5-538	

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### INTENDED USE

Diagnostic kit for determination of iron concentration used both for manual assay and in several automatic analysers. The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

### INTRODUCTION

Iron is the most abundant trace element in the organism. Most of the iron in humans is located within heme molecule which is incorporated into hemoglobin, myoglobin, catalase, peroxidase and cytochromes. Iron is stored bound to ferritin or hemosiderin and is transported by transferrin. Measurement of iron level is valuable especially in diagnosis of different types of anemia.

### METHOD PRINCIPLE

Colorimetric method with ferrozine, without deproteinization. Iron ions ( $Fe^{3+}$ ), bounded in blood to transferrin are released in acid solution and in the presence of detergents and reduced to  $Fe^{2+}$  by ascorbate.  $Fe^{2+}$  forms with 3-(2-pyridyl)-5,6-bis(2-[4-phenyl sulfonic acid])-1,2,4-triazine sodium salt (ferrozine) coloured complex,  $Cu^{2+}$  ions are bound by thiourea. The colour intensity is directly related to the iron concentration.

### REAGENTS

#### Package

	Liquick Cor-FERRUM 30	Liquick Cor-FERRUM 60	Liquick Cor-FERRUM 120
1-REAGENT	5 x 25 ml	5 x 50 ml	5x 100 ml
2-REAGENT	1 x 25 ml	1 x 50 ml	1 x 100 ml
3-STANDARD	1 x 2 ml		

	HC-FERRUM	OS-FERRUM	B50-FERRUM
1-REAGENT	6 x 88 ml	2 x 53 ml	2 x 58 ml
2-REAGENT	6 x 18.5 ml	2 x 13 ml	2 x 14 ml

3-STANDARD is a standard solution of iron ions with a concentration within the 50.4 – 61.6  $\mu\text{g/dl}$  (9.02 – 11.02  $\mu\text{mol/l}$ ). The exact concentration is printed on the label of the each vial.

The reagents are stable up to expiry date printed on the package, if stored at 2-8°C. The reagents are stable for 11 weeks on board the analyser at 2-10°C.

### Concentrations in the reagent

1-REAGENT	
citric acid (pH 1.9)	$\leq 240 \text{ mmol/l}$
thiourea	$\leq 108 \text{ mmol/l}$
detergent	$\leq 7 \%$

FERRUM (II GENERACJA / II GENERATION / II ПОКОЛЕНИЕ)

### 2-REAGENT

sodium ascorbate  $\leq 150 \text{ mmol/l}$   
3-(2-pyridyl)-5,6-bis(2-[5-furyl sulfonic acid])-1,2,4-triazine sodium salt (ferrozine)  $\leq 6 \text{ mmol/l}$   
preservative  
stabilizer

### WARNINGS AND NOTES

- Protect from direct sunlight and avoid contamination!
- When using the manual method, use disposable cuvettes.
- 1-REAGENT and 3-STANDARD meeting the criteria for classification in accordance with Regulation (EC) No 1272/2008.

#### Ingredients:

1-REAGENT contains citric acid.  
3-STANDARD contains hydrochloric acid.

#### Danger

H314 Causes severe skin burns and eye damage.  
P280 Wear protective gloves, protective clothing, eye protection or face protection.  
P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P310 Immediately call a POISON CENTER or doctor.

### ADDITIONAL EQUIPMENT

- automatic analyzer or photometer able to read at 550 nm;
- thermostat at 37°C;
- general laboratory equipment.

### SPECIMEN

Serum free from hemolysis, collected in plastic tubes. Serum should be separated from red blood cells as soon as possible after blood collection. Serum can be stored up to 7 days at 15-25°C or up to 3 weeks at 2-8°C. Nevertheless it is recommended to perform the assay with freshly collected samples!

### PROCEDURE

Applications for analyzers are available on request.

#### Manual procedure

wavelength 550 nm  
temperature 37°C  
cuvette 1 cm

#### Reagent Start method

Pipette into the cuvettes:

	blank (B)	test (T)	standard (S)
1-REAGENT	1000 $\mu\text{l}$	1000 $\mu\text{l}$	1000 $\mu\text{l}$

Bring up to the temperature of determination. Then add:

standard	-	-	80 $\mu\text{l}$
serum	-	80 $\mu\text{l}$	-
distilled water	80 $\mu\text{l}$	-	-

Mix well and after 5 minutes of incubation read the absorbance A1 of test (T) and standard (S) against blank (B). Then add:

2-REAGENT	200 $\mu\text{l}$	200 $\mu\text{l}$	200 $\mu\text{l}$
-----------	-------------------	-------------------	-------------------

Mix well and after 10 minutes of incubation read the absorbance A2 of test (T) and standard (S) against blank (B). The intensity of colour is stable for 30 minutes.

Calculate  $\Delta A$  for the test and standard:

$$\Delta A = (A_2 - 0.84A_1)$$

The 0.84 coefficient compensates the decrease of absorbance after 2-REAGENT addition.

### Calculation

$$\text{iron concentration} = \frac{\Delta A(T)}{\Delta A(S)} \times \text{standard concentration}$$

### REFERENCE VALUES

serum	$\mu\text{g/dl}$	$\mu\text{mol/l}$
newborns	100 – 250	17.9 – 44.8
infants	40 – 100	7.2 – 17.9
children	50 – 120	9.0 – 21.5
adult females	50 – 170	9.0 – 30.4
adult males	65 – 175	11.6 – 31.3

Samples should be taken in the morning from patients in a fasting state, since iron values decrease by 30% during the course of the day.

It is recommended for each laboratory to establish its own reference ranges for local population.

### QUALITY CONTROL

For internal quality control it is recommended to use the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173) with each batch of samples. For calibration when using the manual method IRON STANDARD 56 (Cat. No 5-133) or 3-STANDARD attached to the kit is recommended. For calibration of the automatic analysers systems CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) or CORMAY MULTICALIBRATOR LEVEL 2 (Cat. No 5-175; 5-177) is recommended depending on the calibrator lot number. The calibration curve should be prepared every 11 weeks, with change of reagent lot number or as required e.g. quality control findings outside the specified range.

### PERFORMANCE CHARACTERISTICS

These metrological characteristics have been obtained using automatic analyser Biolis 24i Premium. Results may vary if a different instrument or a manual procedure is used.

- Sensitivity:** 3.6  $\mu\text{g/dl}$  (0.644  $\mu\text{mol/l}$ ).
- Linearity:** up to 1000  $\mu\text{g/dl}$  (179  $\mu\text{mol/l}$ ).

For higher concentrations dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

#### Specificity / Interferences

Ascorbate up to 62 mg/l, bilirubin up to 20 mg/dl, triglycerides up to 1000 mg/dl and copper up to 500  $\mu\text{g/dl}$  do not interfere with the test. Haemoglobin interferes even in small amount with the determination.

#### Precision

Repeatability (run to run) n = 10	Mean [ $\mu\text{g/dl}$ ]	SD [ $\mu\text{g/dl}$ ]	CV [%]
level 1	33.86	0.47	1.39
level 2	317.54	1.76	0.55
Reproducibility (day to day) n = 10	Mean [ $\mu\text{g/dl}$ ]	SD [ $\mu\text{g/dl}$ ]	CV [%]
level 1	243.77	1.97	0.81
level 2	65.91	1.39	2.10

#### Method comparison

A comparison between ferrum values determined at Biolis 24i Premium (y) and at Hitachi 912 (x) using 102 samples gave following results:

$$y = 0.9325x + 7.8482 \mu\text{g/dl};$$

$$R = 0.9925 \quad (R - \text{correlation coefficient})$$

### TRACEABILITY

IRON STANDARD 56 is traceable to the direct spectrophotometric reference method.

### WASTE MANAGEMENT

Please refer to local legal requirements.

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