## ALBUMIN

Nazwa zestawu	(PL) Nr kat.
Liquick Cor-ALBUMIN 30	2-239
Liquick Cor-ALBUMIN 60	2-238
Liquick Cor-ALBUMIN 120	3-326
HC-ALBUMIN	4-538
OS-ALBUMIN	9-403
B50-ALBUMIN	5-500

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### INTENDED USE

Diagnostic kit for determination of albumin concentration, used both for manual assay and in several automatic analysers. The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

### INTRODUCTION

Albumin is the major serum protein, but is present also in other body fluids: cerebrospinal, pleural and peritoneal. Albumin regulates blood oncotic pressure and serves as amino acids reservoir. Beyond of these functions albumin is very important transport protein – binds and keeps dispersed bilirubin, hormones, vitamins, calcium, magnesium, fatty acids and medicines. Decreased albumin blood level is caused usually by liver or kidney disease, malabsorption or malnutrition.

### METHOD PRINCIPLE

Bromocresol green (BCG) forms with albumin, in succinate buffer (acid medium), a coloured complex. The absorbance of this complex is proportional to the albumin concentration in the sample. The colour intensity of the formed complex measured at 630 nm is proportional to albumin concentration in the sample.

### REAGENTS

#### Package

	Liquick Cor- ALBUMIN 30	Liquick Cor- ALBUMIN 60	Liquick Cor- ALBUMIN 120
1-REAGENT	6 x 30 ml	6 x 60 ml	6 x 100 ml
2-STANDARD	1 x 2 ml	1 x 2 ml	-

	HC- ALBUMIN	OS- ALBUMIN	B50- ALBUMIN
1-REAGENT	6 x 100 ml	2 x 53 ml	2 x 58 ml

2-STANDARD is a standard solution of albumin with a concentration within the range 4.50 - 5.50 g/dl. The exact concentration is printed on the label of the each vial.

The reagent is stable up to the kit expiry date printed on the package when stored at 2-8°C. The reagents are stable for 8 weeks on board the analyser at 2-10°C.

Concentrations in the test	
succinate buffer	≤ 108 mmol/l
bromocresol green (BCG)	≤ 0.35 mmol/l
sodium hydroxide	≤ 60 mmol/l
detergent	
preservative	

### WARNINGS AND NOTES

- Protect from direct sunlight and avoid contamination!
- Do not freeze the reagent!
- 2-STANDARD has been tested for the HIV antibody, HBsAg and HCV and found to be non-reactive. However this material should be handled as though capable of transmitting infectious disease.
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.
- 1-REAGENT meets the criteria for classification in accordance with Regulation (EC) No 1272/2008.

### Warning

H319 Causes serious eye irritation.  
P280 Wear protective gloves, protective clothing, eye protection or face protection.  
P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

### ADDITIONAL EQUIPMENT

- automatic analyzer or photometer able to read at 630 nm;
- thermostat at 37°C;
- general laboratory equipment;

### SPECIMEN

Serum free from hemolysis.  
Serum should be separated from red blood cells as soon as possible after blood collection.  
Serum can be stored up to 3 days at 2-8°C or 6 months at -20°C.<sup>2</sup>  
Nevertheless it is recommended to perform the assay with freshly collected samples!

### PROCEDURE

Applications for analysers are available on request.  
The reagent is ready to use.

### Manual procedure

wavelength	630 nm
temperature	37°C
cuvette	1 cm

Pipette into the cuvette:

	blank (B)	test (T)	standard (S)
1-REAGENT	1000 µl	1000 µl	1000 µl

Bring up to the temperature of determination. Then add:

standard	-	-	10 µl
sample	-	10 µl	-
distilled water	10 µl	-	-

Mix well and incubate for 5 minutes. Read the absorbance of test sample A(T) and standard sample A(S) against reagent blank (B).

### Calculation

$$\text{albumin concentration} = \frac{A(T)}{A(S)} \times \text{standard concentration}$$

### REFERENCE VALUES<sup>7</sup>

serum		g/dl	g/l
children	3 – 5 years	3.9 – 5.0	39 - 50
	6 – 15 years	4.1 – 5.1	41 - 51
females	16 – 54 years	3.9 – 5.0	39 – 50
males	16 – 29 years	4.6 – 5.3	46 - 53
	30 – 54 years	4.4 – 5.1	44 - 51
females and males	55 – 79 years	4.2 – 5.0	42 – 50

It is recommended for each laboratory to establish its own reference ranges for local population.

### QUALITY CONTROL

For internal quality control it is recommended to use the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173) with each batch of samples.

For calibration when using the manual method ALBUMIN STANDARD (Cat. No 5-115) is recommended.

For calibration of the automatic analysers systems CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174 and 5-176) and LEVEL 2 (Cat. No 5-175 and 5-177) are recommended.

The calibration curve should be prepared every 7 weeks, with change of reagent lot number or as required e.g. quality control findings outside the specified range.

### PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analyser Biolis 24i Premium. Results may vary if a different instrument or a manual procedure is used.

- Sensitivity:** 1.14 g/dl (11.4 g/l).

- Linearity:** up to 6.5 g/dl (65.0 g/l).

For higher concentration, dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

- Specificity / Interferences**

Haemoglobin up to 2.5 g/dl, ascorbate up to 62 mg/l, bilirubin up to 20 mg/dl and triglycerides up to 1200 mg/dl do not interfere with the test.

- Precision**

Repeatability (run to run) n = 20	Mean [g/dl]	SD [g/dl]	CV [%]
level 1	2.88	0.03	0.96
level 2	4.53	0.03	0.68

Reproducibility (day to day) n = 80	Mean [g/dl]	SD [g/dl]	CV [%]
level 1	4.37	0.04	0.94
level 2	2.86	0.02	0.87

- Method comparison**

A comparison between albumin values determined at Biolis 24i Premium (y) and at ADVIA 1650 (x) using 21 samples gave following results:

$$y = 1.0372 x - 0.0167 \text{ g/dl};$$

$$R = 0.997 \quad (R - \text{correlation coefficient})$$

### TRACEABILITY

ALBUMIN STANDARD is traceable to the ERM-DA470 reference material.

### WASTE MANAGEMENT

Please refer to local legal requirements.

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