PRESTIGE 24i ASO

Nr kat. 4-270, 4-489 (PL)

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Cat. No 4-270, 4-489 (EN)

INTENDED USE

Diagnostic kit for determination of anti-streptolysin O levels, used in automatic analysers Prestige 24i, Biolis 24i and Prestige 24i Premium, Biolis 24i Premium and Biolis 30i. The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Most people infected with hemolytic streptococcus produce anti-streptolysin O (ASO), antibodies against streptolysin O (SLO), an exotoxin of Streptococcus. Measuring the level of ASO is effective for diagnosing, judging the progress of medical treatment, and assessing recovery from diseases caused by hemolytic streptococcus such as rheumatic fever, acute glomerulonephritis, scarlatina and tonsillitis.

METHOD PRINCIPLE

When an antigen-antibody reaction occurs between ASO in a sample and SLO which has been sensitized to latex particles, agglutination results. This agglutination is detected as an absorbance change (572 nm), with the magnitude of the change being proportional to the quantity of ASO in the sample. The actual concentration is then determined by interpolation from a calibration curve prepared from calibrators of known concentration.

REAGENTS

Package

	Cat. No 4-270 (24-TRAY)	Cat. No 4-489 (36-TRAY)
1-Reagent	2 x 20 ml	3 x 13 ml
2-Reagent	2 x 14 ml	3 x 9 ml

The reagents when stored at 2-10°C are stable up to expiry date printed on the package. Stability on board of the analyser at 2-10°C: Prestige 24i – 12 weeks, Biolis 24i Premium – 12 weeks.

Concentrations in the test

suspension of latex particles sensitized with SLO (pH 8.2) 0.17 w/v%
 glycine buffer solution (pH 8.3)
 preservative

Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Reagent bottles should be shaken before use by gently inverting several times.
- After measurements are taken, reagent bottles should be capped and kept at 2-10°C. Care should be taken not to interchange the caps of reagent bottles.
- Reagents with different lot numbers should not be interchanged or mixed.
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.

SPECIMEN

Serum or plasma (Na-EDTA, K-EDTA, Na-Heparin, Li-Heparin, citric acid). After blood has clotted thoroughly, the sample is centrifuged and the serum is separated from blood cells and fibrins.

If the test cannot be done immediately, the sample should be placed in a tightly sealable container and stored at -20°C. Repeated freezing and thawing should be avoided.

Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-Reagent and 2-Reagent are ready to use.

1-Reagent put on basic position in reagent tray.

2-Reagent put on start position in reagent tray.

For reagent blank 0.9% NaCl is recommended.

REFERENCE VALUES ³

serum, plasma	< 160 IU/ml
---------------	-------------

It is recommended for each laboratory to establish its own reference ranges for local population. Diagnosis should only be made after taking clinical symptoms and the results of other tests into consideration.

QUALITY CONTROL

For internal quality control it is recommended to use, with each batch of samples the CORMAY IMMUNO-CONTROL I (Cat. No 4-288).

For the calibration of automatic analysers systems the CORMAY ASO CALIBRATOR kit (Cat. No 4-278) is recommended.

The calibration curve should be prepared every 12 weeks (Prestige 24i, Biolis 24i Premium), with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using the automatic analyser Biolis 30i. Results may vary if a different instrument is used.

- LoB (Limit of Blank):** 0 IU/ml
- LoD (Limit of Detection):** 4 IU/ml
- LoQ (Limit of Quantitation):** 41 IU/ml
- Linearity:** up to 718 IU/ml.

For higher concentrations dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

Specificity / Interferences

Haemoglobin up to 0.5 g/dl, bilirubin up to 20 mg/dl and triglycerides up to 500 mg/dl do not interfere with the test.

Precision

Repeatability (run to run) n = 20	Mean [IU/ml]	SD [IU/ml]	CV [%]
level 1	141	4.01	2.85
level 2	313	3.69	1.18
Reproducibility (day to day) n = 80	Mean [IU/ml]	SD [IU/ml]	CV [%]
level 1	145.4	4.72	3.2
level 2	313.8	6.08	1.9

Method comparison

A comparison between ASO values determined at **Biolis 30i** (y) and at **Beckman Coulter AU 680** (x) using 55 serum samples gave following results:

$$y = 0.9435x + 6.1558 \text{ IU/ml;}$$

$$R = 0.997 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT

Please refer to local legal requirements.

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Item name	51	ASO
Auto Rerun SW		
ON		
Auto Rerun Condition (Absorbance)		
Absorbance Range	Lower	Higher
		OFF
		OFF
Auto Rerun Range (Result)		
	ON	ON
	Lower	Higher
Serum	39	1060
Urine		
Plasma		
CSF		
Dialysis		
Other		
Prozone Range		
	OFF	

• Prestige 24i Premium, Biolis 24i Premium

Item No.	51	Item Name	ASO	Optical
Data information				
Units	IU/ml			
Decimals	0			
Calibration				
Type	Linear1			
Std sample conc.				
Blank	0	#1	*	#2
#3		#4		#5
#6				
Analysis				
Type	END method			
Main Wave Length	570 nm			
Sub Wave Length	800 nm			
Method	Immuno			
Correlation				
Y=	Slope	X+	Intercept	
	1		0	

Item No.	51	Item Name	ASO	Optical
Aspiration				
Kind	Double			
Vol.				
	Kind	Vol.	Add	Units
Sample		3	5	µl
Reagent 1		160	10	µl
Reagent 2		110	10	µl
Data Process				
Read	Main	Start	End	
	Sub	37	38	
Abs.Limit Low High				
	-3		3	
Blank value				
Water Blank				
Correction value				
Blank correction	1			
End Point Limit	2			
Linear Check (%)				
Reaction Monitor				
0 Level Point	1			
Span	3			
Prozone Check				
	Start	End	Limit (%)	
First				
Second				Low
Third				Low

Item No.	51	Item Name	ASO	Optical
Normal Range				
	Male		Female	
	Low	High	Low	High
Serum	0	160	0	160
Urine				
Plasma	0	160	0	160
CSF				
Dialysis				
Other				
Panic Range				
	Male		Female	
	Low	High	Low	High
Serum				
Urine				
Plasma				
CSF				
Dialysis				
Other				

Item No.	51	Item Name	ASO	Optical
Auto Rerun SW				
ON				
Auto Rerun Condition (Absorbance)				
	Lower	Higher		
			OFF	
			OFF	
Auto Rerun Range (Conc.)				
	First Dil	Low	High	
		Re	Value	Dil
Serum			38	1100
Urine				
Plasma				
CSF				
Dialysis				
Other				
Auto Rerun Condition (Prozone)				
	OFF			
Dilution				
	99:Dil1			

PRESTIGE 24i ASO

• **Biolis 30i**

Item no	51	Item name	ASO	Specimen	SERUM / PLASMA	OPTICAL																								
Data information																														
UNITS	IU/mL		Aspiration volume																											
DECIMALS	0		TYPE	Double																										
Analysis		END method		<table border="1"> <tr> <td></td><td>SAMPLE</td><td>REAGENT 1</td><td>REAGENT 2</td></tr> <tr> <td>VOL. (µL)</td><td>3</td><td>160</td><td>110</td></tr> <tr> <td>BOTTLE (ml)</td><td></td><td></td><td></td></tr> </table>				SAMPLE	REAGENT 1	REAGENT 2	VOL. (µL)	3	160	110	BOTTLE (ml)															
	SAMPLE	REAGENT 1	REAGENT 2																											
VOL. (µL)	3	160	110																											
BOTTLE (ml)																														
METHOD			FIRST DIL.																											
Main Wave Length	570 nm		Data processing read																											
Sub Wave Length	800 nm		<table border="1"> <tr> <td></td><td>START</td><td>END</td></tr> <tr> <td>MAIN</td><td>51</td><td>53</td></tr> <tr> <td>SUB</td><td>37</td><td>38</td></tr> </table>					START	END	MAIN	51	53	SUB	37	38															
	START	END																												
MAIN	51	53																												
SUB	37	38																												
CORRELATION (Y= AX + B)			ABS LIMIT																											
A =	1		-3		TO 3																									
B =	0		Collection value																											
Blank value			END POINT																											
	• WATER ° REAGENT		LINEARITY CHECK (%)																											
Calibration			<table border="1"> <tr> <td>TYPE</td><td colspan="2">Linear 1</td></tr> <tr> <td>STABILITY</td><td colspan="2"></td></tr> </table>				TYPE	Linear 1		STABILITY																				
TYPE	Linear 1																													
STABILITY																														
			Prozone check																											
			<table border="1"> <tr> <td></td><td>START</td><td>END</td><td>LIMIT (%)</td></tr> <tr> <td>FIRST</td><td></td><td></td><td></td></tr> <tr> <td>SECOND</td><td></td><td></td><td></td></tr> <tr> <td>°HIGH</td><td colspan="2"></td><td>MINIMUM ABS.</td></tr> <tr> <td>•LOW</td><td colspan="2"></td><td>MEAN</td></tr> <tr> <td></td><td colspan="2"></td><td>VARIATE</td></tr> </table>					START	END	LIMIT (%)	FIRST				SECOND				°HIGH			MINIMUM ABS.	•LOW			MEAN				VARIATE
	START	END	LIMIT (%)																											
FIRST																														
SECOND																														
°HIGH			MINIMUM ABS.																											
•LOW			MEAN																											
			VARIATE																											
Item No	51	Item Name	ASO	Specimen	SERUM / PLASMA	OPTICAL																								
Reference intervals																														
MALE				FEMALE																										
LOW	HIGH	LOW	HIGH																											
0	160	0	160																											
Panic range																														
MALE				FEMALE																										
LOW	HIGH	LOW	HIGH																											
Decision limit																														
MALE				FEMALE																										
Reaction check																														
		°ON		•OFF																										
CHECK																														
LOW																														
HIGH																														
VL CHECK			VH CHECK																											
°ON		•OFF		°ON		•OFF																								

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