## PRESTIGE 24i LQ CHOL

Nr kat. **4-204, 4-404** (PL)

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### LITERATURA

- Siedel J., Hägele E.O., Zigenhorn J., Wahlefeld A.W.: Clin. Chem. 29, 6 (1983).
- Tel R.M., Berends G.T.: J Clin.Chem. Clin. Biochem. 18, 10 (1980).
- Rautela G.S., Liedtke R.J.: Clin. Chem. 24, 1 (1978)
- Schettler G., Nussel E.: Arbeitsmed. Sozialmed. Preventivmed. 10, 25 (1975).
- Richmond W.: Clin. Chem. 19, 1350 (1973).
- Roeschlau P., Bernert E., Gruber W.: J. Clin. Chem. Biochem. 12, 403 (1974).
- Tietz N.W., ed. Clinical Guide to Laboratory Tests, 3rd ed. Philadelphia, PA: WB Saunders, 130 (1995).
- Trinder P.: Ann. Clin. Biochem. 6, 24 (1969).
- Dembńska-Kieć A., Naskalski J.W., Solnica B.: Diagnostyka laboratoryjna z elementami biochemii klinicznej, Wyd. IV, str. 670, 875 (2017).
- Cybulska B., Szostak WB, Filipiak KJ, et al. Polish Forum for Prevention Guidelines on Dyslipidaemia: update 2016. Kardiol Pol. 2017; 75(2): 187–190, doi: 10.5603/KP.2017.0031.
- Piepoli MF, Hoes AW, Agewall S, et al. 2016 European Guidelines on cardiovascular disease prevention in clinical practice: The Sixth Joint Task Force of the European Society of Cardiology and Other Societies on Cardiovascular Disease Prevention in Clinical Practice (constituted by representatives of 10 societies and by invited experts) Developed with the special contribution of the European Association for Cardiovascular Prevention & Rehabilitation (EACPR). Eur Heart J. 2016; 37(29): 2315–2381, doi: 10.1093/eurheartj/ehw106, indexed in Pubmed: 27222591.
- NCEP Expert Panel. Arch Inter Med (148), 36-69, (1988).
- Jacobso D.S., DeMott W.R., Grady H.J., et. al., ed., Laboratory Tests Handbook, 4th ed., Hudson, Lexi-Comp, 143, (1996).

Data wydania: 06. 2023.

## PRESTIGE 24i LQ CHOL

Cat. No **4-204, 4-404** (EN)

### INTENDED USE

Diagnostic kit for determination of total cholesterol concentration intended to use in automatic analyzers: Prestige 24i, Biolis 24i, Prestige 24i Premium, Biolis 24i Premium and Biolis 30i.


The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

### INTRODUCTION

Cholesterol is essential structural component of cell membranes and precursor of bile acids and all steroids hormones. This is why cholesterol has enormous significance for organism normal functioning. But there is also well established association between blood cholesterol concentration and coronary heart disease. Measurement of cholesterol serum level is valuable in prevention and monitoring cardiovascular disease. This determination is useful also for evaluation of intestine absorption, liver and gallbladder function.

### METHOD PRINCIPLE

Enzymatic, colorimetric method with cholesterol esterase and cholesterol oxidase.

The colour intensity is proportional to the cholesterol concentration.

### REAGENTS

#### Package

	Cat. No 4-204 (24-TRAY)	Cat. No 4-404 (36-TRAY)
1-REAGENT	6 x 60 ml	10 x 25 ml

The reagent when stored at 2-8°C is stable up to expiry date printed on the package. Stability on board of the analyzer at 2-10°C: Prestige 24i - 8 weeks, Biolis 24i Premium - 12 weeks.

### Concentrations in the test

Good's buffer (pH 6.4)	< 120 mmol/l
phenol	< 6 mmol/l
4-aminoantipyrene	< 0.4 mmol/l
cholesterol esterase (CHE)	< 4 μkat/l
cholesterol oxidase (CHO)	< 5 μkat/l
peroxidase (POD)	< 24 μkat/l
preservatives, detergents, stabilizers	

### Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.
- 1-REAGENT meeting the criteria for classification in accordance with Regulation (EC) No 1272/2008.

### Warning

H319 Causes serious eye irritation.  
P280 Wear protective gloves, protective clothing, eye protection or face protection.  
P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

### SPECIMEN

Serum, EDTA or heparinized plasma (recommended: heparine lithium, sodium or ammonium salt) free from hemolysis. Blood should be collected only if the patient has been fasting for minimum of 12 hours. Before blood collection patient should stay in rest position for about 30 minutes. Venous blood is recommended for cholesterol measurement.

Plasma cholesterol values have been reported to be 3% to 5% lower than serum cholesterol values.

Serum should be separated from red blood cells as soon as possible after blood collection.

Serum and plasma can be stored up to 3 days at 2-8°C or 6 months at -20°C.

Nevertheless it is recommended to perform the assay with freshly collected samples!

### PROCEDURE

1-REAGENT is ready to use.

1-REAGENT put on basic position in reagent tray.

Deionized water is recommended as a reagent blank.

#### Actions required:

- Biolis 24i Premium:** When performing assays in the analyser, there is a probability of **cross-contamination** affecting the tests results: CHOL – LIPASE, CREATININE - CHOL. To avoid this effect follow the recommendations contained in the instruction 51\_03\_24\_008\_BIOLIS\_24i\_PREMIUM\_CARRYOVER.
- Biolis 30i:** When performing assays in the analyser, there is a probability of **cross-contamination** affecting the tests results: CHOL – LIPASE. To avoid this effect follow the recommendations contained in the instruction 51\_03\_24\_009\_BIOLIS\_30i\_CARRYOVER.

### REFERENCE VALUES

serum / plasma	mg/dl	mmol/l
children <sup>9</sup>		
< 1year	66.1 – 228.5	1.71 – 5.91
1 – 19 years	111.4 – 202.2	2.88 – 5.23
adults <sup>9,10,11</sup>	< 190	< 5.00

It is recommended for each laboratory to establish its own reference ranges for local population.

### QUALITY CONTROL

For internal quality control it is recommended to use, with each batch of samples, the CORMAY SERUM HN (Cat. No 5-172), CORMAY SERUM HP (Cat. No 5-173) and CORMAY LIPID CONTROL LEVEL 1 (Cat. No 5-179), CORMAY LIPID CONTROL LEVEL 2 (Cat. No 5-180).

For the calibration of automatic analyzers the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) are recommended.

The calibration curve should be prepared every 8 weeks (Prestige 24i) or every 12 weeks (Biolis 24i Premium), with change of reagent lot number or as required e.g. quality control findings outside the specified range.

### PERFORMANCE CHARACTERISTICS

The following results have been obtained using the automatic analyser Biolis 30i. Results may vary if a different instrument or a manual procedure is used.

- LoB (Limit of Blank):**  
0.0 mg/dl (0.0 mmol/l)
- LoD (Limit of Detection):**  
1.4 mg/dl (0.036 mmol/l)
- LoQ (Limit of Quantitation):**  
14 mg/dl (0.36 mmol/l)
- Linearity:**  
up to 1000 mg/dl (25.9 mmol/l)

#### Specificity / Interferences

Haemoglobin up to 0.31 g/dl, ascorbate up to 62 mg/l, bilirubin up to 20 mg/dl and triglycerides up to 1000 mg/dl do not interfere with the test.

#### Precision

Repeatability (run to run) n = 20	Mean	SD	CV
	[mg/dl]	[mg/dl]	[%]
level 1	253.9	4.24	1.67
level 2	107.1	2.87	2.68
Reproducibility (day to day) n = 80	Mean	SD	CV
	[mg/dl]	[mg/dl]	[%]
level 1	256.8	6.51	2.5
level 2	109.0	2.54	2.3

#### Method comparison

A comparison between cholesterol values determined at **Biolis 30i** (y) and at **BECKMAN COULTER AU680** (x) using 57 serum samples gave following results:

$$y = 0.9489x + 11.8 \text{ mg/dl};$$

$$R = 0.997 \quad (R - \text{correlation coefficient})$$

### WASTE MANAGEMENT

Please refer to local legal requirements.

### LITERATURE

- Siedel J., Hägele E.O., Zigenhorn J., Wahlefeld A.W.: Clin. Chem. 29, 6 (1983).
- Tel R.M., Berends G.T.: J Clin.Chem. Clin. Biochem. 18, 10 (1980).
- Rautela G.S., Liedtke R.J.: Clin. Chem. 24. 1 (1978)
- Schettler G., Nussel E.: Arbeitsmed. Sozialmed. Preventivmed. 10, 25 (1975).
- Schrimond W.: Clin. Chem. 19, 1350 (1973).
- Roeschlau P., Bernt E., Gruber W.: J. Clin. Chem. Biochem. 12, 403 (1974).
- Tietz N.W., ed. Clinical Guide to Laboratory Tests, 3rd ed. Philadelphia, PA: WB Saunders, 130 (1995).
- Trinder P.: Ann. Clin. Biochem. 6, 24 (1969).
- Dembńska-Kieć A., Naskalski J.W., Solnica B.: Diagnostyka laboratoryjna z elementami biochemii klinicznej. Wyd. IV, str. 670, 875 (2017).
- Cybulska B, Szostak WB, Filipiak KJ, et al. Polish Forum for Prevention Guidelines on Dyslipidaemia: update 2016. Kardiol Pol. 2017; 75(2): 187–190, doi: 10.5603/KP.2017.0031.
- Piepoli MF, Hoes AW, Agewall S, et al. 2016 European Guidelines on cardiovascular disease prevention in clinical practice: The Sixth Joint Task Force of the European Society of Cardiology and Other Societies on Cardiovascular Disease Prevention in Clinical Practice (constituted by representatives of 10 societies and by invited experts)Developed with the special contribution of the European Association for Cardiovascular Prevention & Rehabilitation (EACPR). Eur Heart J. 2016; 37(29): 2315–2381, doi: 10.1093/eurheartj/ehw106, indexed in Pubmed: 27222591.
- NCEP Expert Panel. Arch Inter Med (148), 36-69, (1988).
- Jacobso D.S., DeMott W.R., Grady H.J., et al., ed., Laboratory Tests Handbook, 4th ed., Hudson, Lexi-Comp, 143, (1996).

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- Siedel J., Hägele E.O., Zigenhorn J., Wahlefeld A.W.: Clin. Chem. 29, 6 (1983).
- Tel R.M., Berends G.T.: J Clin.Chem. Clin. Biochem. 18, 10 (1980).
- Rautela G.S., Liedtke R.J.: Clin. Chem. 24. 1 (1978)
- Schettler G., Nussel E.: Arbeitsmed. Sozialmed. Preventivmed. 10, 25 (1975).
- Richmond W.: Clin. Chem. 19, 1350 (1973).
- Roeschlau P., Bernt E., Gruber W.: J. Clin. Chem. Biochem. 12, 403 (1974).
- Tietz N.W., ed. Clinical Guide to Laboratory Tests, 3rd ed. Philadelphia, PA: WB Saunders, 130 (1995).
- Trinder P.: Ann. Clin. Biochem. 6, 24 (1969).
- Dembńska-Kieć A., Naskalski J.W., Solnica B.: Diagnostyka laboratoryjna z elementami biochemii klinicznej, Wyd. IV, str. 670, 875 (2017).
- Cybulska B, Szostak WB, Filipiak KJ, et al. Polish Forum for Prevention Guidelines on Dyslipidaemia: update 2016. Kardiol Pol. 2017; 75(2): 187–190, doi: 10.5603/KP.2017.0031.
- Piepoli MF, Hoes AW, Agewall S, et al. 2016 European Guidelines on cardiovascular disease prevention in clinical practice: The Sixth Joint Task Force of the European Society of Cardiology and Other Societies on Cardiovascular Disease Prevention in Clinical Practice (constituted by representatives of 10 societies and by invited experts)Developed with the special contribution of the European Association for Cardiovascular Prevention & Rehabilitation (EACPR). Eur Heart J. 2016; 37(29): 2315–2381, doi: 10.1093/eurheartj/ehw106, indexed in Pubmed: 27222591.
- NCEP Expert Panel. Arch Inter Med (148), 36-69, (1988).
- Jacobso D.S., DeMott W.R., Grady H.J., et. al., ed., Laboratory Tests Handbook, 4th ed., Hudson, Lexi-Comp, 143, (1996).

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Item name	5	CHOL
<b>Auto Rerun SW</b>		
ON		
<b>Auto Rerun Condition (Absorbance)</b>		
Absorbance Range	Lower	Higher
	OFF	OFF
<b>Auto Rerun Range (Result)</b>		
	ON	ON
	Lower	Higher
Serum	20	750
Urine		
Plasma		
CSF		
Dialysis		
Other		
<b>Prozone Range</b>		
	OFF	

• **Prestige 24i Premium, Biolis 24i Premium**

Item No.	5	Item Name	CHOL	Optical
<b>Data information</b>				
Units	mg/dl			
Decimals	1			
<b>Calibration</b>				
Type	Linear2			
Std sample conc.				
Blank	0	#1	*	#2
#3		#4		#5
#6				
<b>Analysis</b>				
Type	END method			
Main Wave Length	505nm			
Sub Wave Length	700nm			
Method	CHOD-PAP			
<b>Correlation</b>				
Y=	Slope	X+	Intercept	
	1		0	

Item No.	5	Item Name	CHOL	Optical
<b>Aspiration</b>				
Kind	Single			
<b>Data Process</b>				
Read	Start	End		
Main	46	47		
Sub	9	10		
<b>Abs.Limit</b>				
Low	-0.1	High	2	
<b>Correction value</b>				
Blank correction				
End Point Limit				
2				
Linear Check (%)				
<b>Reaction Monitor</b>				
0 Level Point				
1				
Span				
3				
<b>Third mixing</b>				
OFF				

Item No.	5	Item Name	CHOL	Optical
<b>Normal Range</b>				
	Male		Female	
	Low	High	Low	High
Serum	0	200	0	200
Urine				
Plasma				
CSF				
Dialysis				
Other				
<b>Panic Range</b>				
	Male		Female	
	Low	High	Low	High
Serum				
Urine				
Plasma				
CSF				
Dialysis				
Other				

Item No.	5	Item Name	CHOL	Optical
<b>Auto Rerun SW</b>				
ON				
<b>Auto Rerun Condition (Absorbance)</b>				
	Lower	Higher		
	OFF	OFF		
<b>Auto Rerun Range (Conc.)</b>				
	First Dil	Low Value	Dil	High Value
Serum		12.95		730
Urine				
Plasma				
CSF				
Dialysis				
Other				
<b>Auto Rerun Condition (Prozone)</b>				
OFF				
<b>Dilution</b>				
100:Di12				

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