PRESTIGE 24i LQ UREA

Nr kat. **4-206, 4-406** (PL)

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Cat. No **4-206, 4-406** (EN)

INTENDED USE

Diagnostic kit for determination of urea concentration, used in automatic analysers Prestige 24i, Biolis 24i, Prestige 24i Premium, Biolis 24i Premium and Biolis 30i.



The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Urea is a product of amino acids catabolism. It is produced in liver and excreted in urine. Urea in the blood is reported as the blood urea nitrogen (BUN). Increased urea concentration in the serum, called uremia, is observed due to dehydration, renal failure, high-protein diet, increased protein catabolism caused by tissue injury or massive bleeding into the alimentary tract. The reason of reduced urea level could be overhydration, low-protein diet or starvation and severe liver disease.

METHOD PRINCIPLE

Kinetic, enzymatic method with urease and glutamate dehydrogenase.

The rate of absorbance changing at $\lambda=340$ nm is proportional to the urea concentration.

REAGENTS

Package

| | Cat. No 4-206 (24-TRAY) | Cat. No 4-406 (36-TRAY) |
|-----------|--|--|
| 1-REAGENT | 6 x 40 ml | 8 x 23 ml |
| 2-REAGENT | 6 x 12.5 ml | 8 x 7.5 ml |

The reagents when stored at 2-8°C are stable up to expiry date printed on the package. The reagents are stable for 12 weeks (Biolis 30i) on board the analyser at 2-10°C.

Concentrations in the reagent

1-REAGENT

| | |
|--|---------------|
| Tris (pH 7.8) | ≤ 144 mmol/l |
| ADP | ≤ 0.84 mmol/l |
| urease | ≤ 250 μkat/l |
| GLDH | ≤ 10.5 μkat/l |
| stabilizers, detergents, preservatives | |

2-REAGENT

| | |
|----------------------|---------------|
| 2-oxoglutarate | ≤ 48.6 mmol/l |
| NADH | ≤ 1.6 mmol/l |
| buffer, preservative | |

Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.

- 2-REAGENT meeting the criteria for classification in accordance with Regulation (EC) No 1272/2008.

Warning

H319 Causes serious eye irritation
P280 Wear protective gloves, protective clothing, eye protection or face protection.
P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes.

Remove contact lenses, if present and easy to do. Continue rinsing.

SPECIMEN

Serum, EDTA or heparinized plasma free from hemolysis, 24-hours urine.

Do not use heparine ammonium salt and fluoride as anticoagulants.

Specimen can be stored up to 7 days at 2-8°C.

Urine preparation: Samples with visible turbidity or the presence of precipitates should be pre-centrifuged.

Urine sample should be mixed well before analysis, diluted 100-fold with 0.9% NaCl and the results multiplied by 100. Bacterial growth in the specimen may cause erroneously elevated results.

24-hours urine samples should be adjusted to pH < 7 before storage.

Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-REAGENT and 2-REAGENT are ready to use.

1-REAGENT put on basic position in reagent tray.

2-REAGENT put on start position in reagent tray.

For reagent blank deionized water is recommended.

REFERENCE VALUES ⁸

| serum / plasma | mg/dl | mmol/l |
|----------------|---------|-----------|
| | < 50 | < 8.3 |
| 24-hours urine | g/24h | mmol/24h |
| | 20 – 35 | 300 – 550 |

1 mg of urea corresponds to 0.467 mg of urea nitrogen.

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use with each batch of samples, the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173) for determination in serum or CORMAY URINE CONTROL LEVEL 1 (Cat. No 5-161) or LEVEL 2 (Cat. No 5-162) for determination in urine.

For the calibration of automatic analysers systems the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) is recommended.

The calibration curve should be prepared every 9 weeks (Biolis 30i), with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using the automatic analyser Biolis 30i. Results may vary if a different instrument or a manual procedure is used.

- LoB (Limit of Blank):**
1.3 mg/dl (0.22 mmol/l)
- LoD (Limit of Detection):**
2.7 mg/dl (0.45 mmol/l)
- LoQ (Limit of Quantitation):**
4.5 mg/dl (0.75 mmol/l)
- Linearity:**
up to 340 mg/dl (56.44 mmol/l)

Specificity / Interferences

Haemoglobin up to 5 g/dl, ascorbate up to 62 mg/l, bilirubin up to 20 mg/dl and triglycerides up to 1000 mg/dl do not interfere with the test.

Precision

| Repeatability (run to run) n = 20 | Mean [mg/dl] | SD [mg/dl] | CV [%] |
|--|-----------------|---------------|-----------|
| level 1 | 35.8 | 0.98 | 2.7 |
| level 2 | 101.6 | 2.43 | 2.4 |
| Reproducibility (day to day) n = 80 | Mean [mg/dl] | SD [mg/dl] | CV [%] |
| level 1 | 34.7 | 0.87 | 2.5 |
| level 2 | 101.9 | 1.83 | 1.8 |

Method comparison

A comparison between urea values determined at **Biolis 30i** (y) and at **BECKMAN COULTER AU680** (x) using 60 serum samples gave following results:

$$y = 1.073 x - 1.2632 \text{ mg/dl};$$

$$R = 1.000 \quad (R - \text{correlation coefficient})$$

A comparison between urea values determined at **Biolis 30i** (y) and at **BECKMAN COULTER AU680** (x) using 30 plasma samples gave following results:

$$y = 0.9836 x - 0.3576 \text{ mg/dl};$$

$$R = 0.999 \quad (R - \text{correlation coefficient})$$

A comparison between urea values determined at **Biolis 30i** (y) and at **BECKMAN COULTER AU680** (x) using 30 urine samples gave following results:

$$y = 1.0785 x - 16.333 \text{ mg/dl};$$

$$R = 0.999 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT

Please refer to local legal requirements.

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| | | |
|--|--------|--------|
| Item name | 10 | UREA |
| Auto Rerun SW | | |
| ON | | |
| Auto Rerun Range (Result) | | |
| | ON | ON |
| | Lower | Higher |
| Serum | 4.5 | 340 |
| Urine | | |
| Plasma | | |
| CSF | | |
| Dialysis | | |
| Other | | |
| Auto Rerun Condition (Absorbance) | | |
| Absorbance Range | Lower | OFF |
| | Higher | OFF |
| Prozone Range | | |
| | | OFF |

• Prestige 24i Premium, Biolis 24i Premium

| | | | | |
|-------------------------|-------------|-----------|-----------|---------|
| Item No. | 10 | Item Name | UREA | Optical |
| Data information | | | | |
| Units | mg/dl | | | |
| Decimals | 1 | | | |
| Calibration | | | | |
| Type | Linear2 | | | |
| Std sample conc. | | | | |
| Blank | 0 | #1 | * | #2 |
| #3 | | #4 | | #5 |
| #6 | | | | |
| Analysis | | | | |
| Type | RATE method | | | |
| Main Wave Length | 340 nm | | | |
| Sub Wave Length | 450 nm | | | |
| Method | Urease | | | |
| Correlation | | | | |
| Y= | Slope | X+ | Intercept | |
| | 1 | | 0 | |

| | | | | |
|-------------------------|-----------|-----------|-----------|---------|
| Item No. | 10 | Item Name | UREA | Optical |
| Aspiration | | | | |
| Kind | Double | | | |
| Vol. | Kind | Vol. | Add | Units |
| | Sample | 3 | 5 | µl |
| | Reagent 1 | 160 | 10 | µl |
| | Reagent 2 | 40 | 10 | µl |
| Data Process | | | | |
| Read | Main | Start | End | |
| | Sub | | | |
| Abs.Limit | Low | -3 | High | 3 |
| Correction value | | | | |
| Blank correction | 1 | | | |
| End Point Limit | 2 | | | |
| Linear Check (%) | 80 | | | |
| Prozone Check | | | | |
| | Start | End | Limit (%) | |
| First | | | | |
| Second | | | | Low |
| Blank value | | | | |
| Water Blank | | | | |
| Reaction Monitor | | | | |
| 0 Level Point | 1 | | | |
| Span | 3 | | | |
| Third mixing | | | | |
| OFF | | | | |

| | | | | |
|---------------------|------|-----------|------|---------|
| Item No. | 10 | Item Name | UREA | Optical |
| Normal Range | | | | |
| | Male | Female | | |
| | Low | High | Low | High |
| Serum | 0 | 50 | 0 | 50 |
| Urine | | | | |
| Plasma | 0 | 50 | 0 | 50 |
| CSF | | | | |
| Dialysis | | | | |
| Other | | | | |
| Panic Range | | | | |
| | Male | Female | | |
| | Low | High | Low | High |
| Serum | | | | |
| Urine | | | | |
| Plasma | | | | |
| CSF | | | | |
| Dialysis | | | | |
| Other | | | | |

| | | | | |
|--|-----------|-----------|-------|---------|
| Item No. | 10 | Item Name | UREA | Optical |
| Auto Rerun SW | | | | |
| ON | | | | |
| Auto Rerun Range (Conc.) | | | | |
| | First Dil | Low | High | |
| | | Re | Value | Dil |
| Serum | | | 4.5 | |
| Urine | | | | |
| Plasma | | | | |
| CSF | | | | |
| Dialysis | | | | |
| Other | | | | |
| Auto Rerun Condition (Absorbance) | | | | |
| | Lower | OFF | | |
| | Higher | OFF | | |
| Auto Rerun Condition (Prozone) | | | | |
| OFF | | | | |
| Dilution | | | | |
| 100:DiI2 | | | | |

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