PRESTIGE 24i LQ ALP

Nr kat. **4-212, 4-412** (PL)

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
preservative

INTENDED USE

ALP reagent is intended to determine quantitative alkaline phosphatase activity in serum and plasma. It is used to monitoring and as an aid to diagnosis of clinical conditions associated with abnormal ALP activity. ALP reagent is intended to use on automatic analyzers Biolis 24i Premium and Biolis 30i. It is only for *in vitro* diagnosis, for healthcare professional users.

SUMMARY¹⁻⁶

Determinations of alkaline phosphatase activity in serum and plasma are mainly used as an aid in diagnosis, differentiation and monitoring of hepatobiliary disease and bone diseases. Elevated levels of the enzyme in the blood are observed, e.g. in cholestatic liver diseases, osteoporosis or Paget's disease. Decreased ALP levels are associated with Wilson's disease, Hypophosphatasia, and deficiency of magnesium, zinc and vitamin B12. Measurement of ALP activity is also used as a predictive marker in outcomes (death or liver transplantation) of patients with primary sclerosing cholangitis (PSC) and primary biliary cholangitis (PBC). The determination of the level of alkaline phosphatase is also used to monitoring the effectiveness of the PBC and IgG4-related cholangiopathy treatment.

METHOD PRINCIPLE^{7,8}

Kinetic method recommended by International Federation of Clinical Chemistry (IFCC).

2-amino-2-methyl-1-propanol + p-nitrophenylphosphate + H₂O \xrightarrow{ALP} 4-nitrophenol + 2-amino-2-methyl-1-propanol phosphate

The rate of 4-nitrophenol formation is directly proportional to the ALP activity.

REAGENTS

Package

	Cat. No 4-212 (24-TRAY)	Cat. No 4-412 (36-TRAY)
1-REAGENT	4 x 40 ml	8 x 23 ml
2-REAGENT	4 x 12.5 ml	8 x 7.5 ml

Number of tests:

Bolis 24i Premium	710	810
Bolis 30i	720	810

CONCENTRATIONS OF THE ACTIVE INGREDIENTS IN THE REAGENT

1-REAGENT

2-amino-2-methyl-1-propanol (AMP)	425 mmol/l
Mg ²⁺	2.5 mmol/l
Zn ²⁺	1.25 mmol/l

metal ion chelator

pH adjuster

preservative

REAGENT STABILITY

The reagents when stored at 2-8°C are stable up to expiry date printed on the package. The reagents stored on board of the analyzer at 2-10°C are stable for 12 weeks (Bolis 24i Premium, Biolis 30i).

WARNINGS AND NOTES

- Protect from direct sunlight and avoid contamination!
- During the reaction p-nitrophenol is produced. Do not swallow or inhale, avoid contact with skin.
- Do not use reagent beyond expiry date printed on the package.
- Do not mix reagents from different kits or lots.
- Use personal protective equipment to prevent contact with samples, reagents and controls.
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.

SPECIMEN¹

Serum, heparinized plasma free from hemolysis.

Do not use EDTA, citrate and oxalate as anticoagulants because of ALP activity inhibition!

ALP activity remains stable in specimen up to 4 hours at 15-25°C. Freezing of sample causes a decrease of enzyme activity. Frozen specimens should be thawed and kept at room temperature for 18 to 24 hours before measurement to achieve full enzyme reactivation.

Nevertheless it is recommended to perform the assay with freshly collected samples!

Follow tube manufacturers' instructions carefully when using collection tubes.

Human-origin material should be handled as potentially infectious. Standard precautions in normal laboratory work are required.

ASSAY PROCEDURE

Reagents are ready to use.

1-REAGENT put on basic position in reagent tray.

2-REAGENT put on start position in reagent tray.

For reagent blank deionized water is recommended.

Actions required:

- **Bolis 24i Premium:** When performing assays in the analyser there is a probability of **cross-contamination** affecting the tests results: ALP - TG. To avoid this effect follow the recommendations contained in the instruction 51_03_24_008_BIOLIS_24i_PREMIUM_CARRYOVER.

Application guide for Prestige 24i analyser is available on request.

The application guide should be verified and validated by the user prior to testing patient samples.

QUALITY CONTROL

Pay attention to preparation of calibrator and control before ALP determination.

For internal quality control it is recommended to use with each batch of samples, the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173). For the calibration of automatic analysers systems the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) is recommended.

The calibration curve should be prepared every week (Biolis 24i Premium, Biolis 30i).

Calibration is recommended in the following cases:

- after each change of lot,
- after instrument service,
- if controls lie outside the expected range,
- each time a new reagent kit is used.

If the quality control results do not fall within the expected values or within the range determined in the laboratory, despite a successful calibration procedure, do not report results. In this case, please take the following actions:

- verify reagents are not out of expiration date.
- verify that the required maintenance has been carried out.
- verify that the procedure has been performed in accordance with the instructions for use.
- contact the Service Department or distributor for assistance.

REFERENCE VALUES ⁹⁻¹¹

Gender	Age	U/l (37°C)	µkat/l (37°C)
female	1 – 30 days	48 – 406	0.80 – 6.77
	31 days – 1 year	124 – 341	2.07 – 5.68
	1 year – 3 years	108 – 317	1.80 – 5.28
	4 – 15 years	54 – 369	0.91 – 6.23
	16 – 18 years	35 – 124	0.58 – 2.07
	19 – 20 years	39 – 118	0.65 – 1.97
	20 – 50 years	42 – 98	0.71 – 1.67
	50 – 60 years	39 – 118	0.65 – 1.97
≥ 60 years	53 – 141	0.90 – 2.40	
male	1 – 30 days	75 – 316	1.25 – 5.27
	31 days – 1 year	82 – 383	1.37 – 6.38
	1 year – 3 years	104 – 345	1.73 – 5.75
	4 – 15 years	54 – 369	0.91 – 6.23
	16 – 18 years	58 – 331	0.97 – 5.52
	19 – 20 years	41 – 137	0.68 – 2.28
	20 – 50 years	53 – 128	0.90 – 2.18
	50 – 60 years	41 – 137	0.68 – 2.28
≥ 60 years	56 – 119	0.95 – 2.02	

It is recommended for each laboratory to establish its own reference ranges for local population and method principle.

The assay results should be used in conjunction with other data, such as symptoms, results of other tests and clinical history to make clinical decisions. It is not recommended to make clinical diagnosis based on a single result.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analyser Biolis 30i. Results may vary if a different instrument is used.

- LoB (Limit of Blank):**
1.7 U/l (0.028 µkat/l)
- LoD (Limit of Detection):**
3.9 U/l (0.065 µkat/l)
- LoQ (Limit of Quantitation):**
7 U/l (0.12 µkat/l)
- Linearity:**
up to 900 U/l (µkat/l)

For higher activity, dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

- Measurement range:**
7 U/l (0.12 µkat/l) – 900 U/l (15 µkat/l)
- Specificity / Interferences**
Haemoglobin up to 0.625 g/dl, ascorbate up to 62 mg/l, bilirubin up to 20 mg/dl and triglycerides up to 1000 mg/dl do not interfere with the test.

Precision

Repeatability (run to run) n = 20	Mean [U/l]	SD [U/l]	CV [%]
level 1	103.9	0.80	0.77
level 2	445.8	1.65	0.37
Reproducibility (day to day) n = 80	Mean [U/l]	SD [U/l]	CV [%]
level 1	106.7	4.41	4.1
level 2	484.5	15.86	3.3

Method comparison

A comparison between ALP values determined at **Biolis 30i** (y) and at **ADVIA SIEMENS 1800** (x) using 63 serum samples gave following results:
 $y = 0.9829x - 2.1291$ U/l;
 $R = 0.999$ (R – correlation coefficient)

WASTE MANAGEMENT ¹²

After use, the reagents should be handled as potentially infectious and disposed of in accordance with local legal requirements.

- Reagents residues: 18 01 07
- Empty packages: 15 01 02
- Wastewater from the analyzer: 18 01 03*

INCIDENTS ¹³

Any serious incident that has occurred in relations to the device shall be reported to the manufacturer (website address: incidents@cormay.pl) and the competent authority of the Member State in which the user and/or the patient is established.

Serious incident means any incident that directly or indirectly led, might have led or might lead to any of the following:

- the death of a patient, user or other person,

- the temporary or permanent serious deterioration of a patient's, user's or other person's state of health,
- a serious public health threat.

LITERATURE

- Rifai N., Horvath A.R., Wittwer C., ed. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics 6th ed., St. Louis, Missouri, Elsevier, 417, 1393 (2018).
- European Association for the Study of the Liver, EASL Clinical Practice Guidelines: The diagnosis and management of patients with primary biliary cholangitis, Journal of Hepatology, 67, 145-172, (2017).
- Hartleb M. i wsp., Choroby cholestatyczne u dorosłych - wytyczne postępowania Sekcji Hepatologicznej Polskiego Towarzystwa Gastroenterologii, Gastroenterologia Praktyczna, vol. 2, 5-31 (2018).
- Alvarez L, et al., Usefulness of Biochemical Markers of Bone Turnover in Assessing Response to the Treatment of Paget's Disease, Bone, 29 (5), 441-452 (2001).
- Kravets I., Paget's Disease of Bone: Diagnosis and Treatment, The American Journal of Medicine, 131 (11), 1298-1303, (2018)
- Lammers W.J., et al., Levels of Alkaline Phosphatase and Bilirubin Are Surrogate End Points of Outcomes of Patients With Primary Biliary Cirrhosis: An International Follow-up Study, Gastroenterology, 1:47, 1338-1349, (2014)
- Bessey O.A., Lowry O.H. And Brock M.J., A method for the rapid determination of alkaline phosphates with five cubic millimeters of serum, Biol. Chem., 164, 321-9 (1946).
- Bowers G.N.Jr. and McComb R.B., A continuous spectrophotometric method for measuring the activity of serum alkaline phosphatase, Clin Chem., 12(2), 70-89 (1966)
- Burtis C.A., Ashwood E.R., ed. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics 4th ed., PA: WB Saunders, 2290 (2006)
- Soldin S. J., Brugnara C., Wong E. C.: Pediatric reference Ranges, 4th ed, AACC Press, 10-11 (2003).
- Hay, W.W., et al, Current pediatric diagnosis and treatment (15th ed.). New York: Lange MedicalBooks/McGraw Hill. (2000).
- European Commission notice on technical guidance on the classification of waste (2018/C 124/01) of 9 April 2018.
- Regulation (EU) 2017/746 of the European Parliament and of the Council of 5 April 2017 on in vitro diagnostic medical devices.

LIST OF CHANGES

Previous version: 02	Current version: 03
Sections updated: <i>INTENDED USE; SUMMARY; CONCENTRATION OF THE ACTIVE INGREDIENTS IN THE REAGENT; WARNING AND NOTES; SPECIMEN; ASSAY PROCEDURE; QUALITY CONTROL; REFERENCE VALUES; PERFORMANCE CHARACTERISTICS; WASTE MANAGEMENT; LITERATURE</i>	
Sections added: <i>REAGENT STABILITY, INCIDENTS</i>	

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- Rifai N., Horvath A.R., Wittner C., ed. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics 6th ed., St. Louis, Missouri, Elsevier, 417, 1393 (2018).
- European Association for the Study of the Liver, EASL Clinical Practice Guidelines: The diagnosis and management of patients with primary biliary cholangitis, Journal of Hepatology, 67, 145-172 (2017).
- Hartleb M. i wsp., Choroby cholestatyczne u doroslych - wytyczne postepowania Sekcji Hepatologicznej Polskiego Towarzystwa Gastroenterologii, Gastroenterologia Praktyczna, vol. 2, 5-31 (2018).
- Alvarez L, et al., Usefulness of Biochemical Markers of Bone Turnover in Assessing Response to the Treatment of Paget's Disease, Bone, 29 (5), 441-452 (2001).
- Kravets I., Paget's Disease of Bone: Diagnosis and Treatment, The American Journal of Medicine, 131 (11), 1298-1303 (2018).
- Lammers W.J., et al., Levels of Alkaline Phosphatase and Bilirubin Are Surrogate End Points of Outcomes of Patients With Primary Biliary Cirrhosis: An International Follow-up Study, Gastroenterology, 1:47, 1338-1349, (2014)
- Bessey O.A., Lowry O.H. And Brock M.J., A method for the rapid determination of alkaline phosphates with five cubic millimeters of serum, Biol. Chem., 164, 321-9 (1946).
- Bowers G.N.Jr. and McComb R.B., A continuous spectrophotometric method for measuring the activity of serum alkaline phosphatase, Clin Chem., 12(2), 70-89 (1966)
- Burtis C.A., Ashwood E.R., ed. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics 4th ed., PA: WB Saunders, 2290 (2006)
- Soldin S. J., Brugnara C., Wong E. C.: Pediatric reference Ranges, 4th ed, AACCC Press, 10-11 (2003).
- Hay, W.W., et al, Current pediatric diagnosis and treatment (15th ed.). New York: Lange MedicalBooks/McGraw Hill. (2000).
- European Commission notice on technical guidance on the classification of waste (2018/C 124/01) of 9 April 2018.
- Regulation (EU) 2017/746 of the European Parliament and of the Council of 5 April 2017 on in vitro diagnostic medical devices.

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Item No.	17	Item Name	ALP	Optical
Aspiration				
Kind	Double / Single			
Vol.				
Kind	Vol.	Add	Units	
Sample	7	5	µl	
Reagent 1	200	10	µl	
Reagent 2	50	10	µl	
Blank value				
Water Blank				
Reaction Monitor				
0 Level Point	1			
Span	3			
Third mixing				
OFF				
Data Process				
Read				
Main	36	Start	End	
Sub				
Abs.Limit Low ~ High				
	-3		3	
Correction value				
Blank correction				
End Point Limit 1.1				
Linear Check (%) 90				
Prozone Check				
	Start	End	Limit (%)	
First				
Second			Low	

Item No.	17	Item Name	ALP	Optical
Normal Range				
	Male		Female	
	Low	High	Low	High
Serum	53	128	42	98
Urine				
Plasma	53	128	42	98
CSF				
Dialysis				
Other				
Panic Range				
	Male		Female	
	Low	High	Low	High
Serum				
Urine				
Plasma				
CSF				
Dialysis				
Other				

Item No.	17	Item Name	ALP	Optical	
Auto Rerun SW					
ON					
Auto Rerun Range (Conc.)					
	First Dil	Low Value		High Value	
		Re	Dil	Re	Dil
Serum			6.3		760
Urine					
Plasma					
CSF					
Dialysis					
Other					
Auto Rerun Condition (Absorbance)					
Lower OFF					
Higher OFF					
Auto Rerun Condition (Prozone)					
OFF					
Dilution					
100:Dil2					

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