## PRESTIGE 24i LQ ALAT

Nr kat. **4-216, 4-416** (PL)

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
### INTENDED USE

ALAT reagent is intended to determine quantitatively alanine aminotransferase (ALAT) activity in serum. It is used to monitoring and as an aid to diagnosis of clinical conditions associated with abnormal alanine aminotransferase activity. ALAT reagent is intended to use on automatic analysers Biolis 24i Premium and Biolis 30i. It is only for *in vitro* diagnostics, for healthcare professional users.

### SUMMARY <sup>1,2</sup>

Determinations of elevated alanine aminotransferase activity (ALAT, ALT, GPT) in serum are mainly used to aid in the diagnosis and monitor of liver diseases e.g. acute hepatitis, hepatocyte necrosis, hepatic ischemia, cirrhosis, cholestasis, liver tumor, hepatotoxic drugs, mechanical jaundice. Increased levels of ALAT can be observed in severe burns, striated muscle injury, muscle inflammation, pancreatitis, myocardial infarction, infectious mononucleosis, trauma, patients with higher risk of COVID-19 adverse outcome.

### METHOD PRINCIPLE <sup>3</sup>

Optimized, modified method according to International Federation of Clinical Chemistry (IFCC), without pyridoxal phosphate.

L-alanine + 2-oxoglutarate  $\xleftarrow{\text{ALAT}}$  pyruvate + L-glutamate

pyruvate + NADH + H<sup>+</sup>  $\xleftarrow{\text{LDH}}$  lactate + NAD<sup>+</sup>

The rate of absorbance changing at  $\lambda=340$  nm is directly proportional to alanine aminotransferase activity.

### REAGENTS

Package	Cat. No 4-216 (24-TRAY)	Cat. No 4-416 (36-TRAY)
1-REAGENT	6 x 40 ml	8 x 23 ml
2-REAGENT	6 x 12.5 ml	8 x 7.5 ml

### Number of tests:

Biolis 24i Premium	1070	810
Biolis 30i	1080	810

### CONCENTRATION OF THE ACTIVE INGREDIENTS

#### 1-REAGENT

L-alanine 625 mmol/l  
LDH 3.3 U/ml  
buffer Tris  
Stabilizers  
pH adjuster  
detergent  
preservative

### REAGENT STABILITY

The reagents when stored at 2-8°C are stable up to expiry date printed on the package. The reagents are stable for 12 weeks on board the analyser at 2-10°C (Biolis 24i Premium, Biolis 30i).

### WARNINGS AND NOTES

- Protect from direct sunlight and avoid contamination!
- Do not use reagent beyond expiry date printed on the package.
- Do not mix reagents from different kits or lots.
- Use personal protective equipment to prevent contact with samples, reagents and controls.
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.

### SPECIMEN <sup>4,5</sup>

Serum free from hemolysis. Hemolysis should be avoided, since ALAT activity in erythrocytes is 3 to 5 times higher than in normal serum. Do not freeze the samples. ALAT activity remains stable in specimen up to 3 days at 15-25°C or up to 7 days at 2-8°C. Nevertheless it is recommended to perform the assay with freshly collected samples! Follow tube manufacturers' instructions carefully when using collection tubes. Human-origin material should be handled as potentially infectious. Standard precautions in normal laboratory work are required.

### PROCEDURE

1-REAGENT and 2-REAGENT are ready to use.  
1-REAGENT put on basic position in reagent tray.  
2-REAGENT put on start position in reagent tray.  
For reagent blank deionized water is recommended.

#### Actions required:

- Biolis 24i Premium:** When performing assays in the analyser, there is a probability of **cross-contamination** affecting the tests results: ALAT - LDH, ALAT - TG, LACTATE - ALAT. To avoid this effect follow the recommendations contained in the instruction 51\_03\_24\_008\_BIOLIS\_24i\_PREMIUM\_CARRYOVER
- Biolis 30i:** When performing assays in the analyser, there is a probability of **cross-contamination** affecting the tests results: ALAT - LDH, ALAT - TG, LDH - ALAT. To avoid this effect follow the recommendations contained in the instruction 51\_03\_24\_009\_BIOLIS\_30i\_CARRYOVER.

Application guide for Prestige 24i analyser is available on request.

The application guide should be verified and validated by the user prior to testing patient samples.

#### QUALITY CONTROL

For internal quality control it is recommended to use with each batch of samples, the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173). For the calibration of automatic analysers systems the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) is recommended.

The calibration curve should be prepared every 12 weeks (Biolis 24i Premium, Biolis 30i).

Calibration is recommended in the following cases:

- after each change of lot,
- after instrument service,
- if controls lie outside the expected range,
- each time a new reagent kit is used.

If the quality control results do not fall within the expected values or within the range determined in the laboratory, despite a successful calibration procedure, do not report results. In this case, please take the following actions:

- verify reagents are not out of expiration date.
- verify that the required maintenance has been carried out.
- verify that the procedure has been performed in accordance with the instructions for use.
- contact the Service Department or distributor for assistance.

#### REFERENCE VALUES <sup>6</sup>

serum	37°C	
women	up to 34 U/l	up to 0.58 µkat/l
men	up to 45 U/l	up to 0.77 µkat/l

It is recommended for each laboratory to establish its own reference ranges for local population

#### PERFORMANCE CHARACTERISTICS

The following results have been obtained using the automatic analyser Biolis 30i and/or BS-400. Results may vary if a different instrument or a manual procedure is used.

- **LoB (Limit of Blank):**  
2.0 U/l (0.03 µkat/l) – BS-400
- **LoD (Limit of Detection):**  
3.0 U/l (0.05 µkat/l) – BS-400
- **LoQ (Limit of Quantitation):**  
7.0 U/l (0.12 µkat/l) – Biolis 30i
- **Linearity:**  
up to 675 U/l (11.25 µkat/l) – Biolis 30i
- **Measurement blank**  
7.0 U/l (0.12 µkat/l) - 675 U/l (11.25 µkat/l) – Biolis 30i
- **Specificity / Interferences**

Ascorbic acid up to 62 mg/l, triglycerides up to 1000 mg/dl, bilirubin up to 20 mg/dl, haemoglobin up to 0.31 g/dl (in samples containing low activity of ALAT)  
ALAT (II GENERACJA/ II GENERATION/ II ПОКОЛЕНИЕ)

and up to 2.5 g/dl (in samples containing high activity ALAT) do not interfere with the test.

#### ▪ Precision

Repeatability (run to run) n=20	Mean [U/l]	SD [U/l]	CV [%]
level 1	31.1	1.33	4.3
level 2	96.0	1.09	1.1
Reproducibility (day to day) n=80	Mean [U/l]	SD [U/l]	CV [%]
level 1	31.1	1.09	3.5
level 2	106.6	1.25	1.2

#### ▪ Method comparison

A comparison between ALAT values determined at **Biolis 30i** (y) and at **ADVIA SIEMENS 1800** (x) using 59 serum samples gave following results:

$$y = 0.9872 x + 0.442 \text{ U/l};$$

$$R = 0.999 \quad (R - \text{correlation coefficient})$$

#### WASTE MANAGEMENT <sup>7</sup>

After use, the reagents should be handled as potentially infectious and disposed of in accordance with local legal requirements.

- Reagents residues: 18 01 07
- Empty packages: 15 01 02
- Wastewater from the analyzer: 18 01 03\*

#### INCIDENTS <sup>8</sup>

Any serious incident that has occurred in relations to the device shall be reported to the manufacturer (website address: [incidents@cormay.pl](mailto:incidents@cormay.pl)) and the competent authority of the Member State in which the user and/or the patient is established.

**Serious incident** means any incident that directly or indirectly led, might have led or might lead to any of the following:

- the death of a patient, user or other person,
- the temporary or permanent serious deterioration of a patient's, user's or other person's state of health,
- a serious public health threat.

#### LITERATURE

1. Pagana K. Pagana T.J., Mosby's Diagnostic and Laboratory Test Reference, 10<sup>th</sup> ed., Elsevier, St. Louis, USA, 27 (2013).
2. Lippi G., Plebani M., Laboratory abnormalities in patients with COVID-2019 infection, Clin Chem Lab Med, 58(7):1131-1134 (2020).
3. Bergmeyer H.U., Horder M., Rej R.: J. Clin. Chem. Clin. Biochem. 24, 481-495 (1986).
4. Kaplan L.A., Pesce A.J., Clinical Chemistry: Theory, analysis, and correlation, 2nd ed., The C.V. Mosby Company, St. Louis, USA, 896 (1989).
5. World Health Organization (2002). Use of Anticoagulants in Diagnostic Laboratory Investigations (Stability of Blood, Plasma and serum samples); WHO/DIL/LAB/99.1 Rev.2: 1-64
6. Burtis C.A., Ashwood E.R., Bruns D.E. ed. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics, 4th ed., Elsevier Saunders, St. Louis, USA, 2254 (2006).
7. European Commission notice on technical guidance on the classification of waste (2018/C 124/01) of 9 April 2018.
8. Regulation (EU) 2017/746 of the European Parliament and of the Council of 5 April 2017 on in vitro diagnostic medical devices

#### LIST OF CHANGES

Previous version: 05	Current version: 06
<i>Sections updated: WASTE MANAGEMENT; DOCUMENT NUMBER</i>	

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Item No.	12	Item Name	ALAT II GEN	Optical
<b>Aspiration</b>				
Kind	Double			
Vol.				
	Kind	Vol.	Add	Units
	Sample	18	5	µl
	Reagent 1	200	10	µl
	Reagent 2	50	10	µl
<b>Blank value</b>				
Water Blank				
<b>Reaction Monitor</b>				
0 Level Point	1			
Span	3			
<b>Third mixing</b>				
OFF				
<b>Data Process</b>				
Read	Start	End		
	Main	36		
	Sub	46		
Abs.Limit Low High				
	-3	~		3
<b>Correction value</b>				
Blank correction		1		
End Point Limit		2		
Linear Check (%)		80		
<b>Prozone Check</b>				
	Start	End	Limit (%)	
First				
Second			Low	

Item No.	12	Item Name	ALAT II GEN	Optical
<b>Normal Range</b>				
	Male		Female	
	Low	High	Low	High
Serum	0	45	0	34
Urine				
Plasma				
CSF				
Dialysis				
Other				
<b>Panic Range</b>				
	Male		Female	
	Low	High	Low	High
Serum				
Urine				
Plasma				
CSF				
Dialysis				
Other				

Item No.	12	Item Name	ALAT II GEN	Optical			
<b>Auto Rerun SW</b>							
ON							
<b>Auto Rerun Range (Conc.)</b>							
	First Dil	Low		High			
		Re	Value	Dil	Re	Value	Dil
Serum			7			675	
Urine							
Plasma							
CSF							
Dialysis							
Other							
<b>Auto Rerun Condition (Absorbance)</b>							
Lower		OFF					
Higher		OFF					
<b>Auto Rerun Condition (Prozone)</b>							
OFF							
<b>Dilution</b>							
100:Di2							

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