PRESTIGE 24i dTIBC

Nr kat. **4-223, 4-440** (PL)

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Cat. No **4-223, 4-440** (EN)

INTENDED USE

Diagnostic kit for determination of total iron binding capacity intended to use in automatic analyzers: Prestige 24i, Biolis 24i, Prestige 24i Premium, Biolis 24i Premium and Biolis 30i.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Iron is transported as Fe (III) bound to the apotransferrin. The apotransferrin-Fe (III) complex is called transferrin. Normally only about one third of the iron binding sites of transferrin are occupied by Fe (III). TIBC (Total iron-binding capacity) is a measurement for the maximum iron concentration that transferrin can bind.


Serum TIBC levels vary in disorders of iron metabolism, e.g. elevated TIBC level when iron levels are low is useful in the diagnosis and monitoring of iron deficiency anaemia and also anaemia in late pregnancy. Decreased level of TIBC may indicate disorders such as anaemia associated with chronic inflammatory disorders, haemochromatosis, other anaemias, thalassemia, cirrhosis and malignancies.

METHOD PRINCIPLE


Colorimetric method, direct, with chromazurol B.

The assay consists of 2 reaction steps:

- The sample is added to 1-Reagent acidic buffer which contains excess iron and iron-binding dye (Chromazurol B). The low pH causes transferrin to release all iron atoms which then are bound by the iron-binding dye along with the excess iron, producing a coloured dye-iron complex.

- Neutral buffer 2-Reagent addition results in pH change and causes a large increase in affinity of transferrin for iron. Then serum transferrin rapidly binding iron ions separating it from the dye-iron complex. The observed decrease in the absorbance of a coloured dye-iron complex is directly proportional to the total iron-binding capacity TIBC in the serum sample.

REAGENTS

Package	Cat. No 4-223 (24-TRAY)	Cat. No 4-440 (36-TRAY)
1-Reagent	1 x 22 ml	1 x 21 ml
2-Reagent	1 x 7.5 ml	1 x 7 ml

The reagents when stored at 2-8°C are stable up to expiry date printed on the package. The reagents are stable for 7 weeks on board the analyzer.

Concentrations in the test

1-Reagent	
sodium acetate	200 mmol/l
iron chloride	0.02 mmol/l
chromazurol B	< 0.25 mmol/l
hexadecyltrimethylammonium bromide (CTAB)	< 3.0 mmol/l
preservatives, conservants	
2-Reagent	
sodium bicarbonate	325 mmol/l
buffer, preservatives, conservants	

Warnings and notes

- Do not freeze the reagents.
- Protect from light and air, avoid contamination!
- Contaminated glassware is the greatest source of error. The use of disposable plastic ware is recommended.
- Because of the great quantities of iron in the environment, scrupulous care is necessary to ensure that glassware, water and reagents do not become contaminated with iron.
- For internal quality control not frozen control serum is recommended.
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.

SPECIMEN

Serum free from haemolysis. Presence of fibrin in the serum sample caused by incomplete clot formation causes non-specific interference with TIBC determination.

Do not use plasma samples as anticoagulants affect TIBC test results.

Serum can be stored up to 3 days at 4°C or up to 6 months at -20°C or also at -70°C indefinitely.

Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-Reagent and 2-Reagent are ready to use.

For reagent blank deionised water is recommended.

Actions required:

- Biolis 24i Premium:** When performing assays in the analyzer there is a probability of **cross-contamination** affecting the tests results. To avoid this effect, the tests for determination of total iron binding capacity should be performed **in the separate order** (follow the recommendations contained in the instruction 51_03_24_008_BIOLIS_24i_PREMIUM_CARRYOVER).
- Biolis 30i:** When performing assays in the analyser, there is a probability of **cross-contamination** affecting the tests results. To avoid this effect follow the recommendations contained in the instruction 51_03_24_009_BIOLIS_30i_CARRYOVER.

REFERENCE VALUES ²

serum	µg/dl	µmol/l
adults	250 – 425	44,8 – 71,6

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173) with each batch of samples.

For the calibration of automatic analysers Prestige 24i and Biolis 24i Premium the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) or CORMAY MULTICALIBRATOR LEVEL 2 (Cat. No 5-175; 5-177) is recommended depending on the calibrator lot number.

For the calibration of the automatic analyser Biolis 30i, CORMAY MULTICALIBRATOR LEVEL 2 (Cat. No 5-175; 5-177) is recommended.

The calibration curve should be prepared every 7 weeks, with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analyzer Biolis 30i. Results may vary if a different instrument is used.

- LoB (Limit of Blank):** 13 µg/dl (2.33 µmol/l)
- LoD (Limit of Detection):** 21 µg/dl (3.76 µmol/l)
- LoQ (Limit of Quantitation):** 115 µg/dl (20.59 µmol/l)
- Linearity:** up to 636 µg/dl (113.84 µmol/l)

For higher concentration dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

Specificity / Interferences

Haemoglobin up to 1.25 g/dl, ascorbate up to 31 mg/l, bilirubin up to 20 mg/dl, triglycerides up to 1000 mg/dl, copper up to 3 mg/dl, zinc up to 250 µg/dl, iron up to 357.5 µg/dl, Desferal up to 11.5 µg/ml, Cuprimine up to 250 µg/dl and Iron Dextran (Imferon) up to 1430 µg/ml do not interfere with the test.

Precision

Repeatability (run to run) n = 20	Mean [µg/dl]	SD [µg/dl]	CV [%]
level 1	222	3.03	1.37
level 2	367	2.74	0.75

Method comparison

A comparison between TIBC values determined at **Biolis 30i** (y) and at **BS-400** (x) using 63 serum samples gave following results:

$$y = 1.0742x - 25.156 \mu\text{g/dl};$$

$$R = 0.987 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT

Please refer to local legal requirements.

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Item name	62	dTIBC
Auto Rerun SW		
ON		
Auto Rerun Range (Result)		
	ON	ON
	Lower	Higher
Serum	50	600
Urine		
Plasma		
CSF		
Dialysis		
Other		
Auto Rerun Condition (Absorbance)		
Absorbance Range	Lower	OFF
	Higher	OFF
Prozone Range		
OFF		

PROGRAM NA ANALIZATOR / APPLICATION / АДАПТАЦІЯ: Prestige 24i Premium, Biolis 24i Premium

Item No.	62	Item Name	dTIBC	Optical
Data information				
Units	µg/dl			
Decimals	0			
Analysis				
Type	END method			
Main Wave Length	660 nm			
Sub Wave Length				
Method	Direct			
Correlation				
Y=	Slope	X+	Intercept	
	1		0	
Calibration				
Type	Linear1			
Std sample conc.				
Blank	0	#1	*	#2
#3		#4		#5
#6				

Item No.	62	Item Name	dTIBC	Optical
Aspiration				
Kind	Double			
Vol.	Kind	Vol.	Add	Units
	Sample	13.5	5	µl
	Reagent 1	180	10	µl
	Reagent 2	45	10	µl
Data Process				
Read	Main	Start	End	
	Sub	30	31	
Abs.Limit				
Low	-3	High	3	
Blank value				
Water Blank				
Reaction Monitor				
0 Level Point	1			
Span	3			
Third mixing				
OFF				
Correction value				
Blank correction				
End Point Limit				
Linear Check (%)				
Prozone Check				
	Start	End	Limit (%)	
First				
Second				Low

Item No.	62	Item Name	dTIBC	Optical
Normal Range				
	Male		Female	
	Low	High	Low	High
Serum	250	425	250	425
Urine				
Plasma				
CSF				
Dialysis				
Other				
Panic Range				
	Male		Female	
	Low	High	Low	High
Serum				
Urine				
Plasma				
CSF				
Dialysis				
Other				

Item No.	62	Item Name	dTIBC	Optical
Auto Rerun SW				
ON				
Auto Rerun Range (Conc.)				
	First Dil	Low	High	
		Re	Value	Dil
Serum			67	600
Urine				
Plasma				
CSF				
Dialysis				
Other				
Auto Rerun Condition (Absorbance)				
Lower				
Higher				
OFF				
Auto Rerun Condition (Prozone)				
OFF				
Dilution				
100:Dil 2				

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