PRESTIGE 24i LQ MG

Nr kat. 4-229, 4-429, 4-462 (PL)

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(EN)

INTENDED USE

Diagnostic kit for determination of magnesium concentration used in automatic analysers Prestige 24i, Biolis 24i, Prestige 24i Premium, Biolis 24i Premium and Biolis 30i. The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Magnesium in human organism occurs mainly in bone (about 50%) but is present also intracellularly in other tissues. Magnesium serves as a cofactor for multiple enzymatic reactions involved in nucleic acids synthesis, transport and production of energy. Magnesium is important in neuromuscular conduction and activation. Reduced magnesium level generates: concentration disturbances, fatigue, muscle tremor, anxiety state.

METHOD PRINCIPLE

Magnesium forms a purple coloured complex in alkaline solution. In the presence of EGTA, the reaction is specific. The intensity of the purple colour is proportional to the magnesium concentration.

REAGENTS

Package

| | Cat. No 4-229 (24-TRAY) | Cat. No 4-429 (36-TRAY) | Cat. No 4-462 (36-TRAY) |
|-----------|----------------------------|----------------------------|----------------------------|
| 1-REAGENT | 4 x 60 ml | 10 x 25 ml | 5 x 14.5 ml |

The reagent is stable up to the kit expiry date printed on the package when stored at 2-8°C. The reagent stored on board of the analyser at 2-10°C is stable for 3 weeks (Prestige 24i, Biolis 24i Premium) or 12 weeks (Biolis 30i).

Concentrations in the test

| | |
|---------------------|---------------|
| xylydyl blue | ≤ 0.18 mmol/l |
| EGTA | ≤ 0.12 mmol/l |
| CAPS | ≤ 0.06 mol/l |
| potassium hydroxide | ≤ 0.07 mol/l |
| buffer (pH 11.5) | ≤ 0.14 mol/l |
| detergent | |
| preservative | |

WARNINGS AND NOTES

- Protect from direct sunlight and avoid contamination!
- The reagent is air sensitive, to extend reagents stability it is recommended to keep reagent's bottles recapped on the board of analyser.
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.
- 1-REAGENT meeting the criteria for classification in accordance with Regulation (EC) No 1272/2008.

Ingredients:

1-REAGENT contains potassium hydroxide.

Danger

H314 Causes severe skin burns and eye damage.

P280 Wear protective gloves, protective clothing, eye protection or face protection.

P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P310 Immediately call a POISON CENTER or doctor.

SPECIMEN ^{6,7}

Serum, heparinized plasma free from hemolysis, 24-hours urine.

Recommended anticoagulants: heparine lithium, sodium or ammonium salt.

Serum should be separated from red blood cells as soon as possible after blood collection, because erythrocytes contain approximately 3 times the magnesium concentration found in normal serum.

Urine preparation: acidify urine with some drops of concentrated hydrochloride acid to pH 1.0. Then dilute 1 part of acidified urine with 4 parts of distilled water. Multiply the result by 5. Mix well samples before analysis.

Serum and plasma can be stored up to 7 days at 2-8°C. For longer storage samples should be frozen at -20°C.

24-hours urine samples can be stored up to 7 days at 2-8°C. Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-REAGENT is ready to use.

Avoid foaming.

1-REAGENT put on basic position in reagent tray.

For reagent blank deionized water is recommended.

Actions required:

- Biolis 24i Premium:** When performing assays in the analyzer there is a probability of **cross-contamination** affecting the tests results. To avoid this effect, it is recommended to **precede the MG test by wash cycle** with BIOLIS WASHING SOLUTION (Cat. No 3-126, 3-127) according to the guidelines contained in the instruction: 51_03_24_008_BIOLIS_24i_PREMIUM_CARRYOVER.
- Biolis 30i:** When performing assays in the analyzer, there is a probability of **cross-contamination** affecting the tests results. To avoid this effect, the tests for determination of magnesium concentration should be performed **in the separate order** (follow the recommendations contained in the instruction 51_03_24_009_BIOLIS_30i_CARRYOVER).

REFERENCE VALUES

| serum / plasma | mg/dl | mmol/l |
|---------------------|--------------|-------------|
| newborn 2 – 4 d | 1.5 – 2.2 | 0.62 – 0.91 |
| children 5 mo – 6 y | 1.7 – 2.3 | 0.70 – 0.95 |
| 6 – 12 y | 1.7 – 2.1 | 0.70 – 0.86 |
| 12 – 20 y | 1.7 – 2.2 | 0.70 – 0.91 |
| adults | 1.6 – 2.6 | 0.66 – 1.07 |
| 24-hours urine | mg/24h | mmol/24h |
| | 72.9 – 145.8 | 3 – 5 |

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use, with each batch of samples, the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173) for determination in serum, CORMAY URINE CONTROL LEVEL 1 (Cat. No 5-161) and LEVEL 2 (Cat. No 5-162) for determination in urine.

For the calibration of automatic analysers systems the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) is recommended.

The calibration curve should be prepared every week (Prestige 24i, Biolis 24i Premium) or every 2 weeks (Biolis 30i), with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using the automatic analyser Biolis 30i. Results may vary if a different instrument or a manual procedure is used.

- LoB (Limit of Blank):**
0.02 mg/dl (0.008 mmol/l)
- LoD (Limit of Detection):**
0.03 mg/dl (0.012 mmol/l)
- LoQ (Limit of Quantitation):**
0.07 mg/dl (0.029 mmol/l)
- Linearity:**
up to 5.0 mg/dl (2.06 mmol/l)

Specificity / Interferences

Haemoglobin up to 0.313 g/dl, ascorbate up to 62 mg/l, bilirubin up to 15 mg/dl, triglycerides up to 1000 mg/dl and calcium up to 20 mg/dl do not interfere with the test.

Precision

| Repeatability (run to run) n = 20 | Mean [mg/dl] | SD [mg/dl] | CV [%] |
|----------------------------------------|-----------------|---------------|-----------|
| level 1 | 2.02 | 0.02 | 1.04 |
| level 2 | 4.31 | 0.04 | 0.92 |
| Reproducibility (day to day) n = 80 | Mean [mg/dl] | SD [mg/dl] | CV [%] |
| level 1 | 2.0 | 0.08 | 3.7 |
| level 2 | 4.1 | 0.11 | 2.7 |

Method comparison

A comparison between magnesium values determined at **Biolis 30i** (y) and at **ADVIA SIEMENS 1800** (x) using 59 serum samples gave following results:

$$y = 0.9927x + 0.0361 \text{ mg/dl};$$

$$R = 0.975 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT

Please refer to local legal requirements.

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| | | |
|------------------------------------------|--------|--------|
| Item name | 22 | MG |
| Auto Rerun SW | | |
| ON | | |
| Auto Rerun Condition (Absorbance) | | |
| Absorbance Range | Lower | OFF |
| | Higher | OFF |
| Auto Rerun Range (Result) | | |
| | ON | ON |
| | Lower | Higher |
| Serum | 0.14 | 5.0 |
| Urine | | |
| Plasma | | |
| CSF | | |
| Dialysis | | |
| Other | | |
| Prozone Range | | |
| | | OFF |

- Prestige 24i Premium, Biolis 24i Premium

| | | | | | |
|-------------------------|--------------|-----------|----------------|---------|---|
| Item No. | 22 | Item Name | MG | Optical | |
| Data information | | | | | |
| Units | mg/dl | | | | |
| Decimals | 2 | | | | |
| Calibration | | | | | |
| Type | Linear2 | | | | |
| Std sample conc. | | | | | |
| Blank | 0 | #1 | * | #2 | * |
| #3 | | #4 | | #5 | |
| #6 | | | | | |
| Analysis | | | | | |
| Type | END method | | | | |
| Main Wave Length | 505 nm | | | | |
| Sub Wave Length | 700 nm | | | | |
| Method | Xylidyl blue | | | | |
| Correlation | | | | | |
| Y= | Slope 1 | X+ | Intercept 0 | | |

| | | | | |
|-------------------------|--------|-----------|-----------|---------|
| Item No. | 22 | Item Name | MG | Optical |
| Aspiration | | | | |
| Kind | Single | | | |
| Vol. | Kind | Vol. | Add | Units |
| Sample | | 3 | 5 | µl |
| Reagent 1 | | 250 | 10 | µl |
| Data Process | | | | |
| Read | Main | Start | End | |
| | Sub | 9 | 10 | |
| Abs.Limit | | | | |
| Low | -0.1 | High | 2 | |
| Correction value | | | | |
| Blank correction | | | | |
| End Point Limit | | | | |
| 2 | | | | |
| Linear Check (%) | | | | |
| Prozone Check | | | | |
| | Start | End | Limit (%) | |
| First | | | | |
| Second | | | | Low |
| Blank value | | | | |
| Water Blank | | | | |
| Reaction Monitor | | | | |
| 0 Level Point | | | | |
| 1 | | | | |
| Span | | | | |
| 3 | | | | |
| Third mixing | | | | |
| OFF | | | | |

| | | | | |
|---------------------|------|-----------|--------|---------|
| Item No. | 22 | Item Name | MG | Optical |
| Normal Range | | | | |
| | Male | | Female | |
| | Low | High | Low | High |
| Serum | 1.6 | 2.6 | 1.6 | 2.6 |
| Urine | | | | |
| Plasma | 1.6 | 2.6 | 1.6 | 2.6 |
| CSF | | | | |
| Dialysis | | | | |
| Other | | | | |
| Panic Range | | | | |
| | Male | | Female | |
| | Low | High | Low | High |
| Serum | | | | |
| Urine | | | | |
| Plasma | | | | |
| CSF | | | | |
| Dialysis | | | | |
| Other | | | | |

| | | | | | |
|------------------------------------------|-----------|-----------|-----|------------|-----|
| Item No. | 22 | Item Name | MG | Optical | |
| Auto Rerun SW | | | | | |
| ON | | | | | |
| Auto Rerun Condition (Absorbance) | | | | | |
| | Lower | OFF | | | |
| | Higher | OFF | | | |
| Auto Rerun Range (Conc.) | | | | | |
| | First Dil | Low Value | Dil | High Value | Dil |
| Serum | | 0.06 | | 5.0 | |
| Urine | | | | | |
| Plasma | | | | | |
| CSF | | | | | |
| Dialysis | | | | | |
| Other | | | | | |
| Auto Rerun Condition (Prozone) | | | | | |
| OFF | | | | | |
| Dilution | | | | | |
| 100:Dil2 | | | | | |

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