PRESTIGE 24i MICROALBUMIN

Nr kat. 4-244, 4-444

(PL)

2-Reagent

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Cat. No **4-244, 4-444** (EN)

INTENDED USE

Diagnostic kit for determination of albumin concentration in urine and cerebrospinal fluid intended to use in automatic analyzers: Prestige 24i, Biolis 24i, Prestige 24i Premium, Biolis 24i Premium and Biolis 30i.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Albumin is a protein that is formed within the liver and it makes up approximately 60% of the serum protein. Normally only small amounts of albumin are filtered through the renal glomeruli, and that small quantity can be reabsorbed by the renal tubules. In that case there is a low albumin concentration in the urine. When renal disorders appear, level of urine albumin increase but remains still undetectable by routine screening tests (microalbuminuria). The appearance of low but abnormal levels (30-300 mg/24h) of albumin in the urine is an early clinical evidence of nephropathy (mostly diabetic) and cardiovascular disorders.

To avoid the necessity of 24-hour urine collection it is common in clinical practice to measure albumin and creatinine simultaneously and give the result as a albumin/creatinine ratio.

METHOD PRINCIPLE

Immunoturbidimetric method. Albumin in the sample forms with anti-albumin antibodies in the reagent an insoluble complex. The turbidity caused by the complexes is measured spectrophotometrically at 340 nm and is proportional to the amount of albumin in the sample.

REAGENTS

Package

| | Cat. No 4-244 (24-TRAY) | Cat. No 4-444 (36-TRAY) |
|-----------|----------------------------|----------------------------|
| 1-Reagent | 2 x 15.1 ml | 2 x 17.3 ml |
| 2-Reagent | 2 x 4.5 ml | 2 x 4.8 ml |

The reagents are stable up to the expiry date printed on the package when stored at 2-8°C. Stability on board of the analyser at 2-10°C: Prestige 24i – 8 weeks, Biolis 24i Premium – 12 weeks.

Concentrations in the test

| | |
|-------------------------------|--------------|
| 1-Reagent | |
| Tris buffer (pH 7.6) | 18.2 mmol/l |
| sodium chloride | 123.2 mmol/l |
| PEG | < 4% |
| 2-Reagent | |
| sodium chloride | 154 mmol/l |
| anti-human albumin antibodies | |
| preservatives | |

Warnings and notes

- Protect from direct sunlight and avoid contamination.
- Do not freeze the reagents.
- Before use mix reagent by gently inverting each bottle.
- Do not use after expiry date.
- Do not interchange caps.
- Reagents with different lot numbers should not be interchanged or mixed.
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.

SPECIMEN

Urine. Urine used for analysis may come from the first morning sample, random sample or timed collection sample ^[3].

Samples with visible turbidity should be centrifuged before analysis. Determination of uncentrifuged samples may give increased results.

Urine samples are stable for 2 days at room temperature, 14 days at the 8°C ^[7].

Nevertheless it is recommended to perform the assay with freshly collected samples!

Cerebrospinal fluid. CSF should be centrifuged before analysis. If the total protein in CSF is greater than 2000 mg/l, the CSF sample needs to be diluted 1:9 and the result multiplied by 10.

It is recommended to follow NCCLS procedures regarding specimen collecting and handling.

Samples should be stored at 2-4°C and analyzed within 2 hours after collection.

Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-Reagent and 2-Reagent are ready to use.

1-Reagent put on basic position in reagent tray.

2-Reagent put on start position in reagent tray.

For reagent blank 0.9% NaCl is recommended.

Actions required:

- Biolis 30i:** When performing assays in the analyzer, there is a probability of **cross-contamination** affecting the tests results. To avoid this effect follow the recommendations contained in the instruction: 51_03_24_009_BIOLIS_30i_CARRYOVER.

Calculation

For the calculation of albumin 24 hours quantity, multiply the concentration (mg/l) with the volume (l) of the 24 hours urines.

REFERENCE VALUES ³

| urine | mg/24h | µg/min | mg/g creatinine |
|--|----------------|----------|-----------------|
| normal | < 30 | < 20 | < 30 |
| microalbuminuria | 30 – 300 | 20 – 200 | 30 – 300 |
| clinical albuminuria (overt nephropathy) | > 300 | > 200 | > 300 |
| cerebrospinal fluid, lumbar | 177 – 251 mg/l | | |

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use, with each batch of samples, the CORMAY MICROALBUMIN CONTROL (Cat. No 4-461).

For the calibration of automatic analysers systems the CORMAY MICROALBUMIN CALIBRATOR (Cat. No 5-193) is recommended.

The calibration curve should be prepared every 6 weeks (Prestige 24i) or 12 weeks (Biolis 24i Premium), with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results been obtained using the automatic analyser Biolis 30i. Results may vary if a different instrument or a manual procedure is used.

- LoB (Limit of Blank):** 0.0 mg/l
- LoD (Limit of Detection):** 1.21 mg/l
- LoQ (Limit of Quantitation):** 5.2 mg/l
- Linearity:** up to 331 mg/l

For higher concentration dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

Specificity / Interferences

Hemoglobin up to 2.5 g/dl, ascorbate up to 200 mg/dl, creatinine up to 6 g/l, uric acid up to 100 mg/dl, glucose up to 35 g/l, urea up to 50 g/l, bilirubin conjugated up to 60 mg/dl, calcium ion up to 130 mg/dl, magnesium ion up to 1,8 g/l, do not interfere with the test.

Precision

| Repeatability (run to run) n = 20 | Mean [mg/l] | SD [mg/l] | CV [%] |
|--|----------------|--------------|-----------|
| level 1 | 66.28 | 0.57 | 0.86 |
| Reproducibility (day to day) n = 80 | Mean [mg/l] | SD [mg/l] | CV [%] |
| level 1 | 68.9 | 1.51 | 2.2 |

Method comparison

A comparison between microalbumin values determined at **Biolis 30i** (y) and at **BS-400** (x) using 59 urine samples gave following results:

$$y = 0.9979 x - 5.6224 \text{ mg/l};$$

$$R = 0.994 \quad (R - \text{correlation coefficient})$$

A comparison between microalbumin values determined at **Biolis 30i** (y) and at **BS-400** (x) using 57 cerebrospinal fluid samples gave following results:

$$y = 0.9973 x - 2.559 \text{ mg/l};$$

$$R = 0.999 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT

Please refer to local legal requirements.

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| | | |
|--|--------|--------|
| Item name | 56 | MALB |
| Auto Rerun SW | | |
| OFF | | |
| Auto Rerun Range (Result) | | |
| | OFF | OFF |
| | Lower | Higher |
| Serum | | |
| Urine | | |
| Plasma | | |
| CSF | | |
| Dialysis | | |
| Other | | |
| Auto Rerun Condition (Absorbance) | | |
| Absorbance Range | | |
| | Lower | OFF |
| | Higher | OFF |
| Auto Rerun Condition (Prozone) | | |
| Prozone Range | | |
| OFF | | |

▪ **Prestige 24i Premium, Biolis 24i Premium**

| | | | | |
|-------------------------|---------------|-----------|------|---------|
| Item No. | 56 | Item Name | MALB | Optical |
| Data information | | | | |
| Units | mg/l | | | |
| Decimals | 2 | | | |
| Calibration | | | | |
| Type | Spline1 | | | |
| Std sample conc. | | | | |
| Blank | 0 | #1 | * | #2 * |
| #3 | * | #4 | * | #5 * |
| #6 | | | | |
| Analysis | | | | |
| Type | END method | | | |
| Main Wave Length | 340 nm | | | |
| Sub Wave Length | 700 nm | | | |
| Method | Turbidimetric | | | |
| Correlation | | | | |
| Slope | | Intercept | | |
| Y= | 1 | X+ | 0 | |

| | | | | |
|-------------------------|--------|-----------|-----------|---------|
| Item No. | 56 | Item Name | MALB | Optical |
| Aspiration | | | | |
| Kind | Double | | | |
| Vol. | | | | |
| | Kind | Vol. | Add | Units |
| Sample | 12 | 200 | 5 | µl |
| Reagent 1 | 40 | 200 | 10 | µl |
| Reagent 2 | 40 | 200 | 10 | µl |
| Data Process | | | | |
| Read | Main | Start | End | |
| | Sub | 30 | 52 | 31 |
| Abs.Limit Low | | | | |
| | -3.000 | ~ | High | 3.000 |
| Correction value | | | | |
| Blank correction | | | | 0.8413 |
| End Point Limit | | | | 2 |
| Linear Check (%) | | | | |
| Prozone Check | | | | |
| | Start | End | Limit (%) | |
| First | | | | |
| Second | | | Low | |
| Reaction Monitor | | | | |
| 0 Level Point | | | | 1 |
| Span | | | | 3 |
| Third mixing | | | | |
| OFF | | | | |

| | | | | |
|---------------------|------|-----------|--------|---------|
| Item No. | 56 | Item Name | MALB | Optical |
| Normal Range | | | | |
| | Male | | Female | |
| | Low | High | Low | High |
| Serum | | | | |
| Urine | | | | |
| Plasma | | | | |
| CSF | | | | |
| Dialysis | | | | |
| Other | | | | |
| Panic Range | | | | |
| | Male | | Female | |
| | Low | High | Low | High |
| Serum | | | | |
| Urine | | | | |
| Plasma | | | | |
| CSF | | | | |
| Dialysis | | | | |
| Other | | | | |

| | | | | | | | |
|--|-----------|-----------|-------|---------|----|-------|-----|
| Item No. | 56 | Item Name | MALB | Optical | | | |
| Auto Rerun SW | | | | | | | |
| OFF | | | | | | | |
| Auto Rerun Range (Conc.) | | | | | | | |
| | First Dil | Low | | High | | | |
| | | Re | Value | Dil | Re | Value | Dil |
| Serum | | | | | | | |
| Urine | | | | | | | |
| Plasma | | | | | | | |
| CSF | | | | | | | |
| Dialysis | | | | | | | |
| Other | | | | | | | |
| Auto Rerun Condition (Absorbance) | | | | | | | |
| Lower | | | | | | | |
| Higher | | | | | | | |
| OFF | | | | | | | |
| Auto Rerun Condition (Prozone) | | | | | | | |
| OFF | | | | | | | |
| Dilution | | | | | | | |
| 100:Dil 2 | | | | | | | |

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