PRESTIGE 24i CALCIUM ARSENAZO

Nr kat. 4-247, 4-447 (PL)

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Cat. No **4-247, 4-447** (EN)

INTENDED USE

Diagnostic kit for determination of calcium concentration used in automatic analysers Prestige 24i, Biolis 24i, Prestige 24i Premium, Biolis 24i Premium and Biolis 30i.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Calcium plays an essential role in many cell functions: intracellularly in muscle contraction and glycogen metabolism, extracellularly, in bone mineralization, in blood coagulation and in transmission of nerve impulses. Calcium is present in plasma in three forms: free, bound to proteins or complexed with anions as phosphate, citrate and bicarbonate. Decreased total calcium levels can be associated with diseases of the bone apparatus (especially osteoporosis), kidney diseases (especially under dialysis), defective intestinal absorption and hypoparathyroidism. Increased total calcium can be measured in hyperparathyroidism, malignant diseases with metastases and sarcoidosis. Calcium measurements also help in monitoring of calcium supplementation mainly in the prevention of osteoporosis.

METHOD PRINCIPLE

Photometric test using arsenazo III. Calcium with arsenazo III at neutral pH yields a blue colored complex, whose intensity is proportional to the calcium concentration. Interference by magnesium is eliminated by addition of 8-hydroxyquinoline-5-sulfonic acid.

REAGENTS

Package

| | Cat. No 4-247 (24-TRAY) | Cat. No 4-447 (36-TRAY) |
|-----------|----------------------------|----------------------------|
| 1-Reagent | 4 x 52 ml | 6 x 22 ml |

The reagents when stored at 2-8°C are stable up to expiry date printed on the package. The reagent is stable for 12 weeks on board the analyser at 2-10°C (Biolis 30i).

Concentrations in the test

| | |
|------------------------------------|--------------|
| phosphate buffer (pH 7.5) | ≤ 80 mmol/l |
| 8-Hydroxyquinoline-5-sulfonic acid | ≤ 6 mmol/l |
| arsenazo III | ≤ 120 μmol/l |
| detergents, preservative | |

Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Do not freeze the reagent.
- Contaminated glassware is the greatest source of error. The use of disposable plastic ware is recommended. Glassware should be soaked for a few hours in 2M HCl solution and then thoroughly rinsed with distilled water.

SPECIMEN

Serum. Random or 24-hours urine.

Serum can be stored up to 7 days at 20-25°C or up to 3 weeks at 4-8°C. Samples frozen at -20°C can be stored up to 8 months. Discard contaminated specimens.

24-hours urine preparation: To prevent calcium salt precipitation specimens should be collected in 10 ml of 6M HCl. In case of presence of precipitants they can be solved by lowering pH of the urine to below 2.0. Prior to determination dilute the sample with 0.9% NaCl in the ratio of 1 to 1. Multiply the result by the dilution factor. Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-Reagent is ready to use.

1-Reagent put on basic position in reagent tray.

For reagent blank deionized water is recommended.

Actions required:

- Biolis 24i Premium:** When performing assays in the analyser, there is a probability of **cross-contamination** affecting the tests results: CALCIUM ARSENAZO - LIPASE, RF - CALCIUM ARSENAZO. To avoid this effect follow the recommendations contained in the instruction 51_03_24_008_BIOLIS_24i_PREMIUM_CARRYOVER.
- Biolis 30i:** When performing assays in the analyser, there is a probability of **cross-contamination** affecting the tests results: CALCIUM ARSENAZO - UIBC, CALCIUM ARSENAZO - PHOSPHORUS, CALCIUM ARSENAZO - LIPASE, RF - CALCIUM ARSENAZO, TG mono - CALCIUM ARSENAZO. To avoid this effect follow the recommendations contained in the instruction: 51_03_24_009_BIOLIS_30i_CARRYOVER.

REFERENCE VALUES ⁴

| serum | mg/dl | mmol/l |
|----------------|------------|--------------|
| adult | 8.6 – 10.3 | 2.15 – 2.57 |
| random urine | mg/dl | mmol/l |
| male | 0.9 – 37.9 | 0.225 – 9.47 |
| female | 0.5 – 35.7 | 0.125 – 8.92 |
| 24-hours urine | mg/24h | mmol/24h |
| adult | 100 – 300 | 2.5 – 7.5 |

It is recommended for each laboratory to establish its own reference ranges for local population.

Calcium concentration 24-hours urine - calculation

| | | | | |
|--|---|---|---|---|
| calcium concentration in 24-hours urine [mg/24h] | = | calcium concentration in sample of 24-hours urine [mg/dl] | × | urine volume of 24-hours urine [dl/24h] |
|--|---|---|---|---|

QUALITY CONTROL

For internal quality control it is recommended to use, with each batch of samples, the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173) for determination in serum or CORMAY URINE CONTROL LEVEL 1 (Cat. No 5-161) and LEVEL 2 (Cat. No 5-162) for determination in urine.

For the calibration of automatic analysers systems the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) are recommended.

The calibration curve should be prepared every 12 weeks (Biolis 30i), with every change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using the automatic analyser Biolis 30i. Results may vary if a different instrument or a manual procedure is used.

- LoB (Limit of Blank):**
0.07 mg/dl (0.0175 mmol/l)
- LoD (Limit of Detection):**
0.13 mg/dl (0.0325 mmol/l)
- LoQ (Limit of Quantitation):**
0.34 mg/dl (0.085 mmol/l)
- Linearity:**
up to 23 mg/dl (5.75 mmol/l)

For higher concentration, dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

Specificity / Interferences

Haemoglobin up to 5 g/dl, ascorbate up to 62 mg/l, bilirubin up to 40 mg/dl, triglycerides up to 2000 mg/dl and magnesium up to 20 mg/dl do not interfere with the test.

Precision

| Repeatability (run to run) n = 20 | Mean | SD | CV |
|--|---------|---------|------|
| | [mg/dl] | [mg/dl] | [%] |
| level 1 | 9.94 | 0.12 | 1.17 |
| level 2 | 12.9 | 0.12 | 0.92 |
| Reproducibility (day to day) n = 80 | Mean | SD | CV |
| | [mg/dl] | [mg/dl] | [%] |
| level 1 | 10.0 | 0.28 | 2.8 |
| level 2 | 12.8 | 0.26 | 2.1 |

Method comparison

A comparison between calcium values determined at **Biolis 30i** (y) and at **BECKMAN COULTER AU680 0** (x) using 60 serum samples gave following results:

$$y = 0.9629x + 0.6479 \text{ mg/dl};$$

$$R = 0.997 \quad (R - \text{correlation coefficient})$$

A comparison between calcium values determined at **Biolis 30i** (y) and at **BECKMAN COULTER AU680** (x) using 72 urine samples gave following results:

$$y = 0.9016x + 1.14 \text{ mg/dl};$$

R = 0.992

(R – correlation coefficient)

WASTE MANAGEMENT

Please refer to local legal requirements.

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| | | |
|--|--------|--------|
| Item name | 57 | CaArs |
| Auto Rerun SW | | |
| ON | | |
| Auto Rerun Range (Result) | | |
| | ON | ON |
| | Lower | Higher |
| Serum | 0.62 | 23 |
| Urine | | |
| Plasma | | |
| CSF | | |
| Dialysis | | |
| Other | | |
| Auto Rerun Condition (Absorbance) | | |
| Absorbance Range | Lower | OFF |
| | Higher | OFF |
| Auto Rerun Condition (Prozone) | | |
| Prozone Range | OFF | |

• Prestige 24i Premium, Biolis 24i Premium

| | | | | |
|-------------------------|------------|-----------|-----------|---------|
| Item No. | 57 | Item Name | CaArs | Optical |
| Data information | | | | |
| Units | mg/dl | | | |
| Decimals | 2 | | | |
| Calibration | | | | |
| Type | Linear2 | | | |
| Std sample conc. | | | | |
| Blank | 0 | #1 | * | #2 |
| #3 | | #4 | | #5 |
| #6 | | | | |
| Analysis | | | | |
| Type | END method | | | |
| Main Wave Length | 660nm | | | |
| Sub Wave Length | 700nm | | | |
| Method | Arsenazo | | | |
| Correlation | | | | |
| Y= | Slope | X+ | Intercept | |
| | 1 | | 0 | |

| | | | | |
|-------------------------|--------|-----------|-----------|---------|
| Item No. | 57 | Item Name | CaArs | Optical |
| Aspiration | | | | |
| Kind | Single | | | |
| Vol. | | | | |
| | Kind | Vol. | Add | Units |
| Sample | | 3 | 5 | µl |
| Reagent 1 | | 300 | 10 | µl |
| Data Process | | | | |
| Read | Main | Start | End | |
| | Sub | 9 | 10 | |
| Abs.Limit | | | | |
| Low | -3 | | | High |
| | | | | 3 |
| Correction value | | | | |
| Blank correction | | | | |
| End Point Limit | 2 | | | |
| Linear Check (%) | 90 | | | |
| Prozone Check | | | | |
| | Start | End | Limit (%) | |
| First | | | | |
| Second | | | | Low |
| Reaction Monitor | | | | |
| 0 Level Point | 1 | | | |
| Span | 3 | | | |
| Third mixing | | | | |
| OFF | | | | |

| | | | | |
|---------------------|------|-----------|--------|---------|
| Item No. | 57 | Item Name | CaArs | Optical |
| Normal Range | | | | |
| | Male | | Female | |
| | Low | High | Low | High |
| Serum | 8.6 | 10.3 | 8.6 | 10.3 |
| Urine | 0.9 | 37.9 | 0.5 | 35.7 |
| Plasma | | | | |
| CSF | | | | |
| Dialysis | | | | |
| Other | | | | |
| Panic Range | | | | |
| | Male | | Female | |
| | Low | High | Low | High |
| Serum | | | | |
| Urine | | | | |
| Plasma | | | | |
| CSF | | | | |
| Dialysis | | | | |
| Other | | | | |

| | | | | |
|--|-----------|-----------|-------|---------|
| Item No. | 57 | Item Name | CaArs | Optical |
| Auto Rerun SW | | | | |
| ON | | | | |
| Auto Rerun Range (Conc.) | | | | |
| | First Dil | Low | High | |
| | Re | Value | Dil | Re |
| Serum | | 0.3 | | 23 |
| Urine | | | | |
| Plasma | | | | |
| CSF | | | | |
| Dialysis | | | | |
| Other | | | | |
| Auto Rerun Condition (Absorbance) | | | | |
| | Lower | OFF | | |
| | Higher | OFF | | |
| Auto Rerun Condition (Prozone) | | | | |
| | OFF | | | |
| Dilution | | | | |
| | 100:Dil12 | | | |

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