PRESTIGE 24i TOTAL IgE

Nr kat. 4-265, 4-484 (PL)

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Cat. No **4-265, 4-484** (EN)

INTENDED USE

Diagnostic kit for determination of IgE levels intended to use in automatic analysers Prestige 24i, Biolis 24i, Prestige 24i Premium, Biolis 24i Premium and Biolis 30i

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

IgE is an immunoglobulin with a molecular weight of approximately 190 kD normally present in the blood in trace amounts. Continual production of IgE antibodies in response to common naturally occurring allergens, however, often results in elevated serum levels and in the development of such clinically important Type I allergic reactions as asthma, hay fever, dermatitis and food allergies. Elevated IgE levels are also seen in parasitic (helminth) diseases, IgE myeloma, and in hepatitis. The measurement of IgE in human serum is thus considered to be useful in the diagnosis, treatment, assessment of disease progression, or postoperative prognosis for such conditions.

METHOD PRINCIPLE

When an antigen-antibody reaction occurs between IgE in a sample and anti-IgE antibody which has been sensitized to latex particles, agglutination results. This agglutination is detected as an absorbance change (572 nm), with the magnitude of the change being proportional to the quantity of IgE in the sample. The actual concentration is then determined by interpolation from a calibration curve prepared from calibrators of known concentration.

REAGENTS

Package

	Cat. No 4-265 (24-TRAY)	Cat. No 4-484 (36-TRAY)
1-Reagent	1 x 38 ml	2 x 23 ml
2-Reagent	1 x 20 ml	2 x 12.5 ml

The reagents when stored at 2-10°C are stable up to expiry date printed on the package. The reagents are stable for 5 weeks on board the analyser at 2-10°C.

Concentrations in the test

suspension of latex particles sensitized with (mouse) anti-IgE antibodies (pH 7.3) 0.125 w/v%
 glycine buffer solution (pH 8.3)
 preservative

Warnings and notes

- Protect from direct sunlight and avoid contamination!
- After measurements are taken, reagent bottles should be capped and kept at 2-10°C. Care should be taken not to interchange the caps of reagent bottles.

- Reagents with different lot numbers should not be interchanged or mixed.
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.

SPECIMEN

Serum or plasma (Na-EDTA, K-EDTA, Na-Heparin, Li-Heparin, citric acid).

If the test cannot be done immediately, the sample should be placed in a tightly sealable container and stored at -20°C. Repeated freezing and thawing should be avoided. Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

- 1-Reagent and 2-Reagent are ready to use.
- 1-Reagent put on basic position in reagent tray.
- 2-Reagent put on start position in reagent tray.
- For reagent blank 0.9% NaCl is recommended.

REFERENCE VALUES ³

serum, plasma	< 358 IU/ml
---------------	-------------

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use, with each batch of samples the CORMAY IMMUNO-CONTROL II (Cat. No 4-290).

For the calibration of automatic analysers systems the CORMAY IgE CALIBRATORS kit (Cat. No 4-280) is recommended.

The calibration curve should be prepared every 5 weeks, with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using an automatic analyser Biolis 30i. Results may vary if a different instrument is used.

- LoB (Limit of Blank):** 4.5 IU/ml
- LoD (Limit of Detection):** 9.0 IU/ml
- LoQ (Limit of Quantitation):** 25 IU/ml
- Linearity:** up to 1000 IU/ml

For higher concentration of IgE dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

Specificity / Interferences

Haemoglobin up to 0.5 g/dl, bilirubin up to 30 mg/dl, triglycerides up to 1500 mg/dl, RF up to 500 IU/ml do not interfere with the test.

Precision

Repeatability (run to run) n = 20	Mean [IU/ml]	SD [IU/ml]	CV [%]
level 1	71.8	2.70	3.76
level 2	393.8	4.56	1.16
Reproducibility (day to day) n=80	Mean [IU/ml]	SD [IU/ml]	CV [%]
level 1	72.0	4.81	6.7
level 2	406.1	9.76	2.4

Method comparison

A comparison between IgE levels determined at **Biolis 30i** (y) and at **BS-800** (x) using 60 serum samples gave following results:

$$y = 0.98x + 5.2362 \text{ IU/ml;}$$

$$R = 0.997 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT

Please refer to local legal requirements.

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Item name	46	IgE	
Auto Rerun SW			
OFF			
Auto Rerun Range (Result)			
	OFF	OFF	
	Lower	Higher	
Serum			
Urine			
Plasma			
CSF			
Dialysis			
Other			
Auto Rerun Condition (Absorbance)			
Absorbance Range		Lower	OFF
		Higher	OFF
Prozone Range		OFF	

• Prestige 24i Premium, Biolis 24i Premium

Item No.	46	Item Name	IgE	Optical	
Data information					
Units	IU/ml				
Decimals	1				
Analysis					
Type	RATE method				
Main Wave Length	570 nm				
Sub Wave Length					
Method	Immuno				
Correlation					
Y=	Slope	X+	Intercept		
	1		0		
Calibration					
Type	Spline1				
Std sample conc.					
Blank	0	#1	*	#2	*
#3	*	#4	*	#5	*
#6					

Item No.	46	Item Name	IgE	Optical
Aspiration				
Kind	Double			
Data Process				
Read		Start	End	
		Main	36	46
		Sub		
Abs.Limit		Low	High	
		-3	3	
Blank value				
Water Blank				
Reaction Monitor				
0 Level Point	1			
Span	3			
Correction value				
Blank correction		1.00		
End Point Limit		2		
Linear Check (%)		90		
Prozone Check				
	Start	End	Limit (%)	
First				
Second				Low

Item No.	46	Item Name	IgE	Optical
Normal Range				
	Male		Female	
	Low	High	Low	High
Serum	0	358	0	358
Urine				
Plasma				
CSF				
Dialysis				
Other				
Panic Range				
	Male		Female	
	Low	High	Low	High
Serum				
Urine				
Plasma				
CSF				
Dialysis				
Other				

Item No.	46	Item Name	IgE	Optical			
Auto Rerun SW							
OFF							
Auto Rerun Range (Conc.)							
	First Dil	Low		High			
		Re	Value	Dil	Re	Value	Dil
Serum							
Urine							
Plasma							
CSF							
Dialysis							
Other							
Auto Rerun Condition (Absorbance)							
		Lower	OFF				
		Higher	OFF				
Auto Rerun Condition (Prozone)							
OFF							
Dilution							
99:Di1 1							

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