PRESTIGE 24i TRANSFERRIN

Nr kat. 4-274, 4-333 (PL)

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Cat. No **4-274, 4-333** (EN)

INTENDED USE

Diagnostic kit for determination of transferrin concentration intended to use in automatic analysers Prestige 24i, Biolis 24i, Prestige 24i Premium, Biolis 24i Premium and Biolis 30i.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Transferrin (siderophilin) is a glycoprotein synthesised in the liver. Transferrin is the major plasma transport protein for iron. Its concentration correlates with the total iron-binding capacity (TIBC).

Evaluation plasma/serum transferrin levels is useful for the differential diagnosis of anaemia and for monitoring its treatment. Causes of decreased synthesis and low plasma levels may be chronic liver disease and malnutrition. High levels of transferrin occur in pregnancy and during estrogen administration.

METHOD PRINCIPLE

Transferrin, present in a sample, form with the specific antibody an immunological complex. The increase of turbidity after the addition of antiserum is proportional to transferrin concentration in the sample.

REAGENTS

Package

	Cat. No 4-274 (24-TRAY)	Cat. No 4-333 (36-TRAY)
1-Reagent	1 x 42 ml	2 x 24 ml
2-Reagent	1 x 7.5 ml	1 x 7.5 ml

Buffer (1-Reagent) stored at 2-25°C and antiserum (2-Reagent) stored at 2-8°C are stable until expiry date printed on the package.

Reagent components

1-Reagent: imidazole buffer, polymer, inorganic salt, buffer, preservative.

2-Reagent: anti-human transferrin antiserum, inorganic salt, buffer, preservative.

Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Store closed.
- Do not freeze the reagents.
- Nanoparticle-based reagents can settle over time. It may be necessary to delicately mix by repeated turning.
- Human-origin products have been tested for HBsAg and antibodies to HIV 1, HIV 2, HCV and found to be non-reactive. However this material should be handled as potentially infectious.
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.

SPECIMEN ⁶

Serum or plasma. Nonhemolyzed, nonlipaemic, fresh serum is recommended.

Serum can be stored up to 3 days at 4°C or up to 6 months at -20°C.

Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-Reagent and 2-Reagent are ready to use.

1-Reagent put on basic position in reagent tray.

2-Reagent put on start position in reagent tray

0.9% NaCl is recommended as a reagent blank.

REFERENCE VALUES ⁶

newborns	1.30 – 2.75 g/l
children (3 months – 16 years)	2.03 – 3.60 g/l
adults F	2.50 – 3.80 g/l
M	2.15 – 3.65 g/l

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use with each batch of samples the CORMAY IMMUNO-CONTROL III (Cat. No 4-291).

For the calibration of automatic analysers the CORMAY IMMUNO-MULTICAL (Cat. No 4-287) is recommended.

The calibration curve should be prepared with every change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using the automatic analysers Hitachi or Biolis 30i. Results may vary if a different instrument is used.

- LOQ (Limit of Quantitation):** 0.076 g/l

- Linearity:** up to 8 g/l

- Specificity / Interferences**

Haemoglobin up to 0.32 g/dl, bilirubin up to 29.5 mg/dl, triglycerides up to 312 mg/dl, heparin up to 0.5 g/l, sodium fluoride up to 4 g/l, EDTA up to 5 g/l, sodium citrate up to 5 g/l do not interfere with the test.

- Precision**

Repeatability (run to run) n = 30	Mean [g/l]	SD [g/l]	CV [%]
level 1	1.44	0.02	1.12
level 2	2.52	0.04	1.52
level 3	3.47	0.06	1.83

Reproducibility (day to day) n = 66	Mean [g/l]	SD [g/l]	CV [%]
level 1	1.49	0.11	4.56
level 2	2.51	0.12	4.06
level 3	3.51	0.12	1.64

- Method comparison**

A comparison between transferrin values determined at **Biolis 30i** (y) and at **ADVIA SIEMENS 1800** (x) using 70 serum samples, gave following results:

$$y = 1.0081x + 0.0898 \text{ g/l};$$

$$R = 0.979 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT

Please refer to local legal requirements.

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Item name	33	TRF
Auto Rerun SW		
OFF		
Auto Rerun Range (Result)		
	OFF	OFF
	Lower	Higher
Serum		
Urine		
Plasma		
CSF		
Dialysis		
Other		
Auto Rerun Condition (Absorbance)		
Absorbance Range	Lower	OFF
	Higher	OFF
Auto Rerun Condition (Prozone)		
Prozone Range	OFF	

• **Prestige 24i Premium, Biolis 24i Premium**

Item No.	33	Item Name	TRF	Optical
Data information				
Units	g/l			
Decimals	2			
Calibration				
Type	Logit2			
Std sample conc.				
Blank	0	#1 *	#2 *	*
#3	*	#4 *	#5 *	*
#6				
Analysis				
Type	END method			
Main Wave Length	340 nm			
Sub Wave Length	700 nm			
Method	Immuno			
Correlation				
Slope	Intercept			
Y=	1	X+	0	

Item No.	33	Item Name	TRF	Optical
Aspiration				
Kind	Double			
Vol.				
	Kind	Vol.	Add	Units
Sample		2.5	5	µl
Reagent 1		315	10	µl
Reagent 2		35	10	µl
Data Process				
Read	Start	End		
	Main	51	52	
	Sub	21	22	
Abs.Limit	Low	High		
	-3	3		
Blank value				
Water Blank				
Reaction Monitor				
0 Level Point	1			
Span	3			
Third mixing				
ON				
Correction value				
Blank correction	0.9007			
End Point Limit	2			
Linear Check (%)	0			
Prozone Check				
	Start	End	Limit (%)	
First				
Second				Low

Item No.	33	Item Name	TRF	Optical
Normal Range				
	Male		Female	
	Low	High	Low	High
Serum	2	4	2	4
Urine				
Plasma				
CSF				
Dialysis				
Other				
Panic Range				
	Male		Female	
	Low	High	Low	High
Serum				
Urine				
Plasma				
CSF				
Dialysis				
Other				

Item No.	33	Item Name	TRF	Optical			
Auto Rerun SW							
OFF							
Auto Rerun Condition (Absorbance)							
	Lower	OFF					
	Higher	OFF					
Auto Rerun Range (Conc.)							
	First Dil	Low		High			
		Re	Value	Dil	Re	Value	Dil
Serum							
Urine							
Plasma							
CSF							
Dialysis							
Other							
Auto Rerun Condition (Prozone)							
OFF							
Dilution							
99:Di11							

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