## CORMAY CRP ULTRA CALIBRATORS

Nr kat. 4-276 (PL)  
CORMAY CRP ULTRA CALIBRATORS 5 x 2 ml

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Cat. No 4-276 (EN)  
CORMAY CRP ULTRA CALIBRATOR 5 x 2 ml

### INTENDED USE

CORMAY CRP ULTRA CALIBRATORS kit is dedicated for calibration of C-reactive protein determinations on automatic analysers.

The calibrators must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

### INTRODUCTION

CRP (C-reactive protein) is an acute phase protein whose concentration is seen to increase as a result of the inflammatory process, most notably in response to pneumococcal (bacterial) infectious, histolytic disease and a variety of disease states. CRP to be used as a marker or general diagnostic indicator of infections and inflammation, in addition to serving as a monitor of patient response to therapy and surgery. Furthermore, regular measurements of CRP in infants can be a useful aid in the early diagnosis of infectious disease.

CORMAY CRP ULTRA CALIBRATORS kit is prepared by diluting C-reactive protein with normal human serum at various concentrations. Concentrations are shown on each vial and in this kit insert.

### PACKAGE

Cap colour	CRP concentration	Size
WHITE	0,5 mg/dl	1 x 2 ml
YELLOW-GREEN	2 mg/dl	1 x 2 ml
BLUE	4 mg/dl	1 x 2 ml
DARK RED	16 mg/dl	1 x 2 ml
DARK BROWN	32 mg/dl	1 x 2 ml

### PREPARATION AND STORAGE

The calibrators are ready to use.

Unopened calibrator is stable at 2-10°C until the expiration date printed on the vial.

After first opening calibrator is stable for 12 weeks at 2-10°C.

### WARNINGS AND NOTES

- Protect from direct sunlight and avoid contamination!
- As the colours of caps vary according to CRP concentrations in the vial, care should be taken not to interchange the caps.
- The calibrators should not be mixed with ones in the other vial.
- The calibrator should be used according to this kit insert.
- The calibrator is based on human serum that has been tested for HBsAg and antibodies to HIV and HCV and found to be non-reactive. However this material should be handled as though capable of transmitting infectious disease.
- Product contain sodium azide (< 0.1%) as a preservative. Avoid contact with skin and mucous membranes.

### TRACEABILITY

CRP calibrators are prepared referring to the ERM-DA472 reference material.

### WASTE MANAGEMENT

Please refer to local legal requirements.

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