CORMAY AFP CALIBRATORS

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INTENDED USE

Calibrators for determination of α -fetoprotein concentration. The calibrators must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

α -fetoprotein (AFP) is fetoprotein with a molecular weight of approximately 70 kD containing about 3 % sugar. While it is present in high concentrations during fetal growth, its concentration rapidly decrease after birth and is present at extremely low levels in normal human blood.

AFP shows a notable increase in primary hepatic cancer and is considered to be of great diagnostic importance. It is also thought that fluctuations in blood AFP are useful for evaluating the progress, effects of therapy, and postoperative prognosis of hepatoma.

CORMAY AFP CALIBRATORS kit is prepared by adding α -fetoprotein to pooled human serum at various concentrations. Concentrations are shown on each vial and in this kit insert.

CORMAY AFP CALIBRATORS kit is dedicated for calibration of α -fetoprotein determinations on automatic analysers.

PACKAGE

Calibrator	Size
0 ng/ml	1 x 1ml
25 ng/ml	1 x 1ml
100 ng/ml	1 x 1ml
250 ng/ml	1 x 1ml

The calibrators are ready to use.

STABILITY AND STORAGE

Unopened calibrator is stable at 2-10°C until the expiration date printed on the vial.

WARNINGS AND NOTES

- Protect from light.
- Do not change the caps between vials.
- The calibrators should not be mixed with ones in the other vial.
- The calibrators should be used according to this kit insert.
- The calibrators are based on human serum that has been tested for HBsAg and antibodies to HIV and HCV and found to be non-reactive. However, this material should be handled as though capable of transmitting infectious disease.
- The calibrators contain sodium azide (< 0.1%) as a preservative. Avoid contact with skin and mucous membranes.

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