## CORMAY IMMUNO- CONTROL III

(PL)

**CORMAY IMMUNO- CONTROL III**  
Nr kat. 4-291 1 x 2 ml

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### LITERATURA

1. Johnson, A.M.: A new international preparation for proteins in human serum, Arch. Pathol. Lab. Med. 117 (1993) 29-31.

**Data wydania:** 11. 2020.

## CORMAY IMMUNO- CONTROL III

(EN)

**CORMAY IMMUNO- CONTROL III**  
Cat. No 4-291 1 x 2 ml

### INTENDED USE

CORMAY IMMUNO-CONTROL III is designed for use as a multianalyte control sera when measuring various serum proteins in human serum. Use in accordance with instructions enclosed to reagents for the performance of particular proteins.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

### CHARACTERISTICS

The CORMAY IMMUNO-CONTROL III is a human biological fluid diluted in HEPES pH 7.4 buffer containing stabilisers, sodium azide (<1 g/l) and the following human proteins: alpha-1 antitrypsin, alpha-1 glycoprotein acid, antithrombin III, complement C3, complement C4, ceruloplasmin, haptoglobin, IgA, IgG, IgM and transferrin.

### PREPARATION AND STORAGE

Serum is ready for use.

Serum stored at 2-8°C is stable until expiry date printed on the package. Store closed and avoid contamination.

### NOTES

1. Control sera must be used only for the purpose intended by suitably qualified laboratory personnel, under appropriate laboratory conditions.
2. This serum has been tested for the HIV antibody, HBsAg and HCV and found to be non-reactive. However this material should be handled as thought capable of transmitting infectious disease.
3. Control sera contain sodium azide (< 0,1%) as a preservative. Avoid contact with skin and mucous membranes.

### ASSIGNED VALUES

Assigned values for particular proteins vary from lot to lot are given in values sheet.

Values assigned from the reference material-CRM 470.

### WASTE MANAGEMENT

Please refer to local legal requirements.

### LITERATURE

1. Johnson, A.M.: A new international preparation for proteins in human serum, Arch. Pathol. Lab. Med. 117 (1993) 29-31.

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