## PRESTIGE 24i HbA1c DIRECT

Nr kat **4-300, 4-301** (PL)

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Cat. No **4-300, 4-301** (EN)

### INTENDED USE

Diagnostic kit for determination of haemoglobin A1c concentration intended to use in automatic analyzers: Biolis 24i Premium and Biolis 30i.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

### INTRODUCTION

The determination of HbA1c is most commonly performed for the evaluation of glycemic control in diabetes mellitus. HbA1c values provide an indication of glucose levels over the preceding 4-8 weeks.

Throughout the circulatory life of the red cell, hemoglobin A1c is formed continuously by the adduction of glucose to the N-terminal of the hemoglobin beta chain. This process, which is non-enzymatic, reflects the average exposure of hemoglobin to glucose over an extended period. In a classical study, Trivelli et al [1] showed hemoglobin Hb in diabetic subjects to be elevated 2-3 fold over the levels found in normal individuals. Several investigators have recommended that hemoglobin A1c serve as an indicator of metabolic control of the diabetic, since hemoglobin A1c levels approach normal values for diabetics in metabolic control [2,3,4].

Hemoglobin A1c has been defined operationally as the "fast fraction" hemoglobins (HbA1a, HbA1b, HbA1c) that elute first during column chromatography with cation-exchange resins. The non-glycosylated hemoglobin, which consists of the bulk of the hemoglobin has been designated HbA0.

### METHOD PRINCIPLE

Method for hemoglobin A1c determination complies with standardized method certified by the National Glycohemoglobin Standardization Program (NGSP).

The present method utilizes the interaction of antigen and antibody to directly determine the HbA1c concentration in whole blood.

Total hemoglobin and HbA1c have the same unspecific absorption rate to latex particles. When mouse antihuman HbA1c monoclonal antibody is added, latex-HbA1c-mouse anti human HbA1c antibody complex is formed. Agglutination is formed when goat anti-mouse IgG polyclonal antibody interacts with the monoclonal antibody. The amount of agglutination is proportional to the amount of HbA1c absorbed on to the surface of latex particles.

The amount of agglutination is measured as absorbance.

The HbA1c value is obtained from a calibration curve.

### REAGENTS

#### Package

	Cat. No <b>4-300</b> (24-tray)	Cat. No <b>4-301</b> (36-tray)
1-REAGENT	1 x 26 ml	2 x 16 ml
2-REAGENT	1 x 10 ml	2 x 6.3 ml
HEMOLYSING REAGENT	2 x 34 ml	2 x 42 ml

The reagents (1-REAGENT, 2-REAGENT) stored at 2-8°C and HEMOLYSING REAGENT stored at 2-25°C are stable until expiry date printed on the package. The reagents are stable for 12 weeks on board the analyser at 2-10°C (Biolis 24i Premium, Biolis 30i).

### Concentrations in the test

latex	0.13%
mouse anti-human HbA1c monoclonal antibody	0.05 mg/ml
goat anti-mouse IgG polyclonal antibody stabilizers	0.08 mg/dl
buffer	

### Warnings and notes

- Protect from direct sunlight and avoid contamination!
- It has been reported that results may be inconsistent in patients who have the following conditions: opiate addiction, lead-poisoning, alcoholism, ingest large doses of aspirin [6, 7, 8, 9].
- It has been reported that elevated levels of HbF may lead to underestimation of HbA1c and, that uremia does not interfere with HbA1c determination by immunoassay [10].
- This assay should not be used for the diagnosis of diabetes mellitus, but for monitor diabetic patients.
- In using Hemoglobin A1c to monitor diabetic patients, results should be interpreted individually.
- Reagent HEMOLYSING REAGENT (Cat. No 4-398) can be ordered separately.
- Any clinical case with shortened erythrocyte survival (e.g. hemolytic anemia, blood loss, pregnancy) might cause decrease in HbA1c values.

### SPECIMEN

Venous blood collected with EDTA.

Hemoglobin A1c in whole blood collected with EDTA is stable for 7 days at 2-8°C.

Nevertheless it is recommended to perform the assay with freshly collected samples!

#### Sample pretreatment:

- Dispense 500 µl HEMOLYSING REAGENT into tubes labeled: Control, Patients, etc.
- Place 10 µl of well mixed whole blood into the appropriately labeled lyse reagent tube. Mix well and allow to stand for minimum 5 minutes, until complete lysis is evident. Next mix sample for 5 minutes.
- The treated sample may be stored up to 10 days at 2-8°C. Mix sample again for 5 minutes before measurement.
- Note:** calibrators and controls should be also hemolysed according to sample pretreatment.

### PROCEDURE

1-REAGENT and 2-REAGENT are ready to use.

1-REAGENT put on basic position in reagent tray.

2-REAGENT put on start position in reagent tray.

For reagent blank 0.9% NaCl is recommended.

Application guide for Prestige 24i analyser is available on request.

The application guide should be verified and validated by the user prior to testing patient samples.

Test result is read automatically and the value is reported in % of hemoglobin unit in accordance with NGSP standardization. In order to convert the result reported in % haemoglobin (NGSP) to value reported in SI units mmol/mol in accordance with IFCC standardization, the following master equation should be used:

$$\text{HbA1c [mmol/mol IFCC]} = (\text{HbA1c [\% NGSP]} - 2.15) \times 10.929$$

### REFERENCE VALUES <sup>11</sup>

Non-diabetes < 6%

Patients with diabetes, control of glycaemia < 7%

It is recommended for each laboratory to establish its own reference ranges for local population.

### QUALITY CONTROL

For internal quality control it is recommended to use, with each batch of samples, the CORMAY HbA1c DIRECT CONTROLS (Cat. No 4-328).

For the calibration of automatic analysers systems the CORMAY HbA1c DIRECT CALIBRATORS (Cat. No 4-308) are recommended.

Controls and calibrators should be treated with HEMOLYSING REAGENT.

The calibration curve should be prepared every 4 weeks (Biolis 30i) or 12 weeks (Biolis 24i Premium), with change of reagent lot number or as required e.g. quality control findings outside the specified range.

### PERFORMANCE CHARACTERISTICS

The following results have been obtained using the automatic analysers Hitachi 717 and Biolis 30i. Results may vary if a different instrument is used.

- Analytical range:** 2 – 16% (up to 151 mmol/mol).

#### Specificity / Interferences

Bilirubin up to 50 mg/dl, triglycerides up to 2000 mg/dl, ascorbate up to 50 mg/dl, carbamylated Hb up to 7.5 mmol/l and acetylated Hb up to 5.0 mmol/l do not interfere with the test.

#### Precision

Repeatability (run to run) n = 20	Mean [%]	SD [%]	CV [%]
level 1	5.56	0.04	0.8
level 2	12.21	0.17	1.4
Reproducibility (day to day) n = 80	Mean [%]	SD [%]	CV [%]
level 1	6.48	0.09	1.3
level 2	12.42	0.25	2.0

#### Method comparison

A comparison between HbA1c values determined at **Biolis 30i** (y) and at **Bio-Rad Variant II Turbo** (x) using 89 whole blood samples gave following results:

$$y = 0.9081 x + 0.532 \%$$

$$R = 0.982 \quad (R - \text{correlation coefficient})$$

### WASTE MANAGEMENT

Please refer to local legal requirements.

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Item No.	65	Item Name	HbA1c D	Optical
<b>Aspiration</b>				
Kind	Double			
Vol.				
Kind	Vol.	Add	Units	
Sample	4	5	µl	
Reagent 1	180	10	µl	
Reagent 2	60	10	µl	
<b>Blank value</b>				
Water Blank				
<b>Reaction Monitor</b>				
0 Level Point	1			
Span	3			
<b>Third mixing</b>				
ON				
<b>Data Process</b>				
Read	Main	Start	End	
		50	52	
Sub				
<b>Abs.Limit Low High</b>				
	-3			3
<b>Correction value</b>				
Blank correction	1			
End Point Limit	2			
Linear Check (%)				
<b>Prozone Check</b>				
	Start	End	Limit (%)	
First				
Second				Low

Item No.	65	Item Name	HbA1c D	Optical
<b>Normal Range</b>				
	Male		Female	
	Low	High	Low	High
Serum				
Urine				
Plasma				
CSF				
Dialysis				
Other	#	#	#	#
<b>Panic Range</b>				
	Male		Female	
	Low	High	Low	High
Serum				
Urine				
Plasma				
CSF				
Dialysis				
Other	#	#	#	#

Item No.	65	Item Name	HbA1c D	Optical			
<b>Auto Rerun SW</b>							
OFF							
<b>Auto Rerun Range (Conc.)</b>							
	First Dil	Low		High			
		Re	Value	Dil	Re	Value	Dil
Serum							
Urine							
Plasma							
CSF							
Dialysis							
Other							
<b>Auto Rerun Condition (Absorbance)</b>							
Lower	OFF						
Higher	OFF						
<b>Auto Rerun Condition (Prozone)</b>							
OFF							
<b>Dilution</b>							
100:Dil2							

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Item No	65	Item Name	HbA1c	Specimen	OTHER	OPTICAL
<b>Reference intervals</b>						
MALE		FEMALE				
LOW	HIGH	LOW	HIGH			
#	#	#	#			
<b>Panic range</b>						
MALE		FEMALE				
LOW	HIGH	LOW	HIGH			
#	#	#	#			
<b>Decision limit</b>						
	MALE		FEMALE			
<b>Auto rerun</b>						
	°ON		•OFF			
<b>Auto rerun range (conc.)</b>						
Re	Value	Dil.	Re	Value	Dil.	
<b>Auto rerun condition (abs.)</b>						
	°ON		•OFF		DIL.	
LOWER						
HIGH						
<b>Auto rerun condition (prozone)</b>						
	°ON		•OFF			
<b>SAMPLE VOL.</b>						
<b>Dilution</b>						
	•DIL 1		° DIL 2			
<b>Reaction check</b>						
	°ON		•OFF			
	CHECK					
	LOW					
	HIGH					
<b>VL CHECK</b>						
°ON	•OFF		°ON		•OFF	

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