PRESTIGE 24i FIBRINOGEN

Nr kat. 4-313, 4-339 (PL)

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Cat. No **4-313, 4-339** (EN)

INTENDED USE

Diagnostic kit for determination of fibrinogen concentration intended to use in automatic analyzers: Prestige 24i, Biolis 24i, Prestige 24i Premium, Biolis 24i Premium and Biolis 30i.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Fibrinogen is the plasma protein precursor of fibrin, which when cross-linked becomes the principal component of fibrin clot.

Fibrinogen is a sensitive acute phase protein whose concentration rises several fold during inflammation. It can be increased as a response to inflammatory states, with infections, during pregnancy and after trauma.

Evidence has shown that plasma levels above the reference range constitute a significant independent risk factor for both coronary artery and cerebrovascular diseases.

Fibrinogen can be deficient in congenital afibrinogenemia. Levels may also fall significantly as the result of sequestration in extravascular spaces such as the body cavities and sites of recent trauma.

METHOD PRINCIPLE

Fibrinogen presents in a sample form with the specific antibody an immunological complex. The increase of turbidity after the addition of antiserum measured at $\lambda=340$ nm is proportional to fibrinogen concentration in the sample.

REAGENTS

Package

| | Cat. No 4-313 (24-TRAY) | Cat. No 4-339 (36-TRAY) |
|-----------|----------------------------|----------------------------|
| 1-Reagent | 1 x 40 ml | 2 x 23 ml |
| 2-Reagent | 1 x 10 ml | 2 x 6 ml |

Buffer (1-Reagent) stored at 2-8°C and antiserum (2-Reagent) stored at 2-8°C are stable until expiry date printed on the package.

Reagent components

1-Reagent: TRIS buffer, polymer, inorganic salt, preservative.

2-Reagent: anti-human fibrinogen antiserum, buffer, inorganic salt, preservative.

Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Store closed.
- Do not freeze the reagents.

- Nanoparticle-based reagents can settle over time. It may be necessary to delicately mix by repeated turning.
- Human-origin products have been tested for HBsAg and antibodies to HIV 1, HIV 2, HCV and found to be non-reactive. However this material should be handled as potentially infectious.
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.

SPECIMEN ⁵

Plasma collected on sodium citrate.

Before analysis calibrator and samples should be diluted 20-fold with 0.9% NaCl and mixed gently.

Samples may be stored up to 4 hours at room temperature, or up to 18 months, frozen at -70°C.

It is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-Reagent and 2-Reagent are ready to use.

1-Reagent put on basic position in reagent tray.

2-Reagent put on start position in reagent tray.

For reagent blank 0.9% NaCl is recommended.

REFERENCE VALUES ⁴

| | |
|----------|-----------------|
| adults | 200 – 400 mg/dl |
| children | 125 – 300 mg/dl |

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use with each batch of samples the control plasma. For the calibration of automatic analysers systems the FIBRINOGEN CALIBRATOR (Cat. No 4-292) is recommended.

The calibration curve should be prepared with every change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using an automatic analyser Hitachi or Biolis 30i. Results may vary if a different instrument is used.

- Analytical range:** 5,46 mg/dl to 1000 mg/dl.

- Interferences:**

Hemoglobin up to 0.13 g/dl, bilirubin up to 29.5 mg/dl, triglycerides up to 1053 mg/dl do not interfere with the test.

- Precision**

| Repeatability (run to run) n = 30 | Mean [mg/dl] | SD [mg/dl] | CV [%] |
|--------------------------------------|-----------------|---------------|-----------|
| level 1 | 284.7 | 4.83 | 1.7 |
| level 2 | 122.7 | 4.06 | 3.31 |
| level 3 | 80.2 | 2.97 | 3.71 |

| Reproducibility (day to day) n = 64 | Mean [mg/dl] | SD [mg/dl] | CV [%] |
|--|-----------------|---------------|-----------|
| level 1 | 299.6 | 9.1 | 2.9 |
| level 2 | 139.6 | 6.3 | 4.4 |
| level 3 | 91.4 | 5.5 | 5.2 |

- Method comparison**

A comparison between haptoglobin values determined at **Biolis 30i** (y) and at **BS-800** (x) using 60 plasma samples gave following results:

$$y = 1.0315x - 9.0873 \text{ g/l;}$$

$$R = 0.998 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT

Please refer to local legal requirements.

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| | | |
|--|---------------|--------|
| Item name | 38 | FIBR |
| Auto Rerun SW | | |
| OFF | | |
| Auto Rerun Condition (Absorbance) | | |
| Absorbance Range | Lower | OFF |
| | Higher | OFF |
| Auto Rerun Range (Result) | | |
| | OFF | OFF |
| | Lower | Higher |
| Auto Rerun Range (Conc.) | | |
| | Prozone Range | OFF |
| Serum | | |
| Urine | | |
| Plasma | | |
| CSF | | |
| Dialysis | | |
| Other | | |

• **Prestige 24i Premium, Biolis 24i Premium**

| | | | | |
|-------------------------|------------|-----------|------|---------|
| Item No. | 38 | Item Name | FIBR | Optical |
| Data information | | | | |
| Units | mg/dl | | | |
| Decimals | 1 | | | |
| Analysis | | | | |
| Type | END method | | | |
| Main Wave Length | 340nm | | | |
| Sub Wave Length | 700nm | | | |
| Method | Immuno | | | |
| Calibration | | | | |
| Type | Linear1 | | | |
| Std sample conc. | | | | |
| Blank | 0 | #1 | * | #2 |
| #3 | | #4 | | #5 |
| #6 | | | | |
| Correlation | | | | |
| Slope | Intercept | | | |
| Y= | 1 | X+ | 0 | |

| | | | | |
|-------------------------|--------|-----------|-----------|---------|
| Item No. | 38 | Item Name | FIBR | Optical |
| Aspiration | | | | |
| Kind | Double | | | |
| Vol. | | | | |
| Kind | Vol. | Add | Units | |
| Sample | 5 | 5 | µl | |
| Reagent 1 | 250 | 10 | µl | |
| Reagent 2 | 50 | 10 | µl | |
| Data Process | | | | |
| Read | Main | Start | End | |
| | Sub | 30 | 31 | |
| Abs.Limit Low High | | | | |
| | -3 | ~ | 3 | |
| Blank value | | | | |
| Water Blank | | | | |
| Reaction Monitor | | | | |
| 0 Level Point | 1 | | | |
| Span | 3 | | | |
| Correction value | | | | |
| Blank correction | | | | 0.8361 |
| End Point Limit | | | | 2 |
| Linear Check (%) | | | | 0 |
| Prozone Check | | | | |
| | Start | End | Limit (%) | |
| First | | | | |
| Second | | | | Low |

| | | | | |
|---------------------|------|-----------|--------|---------|
| Item No. | 38 | Item Name | FIBR | Optical |
| Normal Range | | | | |
| | Male | | Female | |
| | Low | High | Low | High |
| Serum | 200 | 400 | 200 | 400 |
| Urine | | | | |
| Plasma | | | | |
| CSF | | | | |
| Dialysis | | | | |
| Other | | | | |
| Panic Range | | | | |
| | Male | | Female | |
| | Low | High | Low | High |
| Serum | | | | |
| Urine | | | | |
| Plasma | | | | |
| CSF | | | | |
| Dialysis | | | | |
| Other | | | | |

| | | | | |
|--|-----------|-----------|-------|---------|
| Item No. | 38 | Item Name | FIBR | Optical |
| Auto Rerun SW | | | | |
| OFF | | | | |
| Auto Rerun Condition (Absorbance) | | | | |
| | Lower | OFF | | |
| | Higher | OFF | | |
| Auto Rerun Range (Conc.) | | | | |
| | First Dil | Low | High | |
| | | Re | Value | Dil |
| Serum | | | | |
| Urine | | | | |
| Plasma | | | | |
| CSF | | | | |
| Dialysis | | | | |
| Other | | | | |
| Auto Rerun Condition (Prozone) | | | | |
| OFF | | | | |
| Dilution | | | | |
| 99:Dil1 | | | | |

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