Prestige 24i ANTITHROMBIN III

Nr kat. 4-316, 4-342 (PL)

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Cat. No 4-316, 4-342 (EN)

INTENDED USE

Diagnostic kit for determination of antithrombin III concentration used in automatic analysers Prestige 24i, Biolis 24i, Prestige 24i Premium, Biolis 24i Premium and Biolis 30i.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Antithrombin III (AT III) is an α_2 -glycoprotein of MW 58000 and is made in the liver. AT III is one of the most important regulators of the coagulation system. AT III inactivates thrombin and factors Xa, IXa, XIa and XIIIa. This anticoagulant activity is enhanced by the presence of heparin, which forms a ternary complex with AT III and these procoagulant factors. Reduced concentration of AT III in blood means a great and well established risk for thrombotic complications. Clinical low value are associated with congenital or acquired deficiencies caused by a decreased biosynthesis (liver disease, medical treatment) or an increased loss (gastrointestinal disease, nephrotic syndrome) or an increased consumption (sepsis, major trauma due to surgery and burns wounds).

METHOD PRINCIPLE

The antithrombin III presents in a sample form with the specific antibody an immunological complex. The increase of turbidity after the addition of antiserum measured at $\lambda=340$ nm is proportional to antithrombin III concentration in the sample.

REAGENTS

Package

| | Cat. No 4-316 (24-TRAY) | Cat. No 4-342 (36-TRAY) |
|-----------|----------------------------|----------------------------|
| 1-Reagent | 1 x 40 ml | 2 x 23 ml |
| 2-Reagent | 1 x 10 ml | 2 x 6 ml |

Buffer (1-Reagent) stored at 2-25°C and antiserum (2-Reagent) stored at 2-8°C are stable until expiry date printed on the package.

Reagent components

1-Reagent: tricine buffer, polymer, inorganic salt, preservative.

2-Reagent: anti-human antithrombin III antiserum, buffer, inorganic salt, preservative.

Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Store closed.
- Do not freeze the reagents.
- Nanoparticle-based reagents can settle over time. It may be necessary to delicately mix by repeated turning.

- Human-origin products have been tested for HBsAg and antibodies to HIV 1, HIV 2, HCV and found to be non-reactive. However this material should be handled as potentially infectious.
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.

SPECIMEN ⁷

Citrated plasma.

Samples may be stored up to 4 hours on ice. If the test cannot be done immediately, the sample should be stored frozen.

Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-Reagent and 2-Reagent are ready to use.

1-Reagent put on basic position in reagent tray.

2-Reagent put on start position in reagent tray.

For reagent blank 0.9% NaCl is recommended.

REFERENCE VALUES ⁷

| | |
|--------|-----------------|
| plasma | 0.21 – 0.30 g/l |
|--------|-----------------|

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use, with each batch of samples, the CORMAY IMMUNO-CONTROL III (Cat. No 4-291).

For the calibration of automatic analysers systems the CORMAY IMMUNO-MULTICAL (Cat. No 4-287) is recommended.

The calibration curve should be prepared with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

These metrological characteristics have been obtained using the Hitachi or Biolis 30i automatic analysers. Results may vary if a different instrument is used.

- LoQ (Limit of Quantitation):** 0.028 g/l

- Linearity:** up to 0.806 g/l.

- Interferences:**

Hemoglobin up to 0.06 g/dl, bilirubin up to 22 mg/dl, triglycerides up to 155,5 mg/dl, heparin up to 0.5 g/l, sodium fluoride up to 4 g/l, EDTA up to 5 g/l, sodium citrate up to 5 g/l do not interfere with the test.

- Precision**

| Repeatability (run to run) n = 30 | Mean [g/l] | SD [g/l] | CV [%] |
|--------------------------------------|---------------|-------------|-----------|
| level 1 | 0.16 | 0.002 | 0.011 |
| level 2 | 0.32 | 0.002 | 0.005 |
| level 3 | 0.46 | 0.005 | 0.012 |

| Reproducibility (day to day) n = 56 | Mean [g/l] | SD [g/l] | CV [%] |
|--|---------------|-------------|-----------|
| level 1 | 0.18 | 0.005 | 0.03 |
| level 2 | 0.34 | 0.008 | 0.02 |
| level 3 | 0.50 | 0.010 | 0.02 |

- Method comparison**

A comparison between antithrombin III concentration determined at **Biolis 30i** (y) and at **BS-800** (x) using 59 plasma samples gave following results:

$y = 1.0358x - 0.0103$ g/l;

$R = 0.985$ (R – correlation coefficient)

WASTE MANAGEMENT

Please refer to local legal requirements.

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| | | | |
|--|-----------|--------------|--------------|
| Item name | 41 | AT 3 | |
| Auto Rerun SW | | | |
| OFF | | | |
| Auto Rerun Condition (Absorbance) | | | |
| Absorbance Range | Lower | OFF | |
| | Higher | OFF | |
| Auto Rerun Range (Result) | | | |
| | OFF | OFF | |
| | Lower | Higher | |
| Prozone Range | | | |
| | OFF | | |
| Auto Rerun Range (Conc.) | | | |
| | First Dil | Low | High |
| | | Re Value Dil | Re Value Dil |
| Serum | | | |
| Urine | | | |
| Plasma | | | |
| CSF | | | |
| Dialysis | | | |
| Other | | | |

• **Prestige 24i Premium, Biolis 24i Premium**

| | | | | |
|-------------------------|------------|-----------|------|---------|
| Item No. | 41 | Item Name | AT 3 | Optical |
| Data information | | | | |
| Units | g/L | | | |
| Decimals | 3 | | | |
| Calibration | | | | |
| Type | Logit2 | | | |
| Std sample conc. | | | | |
| Blank | 0 | #1 | * | #2 * |
| #3 | * | #4 | * | #5 * |
| #6 | | | | |
| Analysis | | | | |
| Type | END method | | | |
| Main Wave Length | 340 nm | | | |
| Sub Wave Length | 700 nm | | | |
| Method | Immuno | | | |
| Correlation | | | | |
| | Slope | Intercept | | |
| Y= | 1 | X+ | 0 | |

| | | | | |
|-------------------------|--------|-----------|-----------|---------|
| Item No. | 41 | Item Name | AT 3 | Optical |
| Aspiration | | | | |
| Kind | Double | | | |
| Vol. | | | | |
| | Kind | Vol. | Add | Units |
| Sample | | 10 | 5 | µl |
| Reagent 1 | | 250 | 10 | µl |
| Reagent 2 | | 50 | 10 | µl |
| Data Process | | | | |
| Read | Start | End | | |
| | Main | 51 | 52 | |
| | Sub | 30 | 31 | |
| Abs.Limit | | | | |
| Low | | | | High |
| | -3 | ~ | | 3 |
| Blank value | | | | |
| Water Blank | | | | |
| Reaction Monitor | | | | |
| 0 Level Point | 1 | | | |
| Span | 3 | | | |
| Correction value | | | | |
| Blank correction | 0.8387 | | | |
| End Point Limit | 2 | | | |
| Linear Check (%) | 0 | | | |
| Prozone Check | | | | |
| | Start | End | Limit (%) | |
| First | | | | |
| Second | | | Low | |
| Third mixing | | | | |
| ON | | | | |

| | | | | |
|---------------------|------|-----------|--------|---------|
| Item No. | 41 | Item Name | AT 3 | Optical |
| Normal Range | | | | |
| | Male | | Female | |
| | Low | High | Low | High |
| Serum | | | | |
| Urine | | | | |
| Plasma | 0.25 | 0.45 | 0.25 | 0.45 |
| CSF | | | | |
| Dialysis | | | | |
| Other | | | | |
| Panic Range | | | | |
| | Male | | Female | |
| | Low | High | Low | High |
| Serum | | | | |
| Urine | | | | |
| Plasma | | | | |
| CSF | | | | |
| Dialysis | | | | |
| Other | | | | |

| | | | | |
|--|-----------|--------------|--------------|---------|
| Item No. | 41 | Item Name | AT 3 | Optical |
| Auto Rerun SW | | | | |
| OFF | | | | |
| Auto Rerun Condition (Absorbance) | | | | |
| | Lower | OFF | | |
| | Higher | OFF | | |
| Auto Rerun Range (Conc.) | | | | |
| | First Dil | Low | High | |
| | | Re Value Dil | Re Value Dil | |
| Serum | | | | |
| Urine | | | | |
| Plasma | | | | |
| CSF | | | | |
| Dialysis | | | | |
| Other | | | | |
| Auto Rerun Condition (Prozone) | | | | |
| OFF | | | | |
| Dilution | | | | |
| 99:Dil 1 | | | | |

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