PRESTIGE 24i IgM

Nr kat. 4-273, 4-332

(PL)

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Cat. No 4-273, 4-332 (EN)

INTENDED USE

Diagnostic kit for determination of IgM concentration used in automatic analysers Prestige 24i, Biolis 24i, Prestige 24i Premium, Biolis 24i Premium and Biolis 30i.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Immunoglobulins (Igs) are the instrumental proteins of immunity. Immunity is a property of the lymphoid system which is made of organs (spleen, thymus, bone marrow) and of cells (lymphocytes). Circulating immunoglobulins are secreted in the blood by B lymphocytes and they thereby export far-away the specific biological functions of humoral immunity. Immunoglobulin M (IgM) is the first Ig to appear in response to an antigenic stimulus such as an infectious agent. In many cases, the antigen-specific IgM level subsequently falls and remains low as the IgG response appears.

METHOD PRINCIPLE

The IgM present in a sample form with the specific antibody an immunological complex. The increase of turbidity after the addition of antiserum measured at $\lambda=340$ nm is proportional to IgM concentration in the sample.

REAGENTS

Package

	Cat. No 4-273 (24-TRAY)	Cat. No 4-332 (36-TRAY)
1-Reagent	1 x 40 ml	2 x 23 ml
2-Reagent	1 x 10 ml	2 x 6 ml

Buffer (1-Reagent) stored at 2-25°C and antiserum (2-Reagent) stored at 2-8°C are stable until expiry date printed on the package.

Reagent components

1-Reagent: tricine buffer, polymer, inorganic salt, preservative.

2-Reagent: anti-human IgM antiserum, buffer, inorganic salt, preservative.

Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Do not use after expiry date.
- Do not interchange caps.
- Products from human source have been tested for HBsAg and antibodies to HIV and HCV and found to be non-reactive. However this material should be handled as thought capable of transmitting infectious disease.
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.

SPECIMEN

Serum. Nonhemolyzed, nonlipaemic, fresh serum is recommended.

Serum can be stored up to 3 days at 2-8°C or up to 6 months at -20°C.

Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-Reagent and 2-Reagent are ready to use.

1-Reagent put on basic position in reagent tray.

2-Reagent put on start position in reagent tray.

For reagent blank 0.9% NaCl is recommended.

REFERENCE VALUES

adults	0.50 – 3.00 g/l
children (1 year – 12 years)	0.45 – 2.50 g/l
children (1 month – 12 months)	0.20 – 1.50 g/l

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use, with each batch of samples, the CORMAY IMMUNO-CONTROL III (Cat. No 4-291)

For the calibration of automatic analysers systems the CORMAY IMMUNO-MULTICAL (Cat. No 4-287) is recommended.

The calibration curve should be prepared with every change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using an automatic analysers Hitachi or Biolis 30i. Results may vary if a different instrument is used.

- Measurement range:** 0.038 g/l to 7 g/l.

- Specificity / Interferences:**

Hemoglobin up to 0.32 g/dl, bilirubin up to 22 mg/dl, triglycerides up to 312 mg/dl, heparin up to 0.5 g/l, sodium fluoride up to 4 g/l, EDTA up to 5 g/l, sodium citrate up to 5 g/l do not interfere with the test.

- Precision**

Repeatability (run to run) n = 30	Mean [g/l]	SD [g/l]	CV [%]
level 1	0.60	0.004	0.73
level 2	1.25	0.012	0.87
level 3	18.92	0.014	0.72
Reproducibility (day to day) n = 60	Mean [g/l]	SD [g/l]	CV [%]
level 1	0.56	0.024	3.99
level 2	1.30	0.037	2.81
level 3	19.95	0.050	2.49

- Method comparison**

A comparison between IgM values determined at Biolis 30i (y) and at ADVIA SIEMENS 1800 (x) using 62 serum samples gave following results:

$$y = 1.0933x - 0.0387 \text{ g/l};$$

$$R = 0.995 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT

Please refer to local legal requirements.

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Item name	30	IgM
Auto Rerun SW		
OFF		
Auto Rerun Condition (Absorbance)		
Absorbance Range		
	Lower	OFF
	Higher	OFF
Auto Rerun Range (Result)		
	OFF	OFF
	Lower	Higher
Serum		
Urine		
Plasma		
CSF		
Dialysis		
Other		
Prozone Range		
OFF		

• **Prestige 24i Premium, Biolis 24i Premium**

Item No.	30	Item Name	IgM	Optical
Data information				
Units	g/l			
Decimals	2			
Analysis				
Type	END method			
Main Wave Length	340nm			
Sub Wave Length	700nm			
Method	Immuno			
Calibration				
Type	Logit2			
Std sample conc.				
Blank	0	#1 *	#2 *	#3 *
#3	*	#4 *	#5 *	#6 *
#6				
Correlation				
Slope		Intercept		
Y=	1	X+	0	

Item No.	30	Item Name	IgM	Optical
Aspiration				
Kind	Double			
Data Process				
Read	Main	Start	End	
	Sub	29	31	
Abs.Limit				
Low	-3	High	3	
Blank value				
Water Blank				
Reaction Monitor				
0 Level Point	1			
Span	3			
Third mixing				
ON				
Correction value				
Blank correction				0.8350
End Point Limit				2
Linear Check (%)				0
Prozone Check				
	Start	End	Limit (%)	
First				
Second			Low	

Item No.	30	Item Name	IgM	Optical
Normal Range				
	Male		Female	
	Low	High	Low	High
Serum	0.5	3.0	0.5	3.0
Urine				
Plasma				
CSF				
Dialysis				
Other				
Panic Range				
	Male		Female	
	Low	High	Low	High
Serum				
Urine				
Plasma				
CSF				
Dialysis				
Other				

Item No.	30	Item Name	IgM	Optical			
Auto Rerun SW							
OFF							
Auto Rerun Condition (Absorbance)							
Absorbance Range							
	Lower	OFF					
	Higher	OFF					
Auto Rerun Range (Conc.)							
	First Dil	Low			High		
		Re	Value	Dil	Re	Value	Dil
Serum							
Urine							
Plasma							
CSF							
Dialysis							
Other							
Auto Rerun Condition (Prozone)							
OFF							
Dilution							
99:Dil1							

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