## PRESTIGE 24i COMPLEMENT C3

Nr kat. 4-309, 4-335 (PL)

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Cat. No **4-309, 4-335** (EN)

### INTENDED USE

Diagnostic kit for determination of complement C3 concentration intended to use in automatic analysers Prestige 24i, Biolis 24i, Prestige 24i Premium, Biolis 24i Premium and Biolis 30i.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

### INTRODUCTION

Complement is a group of 20 immunologically distinct proteins present in blood and tissues. They are able to interact with antigen-antibody complexes, with each other and with cell membranes, in a complex way intended at destroying viruses and bacteria. They are synthesised in liver and are present in serum as functionally inactive molecules. They are activated by antigen-antibody complexes. C3 complement is a  $\alpha$ -glycoprotein of 2 subunits. It is an acute phase reactant whose levels are increased during the acute phase. Low levels are found in immune complex diseases and in inherited deficiency which results in recurrent infections.

### METHOD PRINCIPLE

The complement C3, present in a sample, form with the specific antibody an immunological complex. The increase of turbidity after the addition of antiserum measured at  $\lambda=340$  nm is proportional to complement C3 concentration in the sample.

### REAGENTS

#### Package

	Cat. No 4-309 (24-TRAY)	Cat. No 4-335 (36-TRAY)
1-Reagent	1 x 40 ml	2 x 23 ml
2-Reagent	1 x 10 ml	2 x 6 ml

Buffer (1-Reagent) stored at 2-25°C and antiserum (2-Reagent) stored at 2-8°C are stable until expiry date printed on the package.

#### Reagent components

**1-Reagent:** Imidazole buffer, polymer, inorganic salt, buffer, preservative.

**2-Reagent:** anti-human complement C3 antiserum, inorganic salt, buffer, preservative.

#### Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Store closed.
- Do not freeze the reagents.
- Nanoparticle-based reagents can settle over time. It may be necessary to delicately mix by repeated turning.

- Human-origin products have been tested for HBsAg and antibodies to HIV 1, HIV 2, HCV and found to be non-reactive. However this material should be handled as potentially infectious.
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.

### SPECIMEN

Serum or plasma.

Serum should be separated from red blood cells as soon as possible after blood collection. If the test cannot be done immediately, the sample should be stored at -70°C. Nevertheless it is recommended to perform the assay with freshly collected samples!

### PROCEDURE

1-Reagent and 2-Reagent are ready to use.

For reagent blank 0.9% NaCl is recommended.

### REFERENCE VALUES <sup>2</sup>

adults	0.9 – 1.8 g/l
--------	---------------

It is recommended for each laboratory to establish its own reference ranges for local population.

### QUALITY CONTROL

For internal quality control it is recommended to use with each batch of samples, the CORMAY IMMUNO-CONTROL III (Cat. No 4-291).

For the calibration of automatic analysers the CORMAY IMMUNO-MULTICAL (Cat. No 4-287) is recommended.

The calibration curve should be prepared with change of reagent lot number or as required e.g. quality control findings outside the specified range.

### PERFORMANCE CHARACTERISTICS

The following results have been obtained using the automatic analyzer Hitachi and Biolis 30i. Results may vary if a different instrument is used.

- LOQ:** 0.015 g/l
- Linearity:** up to 4.31 g/l
- Specificity / Interferences**  
 Haemoglobin up to 0.32 g/dl, bilirubin up to 22 mg/dl, triglycerides up to 312 mg/dl, heparin up to 0.5 g/l, sodium fluoride up to 4 g/l, EDTA up to 5 g/l, sodium citrate up to 5 g/l do not interfere with the test.

#### Precision

Repeatability (run to run) n = 30	Mean [g/l]	SD [g/l]	CV [%]
level 1	0.759	0.005	0.66
level 2	1.615	0.013	0.80
level 3	2.367	0.021	0.87

Reproducibility (day to day) n = 59	Mean [g/l]	SD [g/l]	CV [%]
level 1	0.778	0.019	2.42
level 2	1.634	0.026	1.59
level 3	2.428	0.045	1.87

#### Method comparison

A comparison between complement C3 values determined at **Biolis 30i** (y) and at **SIEMENS ADVIA 1800** (x) using 49 serum samples gave following results:

$$y = 0.9621x + 0.0601 \text{ g/l};$$

$$R = 0.958 \quad (R - \text{correlation coefficient})$$

### WASTE MANAGEMENT

Please refer to local legal requirements.

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Item name	34	C3
<b>Auto Rerun SW</b>		
OFF		
<b>Auto Rerun Range (Result)</b>		
	OFF	OFF
	Lower	Higher
Serum		
Urine		
Plasma		
CSF		
Dialysis		
Other		
<b>Auto Rerun Condition (Absorbance)</b>		
Absorbance Range	Lower	OFF
	Higher	OFF
<b>Prozone Range</b>		
	OFF	

#### • Prestige 24i Premium, Biolis 24i Premium

Item No.	34	Item Name	C3	Optical
<b>Data information</b>				
Units	g/l			
Decimals	2			
<b>Calibration</b>				
Type	Logit2			
Std sample conc.				
Blank	0	#1	*	#2
#3	*	#4	*	#5
#6				*
<b>Analysis</b>				
Type	END method			
Main Wave Length	340 nm			
Sub Wave Length	700 nm			
Method	Immuno			
<b>Correlation</b>				
	Slope	Intercept		
Y=	1	X+	0	

Item No.	34	Item Name	C3	Optical
<b>Aspiration</b>				
Kind	Double			
Vol.				
	Kind	Vol.	Add	Units
Sample		4	5	µl
Reagent 1		250	10	µl
Reagent 2		50	10	µl
<b>Data Process</b>				
Read	Start	End		
	Main	51		52
	Sub	30		31
<b>Abs.Limit</b>				
Low	-3		High	3
<b>Correction value</b>				
Blank correction	0.8355			
End Point Limit	2			
Linear Check (%)	0			
<b>Prozone Check</b>				
	Start	End	Limit (%)	
First				
Second			Low	
<b>Reaction Monitor</b>				
0 Level Point	1			
Span	3			
<b>Third mixing</b>				
ON				

Item No.	34	Item Name	C3	Optical
<b>Normal Range</b>				
	Male	Female		
	Low	High	Low	High
Serum	0.90	1.80	0.90	1.80
Urine				
Plasma				
CSF				
Dialysis				
Other				
<b>Panic Range</b>				
	Male	Female		
	Low	High	Low	High
Serum				
Urine				
Plasma				
CSF				
Dialysis				
Other				

Item No.	34	Item Name	C3	Optical	
<b>Auto Rerun SW</b>					
OFF					
<b>Auto Rerun Range (Conc.)</b>					
	First Dil	Low		High	
		Re	Value	Dil	Re
Serum					
Urine					
Plasma					
CSF					
Dialysis					
Other					
<b>Auto Rerun Condition (Absorbance)</b>					
	Lower	OFF			
	Higher	OFF			
<b>Auto Rerun Condition (Prozone)</b>					
OFF					
<b>Dilution</b>					
99:Dil1					

PRESTIGE 24i COMPLEMENT C3

• **Biolis 30i**

Item no	34	Item name	C3	Specimen	SERUM / PLASMA	OPTICAL												
<b>Data information</b>																		
UNITS	g/L		<b>Aspiration volume</b>															
DECIMALS	2		TYPE Double															
<b>Analysis</b>																		
METHOD	END method		<table border="1"> <tr> <td>VOL. (µL)</td> <td>SAMPLE</td> <td>REAGENT 1</td> <td>REAGENT 2</td> </tr> <tr> <td>BOTTLE (ml)</td> <td>4</td> <td>250</td> <td>50</td> </tr> </table>				VOL. (µL)	SAMPLE	REAGENT 1	REAGENT 2	BOTTLE (ml)	4	250	50				
VOL. (µL)	SAMPLE	REAGENT 1	REAGENT 2															
BOTTLE (ml)	4	250	50															
Main Wave Length	340 nm		FIRST DIL.															
Sub Wave Length	700 nm																	
<b>CORRELATION (Y= AX + B)</b>																		
A =	1		<b>Data processing read</b>															
B =	0		<table border="1"> <tr> <td></td> <td>START</td> <td>END</td> </tr> <tr> <td>MAIN</td> <td>51</td> <td>52</td> </tr> <tr> <td>SUB</td> <td>30</td> <td>31</td> </tr> </table>					START	END	MAIN	51	52	SUB	30	31			
	START	END																
MAIN	51	52																
SUB	30	31																
<b>Blank value</b>																		
• WATER ° REAGENT																		
<b>Calibration</b>																		
TYPE	Logit 2		ABS LIMIT															
STABILITY			-3 TO 3															
<b>Collection value</b>																		
END POINT 2.5																		
LINEARITY CHECK (%) 0																		
<b>Prozone check</b>																		
<table border="1"> <tr> <td></td> <td>START</td> <td>END</td> <td>LIMIT (%)</td> </tr> <tr> <td>FIRST</td> <td></td> <td></td> <td></td> </tr> <tr> <td>SECOND</td> <td></td> <td></td> <td></td> </tr> </table>								START	END	LIMIT (%)	FIRST				SECOND			
	START	END	LIMIT (%)															
FIRST																		
SECOND																		
<table border="1"> <tr> <td>°HIGH</td> <td colspan="2">MINIMUM ABS.</td> </tr> <tr> <td>•LOW</td> <td>MEAN</td> <td></td> </tr> <tr> <td></td> <td>VARIATE</td> <td></td> </tr> </table>							°HIGH	MINIMUM ABS.		•LOW	MEAN			VARIATE				
°HIGH	MINIMUM ABS.																	
•LOW	MEAN																	
	VARIATE																	
Item No	34	Item Name	C3	Specimen	SERUM / PLASMA	OPTICAL												
<b>Reference intervals</b>																		
MALE				FEMALE														
LOW	HIGH	LOW	HIGH															
0.9	1.8	0.9	1.8															
<b>Panic range</b>																		
MALE				FEMALE														
LOW	HIGH	LOW	HIGH															
<b>Decision limit</b>																		
MALE				FEMALE														
<b>Reaction check</b>																		
<table border="1"> <tr> <td></td> <td>°ON</td> <td>•OFF</td> </tr> <tr> <td>CHECK</td> <td></td> <td></td> </tr> <tr> <td>LOW</td> <td></td> <td></td> </tr> <tr> <td>HIGH</td> <td></td> <td></td> </tr> </table>								°ON	•OFF	CHECK			LOW			HIGH		
	°ON	•OFF																
CHECK																		
LOW																		
HIGH																		
<b>VL CHECK</b>																		
°ON •OFF																		
<b>VH CHECK</b>																		
°ON •OFF																		
<b>Auto rerun</b>																		
•ON °OFF																		
<b>Auto rerun range (conc.)</b>																		
Re	Value	Dil.	Re	Value	Dil.													
	0.06			3														
<b>Auto rerun condition (abs.)</b>																		
						DIL.												
LOWER	°ON	•OFF																
HIGH	°ON	•OFF																
<b>Auto rerun condition (prozone)</b>																		
°ON •OFF																		
SAMPLE VOL.																		
<b>Dilution</b>																		
•DIL 1 ° DIL 2																		