PRESTIGE 24i - CERULOPLASMIN

Nr kat. 4-312, 4-338 (PL)

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Cat. No **4-312, 4-338** (EN)

INTENDED USE

Diagnostic kit for determination of ceruloplasmin concentration used in automatic analysers Prestige 24i, Biolis 24i, Prestige 24i Premium, Biolis 24i Premium and Biolis 30i.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Ceruloplasmin is an α_2 -glycoprotein containing 6-7 atoms of copper per molecule. It has been considered for a long time as a copper transporter but since recently, it has been shown to be a serum ferroxidase playing a major role in oxidising iron (II) to iron (III) in serum and at the cell surface, thereby regulating its binding by transferrin.

Ceruloplasmin is a late acute phase reactant. Low levels are found in malnutrition, nephrosis, severe liver diseases such as primary biliary cirrhosis and in Wilson's disease, an autosomal recessive defect in the regulation of copper metabolism.

METHOD PRINCIPLE

The ceruloplasmin presents in a sample form with the specific antibody an immunological complex. The increase of turbidity after the addition of antiserum measured at $\lambda=340$ nm is proportional to ceruloplasmin concentration in the sample.

REAGENTS

Package

	Cat. No 4-312 (24-TRAY)	Cat. No 4-338 (36-TRAY)
1-Reagent	1 x 40 ml	2 x 23 ml
2-Reagent	1 x 10 ml	2 x 6 ml

Buffer (1-Reagent) stored at 2-25°C and antiserum (2-Reagent) stored at 2-8°C are stable until expiry date printed on the package.

Reagent components

1-Reagent: MES buffer, polymer, inorganic salt, preservative.

2-Reagent: anti-human ceruloplasmin antiserum, buffer, inorganic salt, preservative.

Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Store reagents closed.
- Do not freeze the reagents.
- Nanoparticle-based reagents can settle over time. It may be necessary to delicately mix by repeated turning.
- Human-origin products have been tested for HBsAg and antibodies to HIV 1, HIV 2, HCV and found to be non-reactive. However this material should be handled as potentially infectious.

- Please refer to the MSDS for detailed information concerning safe storage and use of the product.

SPECIMEN ⁷

Serum.

Specimen without lipemia or hemolysis is recommended.

Samples remain stable for 3 days at 2-8°C or 4 weeks at -20°C.

Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-Reagent and 2-Reagent are ready to use.

1-Reagent put on basic position in reagent tray.

2-Reagent put on start position in reagent tray.

For reagent blank 0.9% NaCl is recommended.

REFERENCE VALUES ⁷

age	g/l	mg/dl
1 day – 3 months	0.05 – 0.18	5 – 18
6 – 12 months	0.33 – 0.43	33 – 43
1 – 7 years	0.24 – 0.56	24 – 56
> 7 years	0.18 – 0.45	18 – 45

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use, with each batch of samples, the CORMAY IMMUNO-CONTROL III (Cat. No 4-291)

For the calibration of automatic analysers systems the CORMAY IMMUNO-MULTICAL (Cat. No 4-287) is recommended.

The calibration curve should be prepared with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using the automatic analysers Hitachi or Biolis 30i. Results may vary if a different instrument or a manual procedure is used.

- LoQ (limit of quantitation):**
0.948 mg/dl (0.00948 g/l)
- Linearity:**
up to 105.12 mg/dl (1.0512 g/l).
- Specificity / Interferences**
Hemoglobin up to 0.19 g/dl, bilirubin up to 22 mg/dl triglycerides up to 1669,5 mg/dl, heparin up to 0.5 g/l, sodium fluoride up to 4 g/l, EDTA up to 5 g/l, sodium citrate up to 5 g/l do not interfere with the test.

Precision

Repeatability (run to run) n = 30	Mean	SD	CV
	[g/l]	[g/l]	[%]
level 1	0.25	0.002	0.6
level 2	0.40	0.003	0.6
level 3	0.52	0.004	0.7
Reproducibility (day to day) n = 60	Mean	SD	CV
	[g/l]	[g/l]	[%]
level 1	0.19	0.006	3.0
level 2	0.36	0.006	1.6
level 3	0.49	0.006	1.3

Method comparison

A comparison between ceruloplasmin values determined at **Biolis 30i** (y) and at **BS-800** (x) using 61 serum samples gave following results:

$$y = 1.0481 x - 0.0041 \text{ g/l;}$$

$$R = 0.994 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT

Please refer to local legal requirements.

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Item No.	37	Item Name	CERU	Optical
Aspiration				
Kind	Double			
Vol.	Kind	Vol.	Add	Units
Sample	5	5	5	µl
Reagent 1	250	250	10	µl
Reagent 2	50	50	10	µl
Data Process				
Read	Main	Start	End	
	Sub	30	31	
Abs.Limit	Low	High		
	-3	3		
Correction value				
Blank correction	0.8366			
End Point Limit	2			
Linear Check (%)	0			
Prozone Check				
	Start	End	Limit (%)	
First				
Second			Low	
Blank value				
Water Blank				
Reaction Monitor				
0 Level Point	1			
Span	3			
Third mixing				
ON				

Item No.	37	Item Name	CERU	Optical
Normal Range				
	Male		Female	
	Low	High	Low	High
Serum	0.18	0.45	0.18	0.45
Urine				
Plasma				
CSF				
Dialysis				
Other				
Panic Range				
	Male		Female	
	Low	High	Low	High
Serum				
Urine				
Plasma				
CSF				
Dialysis				
Other				

Item No.	37	Item Name	CERU	Optical
Auto Rerun SW				
OFF				
Auto Rerun Range (Conc.)				
	First Dil	Low Value	High Value	
	Re	Dil	Re	Dil
Serum				
Urine				
Plasma				
CSF				
Dialysis				
Other				
Auto Rerun Condition (Absorbance)				
	Lower	Higher	OFF	
			OFF	
Auto Rerun Condition (Prozone)				
OFF				
Dilution				
99:DilI				

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