PRESTIGE 24i ALPHA 1-GLYCOPROTEIN ACID

Nr kat. 4-314, 4-340 PL)

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Cat. No **4-314, 4-340** (EN)

INTENDED USE

Diagnostic kit for determination of a α 1-glycoprotein acid concentration intended to use in automatic analysers Prestige 24i, Biolis 24i and Prestige 24i Premium, Biolis 24i Premium.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

α 1-glycoprotein acid (orosomucoid) is unique amongst plasma proteins because of its low pH and high carbohydrate content. The protein is a serum transporter for steroid hormone and for many drugs. Its physiological role remains unknown but it is an acute phase reactant. The concentration of this protein in serum is used clinically to monitor acute phase responses and tumour recurrence.

METHOD PRINCIPLE

The α 1-glycoprotein acid presents in a sample form with the specific antibody an immunological complex. The increase of turbidity after the addition of antiserum measured at $\lambda=340$ nm is proportional to α -1-glycoprotein acid concentration in the sample.

REAGENTS

Package

	Cat. No 4-314 (24-TRAY)	Cat. No 4-340 (36-TRAY)
1-Reagent	1 x 40 ml	2 x 23 ml
2-Reagent	1 x 10 ml	2 x 6 ml

Buffer (1-Reagent) stored at 2-25°C and antiserum (2-Reagent) stored at 2-8°C are stable until expiry date printed on the package.

Reagent components

1-Reagent: Glycylglycin buffer, polymer, inorganic salt, preservative.

2-Reagent: anti-human α 1-glycoprotein acid antiserum, buffer, inorganic salt, preservative.

Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Store closed.
- Do not freeze the reagents.
- Nanoparticle-based reagents can settle over time. It may be necessary to delicately mix by repeated turning.
- Human-origin products have been tested for HBsAg and antibodies to HIV 1, HIV 2, HCV and found to be non-reactive. However this material should be handled as potentially infectious.
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.

SPECIMEN ⁴

Serum. Nonhemolyzed, nonlipaemic, fresh serum is recommended.

Serum can be stored up to 6 hours at room temperature, up to 7 days at 2-8°C or up to 3 months at -20°C.

Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-Reagent and 2-Reagent are ready to use.

1-Reagent put on basic position in reagent tray.

2-Reagent put on start position in reagent tray.

For reagent blank 0.9% NaCl is recommended.

Application guide for Biolis 30i analyser is available on request.

REFERENCE VALUES ⁴

serum	0.39 – 1.15 g/l
-------	-----------------

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use with each batch of samples, the CORMAY IMMUNO-CONTROL III (Cat. No 4-291).

For the calibration of automatic analysers the CORMAY IMMUNO-MULTICAL (Cat. No 4-287) is recommended. The calibration curve should be prepared with every change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using an automatic analyser Hitachi. Results may vary if a different instrument is used.

- LOQ (limit of quantitation):** 0.143 g/l

- Linearity:** up to 2.42 g/l.

- Interferences:**

Haemoglobin up to 0.32 g/dl, bilirubin up to 22 mg/dl, triglycerides up to 1669.5 mg/dl, heparin up to 0.5 g/l, sodium fluoride up to 4 g/l, EDTA up to 5 g/l, sodium citrate up to 5 g/l do not interfere with the test.

- Precision**

Repeatability (run to run) n = 30	Mean [g/l]	SD [g/l]	CV [%]
level 1	0.34	0.005	1.36
level 2	0.66	0.005	0.69
level 3	0.92	0.004	0.41

Reproducibility (day to day) n = 60	Mean [g/l]	SD [g/l]	CV [%]
level 1	0.42	0.012	2.48
level 2	0.72	0.011	1.57
level 3	0.96	0.012	1.22

- Method comparison**

A comparison between CORMAY reagent (y) and another commercially available assay (x) using at least 30 samples gave following results:

$$y = 1.0751x - 0.0527 \text{ g/l};$$

$$R = 0.990 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT

Please refer to local legal requirements.

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Item name	39	AGP
Auto Rerun SW		
OFF		
Auto Rerun Range (Result)		
	OFF	OFF
	Lower	Higher
Serum		
Urine		
Plasma		
CSF		
Dialysis		
Other		
Auto Rerun Condition (Absorbance)		
Absorbance Range	Lower	OFF
	Higher	OFF
Prozone Range		
		OFF

• **Prestige 24i Premium, Biolis 24i Premium**

Item No.	39	Item Name	AGP	Optical	
Data information					
Units	g/l				
Decimals	2				
Calibration					
Type	Logit2				
Std sample conc.					
Blank	0	#1	*	#2	*
#3	*	#4	*	#5	*
#6					
Analysis					
Type	END method				
Main Wave Length	340nm				
Sub Wave Length	700nm				
Method	Immuno				
Correlation					
Y=	Slope	Intercept			
	1	0			

Item No.	39	Item Name	AGP	Optical
Aspiration				
Kind	Double			
Vol.				
Kind	Vol.	Add	Units	
Sample	3	5	µl	
Reagent 1	250	10	µl	
Reagent 2	50	10	µl	
Data Process				
Read	Start	End		
Main	44	45		
Sub	30	31		
Abs.Limit				
Low	-3		High	
			3	
Correction value				
Blank correction	0.8350			
End Point Limit	2			
Linear Check (%)	0			
Prozone Check				
First	Start	End	Limit (%)	
Second			Low	
Third mixing				
ON				

Item No.	39	Item Name	AGP	Optical
Normal Range				
	Male		Female	
	Low	High	Low	High
Serum	0.39	1.15	0.39	1.15
Urine				
Plasma				
CSF				
Dialysis				
Other				
Panic Range				
	Male		Female	
	Low	High	Low	High
Serum				
Urine				
Plasma				
CSF				
Dialysis				
Other				

Item No.	39	Item Name	AGP	Optical
Auto Rerun SW				
OFF				
Auto Rerun Range (Conc.)				
	First Dil	Low	High	
		Re	Value	Dil
Serum				
Urine				
Plasma				
CSF				
Dialysis				
Other				
Auto Rerun Condition (Absorbance)				
	Lower	OFF		
	Higher	OFF		
Auto Rerun Condition (Prozone)				
	OFF			
Dilution				
	99:Dil1			

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