Prestige 24i ANTITHROMBIN III

Nr kat. 4-316, 4-342 (PL)

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Cat. No 4-316, 4-342 (EN)

INTENDED USE

Diagnostic kit for determination of antithrombin III concentration used in automatic analysers Prestige 24i, Biolis 24i, Prestige 24i Premium, Biolis 24i Premium and Biolis 30i.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Antithrombin III (AT III) is an α 2-glycoprotein of MW 58000 and is made in the liver. AT III is one of the most important regulators of the coagulation system. AT III inactivates thrombin and factors Xa, IXa, XIa and XIIa. This anticoagulant activity is enhanced by the presence of heparin, which forms a ternary complex with AT III and these procoagulant factors. Reduced concentration of AT III in blood means a great and well established risk for thrombotic complications. Clinical low value are associated with congenital or acquired deficiencies caused by a decreased biosynthesis (liver disease, medical treatment) or an increased loss (gastrointestinal disease, nephrotic syndrome) or an increased consumption (sepsis, major trauma due to surgery and burns wounds).

METHOD PRINCIPLE

The antithrombin III presents in a sample form with the specific antibody an immunological complex. The increase of turbidity after the addition of antiserum measured at $\lambda=340$ nm is proportional to antithrombin III concentration in the sample.

REAGENTS

Package

	Cat. No 4-316 (24-TRAY)	Cat. No 4-342 (36-TRAY)
1-Reagent	1 x 40 ml	2 x 23 ml
2-Reagent	1 x 10 ml	2 x 6 ml

Buffer (1-Reagent) stored at 2-25°C and antiserum (2-Reagent) stored at 2-8°C are stable until expiry date printed on the package.

Reagent components

1-Reagent: tricine buffer, polymer, inorganic salt, preservative.

2-Reagent: anti-human antithrombin III antiserum, buffer, inorganic salt, preservative.

Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Store closed.
- Do not freeze the reagents.
- Nanoparticle-based reagents can settle over time. It may be necessary to delicately mix by repeated turning.

- Human-origin products have been tested for HBsAg and antibodies to HIV 1, HIV 2, HCV and found to be non-reactive. However this material should be handled as potentially infectious.
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.

SPECIMEN ⁷

Citrated plasma.

Samples may be stored up to 4 hours on ice. If the test cannot be done immediately, the sample should be stored frozen.

Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-Reagent and 2-Reagent are ready to use.

1-Reagent put on basic position in reagent tray.

2-Reagent put on start position in reagent tray.

For reagent blank 0.9% NaCl is recommended.

REFERENCE VALUES ⁷

plasma	0.21 – 0.30 g/l
--------	-----------------

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use, with each batch of samples, the CORMAY IMMUNO-CONTROL III (Cat. No 4-291).

For the calibration of automatic analysers systems the CORMAY IMMUNO-MULTICAL (Cat. No 4-287) is recommended.

The calibration curve should be prepared with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

These metrological characteristics have been obtained using the Hitachi or Biolis 30i automatic analysers. Results may vary if a different instrument is used.

- LoQ (Limit of Quantitation):** 0.028 g/l

- Linearity:** up to 0.806 g/l.

- Interferences:**

Hemoglobin up to 0.06 g/dl, bilirubin up to 22 mg/dl, triglycerides up to 155,5 mg/dl, heparin up to 0.5 g/l, sodium fluoride up to 4 g/l, EDTA up to 5 g/l, sodium citrate up to 5 g/l do not interfere with the test.

- Precision**

Repeatability (run to run) n = 30	Mean [g/l]	SD [g/l]	CV [%]
level 1	0.16	0.002	0.011
level 2	0.32	0.002	0.005
level 3	0.46	0.005	0.012

Reproducibility (day to day) n = 56	Mean [g/l]	SD [g/l]	CV [%]
level 1	0.18	0.005	0.03
level 2	0.34	0.008	0.02
level 3	0.50	0.010	0.02

- Method comparison**

A comparison between antithrombin III concentration determined at **Biolis 30i** (y) and at **BS-800** (x) using 59 plasma samples gave following results:

$y = 1.0358 x - 0.0103$ g/l;

$R = 0.985$ (R – correlation coefficient)

WASTE MANAGEMENT

Please refer to local legal requirements.

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Item name	41	AT 3	
Auto Rerun SW			
OFF			
Auto Rerun Condition (Absorbance)			
Absorbance Range			
	Lower	OFF	
	Higher	OFF	
Auto Rerun Range (Result)			
	OFF	OFF	
	Lower	Higher	
Auto Rerun Range (Conc.)			
	Prozone Range		OFF
Serum			
Urine			
Plasma			
CSF			
Dialysis			
Other			

• **Prestige 24i Premium, Biolis 24i Premium**

Item No.	41	Item Name	AT 3	Optical
Data information				
Units	g/L			
Decimals	3			
Calibration				
Type	Logit2			
Std sample conc.				
Blank	0	#1	*	#2 *
#3	*	#4	*	#5 *
#6				
Analysis				
Type	END method			
Main Wave Length	340 nm			
Sub Wave Length	700 nm			
Method	Immuno			
Correlation				
	Slope	Intercept		
Y=	1	X+	0	

Item No.	41	Item Name	AT 3	Optical
Aspiration				
Kind	Double			
Vol.				
	Kind	Vol.	Add	Units
Sample		10	5	µl
Reagent 1		250	10	µl
Reagent 2		50	10	µl
Data Process				
Read	Start	End		
	Main	51	52	
	Sub	30	31	
Abs.Limit				
Low				High
	-3	~	3	
Blank value				
Water Blank				
Reaction Monitor				
0 Level Point	1			
Span	3			
Third mixing				
ON				
Correction value				
Blank correction				0.8387
End Point Limit				2
Linear Check (%)				0
Prozone Check				
	Start	End	Limit (%)	
First				
Second			Low	

Item No.	41	Item Name	AT 3	Optical
Normal Range				
	Male		Female	
	Low	High	Low	High
Serum				
Urine				
Plasma	0.25	0.45	0.25	0.45
CSF				
Dialysis				
Other				
Panic Range				
	Male		Female	
	Low	High	Low	High
Serum				
Urine				
Plasma				
CSF				
Dialysis				
Other				

Item No.	41	Item Name	AT 3	Optical			
Auto Rerun SW							
OFF							
Auto Rerun Condition (Absorbance)							
Absorbance Range							
	Lower	OFF					
	Higher	OFF					
Auto Rerun Range (Conc.)							
	First Dil	Low		High			
		Re	Value	Dil	Re	Value	Dil
Serum							
Urine							
Plasma							
CSF							
Dialysis							
Other							
Auto Rerun Condition (Prozone)							
OFF							
Dilution							
99:Dil 1							

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