PRESTIGE 24i ETHANOL

Nr kat. 4-378, 4-478 (PL)

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LITERATURA

- Burtis C.A., Ashwood E.R., ed. Tietz Textbook of Clinical Chemistry, Toronto: WB Saunders Co., 1171-1175 (1994).
- William H., Porter Ph.D. Clinical Toxicology. In: Burtis CA, Ashwood ER, editors. Tietz Textbook of Clinical Chemistry, 3rd ed. Philadelphia: W.B Saunders Company; 1999. p. 922-923.
- Bonnichsen, R.K., Theorell, H., An Enzymatic Method for the Micro Determination of Ethanol, Scand. J. Clin. Lab. Invest. 3, 58 (1951).
- Lundquist, F., The Determination of Ethyl Alcohol in Blood and Tissues., Methods of Biochemical Analysis, Vol III, Editor, D. Glick, Interscience Publishers, New York, 1957 p. 217-251.
- CLSI Guidelines and Standards, Clinical Laboratory Tests, AACC Press, Washington, 1990.
- Preeti Dalawari: medicine.medscape.com/article/2090019/overview#3.
- Burtis C.A., Ashwood E.R., ed. Tietz Textbook of Clinical Chemistry, 3rd ed. Philadelphia, PA: WB Saunders, 923, (1999).
- Krajowa Izba Diagnostów Laboratoryjnych, Oznaczanie etanolu w krwi i moczu pacjentów w laboratorium medycznym, Warszawa 2016. p.18-19.

Data wydania: 04. 2021.

PRESTIGE 24i ETHANOL

Cat. No **4-378, 4-478** (EN)

INTENDED USE

Diagnostic kit for determination of ethyl alcohol concentration intended to use in automatic analyzers: Prestige 24i, Biolis 24i, Prestige 24i Premium, Biolis 24i Premium and Biolis 30i.


The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Ethanol after consumption is absorbed in alimentary tract and transported with blood to whole organism. A little amount of ethanol is excreted with urine or exhaled with air, but majority is metabolized by liver to acetic aldehyde, then to acetate, finally to carbon dioxide and water. Ethanol consumption is often cause of various types of accidents, and ethanol poisoning may be cause of decease.

METHOD PRINCIPLE

Enzymatic method with alcohol dehydrogenase.

The rate of absorbance changes measured at $\lambda=340$ nm or $\lambda=380$ nm is proportional to the quantity of ethyl alcohol in the sample.

REAGENTS

Package	Cat. No 4-378 (24-TRAY)	Cat. No 4-478 (36-TRAY)
1-Reagent	2 x 30 ml	2 x 30 ml

The reagent when stored tightly capped at 2-8°C is stable up to expiry date printed on the package. Stability on board of the analyzer at 2-10°C: Prestige 24i – 5 weeks, Biolis 24i Premium – 4 weeks.

Concentrations in the test

buffer (pH 7,3)	
NADanalog	1.35 mmol/l
alcohol dehydrogenase	>203 kU/l
stabilizers, preservatives	

Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Avoid contact with skin and mucous membranes.
- Concentration of ethyl alcohol in serum or plasma is approximately 1.2 times higher than in whole blood.⁶
- Turbidity would indicate contamination or deterioration of reagent.

SPECIMEN

Serum or plasma (heparin, EDTA or citrate) free from haemolysis; urine.

Samples may be stored tightly closed in temp. 2-8°C up to 3 days. Urine samples should be stored with as little dead air space as possible.

Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-Reagent is ready to use.

1-Reagent put on basic position in reagent tray.

Deionised water or 0.9% NaCl is recommended as a reagent blank.

REFERENCE VALUES 7

Ethanol consumption cause following symptoms:

blood ethanol concentration (mg / dl)	clinical symptoms*
30-120	euphoria; sociability; talkativeness; diminution of attention, judgment and control; some sensory-motor impairment
90-250	emotional instability; loss of critical judgment; impairment of perception, memory and comprehension; increased reaction time; sensory-motor incoordination; drowsiness;
180-300	disorientation; mental confusion; dizziness; staggering gait; slurred speech; disturbances of vision; increased pain threshold; apathy; lethargy;
250-400	general inertia; markedly decreased response to stimuli; vomiting; incontinence of urine and feces; sleep or stupor;
350-500	complete unconsciousness; coma; depressed or abolished reflexes; impairment of circulation and respiration; possible death.

*Effects are dependent on individual use, chronicity, tolerance and other factors.⁸

To calculate the concentration of ethyl alcohol in blood the result received for serum or plasma should be divided by 1.2 factor. Otherwise, doctor should be informed that the result is overstated and requires a proper interpretation.⁸

Please refer to the local regulations concerning legal ethanol limits for drivers.

Conversion factor of the blood ethanol concentration:
100 mg/dl = 1 %.

QUALITY CONTROL

For internal quality control it is recommended to use, with each batch of samples, the CORMAY AMMONIA/ETHANOL CONTROLS (Cat. No 5-163). For the calibration of automatic analyzers the CORMAY ETHANOL CALIBRATORS (Cat. No 5-105) is recommended. Two levels of calibrator and deionised water or 0.9% NaCl (as a calibrator 0) should be used for calibration. The calibration curve should be prepared every 2 weeks (Prestige 24i) or every 4 weeks (Biolis 24i Premium), with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analyzers Biolis 30i or Beckman CX. Results may vary if a different instrument or a manual procedure is used.

- LoB (Limit of Blank):**
0.4 mg/dl (0.087 mmol/l)
- LoD (Limit of Detection):**
2.1 mg/dl (0.46 mmol/l)
- LoQ (Limit of Quantitation):**
7.0 mg/dl (1.52 mmol/l)
- Linearity:**
up to 800 mg/dl (173.6 mmol/l)

Specificity / Interferences

Haemoglobin up to 0.6 g/dl, bilirubin up to 40 mg/dl, intralipid up to 1000 mg/dl, LDH up to 10 000 U/l (in the presence of 200 mg/dl lactate) do not interfere with the test.

Acetone up to 2000 mg/dl, butanol up to 200 mg/dl, ethylene glycol up to 2000 mg/dl, methanol up to 2000 mg/dl and isopropanol up to 2000 mg/dl do not interfere with the test, which confirms that assay is specific to ethanol.

Precision

Repeatability (run to run) n = 20	Mean [mg/dl]	SD [mg/dl]	CV [%]
level 1	42.4	1.13	2.65
level 2	97.1	2.51	2.58
level 3	277.3	5.74	2.07

Method comparison

A comparison between CORMAY reagent (y) and commercially available assay (x) using 40 serum samples, gave following results:

$$y = 1.064 x - 2.33 \text{ mg/dl;}$$

$$R = 0.9994 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT

Please refer to local legal requirements.

LITERATURE

- Burtis C.A., Ashwood E.R., ed. Tietz Textbook of Clinical Chemistry, Toronto: WB Saunders Co., 1171-1175 (1994).
- William H., Porter Ph.D. Clinical Toxicology. In: Burtis CA, Ashwood ER, editors. Tietz Textbook of Clinical Chemistry, 3rd ed. Philadelphia: W.B Saunders Company; 1999. p. 922-923.
- Bonnichsen, R.K., Theorell, H., An Enzymatic Method for the Micro Determination of Ethanol, Scand. J. Clin. Lab. Invest. 3, 58 (1951).
- Lundquist, F., The Determination of Ethyl Alcohol in Blood and Tissues., Methods of Biochemical Analysis, Vol III, Editor, D. Glick, Interscience Publishers, New York, 1957 p. 217-251.
- CLSI Guidelines and Standards, Clinical Laboratory Tests, AACC Press, Washington, 1990.
- Preeti Dalawari:
emedicine.medscape.com/article/2090019/overview#a3.
- Burtis C.A., Ashwood E.R., ed. Tietz Textbook of Clinical Chemistry, 3rd ed. Philadelphia, PA: WB Saunders, 923, (1999).
- Polish National Chamber of Laboratory Diagnosticians, Oznaczanie etanolu w krwi i moczu pacjentów w laboratorium medycznym, Warszawa 2016. p.18-19.

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- Burtis C.A., Ashwood E.R., ed. Tietz Textbook of Clinical Chemistry, Toronto: WB Saunders Co., 1171-1175 (1994).
- William H., Porter Ph.D. Clinical Toxicology. In: Burtis CA, Ashwood ER, editors. Tietz Textbook of Clinical Chemistry, 3rd ed. Philadelphia: W.B Saunders Company; 1999. p. 922-923.
- Bonnichsen, R.K., Theorell, H., An Enzymatic Method for the Micro Determination of Ethanol, Scand. J. Clin. Lab. Invest. 3, 58 (1951).
- Lundquist, F., The Determination of Ethyl Alcohol in Blood and Tissues, Methods of Biochemical Analysis, Vol III, Editor, D. Glick, Interscience Publishers, New York, 1957 p. 217-251.
- CLSI Guidelines and Standards, Clinical Laboratory Tests, AACC Press, Washington, 1990.
- Preeti Dalawari: medicine.medscape.com/article/2090019-overview#a3.
- Burtis C.A., Ashwood E.R., ed. Tietz Textbook of Clinical Chemistry, 3rd ed. Philadelphia, PA: WB Saunders, 923, (1999).
- Polish National Chamber of Laboratory Diagnosticians, Oznaczenie etanolu w krwi i moczu pacjentów w laboratorium medycznym, Warszawa 2016. p.18-19.

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Item name	67	ETHANOL
Auto Rerun SW		
ON		
Auto Rerun Range (Result)		
	ON	ON
	Lower	Higher
Serum	7	800
Urine		
Plasma		
CSF		
Dialysis		
Other		
Auto Rerun Condition (Absorbance)		
Absorbance Range	Lower	ON
	Higher	OFF
Prozone Range		
Prozone Range	OFF	

• Prestige 24i Premium, Biolis 24i Premium

Item No.	67	Item Name	ETH	Optical
Data information				
Units	mg/dl			
Decimals	1			
Calibration				
Type	Linear2			
Std sample conc.				
Blank	0	#1	*	#2
#3		#4		#5
#6				*
Analysis				
Type	RATE method			
Main Wave Length	380 nm			
Sub Wave Length	505 nm			
Method	UV			
Correlation				
Y=	Slope	X+	Intercept	
	1		0	

Item No.	67	Item Name	ETH	Optical
Aspiration				
Kind	Single			
Data Process				
Read	Main	Start	End	
	Sub	18	25	
Abs.Limit				
Low	-3	High	3	
Correction value				
Blank correction	1			
End Point Limit	2			
Linear Check (%)	90			
Prozone Check				
	Start	End	Limit (%)	
First				
Second			Low	
Reaction Monitor				
0 Level Point	1			
Span	3			
Third mixing				
OFF				

Item No.	67	Item Name	ETH	Optical
Normal Range				
	Male	Female		
	Low	High	Low	High
Serum				
Urine				
Plasma				
CSF				
Dialysis				
Other				
Panic Range				
	Male	Female		
	Low	High	Low	High
Serum				
Urine				
Plasma				
CSF				
Dialysis				
Other				

Item No.	67	Item Name	ETH	Optical
Auto Rerun SW				
ON				
Auto Rerun Range (Conc.)				
	First Dil	Low	High	
		Re	Value	Dil
Serum			7	800
Urine				
Plasma				
CSF				
Dialysis				
Other				
Auto Rerun Condition (Absorbance)				
	Lower	OFF		
	Higher	OFF		
Auto Rerun Condition (Prozone)				
OFF				
Dilution				
100:DiI2				

PRESTIGE 24i ETHANOL

Biolis 30i

Item no	67	Item name	ETH	Specimen	SERUM/PLASMA	OPTICAL												
Data information																		
UNITS	mg/dL		Aspiration volume															
DECIMALS	1		TYPE		Single													
Analysis			<table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td></td> <td>SAMPLE</td> <td>REAGENT 1</td> <td>REAGENT 2</td> </tr> <tr> <td>VOL. (µL)</td> <td>2</td> <td>160</td> <td></td> </tr> <tr> <td>BOTTLE (ml)</td> <td></td> <td></td> <td></td> </tr> </table>					SAMPLE	REAGENT 1	REAGENT 2	VOL. (µL)	2	160		BOTTLE (ml)			
	SAMPLE	REAGENT 1	REAGENT 2															
VOL. (µL)	2	160																
BOTTLE (ml)																		
METHOD	RATE method		FIRST DIL.															
Main Wave Length	380 nm																	
Sub Wave Length	505 nm																	
CORRELATION (Y= AX + B)																		
A =	1		Data processing read															
B =	0		START		END													
Blank value			MAIN		20													
° WATER ° REAGENT			SUB															
Calibration			ABS LIMIT															
TYPE	Linear 2		-3		TO 3													
STABILITY			Collection value															
			END POINT		2.5													
			LINEARITY CHECK (%)		90													
Prozone check																		
		START	END	LIMIT (%)														
FIRST																		
SECOND																		
		MINIMUM ABS.																
°HIGH		MEAN																
•LOW		VARIATE																

Item No	67	Item Name	ETH	Specimen	SERUM/PLASMA	OPTICAL
Reference intervals						
MALE		FEMALE		Auto rerun		
LOW	HIGH	LOW	HIGH	°ON		°OFF
				Auto rerun range (conc.)		
Panic range			Re	Value	Dil.	
MALE		FEMALE		7		800
LOW	HIGH	LOW	HIGH	Auto rerun condition (abs.)		
						DIL.
Decision limit			LOWER	°ON	•OFF	
			HIGH	°ON	•OFF	
			Auto rerun condition (prozone)			
			°ON		•OFF	
			SAMPLE VOL.			
			Dilution			
			•DIL 1		° DIL 2	
Reaction check						
		°ON		•OFF		
CHECK						
LOW						
HIGH						
VL CHECK			VH CHECK			
°ON		•OFF		°ON		•OFF

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