## PRESTIGE 24i LACTATE

Nr kat. **4-386, 4-466** (PL)

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### LITERATURA

- Tietz Textbook of Clinical Chemistry (Edited by Burtis CA and Ashwood ER Eds): Third Edition WB Saunders Company 787-8, (1999).
- Alan H. B. Wu, Tietz Clinical Guide to Laboratory Tests, W.B. Saunders Company, 4th edition, 650-652, (2006).

Data wydania: 04.2021.

## PRESTIGE 24i LACTATE

Cat. No **4-386, 4-466**

(EN)

### INTENDED USE

Diagnostic kit for determination of lactate concentration intended to use in automatic analyzers: Prestige 24i, Biolis 24i, Prestige 24i Premium, Biolis 24i Premium and Biolis 30i. The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

### INTRODUCTION

Lactate is produced in Cori cycle, by anaerobic conversion of glucose, mainly in skeletal muscle. Its determination, frequently done together with pyruvate, is useful in discovering lactic acidosis due to i.a. reduced tissue oxygenation, enzymatic deficiencies, diabetes mellitus, liver and kidneys diseases.

### METHOD PRINCIPLE

Lactate is oxidized by lactate oxidase to pyruvate and hydrogen peroxide, which, in presence of peroxidase (POD), reacts with 4-aminoantipyrine and phenol forming a compound, which colour intensity is proportional to the concentration of lactate in the examined sample.

### REAGENTS

#### Package

	Cat. No <b>4-386</b> <b>(24-TRAY)</b>	Cat. No <b>4-466</b> <b>(36-TRAY)</b>
1-Reagent	2 x 19,5 ml	2 x 17,5 ml

Unopened reagent is stable up to the kit expiry date printed on the package when stored at 2-8°C. Stability on board of the analyser at 2-10°C: Prestige 24i - 12 weeks, Biolis 24i Premium - 12 weeks.

### Concentrations in the test

Tris buffer (pH 7.5)	≥ 50 mmol/l
lactate oxidase	≥ 0.2 kU/l
peroxidase	≥ 2 kU/l
4-aminoantipyrine	≥ 0.4 mmol/l

### Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Do not freeze the reagent!
- Do not use after expiry date.
- The appearance of turbidity or control sera values outside the manufacturer's acceptable range may indicate of the reagents instability.
- Lactate concentration rapidly increases during physical activities. Normal levels are reached again after usually 30 minutes but it may vary according to individuals.
- Draw blood with lowest venous stasis as possible (max. 30 seconds) from fasting and completely resting patient and avoid using a tourniquet.

### SPECIMEN

Plasma. Avoid haemolysis.

Collect samples in tubes containing sodium fluoride and potassium oxalate. Keep samples on ice. Centrifuge within 15 minutes after collection and separate from cells. Analyze promptly. Note whether sample is venous or arterial. It is recommended to follow NCCLS procedures regarding specimen collecting and handling.

Lactate in plasma is stable up to 8 hours at room temperature or up to 14 days at 2-8°C.

Nevertheless it is recommended to perform the assay with freshly collected samples!

### PROCEDURE

1-Reagent is ready to use.

1-Reagent put on basic position in reagent tray.

For reagent blank deionized water is recommended.

#### Actions required:

- Biolis 24i Premium:** When performing assays in the analyser, there is a probability of **cross-contamination** affecting the tests results: LACTATE – ALAT (II GEN). To avoid this effect follow the recommendations contained in the instruction 51\_03\_24\_008\_BIOLIS\_24i\_PREMIUM\_CARRYOVER.
- Biolis 30i:** When performing assays in the analyser, there is a probability of **cross-contamination** affecting the tests results: LACTATE – CRP ULTRA. To avoid this effect follow the recommendations contained in the instruction 51\_03\_24\_009\_BIOLIS\_30i\_CARRYOVER.

### REFERENCE VALUES <sup>2</sup>

plasma (venous)	4.5 – 19.8 mg/dl	0.5 – 2.2 mmol/l
plasma (arterial)	4.5 – 14.4 mg/dl	0.5 – 1.6 mmol/l

It is recommended for each laboratory to establish its own reference ranges for local population.

### QUALITY CONTROL

For internal quality control it is recommended to use, with each batch of samples, the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173). For the calibration of automatic analysers systems the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) is recommended.

The calibration curve should be prepared every 12 weeks (Prestige 24i, Biolis 24i Premium), with change of reagent lot number or as required e.g. quality control findings outside the specified range.

### PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analyser Biolis 30i. Results may vary if a different instrument or a manual procedure is used.

- LoB (Limit of Blank):**  
0.07 mg/dl (0.0078 mmol/l)

- LoD (Limit of Detection):**  
0.17 mg/dl (0.0189 mmol/l)

- LoQ (Limit of Quantitation):**  
1.40 mg/dl (0.1554 mmol/l)

- Linearity:**  
up to 139 mg/dl (15.43 mmol/l)

For higher concentration dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

#### Specificity / Interferences

- In plasma samples containing approximately 12 mg lactate/dl, there is no interference up to: 0.23 g/dl haemoglobin, 8 mg/dl bilirubin, 337 mg/dl triglycerides, 15.5 mg/l ascorbic acid.
- In plasma samples containing approximately 40 mg lactate/dl, there is no interference up to: 1.25 g/dl haemoglobin, 10 mg/dl bilirubin, 1000 mg/dl triglycerides, 62 mg/l ascorbic acid.

#### Precision

Repeatability (run to run) n = 20	Mean [mg/dl]	SD [mg/dl]	CV [%]
level 1	8.8	0.08	0.92
level 2	42.4	0.30	0.70
Reproducibility (day to day) n = 80	Mean [mg/dl]	SD [mg/dl]	CV [%]
level 1	9.1	0.37	4.1
level 2	40.4	0.93	2.3

#### Method comparison

A comparison between lactate values determined at **Biolis 30i** (y) and at **BS-400** (x) using 58 plasma samples gave following results:

$$y = 1.035 x - 0.3969 \text{ mg/dl;}$$

$$R = 0.999 \quad (R - \text{correlation coefficient})$$

### WASTE MANAGEMENT

Please refer to local legal requirements.

### LITERATURE

- Tietz Textbook of Clinical Chemistry (Edited by Burtis CA and Ashwood ER Eds): Third Edition WB Saunders Company 787-8, (1999).
- Alan H. B. Wu, Tietz Clinical Guide to Laboratory Tests, W.B. Saunders Company, 4th edition, 650-652, (2006).

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Item name	66	LACT
<b>Auto Rerun SW</b>		
ON		
<b>Auto Rerun Condition (Absorbance)</b>		
Absorbance Range	Lower	Higher
		OFF
		OFF
<b>Auto Rerun Range (Result)</b>		
	ON	ON
	Lower	Higher
Serum		
Urine		
Plasma	2.40	90
CSF		
Dialysis		
Other		
<b>Prozone Range</b>		
		OFF

## ▪ Prestige 24i Premium, Biolis 24i Premium

Item No.	66	Item Name	LACT	Optical
<b>Data information</b>				
Units	mg/dl			
Decimals	2			
<b>Analysis</b>				
Type	END method			
Main Wave Length	546nm			
Sub Wave Length	700nm			
Method	Colorimetric			
<b>Calibration</b>				
Type	Linear2			
Std sample conc.				
Blank	0	#1	*	#2
#3		#4		#5
#6				
<b>Correlation</b>				
Y=	Slope	X+	Intercept	
	1		0	

Item No.	66	Item Name	LACT	Optical
<b>Aspiration</b>				
Kind	Single			
<b>Data Process</b>				
Read	Main	Start	End	
	Sub	9	10	
<b>Abs.Limit</b>				
Low				High
	-3			3
<b>Blank value</b>				
Water Blank				
<b>Reaction Monitor</b>				
0 Level Point	1			
Span	3			
<b>Third mixing</b>				
OFF				
<b>Correction value</b>				
Blank correction				
End Point Limit				
2				
Linear Check (%)				
<b>Prozone Check</b>				
	Start	End	Limit (%)	
First				
Second				Low

Item No.	66	Item Name	LACT	Optical
<b>Normal Range</b>				
	Male		Female	
	Low	High	Low	High
Serum				
Urine				
Plasma	4.5	19.8	4.5	19.8
CSF				
Dialysis				
Other				
<b>Panic Range</b>				
	Male		Female	
	Low	High	Low	High
Serum				
Urine				
Plasma				
CSF				
Dialysis				
Other				

Item No.	66	Item Name	LACT	Optical
<b>Auto Rerun SW</b>				
ON				
<b>Auto Rerun Condition (Absorbance)</b>				
	Lower	Higher		
			OFF	
			OFF	
<b>Auto Rerun Range (Conc.)</b>				
	First Dil	Low	High	
		Re	Value	Dil
Serum				
Urine				
Plasma		1.4	90	
CSF				
Dialysis				
Other				
<b>Auto Rerun Condition (Prozone)</b>				
OFF				
<b>Dilution</b>				
100:DiI2				

## PRESTIGE 24i LACTATE

• **Biolis 30i**

Item no	66	Item name	LACT	Specimen	PLASMA	OPTICAL
<b>Data information</b>			<b>Aspiration volume</b>			
UNITS	mg/dL			TYPE	Single	
DECIMALS	2				SAMPLE	REAGENT 1
				VOL. (µL)	3	300
				BOTTLE (ml)		
				FIRST DIL.		
<b>Analysis</b>			<b>Data processing read</b>			
METHOD	END method			START	END	
Main Wave Length	546 nm			MAIN	52	53
Sub Wave Length	700 nm			SUB		
<b>CORRELATION (Y= AX + B)</b>			<b>ABS LIMIT</b>			
A =	1			-3	TO	3
B =	0					
<b>Blank value</b>			<b>Collection value</b>			
		• WATER ° REAGENT		END POINT	2.5	
				LINEARITY CHECK (%)	0	
<b>Calibration</b>			<b>Prozone check</b>			
TYPE	Linear 2			START	END	LIMIT (%)
STABILITY				FIRST		
				SECOND		
				MINIMUM ABS.		
				°HIGH	MEAN	
				•LOW	VARIATE	

Item No	66	Item Name	LACT	Specimen	PLASMA	OPTICAL
<b>Reference intervals</b>			<b>Auto rerun</b>			
MALE		FEMALE		•ON		°OFF
LOW	HIGH	LOW	HIGH			
4.5	19.8	4.5	19.8			
<b>Panic range</b>			<b>Auto rerun range (conc.)</b>			
MALE		FEMALE		Re	Value	Dil.
LOW	HIGH	LOW	HIGH		1.4	139
<b>Decision limit</b>			<b>Auto rerun condition (abs.)</b>			
MALE		FEMALE		DIL.		
LOW	HIGH	LOW	HIGH	LOWER	°ON	•OFF
				HIGH	°ON	•OFF
<b>Reaction check</b>			<b>Auto rerun condition (prozone)</b>			
				°ON •OFF		
				SAMPLE VOL.		
<b>VL CHECK</b>			<b>Dilution</b>			
°ON •OFF		°ON •OFF		•DIL 1 ° DIL 2		

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