## PRESTIGE 24i LQ UA

Nr kat. 4-208, 4-408 (PL)

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Cat. No **4-208, 4-408** (EN)

### INTENDED USE

Diagnostic kit with ascorbate oxidase for determination of uric acid concentration used in automatic analysers analysers Prestige 24i, Biolis 24i and Prestige 24i Premium, Biolis 24i Premium, Biolis 30i.



The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

### INTRODUCTION

Uric acid is a product of purine catabolism. It is produced in the liver and excreted in the urine. Both, the amount of uric acid production and the efficiency of renal excretion, affect serum urate level. Elevated serum uric acid level is caused usually by gout, leukemia, diabetes mellitus, hyperfunction of parathyroid and thyroid, renal failure, renal calculus. Urate concentration in serum and in urine depends on glomerular filtration, thus is useful for renal function monitoring.

### METHOD PRINCIPLE

Enzymatic, colorimetric method with uricase and peroxidase.

The colour intensity is proportional to the uric acid concentration.

### REAGENTS

#### Package

	Cat. No <b>4-208</b> <b>(24-TRAY)</b>	Cat. No <b>4-408</b> <b>(36-TRAY)</b>
1-Reagent	6 x 40 ml	8 x 23 ml
2-Reagent	6 x 12.5 ml	8 x 7.5 ml

The reagents when stored at 2-8°C are stable up to expiry date printed on the package.

### Concentrations in the reagent

<b>1-Reagent</b>	
ascorbate oxidase	≤ 104 µkat/l
peroxidase (POD)	≤ 22.4 µkat/l
4-aminoantipyrine	≤ 1.2 mmol/l
sodium hydroxide	≤ 0.8 %
buffer PIPES (pH 7.0)	≤ 120 mmol/l
stabilizers, preservatives, detergent	
<b>2-Reagent</b>	
buffer PIPES (pH 7.0)	≤ 60 mmol/l
ADPS	≤ 2 mmol/l
uricase	≤ 9.9 µkat/l
ferricyanide potassium	≤ 22.8 µmol/l
sodium hydroxide	≤ 0.4 %
stabilizers, preservatives, detergent	

### Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.
- 1-Reagent meeting the criteria for classification in accordance with Regulation (EC) No 1272/2008.

### Warning

H315 Causes skin irritation.

H319 Causes serious eye irritation.

P280 Wear protective gloves/protective clothing/eye protection/face protection.

P302+P352 IF ON SKIN: Wash with plenty of soap and water.

P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

### SPECIMEN

24- hours urine, serum, heparinized plasma free from hemolysis.

Do not use EDTA, fluoride and oxalate as anticoagulants!

**Urine preparation:** To prevent precipitation of salts of uric acid, 10 ml of NaOH (500 g/L) should be added to the collection bottle before collection of a 24-hour specimen. Urine should be diluted with distilled water in the ratio of 1 to 4 (multiply the result by 5).

Serum and plasma can be stored 3-5 days at 2-8°C or 6 months at -20°C. 24-hours urine samples can be stored approximately 3 days at room temperature.

Nevertheless it is recommended to perform the assay with freshly collected samples!

### PROCEDURE

1-Reagent and 2-Reagent are ready to use.

1-Reagent put on basic position in reagent tray.

2-Reagent put on start position in reagent tray.

For reagent blank deionized water is recommended.

#### Actions required:

- Biolis 24i Premium:** When performing assays in the analyser, there is a probability of **cross-contamination** affecting the tests results: CREATININE – UA, GGT – UA, HDL DIRECT – UA, LDL DIRECT – UA, URINE PROTEINS – UA. To avoid this effect follow the recommendations contained in the instruction 51\_03\_24\_008\_BIOLIS\_24i PREMIUM CARRYOVER.
- Biolis 30i:** When performing assays in the analyser, there is a probability of **cross-contamination** affecting the tests results: UA – URINE PROTEINS, CREATININE – UA, HDL DIRECT – UA, LDL DIRECT – UA. To avoid this effect follow the recommendations contained in the instruction 51\_03\_24\_009\_BIOLIS\_30i CARRYOVER.

### REFERENCE VALUES <sup>5</sup>

serum / plasma	mg/dl	µmol/l
females	2.5 – 6.8	149 – 405
males	3.6 – 7.7	214 – 458

24-hours urine	mg/24h	mmol/24h
	250 – 750	1.49 – 4.46

It is recommended for each laboratory to establish its own reference ranges for local population.

### QUALITY CONTROL

For internal quality control it is recommended to use, with each batch of samples, the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173) for determination in serum or CORMAY URINE CONTROL LEVEL 1 (Cat. No 5-161) and LEVEL 2 (Cat. No 5-162) for determination in urine.

For the calibration of automatic analysers systems the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) is recommended.

The calibration curve should be prepared, with every change of reagent lot number or as required e.g. quality control findings outside the specified range.

### PERFORMANCE CHARACTERISTICS

The following results have been obtained using the automatic analyser Biolis 30i. Results may vary if a different instrument or a manual procedure is used.

- LoB (Limit of Blank):**  
0.04 mg/dl (2.38 µmol/l)
- LoD (Limit of Detection):**  
0.09 mg/dl (5.35 µmol/l)
- LoQ (Limit of Quantitation):**  
0.3 mg/dl (17.84 µmol/l) – serum/plasma  
0.14 mg/dl (8.33 µmol/l) – urine
- Linearity:**  
up to 40 mg/dl (2379.2 µmol/l) – serum/plasma  
up to 56 mg/dl (3330.88 µmol/l) – urine

For higher concentration, in serum or plasma, dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

#### Specificity / Interferences

Haemoglobin up to 5 g/dl, ascorbate up to 30 mg/dl for determinations in serum, ascorbate up to 50 mg/dl for determinations in urine, bilirubin up to 20 mg/dl and triglycerides up to 1000 mg/dl do not interfere with the test.

#### Precision

Repeatability (run to run) n = 20	Mean	SD	CV
	[mg/dl]	[mg/dl]	[%]
level 1	4.87	0.06	1.27
level 2	9.35	0.06	0.67
Reproducibility (day to day) n = 80	Mean	SD	CV
	[mg/dl]	[mg/dl]	[%]
level 1	4.86	0.07	1.5
level 2	9.19	0.12	1.3

#### Method comparison

A comparison between uric acid values determined at **Biolis 30i** (y) and at **BECKMAN COULTER AU680** (x) using 60 serum samples gave following results:

$$y = 0.973 x + 0.386 \text{ mg/dl;}$$

$$R = 0.999 \quad (R - \text{correlation coefficient})$$

A comparison between uric acid values determined at **Biolis 30i** (y) and at **BECKMAN COULTER AU680** (x) using 34 plasma samples gave following results:

$$y = 0.9416 x + 0.4715 \text{ mg/dl;}$$

$$R = 0.995 \quad (R - \text{correlation coefficient})$$

A comparison between uric acid values determined at **Biolis 30i** (y) and at **BECKMAN COULTER AU680** (x) using 30 urine samples gave following results:

$$y = 1.0126 x - 0.1377 \text{ mg/dl;}$$

$$R = 0.997 \quad (R - \text{correlation coefficient})$$

### WASTE MANAGEMENT

Please refer to local legal requirements.

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Item name	9	UA
<b>Auto Rerun SW</b>		
ON		
<b>Auto Rerun Range (Result)</b>		
	ON	ON
	Lower	Higher
Serum	0.16	32
Urine	0.16	60
Plasma	0.16	32
CSF		
Dialysis		
Other		
<b>Auto Rerun Condition (Absorbance)</b>		
Absorbance Range		
Lower	OFF	
Higher	OFF	
Prozone Range		OFF

Item No.	9	Item Name	UA	Optical	
<b>Data information</b>					
Units	mg/dl				
Decimals	2				
<b>Calibration</b>					
Type	Linear2				
Std sample conc.					
Blank	0	#1	*	#2	*
#3		#4		#5	
#6					
<b>Analysis</b>					
Type	END method				
Main Wave Length	546nm				
Sub Wave Length	700nm				
Method	Uricase				
<b>Correlation</b>					
Y=	Slope	X+	Intercept		
	1		0		

Item No.	9	Item Name	UA	Optical
<b>Aspiration</b>				
Kind	Double			
Vol.				
	Kind	Vol.	Add	Units
Sample	4	5		µl
Reagent 1	160	10		µl
Reagent 2	40	10		µl
<b>Data Process</b>				
Read	Start	End	<b>Abs.Limit</b>	
Main	51	52	Low	-0.1
Sub	30	31	High	3
<b>Blank value</b>				
Water Blank				
<b>Reaction Monitor</b>				
0 Level Point	1			
Span	3			
<b>Third mixing</b>				
OFF				
<b>Correction value</b>				
Blank correction				
End Point Limit				2
Linear Check (%)				0
<b>Prozone Check</b>				
	Start	End	Limit (%)	
First				
Second			Low	

Item No.	9	Item Name	UA	Optical
<b>Normal Range</b>				
	Male		Female	
	Low	High	Low	High
Serum	3.6	7.7	2.5	6.8
Urine				
Plasma	3.6	7.7	2.5	6.8
CSF				
Dialysis				
Other				
<b>Panic Range</b>				
	Male		Female	
	Low	High	Low	High
Serum				
Urine				
Plasma				
CSF				
Dialysis				
Other				

Item No.	9	Item Name	UA	Optical
<b>Auto Rerun SW</b>				
ON				
<b>Auto Rerun Range (Conc.)</b>				
	First Dil	Low	High	
		Re	Value	Dil
Serum			0.15	32
Urine			0.18	56
Plasma			0.15	32
CSF				
Dialysis				
Other				
<b>Auto Rerun Condition (Absorbance)</b>				
Absorbance Range		Lower	OFF	
		Higher	OFF	
<b>Auto Rerun Condition (Prozone)</b>				
OFF				
<b>Dilution</b>				
100:Dil2				

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