PRESTIGE 24i LQ CREATININE

Nr kat. 4-233, 4-433 (PL)

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Cat. No **4-233, 4-433** (EN)

INTENDED USE

Diagnostic kit for determination of creatinine concentration used in automatic analysers Prestige 24i, Biolis 24i, Prestige 24i Premium, Biolis 24i Premium and Biolis 30i.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Creatinine is a product of creatine nonenzymatic dehydration in skeletal muscle. The amount of creatinine generated and excreted by kidney is proportional to muscle mass and usually is higher in men than women. Daily creatinine generation remains fairly constant, with the exception of crushing injury or degenerative diseases that cause massive damage to muscle. Creatinine blood and urine level depends on glomerular filtration so creatinine clearance is excellent index of renal function.

METHOD PRINCIPLE

Modified Jaffe's method, without deproteinization. In alkaline solution picrate reacts with creatinine to form a yellow-red 2,4,6-trinitrocyclohexadienate. The colour intensity is proportional to the creatinine concentration.

REAGENTS

Package

	Cat. No 4-233 (24-TRAY)	Cat. No 4-433 (36-TRAY)
1-Reagent	6 x 40 ml	8 x 23 ml
2-Reagent	6 x 12.5 ml	8 x 7.5 ml

The reagents when stored at 15-25°C are stable up to expiry date printed on the package. Stability on board of the analyser at 2-10°C: Prestige 24i – 3 weeks, Biolis 24i Premium – 2 weeks.

Concentrations in the test

sodium hydroxide	≤ 450 mmol/l
carbonate buffer	≤ 150 mmol/l
picric acid	≤ 38.8 mmol/l

Warnings and notes

- Protect from direct sunlight and avoid contamination!
- 1-Reagent and 2-Reagent meeting the criteria for classification in accordance with Regulation (EC) No 1272/2008.

Ingredients:

1-REAGENT contains sodium hydroxide.
2-REAGENT contains picric acid.

Danger

H314 Causes severe skin burns and eye damage.
P280 Wear protective gloves/protective clothing/eye protection/face protection.
P301+P330+P331 IF SWALLOWED: rinse mouth. Do NOT induce vomiting.

P303+P361+P353 IF ON SKIN (or hair): Remove/Take off immediately all contaminated clothing. Rinse skin with water/shower.

P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P310 Immediately call a POISON CENTER or doctor/physician.

SPECIMEN

Serum, EDTA or heparinized plasma free from hemolysis, 24-hours urine collected without preservatives.⁹

Urine preparation: before analysis urine sample should be diluted 100-fold with 0.9% NaCl and the results multiplied by 100. Mix well probes before measurement.

Specimen can be stored up to 7 days at 2-8°C. For longer storage samples should be frozen at -20°C.

Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-Reagent and 2-Reagent are ready to use.

1-Reagent put on basic position in reagent tray.

2-Reagent put on start position in reagent tray.

For reagent blank deionized water is recommended.

Actions required:

- Biolis 24i Premium:** When performing assays in the analyser, there is a probability of **cross-contamination** affecting the tests results: CREATININE – BIL TOTAL, CREATININE – CHOL, CREATININE – CK-MB, CREATININE - URINE PROTEINS, CREATININE – UA, CREATININE – BIL DIRECT, CK – CREATININE, CK-MB – CREATININE, URINE PROTEINS - CREATININE. To avoid this effect follow the recommendations contained in the instruction 51_03_24_008_BIOLIS_24i_PREMIUM_CARRYOVER.
- Biolis 30i:** When performing assays in the analyser, there is a probability of **cross-contamination** affecting the tests results: CREATININE – UA, CREATININE – BIL DIRECT, CK – CREATININE, CK-MB – CREATININE, GLUCOSE HEX – CREATININE, URINE PROTEINS - CREATININE. To avoid this effect follow the recommendations contained in the instruction 51_03_24_009_BIOLIS_30i_CARRYOVER.

REFERENCE VALUES ⁷

serum / plasma	mg/dl	µmol/l
female	0.6 – 1.1	53 – 97
male	0.7 – 1.3	62 – 115
24-hours urine	mg/kg/24h	µmol/kg/24h
female	11 – 20	97 – 177
male	14 – 26	124 – 230

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use, with each batch of samples, the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173) for determination in serum or CORMAY URINE CONTROL LEVEL 1 (Cat. No 5-161) or LEVEL 2 (Cat. No 5-162) for determination in urine.

For the calibration of automatic analysers systems the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) is recommended.

The calibration curve should be prepared every week (Prestige 24i, Biolis 24i Premium), with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using the automatic analyser Biolis 30i. Results may vary if a different instrument or a manual procedure is used.

- LoB (Limit of Blank):** 0.01 mg/dl (0.88 µmol/l)
- LoD (Limit of Detection):** 0.04 mg/dl (3.54 µmol/l)
- LoQ (Limit of Quantitation):** 0.6 mg/dl (53.04 µmol/l)
- Linearity:** up to 17.5 mg/dl (1547 µmol/l)

For higher concentration dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

Specificity / Interferences

Haemoglobin up to 2.5 g/dl, triglycerides up to 500 mg/dl, ascorbate up to 62 mg/l and bilirubin up to 20 mg/dl do not interfere with the test.

Precision

Repeatability (run to run) n = 20	Mean [mg/dl]	SD [mg/dl]	CV [%]
level 1	2.03	0.03	1.51
level 2	4.89	0.04	0.89
Reproducibility (day to day) n = 80	Mean [mg/dl]	SD [mg/dl]	CV [%]
level 1	2.1	0.07	3.6
level 2	5.0	0.14	2.7

Method comparison

A comparison between creatinine values determined at **Biolis 30i** (y) and at **BECKMAN COULTER AU680** (x) using 62 serum samples gave following results:

$$y = 0.948x + 0.0411 \text{ mg/dl;}$$

$$R = 0.998 \quad (R - \text{correlation coefficient})$$

A comparison between creatinine values determined at **Biolis 30i** (y) and at **BECKMAN COULTER AU680** (x) using 61 urine samples gave following results:

$$y = 0.9441x - 0.0069 \text{ mg/dl;}$$

$$R = 0.999 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT

Please refer to local legal requirements.

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Item name	8	CREA
Auto Rerun SW		
ON		
Auto Rerun Condition (Absorbance)		
Absorbance Range	Lower	OFF
	Higher	OFF
Auto Rerun Range (Result)		
	ON	ON
	Lower	Higher
Serum	0.04	20
Urine		
Plasma		
CSF		
Dialysis		
Other		
Prozone Range		
		OFF

• Prestige 24i Premium, Biolis 24i Premium

Item No.	8	Item Name	CREA	Optical
Data information				
Units	mg/dl			
Decimals	2			
Calibration				
Type	Linear2			
Std sample conc.				
Blank	0	#1	*	#2
#3		#4		#5
#6				
Analysis				
Type	RATE method			
Main Wave Length	505 nm			
Sub Wave Length	700 nm			
Method	Jaffe			
Correlation				
Y=	Slope	X+	Intercept	
	1	0		

Item No.	8	Item Name	CREA	Optical
Aspiration				
Kind	Double			
Data Process				
Read	Start	End		
Main	37	39		
Sub	35	37		
Abs.Limit Low ~ High				
	-0.1			1.350
Blank value				
Water Blank				
Correction value				
Blank correction				1
End Point Limit				1.350
Linear Check (%)				20
Reaction Monitor				
0 Level Point				1
Span				1.350
Third mixing				
OFF				
Prozone Check				
	Start	End	Limit (%)	
First				
Second				Low

Item No.	8	Item Name	CREA	Optical
Normal Range				
	Male		Female	
	Low	High	Low	High
Serum	0.70	1.30	0.60	1.10
Urine				
Plasma	0.70	1.30	0.60	1.10
CSF				
Dialysis				
Other				
Panic Range				
	Male		Female	
	Low	High	Low	High
Serum				
Urine				
Plasma				
CSF				
Dialysis				
Other				

Item No.	8	Item Name	CREA	Optical
Auto Rerun SW				
ON				
Auto Rerun Condition (Absorbance)				
	Lower	OFF		
	Higher	OFF		
Auto Rerun Range (Conc.)				
	First Dil	Low	High	
		Re	Value	Dil
Serum			0.2	25
Urine				
Plasma				
CSF				
Dialysis				
Other				
Auto Rerun Condition (Prozone)				
OFF				
Dilution				
100:Dil2				

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