## PRESTIGE 24i LQ HBDH

Nr kat. 4-241, 4-441 (PL)

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Cat. No **4-241, 4-441** (EN)

### INTENDED USE

Diagnostic kit for determination of  $\alpha$ -hydroxybutyrate dehydrogenase activity used in automatic analysers Prestige 24i, Biolis 24i and Prestige 24i Premium, Biolis 24i Premium.


The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

### INTRODUCTION

Lactate dehydrogenase (LDH, LD) is a tetrameric molecule containing two possible forms of subunits (H and M). The result is five isoenzymes, one of which is hydroxybutyrate dehydrogenase (HBDH, LD-1) formed by four H subunits. HBDH is present mainly in heart muscle but occur also in kidney and erythrocytes. Normal serum contains mostly LD-2 with lesser amount of LD-1. Changes in the ratio of LD-1 to LD-2 indicate myocardial infarction or hemolysis.

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The rate of absorbance changing at  $\lambda=340$  nm is directly proportional to  $\alpha$ -hydroxybutyrate dehydrogenase activity.

### REAGENTS

#### Package

	Cat. No 4-241 (24-TRAY)	Cat. No 4-441 (36-TRAY)
1-Reagent	4 x 40 ml	4 x 23 ml
2-Reagent	4 x 12.5 ml	4 x 7.5 ml

The reagents when stored at 2-8°C are stable up to expiry date printed on the package. Stability on board of the analyser at 2-10°C: Prestige 24i - 8 weeks, Biolis 24i Premium - 8 weeks.

#### Concentrations in the test

phosphate buffer (pH 7.5)	50 mmol/l
2-oxybutyrate	3 mmol/l
NADH	0.25 mmol/l
preservative	

#### Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.

### SPECIMEN

Serum.

Do not use hemolyzed blood because erythrocytes contain very high HBDH activity. Do not chill or freeze samples. HBDH activity is unstable and is rapidly lost during storage. Specimens can be stored up to 6 hours at 15-25°C. Nevertheless it is recommended to perform the assay with freshly collected samples!

### PROCEDURE

1-Reagent and 2-Reagent are ready to use.

1-Reagent put on basic position in reagent tray.

2-Reagent put on start position in reagent tray.

For reagent blank deionized water is recommended.

Application guide for Biolis 30i analyser is available on request.

### REFERENCE VALUES <sup>7</sup>

serum	37°C
adults	< 182 U/l (< 3.04 $\mu$ kat/l)

It is recommended for each laboratory to establish its own reference ranges for local population.

### QUALITY CONTROL

For internal quality control it is recommended to use, with each batch of samples, the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173).

For the calibration of automatic analysers systems the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) is recommended.

The calibration curve should be prepared every 8 weeks (Prestige 24i, Biolis 24i Premium), with change of reagent lot number or as required e.g. quality control findings outside the specified range.

### PERFORMANCE CHARACTERISTICS

The following results have been obtained using the automatic analysers Prestige 24i and Biolis 24i Premium. Results may vary if a different instrument or a manual procedure is used.

- Sensitivity:**  
25 U/l (5.43  $\mu$ kat/l) - Prestige 24i  
9.21 U/l (0.153  $\mu$ kat/l) - Biolis 24i Premium
- Linearity:**  
up to 500 U/l (8.33  $\mu$ kat/l) - Prestige 24i  
up to 500 U/l (8.33  $\mu$ kat/l) - Biolis 24i Premium

For higher activity dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

### Specificity / Interferences

Haemoglobin up to 2.5 g/dl, ascorbate up to 62 mg/l, bilirubin up to 20 mg/dl and triglycerides up to 1000 mg/dl do not interfere with the test.

### Precision

Repeatability (run to run)		Mean [U/l]	SD [U/l]	CV [%]
<b>Prestige 24i</b> n=20	level 1	117.19	2.45	2.09
	level 2	343.09	3.43	1.00
<b>Biolis 24i Premium</b> n=20	level 1	141.63	2.12	1.50
	level 2	362.68	4.52	1.25
Reproducibility (day to day)		Mean [U/l]	SD [U/l]	CV [%]
<b>Prestige 24i</b> n=80	level 1	131.44	3.24	2.47
	level 2	383.35	4.78	1.25
<b>Biolis 24i Premium</b> n=80	level 1	148.94	2.03	1.37
	level 2	380.33	3.28	0.86

### Method comparison

A comparison between HBDH values determined at **Prestige 24i** (y) and at **ACCENT-300** (x) using 63 samples gave following results:

$$y = 1.0514x - 10.717 \text{ U/l;}$$

$$R = 0.9613 \quad (R - \text{correlation coefficient})$$

A comparison between HBDH values determined at **Biolis 24i Premium** (y) and at **Prestige 24i** (x) using 101 samples gave following results:

$$y = 1.0432x - 7.6087 \text{ U/l;}$$

$$R = 0.9945 \quad (R - \text{correlation coefficient})$$

### WASTE MANAGEMENT

Please refer to local legal requirements.

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Item name	16	HBDH
<b>Auto Rerun SW</b>		
ON		
<b>Auto Rerun Range (Result)</b>		
	ON	ON
	Lower	Higher
Serum	25	500
Urine		
Plasma		
CSF		
Dialysis		
Other		
<b>Auto Rerun Condition (Absorbance)</b>		
Absorbance Range	Lower	OFF
	Higher	OFF
<b>Auto Rerun Condition (Prozone)</b>		
Prozone Range	OFF	

Item name	16	HBDH	
<b>Auto Rerun SW</b>			
ON			
<b>Auto Rerun Range (Conc.)</b>			
	First Dil	Low Value	High Value
		Re	Dil
Serum		9.2	500
Urine			
Plasma			
CSF			
Dialysis			
Other			
<b>Auto Rerun Condition (Absorbance)</b>			
Lower	OFF		
Higher	OFF		
<b>Auto Rerun Condition (Prozone)</b>			
OFF			
<b>Dilution</b>			
100:Dil2			

• **Prestige 24i Premium, Biolis 24i Premium**

Item No.	16	Item Name	HBDH	Optical	
<b>Data information</b>					
Units	U/L				
Decimals	1				
<b>Calibration</b>					
Type	Linear2				
Std sample conc.					
Blank	0	#1	*	#2	*
#3		#4		#5	
#6					
<b>Analysis</b>					
Type	RATE method				
Main Wave Length	340nm				
Sub Wave Length	450nm				
Method	DGKC				

<b>Correlation</b>				
Slope		Intercept		
Y=	1	X+	0	

Item No.	16	Item Name	HBDH	Optical
<b>Aspiration</b>				
Kind	Double			
Vol.				
	Kind	Vol.	Add	Units
Sample	3	5		µl
Reagent 1	200	10		µl
Reagent 2	50	10		µl
<b>Data Process</b>				
Read	Main	Start	End	
	Sub	35	52	
Abs.Limit Low High				
	1	~		2

<b>Blank value</b>				
Water Blank				
<b>Reaction Monitor</b>				
0 Level Point	1			
Span	3			
<b>Correction value</b>				
Blank correction				
End Point Limit 2				
Linear Check (%) 90				

<b>Third mixing</b>				
OFF				
<b>Prozone Check</b>				
	Start	End	Limit (%)	
First				
Second			Low	

Item No.	16	Item Name	HBDH	Optical
<b>Normal Range</b>				
	Male		Female	
	Low	High	Low	High
Serum	0	182	0	182
Urine				
Plasma				
CSF				
Dialysis				
Other				
<b>Panic Range</b>				
	Male		Female	
	Low	High	Low	High
Serum				
Urine				
Plasma				
CSF				
Dialysis				
Other				

Item No.	16	Item Name	HBDH	Optical
<b>Auto Rerun SW</b>				
ON				
<b>Auto Rerun Range (Conc.)</b>				
	First Dil	Low Value	High Value	
		Re	Dil	
Serum		9.2	500	
Urine				
Plasma				
CSF				
Dialysis				
Other				
<b>Auto Rerun Condition (Absorbance)</b>				
Lower	OFF			
Higher	OFF			
<b>Auto Rerun Condition (Prozone)</b>				
OFF				
<b>Dilution</b>				
100:Dil2				