PRESTIGE 24i URINE PROTEINS

Nr kat. 4-242, 4-442 (PL)

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Cat. No **4-242, 4-442** (EN)

INTENDED USE

Diagnostic kit for determination of total proteins in urine and cerebrospinal fluid intended to use in automatic analysers Prestige 24i, Biolis 24i, Prestige 24i Premium, Biolis 24i Premium and Biolis 30i.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

In healthy people with properly functioning kidneys, proteins are actively reabsorbed in the proximal tubules and only small amounts of proteins (several mg per day) are excreted in urine. The measurement of total proteins concentration in urine is used in the diagnosis and treatment of heart and thyroid diseases, which are characterized by proteinuria or albuminuria.

The measurement of total proteins concentration in cerebrospinal fluid (CSF) is especially useful in detecting increased permeability of the blood-brain barrier and in detecting increased intrathecal synthesis of immunoglobulins. Increased concentration of protein in CSF may indicate brain tumors, intracerebral hemorrhage, brain injury, bacterial and viral encephalitis and multiple sclerosis.

METHOD PRINCIPLE

Direct, colorimetric method with pyrogallol red.

At an acidic pH the protein amino acid groups, with the pyrogallol red-molybdate complex, form a coloured compound which shows a maximum absorbance at a wavelength of 600 nm. Colour intensity is proportional to the concentration of proteins in the sample.

REAGENTS

Package

	Cat. No 4-242 (24-TRAY)	Cat. No 4-442 (36-TRAY)
1-Reagent	2 x 44.5 ml	3 x 23.2 ml

The reagent when stored at 15-25°C is stable up to expiry date printed on the package. The reagents stored on board of the analyser at 2-10°C are stable for 12 weeks (Prestige 24i, Biolis 24i Premium)..

Concentrations in the reagent

succinate buffer	≤ 60 mmol/l
pyrogallol red	≤ 0.07 mmol/l
sodium molybdate	≤ 0.05 mmol/l
stabilizers, preservative	

Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Do not freeze the reagent.
- Reagent should be mixed before use by gentle inverting the bottle several times.
- The appearance of turbidity or control urine values outside the manufacturer's acceptable range may indicate of reagent instability.

SPECIMEN

Urine: Urine used for analysis may come from the first morning sample, random sample or timed collection sample according to the classification by ECLM².

In order to collect and prepare the samples only dedicated tubes and containers should be used.

Do not use preservatives.

Freshly collected urine should be kept at room temperature for about an hour and then cooled to 4°C. Direct reduction of the temperature on freshly collected urine can cause precipitation of minerals.

Samples with visible turbidity should be centrifuged before analysis. Determination of uncentrifuged samples may give increased results.

Urine samples are stable for 2 days at 2-8°C⁶.

Nevertheless it is recommended to perform the assay with freshly collected samples!

Cerebrospinal fluid: Cerebrospinal fluid should be centrifuged before analysis. The presence of blood in samples of cerebrospinal fluid may result in false results of protein determination. For proper interpretation of the results, cerebrospinal fluid must be determined simultaneously with a sample of blood taken from a patient at the same time.

The stability of the cerebrospinal fluid is 3 days at 2-8°C, 6 months at the temperature -20°C⁶.

Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-Reagent is ready to use.

1-Reagent put on basic position in reagent tray.

For reagent blank 0.9% NaCl is recommended.

Reagent blank is required each day.

Actions required:

- Biolis 24i Premium:** When performing assays in the analyser, there is a probability of **cross-contamination** affecting the tests results: URINE PROTEINS – CREATININE, URINE PROTEINS – CREA ENZYMATYC, URINE PROTEINS – GGT, URINE PROTEINS – UA, CREATININE - URINE PROTEINS, CREA ENZYMATYC – URINE PROTEINS. To avoid this effect follow the recommendations contained in the instruction: 51_03_24_008_BIOLIS_24i_PREMIUM_CARRYOVER.
- Biolis 30i:** When performing assays in the analyser, there is a probability of **cross-contamination** affecting the tests results: URINE PROTEINS – CREATININE, CREA ENZYMATYC – URINE PROTEINS, FERRUM – URINE PROTEINS, GGT – URINE PROTEINS, HDL DIRECT – URINE PROTEINS, LDL DIRECT – URINE PROTEINS, UA – URINE PROTEINS. To avoid this effect follow the recommendations contained in the instruction 51_03_24_009_BIOLIS_30i_CARRYOVER.

Calculation

For the calculation of proteins excreted over 24 hours, multiply the concentration (mg/dl) by the volume (dl) of the 24 hours urine.

REFERENCE VALUES^{7,10}

urine (adults)	< 15 mg/dl (0.15 g/l)
urine 24-h (adults)	< 100 mg (0.10 g)

cerebrospinal fluid	mg/dl	g/l
children 0 - 4 weeks	20 – 80	0.20 – 0.80
children >4 weeks, adults	15 – 45	0.15 – 0.45

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use with each batch of samples the CORMAY URINE CONTROL LEVEL 1 (Cat. No 5-161) and LEVEL 2 (Cat. No 5-162).

For the calibration of automatic analysers the CORMAY URINE PROTEINS CALIBRATORS (Cat. No 5-181) is recommended.

The calibration curve should be prepared every 12 weeks (Prestige 24i, Biolis 24i Premium), with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analyser Biolis 30i. Results may vary if a different instrument or a manual procedure is used.

- LoB (Limit of Blank):**
1.6 mg/dl (0.016 g/l)
- LoD (Limit of Detection):**
2.0 mg/dl (0.02 g/l)
- LoQ (Limit of Quantitation):**
10 mg/dl (0.1 g/l)
- Linearity:**
up to 137 mg/dl (1.37 g/l)

For higher concentration dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

Specificity / Interferences

Hemoglobin up to 0.004 g/dl, ascorbic acid up to 20 mg/dl, creatinine up to 6 g/l, bilirubin up to 5 mg/dl, conjugated bilirubin up to 60 mg/dl, uric acid up to 85 mg/dl, glucose up to 35 g/l, citrates up to 250 mg/dl, oxalates up to 90 mg/dl, calcium ions up to 130 mg/dl, magnesium ions up to 1.8 g/l, phosphate ions up to 1.2 g/l, urea up to 50 g/l do not affect the results of the determination.

High concentration of iron (II) ions in the test sample may cause interferences¹².

Acetaminophen and some antibiotics from penicillins and aminoglycosides can interfere^{4,11-13}.

Precision

Repeatability (run to run) n=20	Mean [mg/dl]	SD [mg/dl]	CV [%]
level 1	19.0	0.59	3.11
level 2	62.9	1.11	1.77
Reproducibility (day to day) n=80	Mean [mg/dl]	SD [mg/dl]	CV [%]
level 1	22.2	1.43	6.5
level 2	69.5	2.25	3.2

Method comparison

A comparison between total proteins values determined at **Biolis 30i** (y) and at **ADVIA SIEMENS 1800** (x) using 54 urine samples gave following results:

$y = 1.0593x - 0.1592$ mg/dl;

$R = 0.995$

(R – correlation coefficient)

A comparison between total proteins values determined at **Biolis 30i** (y) and at **BECKMAN COULTER AU680** (x) using 67 CSF samples gave following results:

$y = 0.8823x + 6.2421$ mg/dl;

$R = 0.994$

(R – correlation coefficient)

WASTE MANAGEMENT

Please refer to local legal requirements.

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Item name	54	PROT U
Auto Rerun SW		
ON		
Auto Rerun Range (Result)		
	ON	ON
	Lower	Higher
Serum		
Urine	7.3	143
Plasma		
CSF		
Dialysis		
Other		
Auto Rerun Condition (Absorbance)		
Absorbance Range	Lower	OFF
	Higher	OFF
Auto Rerun Condition (Prozone)		
Prozone Range		OFF

• **Prestige 24i Premium, Biolis 24i Premium**

Item No.	54	Item Name	PROT U	Optical
Data information				
Units	mg/dl			
Decimals	1			
Calibration				
Type	Linear1			
Std sample conc.				
Blank	0	#1	*	#2
#3		#4		#5
#6				
Analysis				
Type	END method			
Main Wave Length	600nm			
Sub Wave Length				
Method	Pyrogallol red			
Correlation				
Slope	Intercept			
Y=	1	X+	0	

Item No.	54	Item Name	PROT U	Optical
Aspiration				
Kind	Single			
Data Process				
Read	Start	End		
Main	29	30		
Sub	9	10		
Abs.Limit				
Low	-3	High	3	
Blank value				
Water Blank				
Reaction Monitor				
0 Level Point	1			
Span	3			
Third mixing				
OFF				
Correction value				
Blank correction	1			
End Point Limit	2			
Linear Check (%)	0			
Prozone Check				
First	Start	End	Limit (%)	
Second				Low

Item No.	54	Item Name	PROT U	Optical
Normal Range				
	Male	Female		
	Low	High	Low	High
Serum				
Urine	0	15	0	15
Plasma				
CSF	15	45	15	45
Dialysis				
Other				
Panic Range				
	Male	Female		
	Low	High	Low	High
Serum				
Urine				
Plasma				
CSF				
Dialysis				
Other				

Item No.	54	Item Name	PROT U	Optical
Auto Rerun SW				
ON				
Auto Rerun Range (Conc.)				
	First Dil	Low	High	
		Re	Value	Dil
Serum				
Urine		1.8	160	
Plasma				
CSF				
Dialysis				
Other				
Auto Rerun Condition (Absorbance)				
Lower	OFF			
Higher	OFF			
Auto Rerun Condition (Prozone)				
OFF				
Dilution				
99:DiI				

• **Biolis 30i**

PRESTIGE 24i URINE PROTEINS

Item no	54	Item name	UP	Specimen	CSF / URINE	OPTICAL
Data information				Aspiration volume		
UNITS	mg/dL		TYPE			
DECIMALS	1		Single			
Analysis				Data processing read		
METHOD	END method		SAMPLE	REAGENT 1	REAGENT 2	
Main Wave Length	600 nm		VOL. (µL)	6	180	
Sub Wave Length			BOTTLE (ml)			
CORRELATION (Y= AX + B)				FIRST DIL.		
A =	1		Data processing read			
B =	0		START	END		
Blank value				MAIN	27	
• WATER ° REAGENT				SUB		
Calibration				ABS LIMIT		
TYPE	Linear 1		-3	TO		3
STABILITY			Collection value			
Prozone check				END POINT	2.5	
				LINEARITY CHECK (%)	0	
				START	END	LIMIT (%)
				FIRST		
				SECOND		
				MINIMUM ABS.		
				MEAN		
				VARIATE		

Item No	54	Item Name	UP	Specimen	CSF / URINE	OPTICAL
Reference intervals				Auto rerun		
MALE		FEMALE		°ON °OFF		
LOW	HIGH	LOW	HIGH			
0	15	0	15			
Panic range				Auto rerun range (conc.)		
MALE		FEMALE		Re	Value	Dil.
LOW	HIGH	LOW	HIGH		10	
				Re	Value	Dil.
					137	
Decision limit				Auto rerun condition (abs.)		
MALE		FEMALE		DIL.		
				LOWER	°ON	°OFF
				HIGH	°ON	°OFF
				Auto rerun condition (prozone)		
				°ON °OFF		
				SAMPLE VOL.		
				DILUTION		
				°DIL 1 ° DIL 2		
Reaction check						
				°ON °OFF		
CHECK						
LOW						
HIGH						
VL CHECK				VH CHECK		
°ON °OFF		°ON °OFF				

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