PRESTIGE 24i CALCIUM ARSENAZO

Nr kat. 4-247, 4-447 (PL)

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Cat. No **4-247, 4-447** (EN)

INTENDED USE

Diagnostic kit for determination of calcium concentration used in automatic analysers Prestige 24i, Biolis 24i, Prestige 24i Premium, Biolis 24i Premium and Biolis 30i.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Calcium plays an essential role in many cell functions: intracellularly in muscle contraction and glycogen metabolism, extracellularly, in bone mineralization, in blood coagulation and in transmission of nerve impulses. Calcium is present in plasma in three forms: free, bound to proteins or complexed with anions as phosphate, citrate and bicarbonate. Decreased total calcium levels can be associated with diseases of the bone apparatus (especially osteoporosis), kidney diseases (especially under dialysis), defective intestinal absorption and hypoparathyroidism. Increased total calcium can be measured in hyperparathyroidism, malignant diseases with metastases and sarcoidosis. Calcium measurements also help in monitoring of calcium supplementation mainly in the prevention of osteoporosis.

METHOD PRINCIPLE

Photometric test using arsenazo III. Calcium with arsenazo III at neutral pH yields a blue colored complex, whose intensity is proportional to the calcium concentration. Interference by magnesium is eliminated by addition of 8-hydroxyquinoline-5-sulfonic acid.

REAGENTS

Package

	Cat. No 4-247 (24-TRAY)	Cat. No 4-447 (36-TRAY)
1-Reagent	4 x 52 ml	6 x 22 ml

The reagents when stored at 2-8°C are stable up to expiry date printed on the package. The reagent is stable for 12 weeks on board the analyser at 2-10°C (Biolis 30i).

Concentrations in the test

phosphate buffer (pH 7.5)	≤ 80 mmol/l
8-Hydroxyquinoline-5-sulfonic acid	≤ 6 mmol/l
arsenazo III	≤ 120 μmol/l
detergents, preservative	

Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Do not freeze the reagent.
- Contaminated glassware is the greatest source of error. The use of disposable plastic ware is recommended. Glassware should be soaked for a few hours in 2M HCl solution and then thoroughly rinsed with distilled water.

SPECIMEN

Serum. Random or 24-hours urine.

Serum can be stored up to 7 days at 20-25°C or up to 3 weeks at 4-8°C. Samples frozen at -20°C can be stored up to 8 months. Discard contaminated specimens.

24-hours urine preparation: To prevent calcium salt precipitation specimens should be collected in 10 ml of 6M HCl. In case of presence of precipitants they can be solved by lowering pH of the urine to below 2.0. Prior to determination dilute the sample with 0.9% NaCl in the ratio of 1 to 1. Multiply the result by the dilution factor. Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-Reagent is ready to use.

1-Reagent put on basic position in reagent tray.

For reagent blank deionized water is recommended.

Actions required:

- Biolis 24i Premium:** When performing assays in the analyser, there is a probability of **cross-contamination** affecting the tests results: CALCIUM ARSENAZO - LIPASE, RF - CALCIUM ARSENAZO. To avoid this effect follow the recommendations contained in the instruction 51_03_24_008_BIOLIS_24i_PREMIUM_CARRYOVER.
- Biolis 30i:** When performing assays in the analyser, there is a probability of **cross-contamination** affecting the tests results: CALCIUM ARSENAZO - UIBC, CALCIUM ARSENAZO - PHOSPHORUS, CALCIUM ARSENAZO - LIPASE, RF - CALCIUM ARSENAZO, TG mono - CALCIUM ARSENAZO. To avoid this effect follow the recommendations contained in the instruction: 51_03_24_009_BIOLIS_30i_CARRYOVER.

REFERENCE VALUES ⁴

serum	mg/dl	mmol/l
adult	8.6 – 10.3	2.15 – 2.57
random urine	mg/dl	mmol/l
male	0.9 – 37.9	0.225 – 9.47
female	0.5 – 35.7	0.125 – 8.92
24-hours urine	mg/24h	mmol/24h
adult	100 – 300	2.5 – 7.5

It is recommended for each laboratory to establish its own reference ranges for local population.

Calcium concentration 24-hours urine - calculation

calcium concentration in 24-hours urine [mg/24h]	=	calcium concentration in sample of 24-hours urine [mg/dl]	×	urine volume of 24-hours urine [dl/24h]
--	---	---	---	---

QUALITY CONTROL

For internal quality control it is recommended to use, with each batch of samples, the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173) for determination in serum or CORMAY URINE CONTROL LEVEL 1 (Cat. No 5-161) and LEVEL 2 (Cat. No 5-162) for determination in urine.

For the calibration of automatic analysers systems the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) are recommended.

The calibration curve should be prepared every 12 weeks (Biolis 30i), with every change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using the automatic analyser Biolis 30i. Results may vary if a different instrument or a manual procedure is used.

- LoB (Limit of Blank):**
0.07 mg/dl (0.0175 mmol/l)
- LoD (Limit of Detection):**
0.13 mg/dl (0.0325 mmol/l)
- LoQ (Limit of Quantitation):**
0.34 mg/dl (0.085 mmol/l)
- Linearity:**
up to 23 mg/dl (5.75 mmol/l)

For higher concentration, dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

Specificity / Interferences

Haemoglobin up to 5 g/dl, ascorbate up to 62 mg/l, bilirubin up to 40 mg/dl, triglycerides up to 2000 mg/dl and magnesium up to 20 mg/dl do not interfere with the test.

Precision

Repeatability (run to run) n = 20	Mean	SD	CV
	[mg/dl]	[mg/dl]	[%]
level 1	9.94	0.12	1.17
level 2	12.9	0.12	0.92
Reproducibility (day to day) n = 80	Mean	SD	CV
	[mg/dl]	[mg/dl]	[%]
level 1	10.0	0.28	2.8
level 2	12.8	0.26	2.1

Method comparison

A comparison between calcium values determined at **Biolis 30i** (y) and at **BECKMAN COULTER AU680 0** (x) using 60 serum samples gave following results:

$$y = 0.9629x + 0.6479 \text{ mg/dl};$$

$$R = 0.997 \quad (R - \text{correlation coefficient})$$

A comparison between calcium values determined at **Biolis 30i** (y) and at **BECKMAN COULTER AU680** (x) using 72 urine samples gave following results:

$$y = 0.9016x + 1.14 \text{ mg/dl};$$

R = 0.992

(R – correlation coefficient)

WASTE MANAGEMENT

Please refer to local legal requirements.

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Item name	57	CaArs
Auto Rerun SW		
ON		
Auto Rerun Range (Result)		
	ON	ON
	Lower	Higher
Serum	0.62	23
Urine		
Plasma		
CSF		
Dialysis		
Other		
Auto Rerun Condition (Absorbance)		
Absorbance Range	Lower	OFF
	Higher	OFF
Auto Rerun Condition (Prozone)		
Prozone Range	OFF	

• Prestige 24i Premium, Biolis 24i Premium

Item No.	57	Item Name	CaArs	Optical
Data information				
Units	mg/dl			
Decimals	2			
Calibration				
Type	Linear2			
Std sample conc.				
Blank	0	#1	*	#2
#3		#4		#5
#6				
Analysis				
Type	END method			
Main Wave Length	660nm			
Sub Wave Length	700nm			
Method	Arsenazo			
Correlation				
Y=	Slope	X+	Intercept	
	1		0	

Item No.	57	Item Name	CaArs	Optical
Aspiration				
Kind	Single			
Vol.				
	Kind	Vol.	Add	Units
Sample		3	5	µl
Reagent 1		300	10	µl
Data Process				
Read	Main	Start	End	
	Sub	9	10	
Abs.Limit				
Low	-3			High
				3
Correction value				
Blank correction				
End Point Limit	2			
Linear Check (%)	90			
Prozone Check				
	Start	End	Limit (%)	
First				
Second				Low
Reaction Monitor				
0 Level Point	1			
Span	3			
Third mixing				
OFF				

Item No.	57	Item Name	CaArs	Optical
Normal Range				
	Male		Female	
	Low	High	Low	High
Serum	8.6	10.3	8.6	10.3
Urine	0.9	37.9	0.5	35.7
Plasma				
CSF				
Dialysis				
Other				
Panic Range				
	Male		Female	
	Low	High	Low	High
Serum				
Urine				
Plasma				
CSF				
Dialysis				
Other				

Item No.	57	Item Name	CaArs	Optical
Auto Rerun SW				
ON				
Auto Rerun Range (Conc.)				
	First Dil	Low	High	
	Re	Value	Dil	Re
Serum		0.3		23
Urine				
Plasma				
CSF				
Dialysis				
Other				
Auto Rerun Condition (Absorbance)				
	Lower	OFF		
	Higher	OFF		
Auto Rerun Condition (Prozone)				
OFF				
Dilution				
100:Dil12				

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