## PRESTIGE 24i ACP

Nr kat. 4-249, 4-449 (PL)

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(EN)

### INTENDED USE

Diagnostic kit for determination of total, non prostatic and prostatic acid phosphatase activity, may be used in automatic analysers Prestige 24i, Biolis 24i and Prestige 24i Premium, Biolis 24i Premium.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

### INTRODUCTION

Acid phosphatases (ACP) are a group of enzymes that can be found in liver, spleen, erythrocytes, platelets, bone marrow and prostate gland. The ratio between total acid phosphatase and prostatic phosphatase is healthy males is about 1:1. Activity of total acid phosphatase increases in case of Paget's disease, hyperparathyroidism, bone cancer, Gaucher's disease, Niemann-Pick disease and myelocytic leukaemia. Prostatic acid phosphatase levels increase in case of prostatic cancer.

### METHOD PRINCIPLE

Acid phosphatases catalyses the hydrolysis, in acid medium, of 1-naphthylphosphate (1-NP) into 1-naphtol and phosphate. 1-naphtol reacts with diazo-2-chloro-5-toluene (Fast Red TR salt), forming an azo dye compound with absorbance increase is proportional to the total acid phosphatase activity (ACP). Prostatic acid phosphatase (ACP-P) is inhibited by tartrate and is detected with an indirect method by subtraction between ACP and non-prostatic acid phosphatase (ACP-NP).

### REAGENTS

#### Package

	Cat. No 4-249 (24-TRAY)	Cat. No 4-449 (36-TRAY)
REAGENT 1a	6 x 20 ml	6 x 20 ml
REAGENT 1b	6 vials (powder)	6 vials (powder)
TARTRATE REAGENT	1 vial (powder)	1 vial (powder)
ACETIC ACID SOLUTION (acetic acid 0.1 mol/l)	1 x 3 ml	1 x 3 ml
small spoon	1 pcs.	1 pcs.

The reagents when stored at 2-8°C are stable up to expiry date printed on the package.

#### Working reagents stability

Stability of working reagents ACP and ACP-NP after reconstitution: 7 days at 2-8°C, if contamination avoided and bottle recapped immediately after use.

#### Concentrations in the working reagent ACP

citrate buffer (pH 5.4)	0.1 mol/l
1-NP	10 mmol/l
Fast Red TR salt	> 1.0 mmol/l

#### Concentrations in the working reagent ACP-NP

citrate buffer (pH 5.4)	0.1 mol/l
1-NP	10 mmol/l
Fast Red TR salt	> 1.0 mmol/l
sodium tartrate	0.18 mol/l

### Warnings and notes

- Protect from direct sunlight and avoid contamination!
- The reagents must be limpid; do not use if turbid.
- It is recommended to carry out the prostatic acid phosphatase test only on samples with total acid phosphatase activity > 3.5 U/l as it is possible to obtain not reliable with a total acid phosphatase activity lower than the above-mentioned limit.
- Controls and calibrators used in determination of ACP should be used according to the guidelines contained in their instruction for use
- A slight pink coloration of working reagents ACP and ACP-NP does not affect the reagents performance.
- ACETIC ACID SOLUTION meeting the criteria for classification in accordance with Regulation (EC) No 1272/2008.

#### Ingredients:

ACETIC ACID SOLUTION contains acetic acid.

#### Danger

H314 Causes severe skin burns and eye damage.  
H318 Causes eye damage.  
P280: Wear protective gloves/protective clothing/eye protection/face protection  
P302+P352: IF ON SKIN: Wash with plenty

of water

P305+P351+P338: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P321: Specific treatment is urgently needed (go to see a doctor with the Safety data sheet for this product).

P501: Dispose of the contents/containers in accordance with the current legislation on waste treatment.

### SPECIMEN

Fresh serum. Do not use haemolized, lipemic and icteric samples. Since the sample activity decreases of about 50% within an hour at room temperature, acidify the specimen immediately after collection by adding 20µl ACETIC ACID SOLUTION for each ml of fresh serum.

Stability of acidified sample: 5 days at 2-8°C.

Nevertheless it is recommended to perform the assay with freshly collected samples!

### PROCEDURE

#### ACP ASSAY

**Working reagent ACP:** Dissolve the contents of one bottle of Reagent 1b with the contents of one bottle of Reagent 1a. Let stand for about 10 minutes and gently mix. The contents pour into bottle of Reagent 1a, put in analyser and use application for ACP.

#### ACP-NP ASSAY

**Working reagent ACP-NP:** Dissolve the contents of one bottle of Reagent 1b with the contents of one bottle of Reagent 1a. Let stand for about 10 minutes and gently mix. In order to determine ACP-NP activity add one small spoon contents (corresponding to 35-40 mg) of Tartrate Reagent for each 5 ml working reagent. Gently mix until complete solution. Such solution put in analyser and use application for ACP-NP.

### Actions required:

- Biolis 24i Premium:** When performing assays in the analyzer there is a probability of **cross-contamination** affecting the tests results. To avoid this effect, the tests for determination of acid phosphatase activity should be performed **in the separate order** (follow the recommendations contained in the instruction: 51\_03\_24\_008\_BIOLIS\_24i\_PREMIUM\_CARRYOVER).

Working reagent put on basic position in reagent tray.

For reagent blank deionized water is recommended.

Application guide for Biolis 30i analyser is available on request.

### CALCULATION OF ACP-P (prostatic acid phosphatase):

$$\text{ACP-P [U/l]} = \text{ACP} - \text{ACP-NP}$$

### REFERENCE VALUES <sup>4</sup>

ACP		
serum	U/l	µkat/l
male	2.5 – 11.7	0.042 – 0.195
female	0.3 – 9.2	0.005 – 0.154
ACP-P		
serum (tartrate inhibition)	U/l	µkat/l
male	0.2 – 3.5	0.003 – 0.058
female	0.0 – 0.8	0.000 – 0.013

It is recommended for each laboratory to establish its own reference ranges for local population.

### QUALITY CONTROL

For internal quality control of ACP it is recommended to use, with each batch of samples, the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173). For internal quality control of ACP-NP it is recommended to use, with each batch of samples, the control sera.

For the calibration of automatic analysers systems Prestige 24i, Biolis 24i the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) or LEVEL 2 (Cat. No 5-175; 5-177) are recommended.

For the calibration of automatic analysers systems Prestige 24i Premium, Biolis 24i Premium the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) are recommended.

The calibration curve should be prepared with change of reagent lot number or as required e.g. quality control findings outside the specified range.

### PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analyser Hitachi 704. Results may vary if a different instrument is used.

- Sensitivity:**  
ACP: 0.35 U/l (0.006 µkat/l).  
ACP-NP: 0.39 U/l (0.007 µkat/l).
- Linearity:**  
ACP/ACP-NP: up to 75 U/l (1.25 µkat/l).

For higher activity dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

#### Specificity / Interferences

Haemoglobin up to 0.15 g/dl, triglycerides up to 400 mg/dl, bilirubin up to 20 mg/dl and do not interfere with the test.

### Precision

ACP			
Repeatability (run to run) n = 20	Mean [U/l]	SD [U/l]	CV [%]
level 1	6.7	0.26	3.87
level 2	13.8	0.37	2.70
level 3	28.0	0.64	2.30
Reproducibility (day to day) n = 20	Mean [U/l]	SD [U/l]	CV [%]
level 1	6.64	0.27	4.06
level 2	13.84	0.37	2.64
level 3	28.10	0.62	2.22
ACP-NP			
Repeatability (run to run) n = 20	Mean [U/l]	SD [U/l]	CV [%]
level 1	2.8	0.19	6.73
level 2	3.9	0.18	4.50
level 3	6.0	0.23	3.86
Reproducibility (day to day) n = 20	Mean [U/l]	SD [U/l]	CV [%]
level 1	2.93	0.17	5.88
level 2	3.90	0.19	4.95
level 3	5.96	0.23	3.82

### Method comparison

A comparison between CORMAY reagent (y) and another commercially available assay (x) using 60 samples gave following results:

ACP

$$y = 1.0162x + 0.10151 \text{ U/l;}$$

$$R = 0.997 \quad (R - \text{correlation coefficient})$$

ACP-NP

$$y = 1.1325x - 0.09887 \text{ U/l;}$$

$$R = 0.998 \quad (R - \text{correlation coefficient})$$

### WASTE MANAGEMENT

Please refer to local legal requirements.

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Item name	55	ACP/ACP-NP	
<b>Auto Rerun SW</b>			
OFF			
<b>Auto Rerun Range (Result)</b>			
	OFF	OFF	
	Lower	Higher	
Serum			
Urine			
Plasma			
CSF			
Dialysis			
Other			
<b>Auto Rerun Condition (Absorbance)</b>			
Absorbance Range	Lower	Higher	OFF
			OFF
<b>Auto Rerun Condition (Prozone)</b>			
Prozone Range	OFF		

• **Prestige 24i Premium, Biolis 24i Premium (ACP i ACP-NP)**

Item No.	55	Item Name	ACP/ACP-NP	Optical	
<b>Data information</b>					
Units	U/l				
Decimals	1				
<b>Calibration</b>					
Type	Linear2				
Std sample conc.					
Blank	0	#1	*	#2	*
#3		#4		#5	
#6					
<b>Analysis</b>					
Type	RATE method				
Main Wave Length	405 nm				
Sub Wave Length	750 nm				
Method	Fast red				
<b>Correlation</b>					
Y=	Slope	Intercept			
	1	0			

Item No.	55	Item Name	ACP/ACP-NP	Optical
<b>Aspiration</b>				
Kind	Single			
<b>Data Process</b>				
Read	Main	Start	End	
	Sub	32	40	
<b>Vol.</b>				
	Kind	Vol.	Add	Units
	Sample	25	5	µl
	Reagent 1	250	10	µl
	Reagent 2			
<b>Blank value</b>				
Water Blank				
<b>Correction value</b>				
Blank correction	1			
End Point Limit	2			
Linear Check (%)	90			
<b>Prozone Check</b>				
	Start	End	Limit (%)	
First				
Second				Low
<b>Reaction Monitor</b>				
0 Level Point	1			
Span	3			
<b>Third mixing</b>				
OFF				

Item No.	55	Item Name	ACP/ACP-NP	Optical
<b>Normal Range</b>				
	Male		Female	
	Low	High	Low	High
Serum				
Urine				
Plasma				
CSF				
Dialysis				
Other				
<b>Panic Range</b>				
	Male		Female	
	Low	High	Low	High
Serum				
Urine				
Plasma				
CSF				
Dialysis				
Other				

Item No.	55	Item Name	ACP/ACP-NP	Optical
<b>Auto Rerun SW</b>				
OFF				
<b>Auto Rerun Range (Conc.)</b>				
	First Dil	Low	High	
		Re	Value	Dil
Serum				
Urine				
Plasma				
CSF				
Dialysis				
Other				
<b>Auto Rerun Condition (Absorbance)</b>				
	Lower	Higher	OFF	
			OFF	
<b>Auto Rerun Condition (Prozone)</b>				
OFF				
<b>Dilution</b>				
100:Di2				

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