PRESTIGE 24i LQ AMYLASE

Nr kat. **4-255, 4-455** (PL)

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Cat. No **4-255, 4-455** (EN)

INTENDED USE

Diagnostic kit for determination of α -amylase activity intended to use in automatic analyzers: Prestige 24i, Biolis 24i, Prestige 24i Premium, Biolis 24i Premium and Biolis 30i.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

α -Amylase is a digestive enzyme secreted by salivary glands and pancreas. Low level of amylase is also found in skeletal muscle, adipose tissue and fallopian tubes. α -Amylase is measured generally in pancreas diseases. Elevation of amylase activity is observed also due to inflammation of abdominal cavity or salivary glands.

METHOD PRINCIPLE

2-Chloro-4-nitrophenyl- α -maltotriose (CNP-G3) is a direct substrate for determination of α -amylase activity, which does not require the presence of ancillary enzymes.

10 CNP-G3 $\xrightarrow{\alpha\text{-amylase}}$ 9 CNP + CNP-G2 + 9 maltotriose + glucose

The rate of 2-chloro-4-nitrophenol formation can be monitored at 405 nm and is proportional to the α -amylase activity.

REAGENTS

Package

	Cat. No 4-255	Cat. No 4-455
	(24-TRAY)	(36-TRAY)
1-REAGENT	4 x 60 ml	10 x 25 ml

The reagent when stored at 2-8°C is stable up to expiry date printed on the package. Stability on board of the analyzer at 2-10°C: Prestige 24i - 8 weeks, Biolis 24i Premium - 12 weeks, Biolis 30i - 12 weeks.

Concentrations in the test

MES buffer	< 120 mmol/l
calcium acetate	< 7 mmol/l
potassium hydroxide	< 40 mmol/l
potassium thiocyanate	< 1100 mmol/l
2-chloro-4-nitrophenyl- α -maltotriose preservative, stabilizer	< 2 mmol/l

Warnings and notes

- Protect from direct sunlight!
- Prevent the reagent from microbiological contamination and from saliva and sweat α -amylase! Saliva and sweat contain α -amylase. Do not pipette by the mouth, avoid skin contact with reagent, specimens, tips, cuvettes. Ensure to use automatic pipettes and laboratory gloves.
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.
- 1-REAGENT meeting the criteria for classification in accordance with Regulation (EC) No 1272/2008.

Warning

H319 Causes serious eye irritation.
P280 Wear protective gloves, protective clothing, eye protection or face protection.
P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

SPECIMEN

Serum or plasma collected on heparin, free from hemolysis, urine.

Do not use anticoagulants: EDTA, citrates and oxalates as they inhibit amylase activity.

Serum / plasma can be stored for 7 days at 15-25°C or for one month at 2-8°C.⁷

Urine can be stored for 2 days at 15-25°C or for 10 days at 2-8°C.⁹ Amylase is very unstable in acid urine. Adjust pH to approximately 7.0 before storage.

Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-REAGENT is ready to use.

1-REAGENT put on basic position in reagent tray.

For reagent blank deionized water is recommended.

REFERENCE VALUES⁸

serum / plasma	U/l	μ kat/l
	20 – 104	0.34 – 1.77
urine	U/l	μ kat/l
	32 – 641	0.54 – 10.90

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use the following controls for each batch of samples: CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173) for determination in serum; CORMAY URINE CONTROL LEVEL 1 (Cat. No 5-161) and LEVEL 2 (Cat. No 5-162) for determination in urine.

For the calibration of automatic analyzers systems the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) is recommended.

The calibration curve should be prepared every 8 weeks (Prestige 24i, Biolis 24i Premium) or 12 weeks (Biolis 30i), with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using the automatic analyzer Biolis 30i. Results may vary if a different instrument or a manual procedure is used.

- LoB (Limit of Blank):**
1.1 U/l (0.018 μ kat/l)
- LoD (Limit of Detection):**
1.7 U/l (0.028 μ kat/l)

- LoQ (Limit of Quantitation):**
4.5 U/l (0.075 μ kat/l)

- Linearity:**
up to 2300 U/l (38.3 μ kat/l)

Specificity / Interferences

Haemoglobin up to 2.5 g/dl, ascorbate up to 62 mg/l, bilirubin up to 20 mg/dl and triglycerides up to 1000 mg/dl do not interfere with the test.

Precision

	Mean [U/l]	SD [U/l]	CV [%]
Repeatability (run to run) n = 20			
level 1	70.4	0.97	1.38
level 2	433.7	3.96	0.91
Reproducibility (day to day) n = 80			
level 1	69.2	2.88	4.2
level 2	440.0	14.70	3.3

Method comparison

A comparison between amylase values determined at **Biolis 30i** (y) and at **BS-800** (x) using 60 serum samples gave following results:

$$y = 1.0142x - 5.8464 \text{ U/l};$$

$$R = 0.999 \quad (R - \text{correlation coefficient})$$

A comparison between amylase values determined at **Biolis 30i** (y) and at **BS-800** (x) using 62 urine samples gave following results:

$$y = 1.0214x - 7.6945 \text{ U/l};$$

$$R = 1.000 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT

Please refer to local legal requirements.

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Item name	14	AMYL
Auto Rerun SW		
ON		
Auto Rerun Range (Result)		
	ON	ON
	Lower	Higher
Serum	5.3	1500
Urine		
Plasma		
CSF		
Dialysis		
Other		
Auto Rerun Condition (Absorbance)		
Absorbance Range		
	Lower	OFF
	Higher	OFF
Prozone Range		
OFF		

■ Prestige 24i Premium, Biolis 24i Premium

Item No.	14	Item Name	AMYL	Optical	
Data information					
Units	U/L				
Decimals	1				
Analysis					
Type	RATE method				
Main Wave Length	405nm				
Sub Wave Length	660nm				
Method	CNP-G3				
Calibration					
Type	Linear2				
Std sample conc.					
Blank	0	#1	*	#2	*
#3		#4		#5	
#6					
Correlation					
Y=	Slope	Intercept			
	1	X+ 0			

Item No.	14	Item Name	AMYL	Optical
Aspiration				
Kind	Single			
Vol.				
Kind	Vol.	Add	Units	
Sample	3	5	µl	
Reagent 1	200	10	µl	
Data Process				
Read	Main	Start	End	
	Sub	20	30	
Abs.Limit Low High				
	-0.1	~	2	
Correction value				
Blank correction				
End Point Limit				2.1
Linear Check (%)				10
Prozone Check				
	Start	End	Limit (%)	
First				
Second				Low

Item No.	14	Item Name	AMYL	Optical
Normal Range				
	Male		Female	
	Low	High	Low	High
Serum	20	104	20	104
Urine	32	641	32	641
Plasma	20	104	20	104
CSF				
Dialysis				
Other				
Panic Range				
	Male		Female	
	Low	High	Low	High
Serum				
Urine				
Plasma				
CSF				
Dialysis				
Other				

Item No.	14	Item Name	AMYL	Optical		
Auto Rerun SW						
ON						
Auto Rerun Range (Conc.)						
	First Dil	Low	High			
	Re	Value	Dil	Re	Value	Dil
Serum		2.5			1500	
Urine						
Plasma						
CSF						
Dialysis						
Other						
Auto Rerun Condition (Absorbance)						
ON						
	Lower	OFF				
	Higher	OFF				
Auto Rerun Condition (Prozone)						
OFF						
Dilution						
100:Di12						

PRESTIGE 24i LQ AMYLASE

Biolis 30i

Item no	14	Item name	AMYL	Specimen	SERUM/PLASMA/URINE	OPTICAL	
Data information				Aspiration volume			
UNITS		U/L		TYPE		Single	
DECIMALS		1					
Analysis							
METHOD		RATE method		SAMPLE		REAGENT 1 REAGENT 2	
Main Wave Length		405 nm		VOL. (µL)		3 200	
Sub Wave Length		660 nm		BOTTLE (ml)			
CORRELATION (Y= AX + B)				FIRST DIL.			
A =		1					
B =		0					
Blank value				Data processing read			
° WATER		° REAGENT		START		END	
				MAIN		17 27	
				SUB			
Calibration				ABS LIMIT			
TYPE		Linear 2		-0.1		TO 2	
STABILITY				Collection value			
				END POINT		2.5	
				LINEARITY CHECK (%)		90	
Prozone check							
				START		END LIMIT (%)	
FIRST							
SECOND							
						MINIMUM ABS.	
°HIGH						MEAN	
°LOW						VARIATE	

Item No	14	Item Name	AMYL	Specimen	SERUM/PLASMA/URINE	OPTICAL	
Reference intervals				Auto rerun			
MALE		FEMALE		°ON		°OFF	
LOW	HIGH	LOW	HIGH				
20	104	20	104	Auto rerun range (conc.)			
				Re	Value	Dil.	Re Value Dil.
					4.5		2300
Panic range				Auto rerun condition (abs.)			
MALE		FEMALE					DIL.
LOW	HIGH	LOW	HIGH	LOWER	°ON	°OFF	
				HIGH	°ON	°OFF	
Decision limit				Auto rerun condition (prozone)			
				°ON		°OFF	
				SAMPLE VOL.			
Reaction check				Dilution			
				°DIL 1		°DIL 2	
VL CHECK				VH CHECK			
°ON		°OFF		°ON		°OFF	

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